

# Protocol for Thawing Cryopreserved Cells

(December 2011)

## I. Introduction

TZM-bl and 293T/17 cell lines are integral reagents for many aspects of the neutralizing antibody assay. Thawing the cryopreserved cells properly is crucial to ensure the viability and functionality of the cells throughout usage in the assays. When thawing cells it is important to remember to wear appropriate personal protective equipment (PPE) to minimize the risk of injury.

## II. Definitions

**GM: Growth Media**

**DMSO: Dimethyl Sulfoxide**

**CO<sub>2</sub>: Carbon Dioxide**

## III. Reagents and Materials

Recommended vendors are listed. Unless otherwise specified, products of equal or better quality than the recommended ones can be used whenever necessary.

### **293T/17 Cells**

*Vendor:* American Tissue Culture Collection (ATCC)

### **TZM-bl Cells**

*Vendor:* NIH AIDS Research and Reference Reagent Program

**Growth Medium** (see Protocol for Reagent Preparation for Use in the Neutralizing Antibody Assay for HIV-1 in TZM-bl Cells)

### **Microliter pipettor tips, sterile**

*Vendor:* Generic

### **Disposable pipettes, sterile, individually wrapped**

*Vendor:* Falcon/VWR

1 ml pipettes

5 ml pipettes

10 ml pipettes

25 ml pipettes

50 ml pipettes

### **Cryogenic vials, 1.5 ml sterile screw cap**

*Vendor:* Sarstedt Brand Products

### **Culture flasks with vented caps, sterile**

*Vendor:* Costar/VWR

T-25 flask  
T-75 flask

**Face shield**

*Vendor:* Generic

**IV. Instrumentation**

Recommended manufacturers are listed. Unless otherwise specified, equipment of equal or better quality than the recommended ones can be used whenever necessary.

**Waterbath**

*Manufacturer:* Precision Scientific

**Biological Safety Cabinet**

*Manufacturer:* Baker Co.

**CO<sub>2</sub> Incubator**

*Manufacturer:* Forma Scientific

**Liquid Nitrogen Freezer / Dewar**

*Manufacturer:* MVE

**Light Microscope**

*Manufacturer:* Olympus

**Pipettor**

*Manufacturer:* Rainin

**V. Protocol**

**1. Thawing cells**

**NOTE 1:** Be sure to wear a full-face shield during the handling of frozen specimens.

**1.1** Transfer cryovials containing frozen cells from liquid nitrogen to a room temperature water bath in the biological safety cabinet.

**1.2** If liquid nitrogen has seeped into the cryovial, loosen the cap slightly to allow the nitrogen to escape during thawing.

**1.3** Hold the cryovial on the surface of the water bath with an occasional gentle “flick” during thawing. Do not leave the cryovial unattended during the thawing process. (It is important for cell viability that the cells are thawed and processed quickly – thawing only takes a few seconds).

**1.4** Dry off the outside of the cryovials and wipe with a 70% ethanol solution before opening the vial to prevent contamination.

**1.5** Transfer the contents of one vial of cells to a T-75 culture flask containing 30 ml of GM.

**NOTE 2:** It is important to dilute a cryoprotectant DMSO present in the cryovial at least 30-fold at this point to avoid cell toxicity.

**1.6** Incubate the cells at 37°C/5% CO<sub>2</sub> overnight keeping the flasks in a horizontal position.

**1.7** The next day, remove the medium and replace with 15 ml of fresh GM. Change the medium every 2-3 days until the cell monolayers are confluent.