

## Step 1. Start NewPyDas


Click this icon on desktop:




You will be your proposal folder.

If you do not see your proposal ID here, type `gohome()` at the prompt.

PyDas (integrated control)



DComClient ( data acquisition)



## Step 2.

Edit a batch run script (python script):

Type:

`edit -x myrun.py`

Or click 'ScriptEditor'

'myrun.py' is a file where data collection commands are stored (see below for details)

The following programs are started by **NewPyDas** as well:

- **Decomclent (Top Right):** where data collection takes place. It can be used to manually start/stop data collection
- **Chopper control**
- **Motor control**
- **Sample environment temperature control**


Most often used instrument parameters can be controlled via the 'control panel' tab:

## Step 3. Start a batch Run:

Type:


`exe myrun.py`

Or click: 'Run' and select a script to run.



# Edit EQSANS Scan Scripts Using the ScriptEditor


(1) Start the ScriptEditor by clicking the 'ScriptEditor' Button in the PyDas window:



(2) Create a script within the script editor.


(3) Save the script to hard disk

(3) Run your script



# Appendix

## A. EQ-SANS Instrument geometry




Note: Frame overlap mirror (FOM) is not installed.

# B. Helper Applications

PostIt Notepad:


PostIt pad can be used to post note for others (e.g. remote users)

It also has buttons for launching other helper apps




## Sample sheet

Fill in the sample sheet for loaded samples



## EQSANS calculator.

Simple calculator for helping experimental setup.



## C. Data Collection Script - Function Reference

---

Highlight:

**Most Frequently Used.**

Frequently Used.

Typically not directly used

### Obtaining instrument info:

1. `calcBandWidth( chopper_speed1 ,  
chopper_phases1 ,  
chopper_phases2 ,  
chopper_phases3 ,  
chopper_phases4 ):`

*Calculate the expected wavelength band give chopper speed and phases.  
Current detector location is used for calculation*

2. `info(),`

3. `getinfo():`

*print the current instrument info*

### Scan/Run:

4. `runSample( sample_id,  
run_condition=ProtonCharge(1e12),  
title=None,  
comments=None):`

*Move sample to position, start a run, stop the run when run condition is met, and save the data*

<code>sample_id</code>	:	<i>Sample cell position on sample changer</i>
<code>run_condition</code>	:	<i>finish run when condition is met</i>
example:		
	<code>ProtonCharge(1e12)</code>	<i>=&gt; run till proton charge &gt;=1e12</i>
	<code>RunTime(1000)</code>	<i>=&gt; run till time &gt;=1000 sec</i>
	<code>DectectorCounts(1e5)</code>	<i>=&gt; run till detector counts &gt;=1e5</i>
	<code>RoiCounts(1e5),</code>	<i>=&gt; run till counts &gt;=1e5</i>
	<code>MonitorCounts(1e6),</code>	<i>=&gt; run till monitor counts&gt;=1e6</i>
<code>title (optional):</code>	<code>title of the run</code>	<i>(a string within quotes)</i>
<code>comments (optional):</code>	<code>comments</code>	<i>(within quotes)</i>

#### **Example:**

`runSample(10,RunTime(3600*2),'My Sample','Run no 1')`      *#=> run sample at position 10 for 2 hours*

**Note:** If decomclient is already collecting data, `runSample()` will not move the sample changer and will wait for the current run to finish (i.e. to meet our new run condition)

5. **runCurrent**(run\_condition=ProtonCharge(1e12), title=None,comments=None):
6. **runNow**(run\_condition=ProtonCharge(1e12), title=None,comments=None):

Same as runSample(), but without moving the sample changer.  
Use these when sample changer is not used.

**Example:**

**runCurrent**(RunTime(3600\*2),'My Sample','Run no 1') #=> run for 2 hours

7. **waitRun**(run\_condition=ProtonCharge(1e12), title=None,comments=None):

Wait for the current run to finish and save it.

## Sample changer:

8. **whereIsSample():**

Print the current sample cell location

9. **setSampleLocation(location) :**

Sample to moderator distance in mm (now overwritten by value from DAS). Used for bandwidth calculation only.

10. **loadSample(wait=None):**

Move the sample changer to loading position.  
If wait == 1, then wait until the sample changer is in position before returns.

11. **selectSample(sample\_id):**

Move sample cell no: sample\_id into the beam

**Example:**

**selectSample(1)** #=> move sample changer to pos. 1

12. **moveSampleChanger(x=0,y=0),**

13. **moveSampleChangerBy(dx=0,dy=0):**

Move the sample change to a absolute coordintate

Sample changer mapping:

```
# dy per row = 24
# dx per column= 41
# row 1 y = 898
# row 15 y = 1964
# column 1 x = 116
# column 3 x = 14
#
#      pos      3          2          1
#(x,y)=      (14,1964)    (65,1964)    (116,1964)
#
#      ...
```

```
#
#   pos   45           44           43
#(x,y)=  (14,989)    (65,989)    (116,989)
```

**14. holdChanger() :**

*Hold the sample changer and prevent it from moving within python control*

*Useful when sample changer is not used for experiment.*

**15. releaseChanger() :**

*Release the sample changer and allow it to move with python*

**16. saveChangerHolderState() :**

*Save the current sample changer holding state to a file*

*Default file: C:\\eqsans\_runs\\logs\\runlog\\ sample\_changer\_holding\_status.txt*

**17. readChangerHolderState() :**

*Read in the current sample changer holding state from the default file*

**18. setTemp(new\_t):**

**19. setTemperature(new\_t):**

*Set the sample changer temperature setpoint in C*

**Example:**

**setTemp(20)**                   #=> set the sample changer temperature set point to 20C

**20. getTemp():**

**21. getTemperature():**

*Get the average sample changer temperature in C*

**Example:**

**getTemp()**                   #=> get the current average sample temperature

**22. waitTemp(timeout=None):**

**23. waitTemperature(timeout=None):**

*Wait sample to reach its temperature of a maximum of 'timeout' seconds.*

*Default: if timeout is not given, wait until temperature is in range*

**24. calcTiCellTemp(cell\_no=None,set\_temp=None) :**

**25. calcAlCellTemp(cell\_no=None,set\_temp=None) :**

**26. calcCellTemp(cell\_type,cell\_no=None,set\_temp=None):**

*Calculate the expected sample (liquid temperature) for a given cell at the set\_temp for Titanium or Aluminum cells*

Note that Al and Ti cells have different behaviors

cell\_no = 1 - 45

If cell\_no < 1 or not given, calculate for all cells

If set\_temp not given, use current set value from DAS

cell\_type = 'Al' for aluminum cells. 'Ti' for titanium cells

**Example:**

`calcCellTemp('Al',2)`

#=> calculate the estimated *sample* temperature in cell 2.

27. `setTiCellTemp(cell_no,set_temp):`

28. `setAlCellTemp(cell_no,set_temp):`

29. `setCellTemp(cell_type,cell_no,set_temp):`

Set the temperate such that cell cell\_no will reach set\_temp for a Titanium or Aluminum cell.

cell\_type = 'Al' for aluminum cells. 'Ti' for titanium cells

**Example:**

`setCellTemp('Al',3,25)`

#=> set the changer temperature such that cell 3 reaches 25C.

## Detector:

30. `whereIsDetector():`

Print detector location

31. `moveDetector( z_value_in_mm,  
y_value_in_mm=None,  
x_value_in_mm=None) :`

Move the detector to new location.

Range of motion:

z\_value\_in\_mm = 1210. to 10100.

x\_value\_in\_mm = -20. to 20.

y\_value\_in\_mm = -20. to 20.

(Hardware range: z = 1205. to 10177.,x = -25. to 25., y = -25. to 25.)

**Note:** the sequence for moving the detector and selecting a bandwidth:

1. `moveDetector(...)`
2. `setWL(...)` or `setFrame(...)`
3. `adjustBeamstop(...)` (optional)

**Example:**

`moveDetector(4000)`

#=> move the detector to 4000mm from sample.

**!!! Note that when detector is moved beyond 5000mm, smaller slits have to be selected to prevent the direct beam hitting the detector. !!!**

32. `setBeamstop(type='c',start_wavelength=None) :`



33. **selectBeamstop**(type='c',start\_wavelength=None) :  
Select a beam stop and move it to beam center according to the current detector location and selected wavelength band

Type ='c' ( default) for circular beamstop, otherwise for square.

**Example:**

`selectBeamstop('c')`

#=> select the circular beamstop, adjust its position according to current detector location and selected wavelength bands.

34. **moveBeamstopX**(x) :

35. **moveBeamstopY**(y) :

*Move the beamstop.*

Approx. range:

-100 <=x <=100 (in mm )

0 <=y <=600 (in mm )

36. **adjustBeamstop**(start\_wavelength=None) :

Adjust the beamstop position according to detector location and selected wavelength.

**Note:** the sequence for moving the detector and selecting a bandwidth:

1. `moveDetector(...)`
2. `setWL(...)` or `setFrame(...)`
3. `adjustBeamstop(...)` (optional)

**Example:**

`adjustBeamstop()`

#=> Adjust the position of the current beamstop according to current detector location and selected wavelength bands.

## Chopper and wavelength selection:

37. **setChopperSpeedsAndPhases**(speed\_for\_all\_choppers,  
chopper\_1\_phase\_in\_microsec,  
chopper\_2\_phase\_in\_microsec,  
chopper\_3\_phase\_in\_microsec,  
chopper\_4\_phase\_in\_microsec):

*Change chopper setting*

38. **setChopperByStartingWavelength**( start\_wavelenth,  
chopper\_speed\_in\_Hz,  
sample\_to\_detector\_in\_mm=None,  
sample\_to\_moderator\_in\_mm=None):

*Set chopper speeds and phases.*

*Detector to sample and sample to detector distance are obtained from DAS.*

39. `setChopperByStartingFrameNumber( start_frame,  
chopper_speed_in_Hz,  
sample_to_detector_in_mm=None,  
sample_to_moderator_in_mm=None):`

*Set chopper speeds and phases to a selected frame.  
Detector to sample and sample to detector distance are obtained from DAS.*

40. `selectFrame( frame_no,  
is_frame_skipping=None,  
sample_to_detector_in_mm=None,  
sample_to_moderator_in_mm=None) :`

41. `setFrame( frame_no,  
is_frame_skipping=None,  
sample_to_detector_in_mm=None,  
sample_to_moderator_in_mm=None) :`

*Select a frame by phasing the choppers.  
Detector to sample and sample to detector distance are obtained from DAS.*

42. `selectStartingWavelength( wavelength,  
is_frame_skipping=None,  
sample_to_detector_in_mm=None,  
sample_to_moderator_in_mm=None) :`

43. `setWL( wavelength,  
is_frame_skipping=None,  
sample_to_detector_in_mm=None,  
sample_to_moderator_in_mm=None) :`

*Set chopper speeds and phases to start at wavelength (in Angstrom).  
Detector to sample and sample to detector distance are obtained from DAS.*

**Note:** the sequence for moving the detector and selecting a bandwidth:

1. `moveDetector(...)`
2. `setWL(...)` or `setFrame(...)`
3. `adjustBeamstop(...)` (optional)

**Example:**

```
selectFrame(2)  
setWL(2.5,'skipping')
```

#=> Select Frame No. 2 at the current detector location. No frame skipping.  
#=> Select start wavelength =2.5 at the current detector location in frame skipping mode.

## Collimation Slits and Shutter

44. `setBeamSlit(slit_no,wheel_no=1):`

45. `selectSlit(slit_no,wheel_no=1):`

46. `setSlit(slit_no,wheel_no=1):`

47. `changeSlit(slit_no,wheel_no=1):`

*Move to slit to selected position on selected slit wheel*

wheel\_no = 1,2,3 (1 being the most upstream one)

Possible value for slit\_no:

For wheel no 1:

slit\_no='closed', 'd10mm', '10x10mm', 'd15mm', '15x15mm', 'd20mm', '20x20mm', 'open'

For wheel no 2:

slit\_no='open1', 'd10mm', '10x10mm', 'd15mm', '15x15mm', 'd20mm', '20x20mm', 'open'

For wheel no 3:

slit\_no='open1', 'd10mm', '10x10mm', 'd15mm', '15x15mm', 'd20mm', '20x20mm', 'open'

**Example:**

`selectSlit('d10mm')`

#=> select slit d10mm on wheel 1

**48. setSecondaryShutter(position) :**

*Set secondary shutter to position 1-8*

**49. openShutter(position=2) :**

*Open the secondary shutter, position =1,2,3,4. Default =2*

**50. closeShutter(position=2) :**

*Close the secondary shutter, position =1,2,3,4. Default =2.*

**Example:**

`openShutter()`

## Plotting:

**51. clearPlotData():**

Clear stored plot data

**52. storePlotData():**

Store data for plotting

**53. plot() :**

plot stored data (default: detector counts vs run number) using PyDas's plotting routine

**54. plotROI() :**

*plot stored ROI data using PyDas's plotting routine*

**55. plotMonitor() :**

*plot monitor data using PyDas's plotting routine*

**56. beginPlot(type=None):**

*Get ready for plotting data. Scans afterwards will be plotted*

Type = 'roi', 'ROI', 'Roi', 'r', 'R' for ROI plot

Type = 'monitor', 'm', 'M', 'Monitor', 'MONITOR' for monitor plot

Default: detector counts

57. **finishPlot():**

*Finish plotting. Scans afterwards will not be plotted*

## High Voltage routines :

58. **HVStatus():**

*The status of the detector High Voltage supply*

59. **isHVOn() :**

*True if all three H.V. are on*

60. **isHVOff() :**

*True if all three H.V. are off*

61. **setHVOn() :**

*Turn H.V. on*

62. **setHVOff() :**

*Turn H.V. off*

## Testing and others:

63. **setTesting(yes=None) :**

*Set a testing flag (1 or 0). No actual command will be sent to the instrument*

64. **checkScript() :**

*Check script for errors*

65. **setVerbose(yes=None) :**

*Verbose 1 or 0. (not yet consistently implemented)*

66. **checkHelperApps ():**

*Check motor and chopper apps are running*

67. **gohome(home=None):**

*Goto home folder. Default : folder that corresponds to current proposal ID.*