



Update on Rapid HIV Test Waivers

CLIAC Meeting
September 22, 2004

Elliot P. Cowan, Ph.D.
Associate Director

Division of Emerging and Transfusion Transmitted Diseases
FDA/CBER/OBRR



CLIA Waived Rapid HIV Tests

- ◆ OraQuick[®] ADVANCE Rapid HIV-1/HIV-2 Antibody Test (OraSure Technologies, Bethlehem, PA)
 - Detection of antibodies to HIV-1 in fingerstick whole blood: 1/31/03
 - Waiver subsequently carried over to detection of antibodies to HIV-1 and HIV-2 in both fingerstick and venipuncture whole blood
 - Detection of antibodies to HIV-1 and HIV-2 in oral fluid: 6/25/04




OraQuick[®]: Oral Fluid Test Procedure

COLLECT SAMPLE



RUN TEST



OraQuick[®]: Test Results

Non-
Reactive



Reactive



Invalid



CLIA Waived Rapid HIV Tests

- ◆ Uni-Gold™
Recombigen® HIV
(Trinity Biotech PLC,
Wicklow, Ireland)
 - Detection of antibodies
to HIV-1 in venipuncture
whole blood: 6/23/04



Uni-Gold™: Test Procedure

1. Add sample



2. Wash



3. Read in 10 min



Uni-Gold™: Test Results

Non-
Reactive



Reactive



Invalid



Rapid HIV Tests: Interpretation

- ◆ Non-reactive = Negative
- ◆ Reactive = Preliminary Positive
 - All reactive results must be confirmed using an appropriate supplemental test
 - Working with HIV supplemental test manufacturers to change labeling to allow use with rapid HIV tests
 - Screening test results are highly accurate, but reactive test result should be confirmed by supplemental testing
 - State law may preclude interpretation



Performance of Rapid HIV Tests

SENSITIVITY*

	Whole Blood	Oral Fluid
--	-------------	------------

OraQuick	99.6%	99.3%
----------	-------	-------

Uni-Gold	100%	-
----------	------	---

SPECIFICITY*

	Whole Blood	Oral Fluid
--	-------------	------------

OraQuick	100%	99.8%
----------	------	-------

Uni-Gold	99.7%	-
----------	-------	---

*Based on data from clinical trials



Rapid HIV Tests are Restricted Devices

- ◆ Sale restricted to clinical laboratories
 - that have an adequate quality assurance program and
 - where there is assurance that operators will receive and use the instructional materials
- ◆ Approved for use only by an agent of a clinical laboratory
 - I.e., not for self-testing
- ◆ Test subjects must receive the “Subject Information” pamphlet and pre-test counseling prior to specimen collection, and appropriate counseling when test results are provided



Rapid HIV Test Restrictions, cont.

- ◆ Not approved for use to screen blood or tissue donors
- ◆ Customer letter included with all kits
 - “By purchasing this device, you are doing so as an agent of a clinical laboratory and agree that you or any of your consignees will abide by the...restrictions on the sale, distribution, and use of the device...”
- ◆ Sales and use restrictions apply to waived rapid HIV tests

Current Rapid HIV Test Issues

- ◆ Reporting a positive test result on the basis of rapid HIV testing
 - Statistical validation of reactive rapid test results by a second rapid test could obviate the need for supplemental testing by Western blot, IFA, or NAT
 - Development of a policy on validation of rapid test results through the use of multiple rapid test results is a joint effort by the PHS Agencies