

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23
24
25

DRAFT BACKGROUND REVIEW DOCUMENT

CURRENT STATUS OF FIVE *IN VITRO* PYROGENICITY TEST METHODS

- **THE HUMAN WHOLE BLOOD (WB)/IL-1 *IN VITRO* PYROGEN TEST: APPLICATION OF CRYOPRESERVED HUMAN WB**
- **AN ALTERNATIVE *IN VITRO* PYROGENICITY TEST USING THE MONOCYTOID CELL LINE MONO MAC 6 (MM6)/IL-6**
- **THE HUMAN PERIPHERAL BLOOD MONONUCLEAR CELL (PBMC)/IL-6 *IN VITRO* PYROGEN TEST**
- **THE HUMAN WB/IL-1 *IN VITRO* PYROGEN TEST**
- **THE HUMAN WB/IL-6 *IN VITRO* PYROGEN TEST**

National Toxicology Program (NTP) Interagency Center for the Evaluation of Alternative Toxicological Methods (NICEATM)

December 1, 2006

25
26
27
28
29
30
31
32
33
34
35
36
37
38
39
40
41
42
43
44
45
46
47
48
49
50
51

[This Page Intentionally Left Blank]

52

TABLE OF CONTENTS

53

Page Number

54

LIST OF TABLES vii

55

LIST OF ACRONYMS AND ABBREVIATIONSix

56

ACKNOWLEDGEMENTS xiii

57

PREFACExvi

58

EXECUTIVE SUMMARYxx

59

60

1.0 Introduction and Rationale for the Proposed Use of *In Vitro*

61

Pyrogenicity Test Methods 1-1

62

1.1 Introduction 1-1

63

1.1.1 Historical Background of *In Vitro* Pyrogenicity Test Methods

64

and the Rationale for Their Development..... 1-1

65

1.1.2 Peer Reviews of *In Vitro* Pyrogenicity Test Method

66

Validation Studies 1-3

67

1.2 Regulatory Rationale and Applicability 1-4

68

1.2.1 Current Regulatory Testing Requirements and ICCVAM

69

Prioritization Criteria..... 1-4

70

1.2.2 Intended Uses of the Proposed *In Vitro* Pyrogenicity Test

71

Methods 1-6

72

1.2.3 Similarities and Differences in the Endpoints Measured

73

by the Proposed Test Methods and the *In Vivo* Reference

74

Test Method 1-7

75

1.2.4 Use of the Proposed Test Methods in an Overall Strategy

76

of Hazard or Safety Assessment..... 1-7

77

1.3 Scientific Basis for the *In Vitro* Pyrogenicity Test Methods..... 1-8

78

1.3.1 Purpose and Mechanistic Basis of the *In Vitro* Pyrogenicity

79	Test Methods.....	1-8
80	1.3.2 Similarities and Differences of Modes of Action Between the	
81	<i>In Vitro</i> Pyrogenicity Test Methods and the Fever Response	
82	in Humans and/or Rabbits	1-8
83	1.3.3 Range of Substances Amenable to the <i>In Vitro</i> Pyrogenicity Test	
84	Methods and Limits of These Methods	1-9
85	1.4 Validation of the <i>In Vitro</i> Pyrogenicity Test Methods.....	1-9
86	1.5 Search Strategies and Selection of Citations for the <i>In Vitro</i>	
87	Pyrogenicity BRD.....	1-10
88		
89	2.0 <i>In Vitro</i> Pyrogenicity Test Method Protocol Components.....	2-1
90	2.1 Overview of How the <i>In Vitro</i> Pyrogenicity Test Methods	
91	Are Conducted.....	2-1
92	2.2 Description and Rationale for the Test Method Components for	
93	Proposed Standardized Protocols	2-1
94	2.2.1 Methods Used to Analyze the Data, Including	
95	Methods to Analyze for Interference With the Assay.....	2-7
96	2.2.2 Decision Criteria and the Basis for the Prediction	
97	Model Used to Identify a Pyrogenic Substance.....	2-7
98	2.2.3 Information and Data to be Included in the Study	
99	Report and Availability of Standard Forms for	
100	Data Collection and Submission.....	2-8
101	2.3 Basis for Selection of the Test Method Systems	2-10
102	2.4 Proprietary Components.....	2-10
103	2.5 Number of Replicates	2-11
104	2.5.1 Number of Donors.....	2-11
105	2.5.2 Number of Assay Replicates.....	2-12
106	2.6 Modifications to the Test Method Protocols Based on ECVAM	
107	Validation Study Results.....	2-12
108	2.7 Differences with Comparable Validated Test Methods with	
109	Established Performance Standards.....	2-12

110	3.0	Substances Used for the Validation of <i>In Vitro</i> Pyrogenicity Test	
111		Methods.....	3-1
112	3.1	Rationale for the Substances or Products Selected for Testing.....	3-1
113	3.2	Number of Substances.....	3-1
114	3.3	Identification and Description of Substances Tested.....	3-2
115	3.4	Sample Coding Procedure.....	3-2
116	3.5	Rationale for the Selection of the Recommended Reference	
117		Substances	3-3
118			
119	4.0	<i>In Vivo</i> Reference Data for the Assessment of Test Method Accuracy.....	4-1
120	4.1	Description of the Protocol Used to Generate <i>In Vivo</i> Data.....	4-1
121		4.1.1 The Rabbit Pyrogen Test (RPT).....	4-1
122		4.1.2 Current <i>In Vivo</i> Pyrogenicity Test Method Protocols	4-4
123	4.2	Reference Data Used to Assess <i>In Vitro</i> Test Method Accuracy	4-5
124	4.3	Availability of Original Records for the <i>In Vivo</i> Reference Data.....	4-6
125	4.4	<i>In Vivo</i> Data Quality.....	4-6
126	4.5	Availability and Use of Toxicity Information from the Species	
127		of Interest.....	4-6
128	4.6	Information on the Accuracy and Reliability of the <i>In Vivo</i>	
129		Test Method.....	4-7
130			
131	5.0	Test Method Data and Results	5-1
132	5.1	Test Method Protocol.....	5-1
133	5.2	Availability of Copies of Original Data Used to Evaluate Test	
134		Method Performance.....	5-1
135	5.3	Description of the Statistical Approaches Used to Evaluate the	
136		Resulting Data	5-1
137	5.4	Summary of Results	5-3
138	5.5	Use of Coded Chemicals and Compliance with GLP Guidelines	5-3
139	5.6	Lot-to-Lot Consistency of Test Substances.....	5-3
140	5.7	Availability of Data for External Audit.....	5-4

141	6.0	Relevance of the <i>In Vitro</i> Pyrogenicity Test Methods	6-1
142	6.1	Accuracy of the <i>In Vitro</i> Pyrogenicity Test Methods	6-1
143	6.1.1	Relevance of the Cryo WB/IL-1 Test Method.....	6-2
144	6.1.2	Relevance of the MM6/IL-6 Test Method.....	6-3
145	6.1.3	Relevance of the PBMC/IL-6 Test Method.....	6-3
146	6.1.4	Relevance of the WB/IL-6 Test Method	6-4
147	6.1.5	Relevance of the WB/IL-1 Test Method	6-4
148	6.2	Summary of the Performance Statistics for <i>In Vitro</i> Pyrogenicity Test	
149		Methods.....	6-4
150	6.2.1	Discordant Results.....	6-5
151	6.2.2	Strengths and Limitations of <i>In Vitro</i> Pyrogenicity Test Methods	6-5
152			
153	7.0	Reliability of the <i>In Vitro</i> Pyrogenicity Test Methods.....	7-1
154	7.1	Selection Rationale for the Substances Used to Evaluate the	
155		Reliability of <i>In Vitro</i> Pyrogenicity Test Methods	7-1
156	7.2	Analysis of Intralaboratory Repeatability and Reproducibility	7-2
157	7.2.1	Intralaboratory Repeatability	7-2
158	7.2.2	Intralaboratory Reproducibility.....	7-3
159	7.2.3	Interlaboratory Reproducibility.....	7-6
160	7.3	Historical Positive and Negative Control Data	7-7
161			
162	8.0	Test Method Data Quality.....	8-1
163	8.1	Adherence to National and International GLP Guidelines	8-1
164	8.2	Data Quality Audits	8-1
165	8.3	Impact of Deviations from GLP Guidelines	8-1
166	8.4	Availability of Laboratory Notebooks or Other Records.....	8-2
167	8.5	Need for Data Quality.....	8-2
168			
169	9.0	Other Scientific Reports and Reviews	9-1
170	9.1	Summaries of <i>In Vitro</i> Pyrogenicity Test Methods and Data	
171		from Published and Unpublished Studies	9-1

172	9.1.1	Andrade et al. (2003).....	9-2
173	9.1.2	Bleeker et al. (1994).....	9-5
174	9.1.3	Carlin and Viitanen (2003).....	9-5
175	9.1.4	Carlin and Viitanen (2005).....	9-9
176	9.1.5	Daneshian et al. (2006).....	9-9
177	9.1.6	Eperon et al. (1996, 1997).....	9-10
178	9.1.7	Pool et al. (1998).....	9-13
179	9.1.8	Taktak et al. (1991).....	9-15
180	9.2	Conclusions from Scientific Literature Based on Independent	
181		Peer-Reviewed Reports and/or Reviews.....	9-17
182	9.2.1	De Groote et al. (1992).....	9-17
183	9.2.2	Fennrich et al. (1999).....	9-18
184	9.2.3	Hansen and Christensen (1990).....	9-19
185	9.2.4	Hartung and Wendel (1996).....	9-20
186	9.2.5	Moesby et al. (1999).....	9-20
187	9.2.6	Nakagawa et al. (2002).....	9-21
188	9.2.7	Pool et al. (1999).....	9-21
189	9.2.8	Poole et al. (2003).....	9-22
190	9.2.9	Schindler et al. (2004).....	9-22
191			
192	10.0	Animal Welfare Considerations (Refinement, Reduction,	
193		and Replacement).....	10-1
194	10.1	How the Five <i>In Vitro</i> Test Methods Will Refine, Reduce, or	
195		Replace Animal Use.....	10-1
196	10.2	Requirement for the Use of Animals.....	10-1
197	10.2.1	Rationale for the Use of Animals.....	10-1
198			
199	11.0	Practical Considerations.....	11-1
200	11.1	Transferability of the <i>In Vitro</i> Pyrogenicity Test Methods.....	11-1
201	11.1.1	Facilities and Major Fixed Equipment.....	11-1
202	11.1.2	General Availability of Other Necessary Equipment	

203		and Supplies.....	11-2
204	11.2	Personnel Training Considerations.....	11-2
205		11.2.1 Required Training and Expertise Needed to Conduct	
206		the <i>In Vitro</i> Pyrogenicity Test Methods.....	11-2
207	11.3	Cost Considerations.....	11-3
208	11.4	Time Considerations.....	11-3
209			
210	12.0	References.....	12-1
211			
212	13.0	Glossary.....	13-1
213			
214	Appendix A	ECVAM BRD Submission.....	A-1
215		The Human Whole Blood (WB)/IL-1 <i>In Vitro</i> Pyrogen Test:	
216		Application of Cryopreserved Human WB.....	TAB A1
217		An Alternative <i>In Vitro</i> Pyrogenicity Test Using the Monocytoid Cell Line	
218		Mono Mac 6 (MM6)/IL-6.....	TAB A2
219		The Human Peripheral Blood Mononuclear Cell (PBMC)/IL-6	
220		<i>In Vitro</i> Pyrogen Test.....	TAB A3
221		The Human WB/IL-1 <i>In Vitro</i> Pyrogen Test.....	TAB A4
222		The Human WB/IL-6 <i>In Vitro</i> Pyrogen Test.....	TAB A5
223			
224	Appendix B	ECVAM Response to ICCVAM Questions.....	B-1
225		ECVAM Information and Additional Unpublished Data.....	TAB B
226			
227	Appendix C	Supplemental Information.....	C-1
228		Key References.....	TAB C1
229		Guidelines for Pyrogenicity Testing.....	TAB C2

230

LIST OF TABLES

231

Page Number

232

Table ES-1 Performance Analysis for Five *In Vitro* Pyrogenicity

233

Test Methods xxiv

234

Table 1-1 Summary of U.S. and European Legislation and Statutory Protocol

235

Requirements for Pyrogenicity Testing..... 1-5

236

Table 2-1 *In Vitro* Pyrogenicity Test Method Components.....2-3

237

Table 2-2 Prediction Model Used for *In vitro* Pyrogenicity Test Methods2-8

238

Table 3-1 Test Substances (Parenteral Drugs) Used in the Validation Studies for

239

Determining Test Method Accuracy.....3-14

240

Table 3-2 Test Substances (Parenteral Drugs) Used in the Validation Studies for

241

Determining Test Method Reproducibility3-14

242

Table 4-1 Test Guidelines for the Rabbit Pyrogen Test.....4-2

243

Table 4-2 Decision Criteria for Determining a Pyrogenic Response in the

244

Rabbit Pyrogen Test.....4-5

245

Table 6-1 Performance Analyses for Five *In Vitro* Pyrogenicity Test Methods.....6-3

246

Table 7-1 Intralaboratory Repeatability Assessed with Saline Spiked with

247

WHO-LPS 94/580.....7-3

248

Table 7-2 Intralaboratory Reproducibility of *In Vitro* Pyrogenicity Test

249

Methods7-5

250

Table 7-3 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity Test

251

Methods7-6

252

Table 7-4 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity Test

253

Methods7-7

254

Table 9-1 Results of Pyrogenicity Testing of Pharmaceutical/Biological

255 Products in the Human PBMC Assay, the BET, and the RPT 9-3

256 Table 9-2 Results of Pyrogenicity Testing of Pharmaceutical/Biological

257 Products by the Human WB Culture Assay, the BET, and the RPT 9-4

258 Table 9-3 IL-6 Production from WB after Exposure to Endotoxin or

259 Five Infanrix[®] Vaccines 9-7

260 Table 9-4 IL-6 Production by MM6 Cells after Exposure to Endotoxin or Five

261 Infanrix[®] Vaccines 9-8

262 Table 9-5 Pyrogenic Activity of Blood Preparations for Parenteral Use 9-12

263 Table 9-6 Pyrogenic Activity of Vaccine Preparations 9-13

264 Table 9-7 Comparison of the WB Test, BET, and the RPT for Detecting the

265 Pyrogenicity of Production Batches of Biological Products 9-14

266 Table 9-8 Results of Pyrogenicity Testing of Batches of Therapeutic HSA

267 Using the MM6/IL-6, BET, and RPT 9-16

268 Table 9-9 Comparison of the Application Spectra of the RPT, the BET,

269 and the Human WB Assay (PyroCheck[®]) 9-19

270 LIST OF ACRONYMS AND ABBREVIATIONS		
271		
272	3 R's	Replacement, Reduction, and Refinement
273	ANOVA	Analysis of variance
274	AWIPT	Adsorb, Wash, <i>In Vitro</i> Pyrogen Test
275	BET	Bacterial Endotoxin Test
276	BP	British Pharmacopoeia
277	BRD	Background Review Document
278	cAMP	Cyclic adenosine monophosphate
279	°C	Degrees centigrade
280	CAS	Chemical Abstract Service
281	CASRN	Chemical Abstracts service Registry Number
282	CBER	Center for Biologics Evaluation and Research
283	CDER	Center for Drug Evaluation and Research
284	CDRH	Center for Devices and Radiological Health
285	CEC	Commission of European Communities
286	CFR	Code of Federal Regulations
287	cm	Centimeter
288	Cryo-	Cryopreserved
289	CSF	Colony stimulating factor
290	CV	Coefficient of variation
291	CVM	Center for Veterinary Medicine
292	EC	European Commission
293	ECVAM	European Centre for the Validation of Alternative Methods
294	EDQM	European Directorate for the Quality of Medicines
295	EEC	European Economic Community
296	ELC	Endotoxin limit concentration
297	ELISA	Enzyme-linked immunosorbent assay
298	EP	European Pharmacopoeia
299	EPA	U.S. Environmental Protection Agency
300	ESAC	ECVAM Scientific Advisory Committee

301	EU	European Union
302	EU/mL	Endotoxin units per milliliter
303	FDA	U.S. Food and Drug Administration
304	Fn	Fibronectin
305	FR	<i>Federal Register</i>
306	g	Gram
307	GLP	Good Laboratory Practice
308	GM-CSF	Granulocyte-macrophage-colony stimulating factor
309	GSK	GlaxoSmithKline
310	Hb	Hemoglobin
311	hGH	Human growth hormone
312	HSA	Human serum albumin
313	IACAC	Institutional Animal Care and Use Committee
314	ICCVAM	Interagency Coordinating Committee on the Validation of
315		Alternative Methods
316	IFN- γ	Interferon- γ
317	IgG	Immunoglobulin G
318	IL-1	Interleukin-1
319	IL-2	Interleukin-2
320	IL-6	Interleukin-6
321	i.m.	Intramuscular
322	ISO	International Standards Organization
323	i.t.	Intrathecal
324	IU	International units (endotoxin)
325	i.v.	Intravenous
326	JP	Japanese Pharmacopoeia
327	K	Threshold pyrogen dose
328	Kg	Kilogram
329	LAL	<i>Limulus</i> Amebocyte Lysate
330	λ	Sensitivity of LAL reagent in EU/mL
331	LPS	Lipopolysaccharide

332	MAPK	Mitogen associated protein kinase
333	M	Maximum human or rabbit pyrogen test dose in kg/hr
334	mg	Milligram
335	min	Minute
336	mL	Milliliter
337	mM	Millimolar
338	mm	Millimeter
339	MM6	Mono Mac 6
340	MVD	Maximum valid dilution
341	NC	Negative control
342	NF	National Formulary
343	ng	Nanogram
344	nm	nanometer
345	NICEATM	National Toxicology Program Interagency Center for the
346		Evaluation of Alternative Toxicological Methods
347	NPC	Negative product control
348	OECD	Organization for Economic Co-Operation and Development
349	OD	Optical density
350	PBMCs	Peripheral blood mononuclear cells
351	PBS	Phosphate buffered saline
352	PEI	Paul Ehrlich Institute
353	pg	Picogram
354	PG	Peptidoglycan
355	PGE ₂	Prostaglandin E ₂
356	PHA	Phytohaemagglutinin
357	P.L.	Public law
358	PMA	Phorbol myristate acetate
359	Poly (I:C)	Polyinosine:polycytidylic acid
360	PPC	Positive product control
361	PWG	Pyrogenicity Working Group
362	QA	Quality assurance

363	rec	Recombinant
364	RPT	Rabbit pyrogen test
365	s.c.	Subcutaneous
366	SD	Standard deviation
367	SEM	Standard error of the mean
368	SOT	Society of Toxicology
369	SMT	Study Management Team
370	SHS	Stabilized human serum
371	SOP	Standard operating procedure
372	TG	Test Guideline
373	THP-1	Acute monocyte leukemia cell line
374	TLR-4	Toll-like receptor-4
375	TNF- α	Tumor necrosis factor- α
376	μ g	Microgram
377	UK	United Kingdom
378	U.S.	United States
379	U.S.C.	United States Code
380	UN	United Nations
381	USP	U.S. Pharmacopeia
382	USPTO	United States Patent Office
383	WB	Whole blood
384	WHO	World Health Organization
385	xg	Times gravity
386		

386

ACKNOWLEDGEMENTS

387

388

*The following individuals are acknowledged for their contributions to the in vitro
pyrogenicity test method review process*

389

390

391

Interagency Coordinating Committee on the Validation of Alternative Methods (ICCVAM) Pyrogenicity Working Group (PWG)

392

393

394

U.S. Environmental Protection Agency (EPA)

Ayaad Assad, D.V.M., Ph.D.
Karen Hamernik, Ph.D.
Louis (Gino) Scarano

U.S. Food and Drug Administration (FDA)

Mustafa Akkoyunlu, M.D., Ph.D.
Pankaj Amin
Christine Anderson
Kimberly Benton, Ph.D.
David Hussong, Ph.D.

Abby Jacobs, Ph.D.
Christopher Joneckis, Ph.D.
Raju Kammula, D.V.M., Ph.D., D.A.B.T.
Richard McFarland, M.D. (Chair)
Ramesh Panguluri,
Penelope Rice, Ph.D.

U.S. Food and Drug Administration (FDA) (Con't)

Leonard Schechtman, Ph.D.
Amy Rosenberg, M.D.
Daniela Verthelyi, M.D., Ph.D.
Jiaqin Yao, Ph.D.

National Institute of Environmental Health Sciences (NIEHS)

William Stokes, D.V.M., D.A.C.L.A.M.
(Director, NICEATM)
Raymond Tice, Ph.D.
(Deputy Director, NICEATM)

U.S. Department of Agriculture (USDA)

Jodie Kulpa-Eddy, D.V.M.

ECVAM Liaison

Marlies Halder, Ph.D.¹

395

396

¹ Dr. Halder was not a participant in the ECVAM validation studies.

396 **National Toxicology Program (NTP) Interagency Center for the Evaluation of**
397 **Alternative Methods (NICEATM)**

398

David Allen, Ph.D.
ILS, Inc.

Deborah McCarley
NIEHS

Bradley Blackard, M.S.P.H.
ILS, Inc.

**William Stokes, D.V.M., Diplomate,
A.C.L.A.M. (Director)**
NIEHS

Sue Brenzel
ILS, Inc.

Raymond Tice, Ph.D. (Deputy Director)
NIEHS

Elizabeth Lipscomb, Ph.D.
ILS, Inc.

James Truax, M.A.
ILS, Inc.

Linda Litchfield
ILS, Inc.

Doug Winters, Ph.D.
ILS, Inc.

399

400

400

401

Additional Reviewers for the *In Vitro* Pyrogenicity Test Methods Background

402

Review Documents

403

404

Joe Haseman, Ph.D.
Consultant
Raleigh, North Carolina

Errol Zeiger, Ph.D., J.D.
Consultant
Chapel Hill, North Carolina

405

406

407

408

PREFACE

409 Endotoxin, a bacterial pyrogen also known as lipopolysaccharide, is an integral component of
410 the Gram-negative bacterial cell membrane. Endotoxin directly interacts with host
411 monocytoïd cells to induce the release of a variety of proinflammatory cytokines (e.g.,
412 interleukin [IL]-1 β , IL-6, tumor necrosis factor [TNF]- α). In addition to an initial febrile
413 reaction, excessive release of these cytokines during Gram-negative bacterial sepsis can lead
414 to multiple organ failure and death. For this reason, it is critical that parenteral
415 pharmaceuticals, fluids for injection, medical devices, and human biological products be
416 properly and accurately evaluated for the presence of endotoxin prior to their clinical or
417 veterinary use. The original pyrogen test, the rabbit pyrogen test (RPT), was developed in
418 1941 to limit to an acceptable level the risks of febrile reaction in the patient to
419 administration of, or contact with, the product of concern. While the RPT continues to serve
420 this purpose well today, an endotoxin test using an extract from the blood cells of the
421 horseshoe crab (i.e., the bacterial endotoxin test [BET]) was developed in the early 1970's as
422 an alternative to the RPT for the detection of this pyrogen. In 1980, the U.S. Food and Drug
423 Administration published guidelines for use of the BET as an end product test for human and
424 animal drug products. The U.S., European, and Japanese Pharmacopoeias currently recognize
425 both test methods for pyrogen testing (i.e., RPT and BET). The BET is recognized for its
426 sensitivity to the presence of Gram-negative endotoxins, but it has well documented
427 limitations, including its inability to respond to non-endotoxin pyrogens, as well as its
428 susceptibility to interference from certain types of materials (e.g., high protein and lipid
429 levels; glucans). In contrast, the RPT is capable of detecting both endotoxin and non-
430 endotoxin pyrogens, but requires the need for interspecies translation from rabbits to humans,
431 and is associated with animal welfare concerns that discourage animal testing.

432 More recent efforts have focused on the development of an *in vitro* test system that combines
433 the sensitivity of the BET with the wide range of pyrogens detectable by the RPT. With this
434 intention, test systems based on the activation of human monocytes *in vitro* have been
435 developed that take advantage of the role of these cells in the fever response. The European
436 Centre for the Validation of Alternative Methods (ECVAM), a unit of the Institute for Health
437 and Consumer Protection (IHCP) at the European Commission's Joint Research Centre,

438 conducted a validation study to independently evaluate the usefulness of six *in vitro*
439 pyrogenicity test methods. The study was financed by the European Commission within the
440 5th Framework Programme of DG Research and recently published (Hoffmann et al, 2005).
441 Since two tests based on THP-1 cell lines did not meet the validation criteria, they are not
442 included in the peer review. During Summer 2004, the University of Konstanz (Germany)
443 carried out catch-up validation studies of two tests using cryo-preserved blood (WB/IL-1) or
444 blood cells (PBMC-IL6), the results of which were recently published (Schindler et al. 2006).
445 Based on these studies, in June 2005, ECVAM submitted background review documents
446 (BRDs) for five methods to the National Toxicology Program (NTP) Interagency Center for
447 the Evaluation of Alternative Toxicological Methods (NICEATM) for consideration as
448 replacements for the RPT. The proposed test methods were:

- 449 • The Human Whole Blood (WB)/IL-1 *In Vitro* Pyrogenicity Test: Application
450 of Cryopreserved Human WB
- 451 • An Alternative *In Vitro* Pyrogenicity Test Using the Monocytoid Cell Line
452 Mono Mac 6 (MM6)/IL-6
- 453 • The Human PBMC/IL-6 *In Vitro* Pyrogen Test
- 454 • The Human WB/IL-1 *In Vitro* Pyrogen Test
- 455 • The Human WB/IL-6 *In Vitro* Pyrogen Test.

456 For simplicity, the submitted studies are referred to collectively as the ECVAM validation
457 study in this document.

458 ICCVAM, which is charged with coordinating the technical evaluations of new, revised, and
459 alternative test methods with regulatory applicability (ICCVAM Authorization Act of 2000,
460 Public Law [P.L.] 106-545), unanimously agreed that the five submitted *in vitro* test methods
461 should have a high priority for evaluation. An ICCVAM Pyrogenicity Working Group
462 (PWG) was established to work with NICEATM to carry out these evaluations. The PWG
463 consists of knowledgeable scientists from ICCVAM member agencies. The PWG functions
464 include the review of draft test method BRDs, recommending proposed performance
465 standards, identifying and recommending scientists for expert scientific review panels,
466 preparing questions for expert or peer review panels, developing draft ICCVAM test method

467 recommendations regarding the usefulness and applicability of the alternative test methods
468 for regulatory testing, and recommending necessary validation studies. ICCVAM and
469 NICEATM also collaborate closely with ECVAM. Accordingly, an ECVAM liaison was
470 designated for the ICCVAM PWG to ensure input and contributions during the evaluation
471 and review process.

472 NICEATM, which administers the ICCVAM and provides scientific support for ICCVAM
473 activities, subsequently prepared a comprehensive background review document (BRD) that
474 provided information and data from the validation studies for each of the five *in vitro* test
475 methods. A request for any other data and information on these test methods was made
476 through a 2005 *Federal Register* (FR) request (Available:
477 <http://iccvam.niehs.nih.gov/methods/pyrogen.htm>), through the ICCVAM electronic mailing
478 list, and through direct requests to over 100 interested stakeholders. No additional data or
479 information was submitted in response to these requests.

480 This BRD is publicly available on the ICCVAM/NICEATM website ([http://](http://iccvam.niehs.nih.gov)
481 iccvam.niehs.nih.gov) or from NICEATM on request. Comments from the public and scientific
482 community are welcome and will be provided to the panel experts and made available on the
483 ICCVAM/NICEATM website (see FR notice TBD, published TBD, available at
484 <http://iccvam.niehs.nih.gov>). The independent review of the usefulness and limitations of the
485 five test methods will take place in a public meeting of the international expert scientific
486 panel on February 6, 2007 at the National Institutes of Health campus in Bethesda, Maryland.
487 The review panel's independent report will be published and made available for public
488 comment following the panel's February meeting. The ICCVAM and the PWG will consider
489 the report and public comments, and prepare final test method recommendations that will be
490 provided to federal agencies and made available to the public. ICCVAM final test method
491 recommendations are forwarded to U.S. Federal agencies for consideration, in accordance
492 with the ICCVAM Authorization Act of 2000 (P.L. 106-545).

493 We want to acknowledge the excellent cooperation and contributions from the many
494 organizations and scientists who provided critical data and information necessary for the
495 original ECVAM BRD and for the subsequent ICCVAM review. The efforts of many
496 individuals who contributed to the preparation of the ICCVAM BRD are also gratefully

497 acknowledged. These include David Allen, Ph.D., Bradley Blackard, M.S.P.H., and James
498 Truax, M.A. of Integrated Laboratory Systems (ILS, Inc.), the NICEATM Support
499 Contractor, as well as the members of the ICCVAM PWG and ICCVAM representatives who
500 subsequently reviewed and provided comments throughout the process leading to this final
501 draft version. We also want to thank Raymond Tice, Ph.D., Deputy Director of NICEATM,
502 for his extensive efforts on this project. Finally, we want to recognize the excellent leadership
503 of the PWG Chair, Dr. Richard McFarland (FDA).

504

505

506 William S. Stokes, D.V.M., Diplomate A.C.L.A.M.

507 Director, NICEATM

508 Executive Director, ICCVAM

509

510 Leonard Schechtman, Ph.D.

511 U.S. Food and Drug Administration

512 National Center for Toxicological Research

513 Chairman, ICCVAM

514

515 *September 15, 2006*

516

517

EXECUTIVE SUMMARY

518 This draft Background Review Document (BRD) provides a comprehensive review of
519 available data and information regarding the usefulness and limitations of five alternative *in*
520 *vitro* pyrogenicity test methods. The test methods are:

- 521 • The Human Whole Blood (WB)/IL-1 *In Vitro* Pyrogen Test: Application of
522 Cryopreserved Human WB
- 523 • An Alternative *In Vitro* Pyrogenicity Test Using the Monocytoid Cell Line
524 Mono Mac 6 (MM6)/IL-6
- 525 • The Human Peripheral Blood Mononuclear Cell (PBMC)/IL-6 *In Vitro*
526 Pyrogen Test
- 527 • The Human WB/IL-1 *In Vitro* Pyrogen Test
- 528 • The Human WB/IL-6 *In Vitro* Pyrogen Test

529 The test methods were reviewed for their ability to detect the presence of Gram-negative
530 endotoxin when spiked into a variety of parenteral pharmaceuticals. The objective of this
531 BRD is to describe the current validation status of the *in vitro* pyrogenicity test methods,
532 including what is known about their relevance² and reliability³, the scope of the substances
533 tested, and the availability of a standardized test method protocol for each test method.

534 The information summarized in this BRD is based on data contained in five individual BRDs
535 submitted by the European Centre for the Validation of Alternative Methods (see **Appendix**
536 **A**) to the National Toxicology Program (NTP) Interagency Center for the Evaluation of
537 Alternative Toxicological Methods (NICEATM). The ECVAM BRDs were prepared
538 according to the Interagency Coordinating Committee on the Validation of Alternative
539 Methods (ICCVAM) Guidelines for the Nomination and Submission of New, Revised, and
540 Alternative Test Methods (ICCVAM 2003) to allow for an independent scientific peer review

² The extent to which a test method correctly predicts or measures the biological effect of interest in humans or another species of interest. Relevance incorporates consideration of the "accuracy" or "concordance" of a test method.

³ A measure of the degree to which a test method can be performed reproducibly within and among laboratories over time. It is assessed by calculating intra- and interlaboratory reproducibility and intralaboratory repeatability.

541 panel to assess the validation status of the proposed test methods, and for agencies to assess
542 the acceptability of the proposed test methods for providing useful information for hazard or
543 risk assessment. Each ECVAM BRD summarizes the validation studies conducted with each
544 *in vitro* pyrogenicity test method, which is in turn compiled into a single document in the
545 current BRD, which compares and contrasts the five test methods.

546 This BRD also summarizes available information obtained from published studies as well as
547 additional unpublished data provided by ECVAM (see **Section 9.0**, Other Scientific Reports
548 and Reviews). This section discusses *in vitro* pyrogenicity test method studies that could not
549 be included in the performance analyses because of the lack of appropriate study details or
550 test method results and/or the lack of appropriate *in vivo* rabbit pyrogen test (RPT) reference
551 data. An online literature search for additional data and information on the proposed *in vitro*
552 pyrogenicity test methods identified eighteen published studies that contained relevant data.
553 ECVAM also provided additional unpublished data in response to a request for additional
554 information related to the validation studies.

555 The proposed *in vitro* pyrogenicity test methods are based on the measurement of
556 proinflammatory cytokines (i.e., interleukin-1 [IL-1] or interleukin-6 [IL-6]), released from
557 monocytoic cells contained in WB, from isolated PBMC, or from a monocytoic cell line in
558 response to exposure to a pyrogen, namely, Gram-negative endotoxin. No data were provided
559 from the validation studies supporting the usefulness of these test method for non-endotoxin-
560 based pyrogens.

561 U.S. regulatory agencies were surveyed to determine whether any of the proposed test
562 methods have been considered for regulatory use where submission of test data is required.
563 General regulatory practice (e.g., in the US and EU) is to accept pyrogenicity test method
564 data when product specific validation of the test method has been demonstrated. It was noted
565 in the ECVAM BRD that the U.S. FDA has accepted data from the PBMC test developed by
566 Novartis and Baxter Healthcare which, in conjunction with RPT and BET data, were used to
567 support the safety testing of a single specific drug product (New Drug Application Number
568 16-267/S-037).

569 Although there are differences among the *in vitro* pyrogenicity test methods based
570 predominantly on the cell type used, there are some basic steps that are consistent across all
571 methods as follows:

- 572 • The test substance is applied to the specific human-derived cells used in the *in*
573 *vitro* test method (i.e., mixed with a suspension of cells).
- 574 • The test substance is incubated with the cells for a specified period of time
- 575 • The concentration of pro-inflammatory cytokines (e.g., IL-1 β , IL-6) is
576 quantified via a cytokine-specific enzyme-linked immunosorbent assay
577 (ELISA) by comparison to a standard curve
- 578 • Using an endotoxin standard curve, the endotoxin content of the product is
579 calculated
- 580 • A product “passes” (i.e., is considered negative for endotoxin) if the endotoxin
581 content is < 0.5 endotoxin units (EU)/mL.

582 A total of ten parenteral pharmaceuticals were used for the determination of the sensitivity
583 and specificity of the five *in vitro* pyrogenicity test methods. The ability of the *in vitro*
584 pyrogenicity test methods to correctly identify the presence of Gram-negative endotoxin was
585 evaluated using parenteral pharmaceuticals spiked endotoxin (WHO-LPS 94/580 [*E. coli*
586 O113:H10:K-]). Each substance was spiked with five concentrations of endotoxin tested once
587 in three different laboratories. As outlined in **Table ES-1**, this analysis indicated that
588 concordance among the test methods ranged from 81% to 93%, sensitivity ranged from 89%
589 to 97%, specificity ranged from 81% to 97%, false negative rates ranged from 3% to 27%,
590 and false positive rates ranged from 3% to 19%.

591

591 **Table ES-1 Performance Analysis for Five *In Vitro* Pyrogenicity Test Methods¹**

Test Method	Concordance ²	Sensitivity	Specificity	False Negative Rate	False Positive Rate
cryo WB/IL-1	91.7% (110/120)	97.4% (75/77)	81.4% (35/43)	2.6% (2/77)	18.6% (8/43)
MM6/IL-6	93.2% (138/148)	95.5% (85/89)	89.8% (53/59)	4.5% (4/89)	10.2% (6/59)
PBMC/IL-6	93.3% (140/150)	92.2% (83/90)	95.0% (57/60)	7.8% (7/90)	5.0% (3/60)
PBMC/IL-6 (cryo) ³	91.9% (136/148)	88.8% (79/89)	96.6% (57/59)	11.2% (10/89)	3.4% (2/59)
WB/IL-6	91.9% (136/148)	88.8% (79/89)	96.6% (57/59)	11.2% (10/89)	3.4% (2/59)
WB/IL-1	81.0% (119/147)	72.7% (64/88)	93.2% (55/59)	27.3% (24/88)	6.8% (4/59)
WB/IL-1 (96-well plate method) ⁴	92.8% (129/139)	98.8% (83/84)	83.6% (46/55)	1.2% (1/84)	16.4% (9/55)

592 Abbreviations: cryo = Cryopreserved; IL-1 = Interleukin-1; IL-6 = Interleukin -6; MM6 = Mono Mac 6; PBMC = Peripheral
593 blood mononuclear cells; WB = Whole blood

594 ¹Based on results of 10 parenteral drugs tested in each of three different laboratories; samples of each drug were tested with
595 or without being spiked with a Gram-negative endotoxin standard (0, 0.25, 0.5, 0.5, or 1.0 EU/mL).

596 ²Percentage (Number of correct runs/total number of runs)

597 ³A modification of the PBMC/IL-6 test method using cryopreserved PBMCs.

598 ⁴A modification of the WB/IL-1 test method using 96-well plates instead of tubes for the test substance incubation.
599

600 It was not possible to make a direct comparison between the RPT and *in vitro* pyrogenicity
601 test results without the availability of parallel testing data (i.e., same test substance tested
602 using the *in vitro* and *in vivo* methods). Therefore, *in vitro* results that are discordant from the
603 RPT could not be identified with these studies. Discordant results reflect a failure of the *in*
604 *vitro* test method(s) to identify Gram-negative endotoxin spiked into a test substance at the
605 threshold concentration (0.5 EU/mL) established based on historical data from the RPT.

606 The limitations of these test methods have not been fully explored and identified. For this
607 reason, pre-testing product specific validation will be necessary to establish if a particular test
608 substance/material is appropriate for evaluation using these *in vitro* test methods. One
609 identified limitation of the *in vitro* methods is the lack of data to determine their responses to,
610 and suitability for, non-endotoxin pyrogens that are known to be detected by the RPT.
611 However, an advantage to these *in vitro* test methods is that they are derived from human
612 tissues, and thus avoid potential uncertainty associated with cross-species extrapolation.

613 Intralaboratory repeatability of each test method was evaluated by testing saline and various
614 endotoxin spikes (0.06 to 0.5 EU/mL) in saline and evaluating the closeness of agreement
615 among optical density readings for cytokine measurements at each concentration. Up to 20
616 replicates per concentration were tested and results indicated that variability in OD
617 measurements increased with increasing endotoxin concentration, but the variability was not
618 so great to interfere with distinguishing the 0.5 EU/mL spike concentration (i.e., the
619 threshold for pyrogenicity) from the lower concentrations.

620 Intralaboratory reproducibility was evaluated using three marketed pharmaceuticals spiked
621 with various concentrations of endotoxin. Three identical, independent runs conducted in
622 each of the three testing laboratories, with the exception of the cryo WB/IL-1 test method⁴.
623 All three possible combinations were compared (i.e., run 1 vs. run 2; run 1 vs. run 3; run 2
624 vs. run 3) and a mean value calculated, intended to provide an overall proportion of inter-run
625 agreement. Inter-run agreement between two runs within each laboratory ranged from 75% to
626 100%, with mean values ranging from 83% to 100%. Agreement across 3 runs within a
627 single laboratory ranged from 75% to 100%. Interlaboratory reproducibility was evaluated in
628 two different studies in which each run from one laboratory was compared to all other runs of
629 another laboratory. The proportion of equally qualified samples provided a measure of
630 reproducibility. In the first study, three marketed pharmaceutical products were spiked with a
631 saline control or various concentrations of endotoxin and tested three times in three different
632 laboratories, except for cryoWB/IL-1, which was tested once in each laboratory in the catch-
633 up validation study. The agreement across the three laboratories for each test method ranged
634 from 58% to 86%⁵, depending on the test method used and 92% for the cryoWB/IL-1 test
635 method.

636 In the second study, reproducibility was determined using the results from the ten substances
637 used in the accuracy analysis. Each substance was spiked with five concentrations of
638 endotoxin and tested once in each of three laboratories. The agreement across three
639 laboratories for each test method ranged from 57% to 88%, depending on the test method

⁴ The ECVAM cryo WB/IL-1 test method BRD states that there was no direct assessment of intralaboratory reproducibility because such an evaluation was performed in the WB IL-1 (fresh blood) test method, and the authors assume that variability is not affected by the change to cryopreserved blood.

640 used. The levels, and order of agreement among laboratories was the same for both studies;
641 the WB/IL-1 test method showed the least agreement (57-58%) and the cryo WB/IL-1 test
642 method showed the most (88-92%).

643 As stated above, this BRD provides a comprehensive summary of the current validation
644 status of five *in vitro* pyrogenicity test methods, including what is known about their
645 relevance and reliability, the scope of the substances tested, and the availability of a
646 standardized test method protocol for each test method. Data for these IVPT methods will be
647 maintained for future use, so that these performance statistics may be updated as additional
648 information becomes available.

649

⁵ However a modification of the WB/IL-1 test method (using 96-well plates for the test substance incubation) tested once in each laboratory resulted in agreement among laboratories of 83% to 92%.

650 **1.0 INTRODUCTION AND RATIONALE FOR THE PROPOSED USE OF *IN***
651 ***VITRO* PYROGENICITY TEST METHODS**

652 **1.1 Introduction**

653 1.1.1 Historical Background of *In Vitro* Pyrogenicity Test Methods and the Rationale for
654 Their Development

655 A brief summary of the historical development the five *in vitro* pyrogenicity test methods
656 was provided in Section 1.1.1 of each European Centre for the Validation of Alternative
657 Methods (ECVAM) Background Review Document (BRD) provided to the Interagency
658 Coordinating Committee on the Validation of Alternative Methods (ICCVAM), and is
659 included in **Appendix A**⁶. This section includes supplementary information and provides a
660 context for U.S. regulatory considerations.

661 Pyrogenic substances increase body temperature by inducing leukocytes (i.e., neutrophils,
662 monocytes/macrophages, lymphocytes) to release proinflammatory cytokines (e.g.,
663 interleukin [IL]-1, IL-6, and tumor necrosis factor- α [TNF- α]) that act as endogenous
664 pyrogens (Dinarello et al. 1999). These pyrogenic substances may originate from a variety of
665 biological or synthetic/manufacturing sources. They may also be released from
666 microbiological organisms such as bacteria, viruses, and fungi during cell death or following
667 immunological attack (i.e., cell damage or death due to a local or systemic immune
668 response). One of the most potent pyrogenic materials is bacterial endotoxin, which is a
669 component of the outer cell wall of Gram-negative bacteria. Pyrogens may also be found in
670 processing and packaging materials, chemicals, raw materials, or equipment used during
671 manufacturing of parenteral drugs or medical devices. The presence of endotoxins in
672 otherwise sterile biological preparations such as parenteral drugs suggests the presence of
673 current or past bacterial contamination.

674 The U.S., European, and Japanese Pharmacopoeias currently recognize two test methods for
675 pyrogen testing, the *in vivo* rabbit pyrogen test (RPT) and the *in vitro* bacterial endotoxin test
676 (BET), commonly referred to as the *Limulus* amoebocyte lysate (LAL) test. The BET is

⁶ References to ECVAM BRD sections are in normal type to distinguish them from references to ICCVAM BRD sections which are in **boldface** type.

677 recognized for its sensitivity to the presence of Gram-negative endotoxins, but it has well
678 documented limitations, including its inability to respond to non-endotoxin pyrogens, as well
679 as its susceptibility to interference from certain types of materials (e.g., high protein and lipid
680 levels; glucans). In contrast, the RPT is capable of detecting both endotoxin and non-
681 endotoxin pyrogens. However, disadvantages of the RPT include the need for interspecies
682 translation from rabbits to humans, along with animal welfare concerns that discourage
683 animal testing.

684 An *in vitro* test system that combines the sensitivity of the BET with the wide range of
685 pyrogens detectable by the RPT would be an obvious improvement for pyrogen testing. With
686 this intention, test systems based on the activation of human monocytes *in vitro* were
687 developed that take advantage of an increased understanding of the biological mechanisms
688 responsible for the human fever reaction (Dinarelo, 1999). Initial efforts focused on
689 peripheral blood mononuclear cells (PBMC), which release proinflammatory cytokines when
690 exposed to endotoxin (Duff and Atkins, 1982; Dinarelo 1984). A number of similar test
691 systems, using either whole blood, PBMCs, or monocytoid cell lines (e.g., MONO MAC 6
692 [MM6], THP-1) were subsequently developed (Poole et al., 1988; Ziegler-Heitbrock et al,
693 1988; Tsuchiya et al, 1980; Hartung and Wendel, 1996; Hartung et al, 2001; Poole et al,
694 2003). Five test systems developed from human monocytoid cells were selected by ECVAM
695 for prevalidation and validation studies with the intent of comparing their effectiveness for
696 replacing the RPT and thereby eliminating the need for using rabbits for pyrogen testing. The
697 results of these studies have been published (Hoffmann et al. 2005). The five tests selected
698 were:

- 699 • The Human Whole Blood/ IL-1 *In Vitro* Pyrogen Test Using Cryopreserved
700 Human Whole Blood (cryo WB/IL-1)
- 701 • An Alternative *In Vitro* Pyrogen Test Using the Human Monocytoid Cell Line
702 MONO MAC6 (MM6/IL-6)
- 703 • *In Vitro* Pyrogen Test Using Human PBMCs (PBMC/IL-6)⁷

⁷ As indicated in the ECVAM BRDs for the WB/IL-1 and PBMC/IL-6 test methods, catch-up validation studies were conducted to evaluate the performance of the WB/IL-1 test method when using 96-well plates, and the PBMC/IL-6 test method when using cryopreserved PBMCs. The plating procedure (WB/IL-1) and the cryopreservation procedure (PBMC/IL-6) are the only differences in the test method protocols (see **Appendix A**). These modifications were not

- 704 • The Human Whole Blood/IL-1 *In Vitro* Pyrogen Test (WB/IL-1)⁵
705 • The Human Whole Blood/IL-6 *In Vitro* Pyrogen Test (WB/IL-6)

706 1.1.2 Peer Reviews of *In Vitro* Pyrogenicity Test Method Validation Studies

707 The ECVAM-sponsored validation studies of each of these *in vitro* methods have been the
708 subject of a recent formal peer review convened by the ECVAM Scientific Advisory
709 Committee (ESAC). Two members of the ESAC served as co-Chairpersons for the review
710 panel, which consisted of five additional U.S. and European reviewers. These reviewers
711 assessed the ability of each test method to serve as a complete replacement for the *in vivo*
712 rabbit pyrogen test (RPT). Subsequent to this review, the ESAC declared that, “these tests
713 have been scientifically validated for the detection of pyrogenicity mediated by Gram-
714 negative endotoxins, and quantification of this pyrogen, in materials currently evaluated and
715 characterized by rabbit pyrogen tests.”

716 However, in their comments to the Study Management Team (SMT), two of the reviewers
717 suggested that, because an adequate link between the *in vitro* pyrogenicity test methods and
718 the RPT was not demonstrated, the validation studies fell short of providing a basis for
719 recommendation of any of these methods as substitutes or replacements for the methods
720 presently indicated in the U.S., European, and Japanese Pharmacopeias (i.e., the RPT and
721 BET).

722 This BRD was prepared for use by an ICCVAM Expert Panel following the review of these
723 *in vitro* pyrogenicity test methods. Because individual BRDs for each method were provided
724 by ECVAM, this ICCVAM BRD serves to combine common information, and references the
725 appropriate sections of the ECVAM BRDs for specifics related to the individual test
726 methods. The results and recommendations of the Expert Panel report, combined with the
727 analyses presented in the ECVAM BRDs, will be used to support ICCVAM
728 recommendations on the usefulness and limitations of each test method, the proposed
729 standardized test method protocols, performance standards, and any additional studies
730 considered necessary to further develop or characterize any or all of these test methods.

submitted by ECVAM as separate test methods, and are therefore not being considered as separate test methods in this BRD. However, where relevant, comparative information is provided (e.g., see **Table 2-1** and **Sections 6.1** and **7.2**).

731 **1.2 Regulatory Rationale and Applicability**

732 1.2.1 Current Regulatory Testing Requirements and ICCVAM Prioritization Criteria

733 This section reviews and summarizes the extent to which the five ICCVAM prioritization
734 criteria (ICCVAM 2003) apply to the *in vitro* pyrogenicity test methods under consideration.

735 **Criteria 1. The extents to which the proposed test methods are (a) applicable to**
736 **regulatory testing needs and (b) applicable to multiple agencies/programs.**

737 Pyrogenicity testing is used by regulatory authorities primarily for end product release of
738 human and animal parenteral drugs, biological products, and medical devices. The results
739 from these assays are used to limit, to an acceptable level, the risks of febrile reaction in the
740 patient exposed to the product of concern by injection and/or implantation. The current U.S.
741 legislation requiring the use of pyrogenicity testing is set forth in the Federal Food, Drug, and
742 Cosmetic Act (U.S.C., Title 21, Chapter 9). In addition, the U.S. Pharmacopeia (USP)
743 maintains sterility requirements for pharmaceuticals that include pyrogenicity testing. As
744 detailed in **Table 1-1**, the U.S. Food and Drug Administration (FDA) is the principal U.S.
745 regulatory agency that requires pyrogenicity testing, with different Centers within the FDA
746 regulating the affected products. The Center for Biologics Evaluation and Research, the
747 Center for Drug Evaluation and Research, the Center for Devices and Radiological Health,
748 and the Center for Veterinary Medicine require that human injectable drugs (including
749 biological products), animal injectable drugs, and medical devices be tested for the presence
750 of pyrogenic substances. **Table 1-1** also shows the statutory protocol requirements used by
751 each FDA Center, along with the comparable enabling legislation and statutory protocol
752 requirements of the European Union (EU) member nations.

753

753 **Table 1-1 Summary of U.S. and European Legislation and Statutory Protocol**
 754 **Requirements for Pyrogenicity Testing**

Agency	Regulated Products	Legislation	Statutory Protocol Requirements	Non-Governmental Standards
United States				
FDA-CBER	Biological products	- Federal Food, Drug, and Cosmetic Act (U.S.C. Title 21, Chapter 9)	- 21 CFR 610.13	- USP28 NF23<85> - USP28 NF23<151> - ISO 10993-11
FDA-CDER	Human parenteral pharmaceuticals			
FDA-CDRH	Medical devices			
FDA-CVM	Veterinary pharmaceuticals			
Europe				
EDQM	Human/veterinary parenteral pharmaceuticals, biological products, medical devices	- Council Regulation (EEC) 230/9/93	- EP5.0 2.6.8 - EP5.0 2.6.14	- ISO 10993-11
EMEA		- Council Directive 93/39/EEC		
Regulatory Authorities for Individual EU Countries		- Council Directive 93/40/EEC		

755 Abbreviations: CBER = Center for Biologics Evaluation and Research; CDER = Center for Drug Evaluation and Research;
 756 CDRH = Center for Devices and Radiological Health; CFR = Code of Federal Regulations; CVM = Center for Veterinary
 757 Medicine; EDQM = European Directorate for the Quality of Medicines; EMEA = European Medicines Agency; EP =
 758 European Pharmacopoeia; EU = European Union; FDA = U.S. Food and Drug Administration; US =: U.S. Pharmacopeia
 759

760 **Criteria 2. Warranted, based on the extent of expected use or application and impact on**
 761 **human, animal, or ecological health.**

762 The proposed test methods are intended to replace a method that is used extensively in
 763 pharmaceutical, biological product, and medical device development and registration (i.e.,
 764 the RPT).

765 **Criteria 3: The potential for the proposed test methods, compared to current test**
 766 **methods accepted by regulatory agencies, to (a) refine animal use (decrease or**
 767 **eliminate pain and distress), (b) reduce animal use, or (c) replace animal use.⁸**

768 The two most common pyrogen tests presently used (i.e., RPT, BET) require the use of
 769 animals. The RPT is performed in rabbits that can be maintained and used for additional
 770 pyrogen tests. Although the *in vitro* BET is performed using haemolymph (the equivalent of

⁸ Refinement alternative is defined as a new or revised test method that refines procedures to lessen or eliminate pain or distress to animals, or enhances animal well-being; Reduction alternative is defined as a new or revised test method that reduces the number of animals required; Replacement alternative is defined as a new or revised test method that replaces animals with nonanimal systems or one animal species with a phylogenetically lower one (e.g., a mammal with an invertebrate) (ICCVAM 1997).

771 blood) drawn from *Limulus polyphemus* (horseshoe crabs), which are subsequently returned
772 to the wild, there is some mortality associated with the procedure (which requires
773 approximately 20% of the animal's total blood volume). The proposed methods will reduce
774 and replace animal use because they rely on human white blood cells that can be obtained
775 with minimal risk from volunteers, or a human white blood cell line that can be maintained
776 by the test laboratory.

777 **Criteria 4: The potential for the proposed test method to provide improved prediction**
778 **of adverse health or environmental effects, compared to current test methods accepted**
779 **by regulatory agencies.**

780 Sufficient data are presented to allow an assessment of the performance of the proposed test
781 methods relative to the RPT (see Section 6.0). Because these methods are conducted using
782 cells of human origin, one might assume that they would better reflect the human
783 physiological response than current, non-human-based methods (i.e., RPT, BET), and thus
784 more effectively predict adverse effects.

785 **Criteria 5: The extent to which the test method provides other advantages (e.g., reduced**
786 **cost and time to perform) compared to current methods.**

787 Specific costs for the various tests have not been provided, and therefore such a definitive
788 determination cannot be made. Because the proposed test methods are reportedly more labor-
789 intensive than the RPT or BET, the costs may be greater. These costs, however, may be
790 offset by the costs of maintaining an animal facility for the rabbits used in the RPT. The
791 possibility of adapting the proposed methods to high throughput procedures may make them
792 more cost effective.

793 1.2.2 Intended Uses of the Proposed *In Vitro* Pyrogenicity Test Methods

794 The proposed test methods are intended as an end product release test for the identification of
795 pyrogens in human and animal parenteral drugs, biological products, and medical devices.
796 Results from pyrogenicity testing are used to limit, to an acceptable level, the risks of febrile
797 reaction to the injection and/or implantation of the product of concern.

798 1.2.3 Similarities and Differences in the Endpoints Measured by the Proposed Test
799 Methods and the *In Vivo* Reference Test Method

800 The endpoint measured in the *in vitro* pyrogenicity test methods is release of
801 proinflammatory cytokines, either IL-1 β or IL-6, in response to a test substance challenge,
802 depending on the specific cell type employed. The RPT involves measuring the rise in body
803 temperature evoked in rabbits by the intravenous injection of a test solution. Although there
804 is no direct association between the endpoints measured in the *in vitro* blood cell assays and
805 the *in vivo* rabbit assay, cytokine release is involved in the development of the inflammatory
806 response, which can lead to an increase in body temperature. Therefore, the *in vitro* release of
807 proinflammatory cytokines, such as IL-1 β and IL-6, is intended to presage the onset of such a
808 cascade of responses resulting in a fever response. The cell types used for the various *in vitro*
809 methods include those that would be directly associated with an inflammatory response (i.e.,
810 monocytes/monocytoid cells) *in vivo*. Both the *in vitro* and *in vivo* tests provide quantitative
811 data that can be applied to specific decision criteria to identify a pyrogenic reaction.

812 1.2.4 Use of the Proposed Test Methods in an Overall Strategy of Hazard or Safety
813 Assessment

814 As detailed in **Table 1-1**, there are current U.S. and European regulatory requirements to test
815 pharmaceutical products, biological products, and medical devices for pyrogenicity. The
816 pyrogenicity tests that are currently acceptable to regulatory authorities require intact animals
817 (rabbits) or an *in vitro* test that requires the use of horseshoe crab haemolymph (BET).
818 According to ECVAM, the *in vitro* human blood cell test methods are intended to replace the
819 RPT for the identification of pyrogens where: (a) the test material is incompatible with the
820 BET; or (b) the test material contains a non-endotoxin mediated pyrogen (although as
821 detailed in **Section 3.0**, only Gram-negative endotoxin was included in the validation study,
822 suggesting that other types of pyrogens have not been adequately validated – see **Section**
823 **1.1.2)**⁹.

⁹ Additional information on testing of other types of pyrogens was also provided by ECVAM and is included in **Appendix B**.

824 1.3 Scientific Basis for the *In Vitro* Pyrogenicity Test Methods

825 1.3.1 Purpose and Mechanistic Basis of the *In Vitro* Pyrogenicity Test Methods

826 The proposed methods are intended to identify pyrogens in parenteral pharmaceuticals,
827 medical devices, and human biological products. These methods provide an *in vitro* model of
828 the initiation of the human fever response by measuring proinflammatory cytokine release
829 (i.e., IL-1 β or IL-6) from human monocytes/monocytoid cells exposed to pyrogens. These
830 proinflammatory cytokines are associated with the initiation of the *in vivo* fever response.

831 1.3.2 Similarities and Differences of Modes of Action Between the *In Vitro* Pyrogenicity 832 Test Methods and the Fever Response in Humans and/or Rabbits

833 As detailed in **Section 1.2.3**, each of the five proposed *in vitro* test methods do not measure
834 actual fever induction, but use proinflammatory cytokine release (i.e., IL-1 β or IL-6) from
835 human monocytoid cells as an indicator of the presence of a pyrogenic substance. By
836 comparison, the RPT involves measuring the effect of concern, i.e., a change in body
837 temperature in rabbits over a specified time period following an intravenous injection of a
838 test substance.

839 Proinflammatory cytokines such as IL-1 β and IL-6 have been characterized as endogenous
840 pyrogens that are released from monocytoid cells in response to the presence of endotoxin
841 and other pyrogenic substances, and have been associated with the fever response induced by
842 pyrogenic substances in both humans and rabbits (Dinarello et al. 1999). The pyrogenicity of
843 endotoxin is believed to be mediated through activation of the Toll-like (TLR-4) receptor on
844 the blood cells. TLR-4 activation induces the mitogen-activated protein kinase (MAPK)
845 signal transduction cascade, which activates proinflammatory cytokine expression. There is
846 no direct evidence that circulating cytokines penetrate the blood-brain barrier, but it is more
847 likely that they produce effects on the rich vascular network or circumventricular organs of
848 the hypothalamus, or induce release of other effectors from endothelial cells (e.g.,
849 prostaglandins such as PGE₂) that increase levels of neurotransmitters (e.g., cAMP) which
850 modulate the thermoregulatory center of the hypothalamus and induce a fever response.
851 (Dinarello 1999; Beutler and Rietschel 2003).

852 As indicated above, the proposed *in vitro* test methods use human cells whereas the RPT and
853 BET procedures use rabbits and horseshoe crab haemolymph, respectively. Although the
854 relative sensitivities of each species to Gram-negative endotoxins vary, the responses of
855 humans, horseshoe crabs (via haemolymph gelation), and rabbits to these pyrogens have been
856 studied extensively, and test methods based on blood products or blood cells from each of
857 these species appear to be capable of responding to pyrogens (Cooper et al 1971; Greisman
858 and Hornick, 1969; Hoffman et al. 2005). However, there are documented cases of medicinal
859 products, interfering substances, and specific pyrogenic substances that yield false-positive or
860 false-negative results in either the RPT or BET tests (Carlin and Viitanen 2003, 2005;
861 Moesby et al. 2000; Hartung et al. 2001; Ochiai et al. 2001).

862 1.3.3 Range of Substances Amenable to the *In Vitro* Pyrogenicity Test Methods and 863 Limits of These Methods

864 The proposed methods are intended to be used for the identification of pyrogenic substances
865 in parenteral pharmaceuticals, biological products, and medical devices. Because they are
866 based on isolated human monocytes/monocytoid cells, they are considered capable of
867 detecting both Gram-negative endotoxin and non-endotoxin-based pyrogens. While **Section**
868 **9.0** summarizes a number of published studies that have measured cytokine release following
869 exposure to non-endotoxin pyrogens (e.g., lipoteichoic acid), the ECVAM validation studies
870 focused specifically on Gram-negative endotoxin due to the unavailability of standardized,
871 non-endotoxin pyrogens (see **Section 3.0**). Because these test methods measure the release
872 of proinflammatory cytokines, test substances that modify this response (e.g., anti-
873 inflammatory or immunosuppressant drugs) are considered inappropriate for testing. Each
874 test protocol includes an interference test for identifying such substances so that they can be
875 excluded from testing.

876 **1.4 Validation of the *In Vitro* Pyrogenicity Test Methods**

877 The ICCVAM Authorization Act (Sec. 4(c) mandates that “[e]ach Federal Agency ... shall
878 ensure that any new or revised ... test method ... is determined to be valid for its proposed
879 use prior to requiring, recommending, or encouraging [its use].” (Public Law [P.L.] 106-
880 545).

881 Validation is the process by which the reliability and relevance of an assay for a specific
882 purpose are established (ICCVAM 1997). Relevance is defined as the extent to which an
883 assay will correctly predict or measure the biological effect of interest (ICCVAM 1997). For
884 the *in vitro* pyrogenicity test methods described in this BRD, relevance is restricted to how
885 well the assays detect the presence of Gram-negative endotoxin. Reliability is defined as the
886 reproducibility of a test method within and among laboratories and should be based on
887 performance with a diverse set of substances that are representative of the types of chemical
888 and product classes that are expected to be tested and the range of responses that needs to be
889 identified. The validation process is designed to provide data and information that will allow
890 U.S. Federal agencies to develop guidance on the development and use of *in vitro*
891 pyrogenicity test methods for the detection of Gram-negative endotoxins.

892 The first stage in the evaluation of a new test procedure is the preparation of a BRD that
893 presents and evaluates the relevant data and information about the method, including its
894 mechanistic basis, proposed uses, reliability, and performance characteristics (ICCVAM
895 1997). This BRD summarizes the available information on each of the five *in vitro*
896 pyrogenicity test methods listed in **Section 1.1.1**.

897 Where adequate data are available, the qualitative and quantitative performances of the
898 assays are evaluated, and the reliability of each new method is compared with the reliability
899 of the other methods and the currently accepted method. This BRD will aid in identifying
900 essential test method components that should be considered during the development and
901 validation of the various methods, and determine whether there are sufficient data to support
902 a recommendation that a standardized protocol be developed for use by the regulatory
903 agencies.

904 **1.5 Search Strategies and Selection of Citations for the *In Vitro* Pyrogenicity BRD**

905 NICEATM conducted an online literature search for relevant information on the proposed
906 test methods using multiple internet databases (i.e., PubMed, SCOPUS, TOXLINE, Web of
907 Science). Specifically, records were sought using various combinations of the terms *in vitro*,
908 WBC, whole blood, PBMN, MONO MAC6, MM6, endotoxin, LPS, pyrogen, LAL, BET,
909 IL-1, and IL-6. This search was conducted to supplement and update the list of peer-reviewed
910 publications related to *in vitro* pyrogen testing that was provided in the ECVAM BRDs. U.S.,

911 EU, and Japanese pyrogenicity test guidelines were obtained from relevant regulatory
912 agencies via the internet or through direct requests. A resulting database of 315 references is
913 currently maintained as an EndNote[®] v8.0 electronic file. This database confirmed that the
914 lists of references included in the ECVAM BRDs are complete and up-to-date.

915 **2.0 IN VITRO PYROGENICITY TEST METHOD PROTOCOL**
916 **COMPONENTS**

917 **2.1 Overview of How the In Vitro Pyrogenicity Test Methods Are Conducted**

918 Although there are differences among the *in vitro* pyrogenicity test methods considered in
919 this BRD, the basic procedural steps are consistent across all five methods:

- 920 • Interference testing is performed to verify that a test substance does not
921 interfere with either the cell system used or with the specific cytokine-specific
922 enzyme-linked immunosorbent assay (ELISA).
- 923 • The test substance is mixed with a suspension of human-derived blood cells
924 and incubated for a specified period of time (i.e., 16 to 24 hr).
- 925 • The concentration of the specific proinflammatory cytokine (e.g., IL-1 β , IL-6)
926 is measured using an ELISA, and is compared to the response curve of an
927 endotoxin standard
- 928 • An internationally accepted endotoxin standard derived from *Escherischia*
929 *coli* (*E. coli*) (i.e., WHO-lipopolysaccharide [LPS] 94/580 [*E. coli*
930 0113:h10:K-]), or an endotoxin standard that has been calibrated against this
931 standard, is used to generate the standard response curve for the assay. The
932 endotoxin activity of a test substance is calculated by comparing the induced
933 cytokine release with that induced by the endotoxin standard.
- 934 • A product “passes” (i.e., is considered negative for endotoxin pyrogen
935 activity) if the cytokine response to the test substance is less than that induced
936 by 0.5 endotoxin units (EU)/mL.

937 **2.2 Description and Rationale for the Test Method Components for Proposed**
938 **Standardized Protocols**

939 The standard operating procedures (SOP) for each method assessed in the ECVAM
940 validation studies are provided as appendices to each BRD (see Appendix A of each BRD).
941 As indicated in **Section 2.1**, the essential principles of each protocol are common among the
942 five methods reviewed. These include isolating and/or culturing human monocytoid cells

943 (either included in whole blood, separated as a fraction [i.e., PBMCs], or as cell line [i.e.,
944 MM6]), performing interference testing with each substance, treating the cells in suspension
945 with a test substance, collecting cytokine release data, and evaluating the data in relation to
946 the proposed prediction model. **Table 2-1** provides a more complete comparison of the
947 similarities and differences among the five test protocols. No rationale was provided by
948 ECVAM for use of the particular primary blood cells in the various test methods; in such
949 cases, their selection is justified by historical use. The use of the MM6 cell line was justified
950 based on mechanistic considerations and its response to endotoxins.

951 Table 2-1 In Vitro Pyrogenicity Test Method Components

	Test Method Component	WB/IL-1 ¹	cryo WB/IL-1	WB/IL-6	PBMC/IL-6 ¹	MM6/IL-6
Materials required	Source of cells	Human whole blood	Human whole blood	Human whole blood	Human whole blood	MM6 cell bank (original cell line maintained by Prof. H. Ziegler-Heitbrock, U. Munich)
	Laboratory equipment	<ul style="list-style-type: none"> Incubator (37°C, 5% CO₂) Adjustable pipetters (2-2000 µL) Multichannel pipetter (8- or 12-channel) Serological pipets Pyrogen-free tips Centrifuge Vortex mixer Heparinized blood tubes Pyrogen-free hypodermic needles Pyrogen-free microfuge tubes OR 96-well plates (depending on tube or plate method) Microtiter plate reader (450 nm and 600-690 nm) Data analysis software 	<ul style="list-style-type: none"> Incubator (37°C, 5% CO₂, Laminar flow hood (recommended)) Adjustable pipetters (2-2000 µL) Multichannel pipetter (8- or 12-channel) Pyrogen-free tips Centrifuge Vortex mixer Pyrogen-free tissue culture plasticware (e.g., flasks, 96-well plates, centrifuge tubes) Microtiter plate reader (450 nm and 600-690 nm) Data analysis software 	<ul style="list-style-type: none"> Incubator (37°C, 5% CO₂, humidified) Laminar flow hood (Class II) Centrifuge Water bath pH meter Microtiter plate reader (450 nm and 600-690 nm capable) Data analysis software Adjustable pipetters (2-2000 µL) Multichannel pipetter (8- or 12-channel) Heparinized blood tubes Pyrogen-free pipet tips Vortex mixer Pyrogen-free hypodermic needles Pyrogen-free tissue culture plasticware (e.g., flasks, 96-well plates, centrifuge tubes) 	<ul style="list-style-type: none"> Incubator (37°C, 5% CO₂, humidified) Inverted microscope Hemocytometer Laminar flow hood (Class II) Centrifuge Water bath pH meter Microtiter plate reader (450 nm and 600-690 nm capable) Data analysis software Adjustable pipetters (2-2000 µL) Multichannel pipetter (8- or 12-channel) Heparinized blood tubes Pyrogen-free pipet tips Vortex mixer Pyrogen-free hypodermic needles Pyrogen-free tissue culture plasticware (e.g., flasks, 96-well plates, centrifuge tubes) 	<ul style="list-style-type: none"> Incubator (37°C, 5% CO₂, humidified) Inverted microscope Hemocytometer Laminar flow hood (Class II) Centrifuge Water bath pH meter Microtiter plate reader (450 nm and 600-690 nm capable) Data analysis software Adjustable pipetters (2-2000 µL) Multichannel pipetter (8- or 12-channel) Pyrogen-free pipet tips Vortex mixer Pyrogen-free tissue culture plasticware (e.g., flasks, 96-well plates, centrifuge tubes)
	Culture medium	None - whole blood is diluted with 0.9% NaCl	RPMI Complete Medium <ul style="list-style-type: none"> RPMI 1640 (part of the Endosafe Kit for cryoblood) no specific additives needed 	None - whole blood is diluted with 0.9% NaCl	RPMI Complete Medium <ul style="list-style-type: none"> RPMI 1640 HSA L-Glutamine (2 mM) Penicillin/streptomycin 	RPMI Complete Medium <ul style="list-style-type: none"> RPMI 1640 HSA L-Glutamine (2 mM) MEM non-essential amino acid solution (0.1 mM) HEPES (20 mM) Bovine insulin (0.23 IU/mL) Oxaloacetic acid (1 mM)

	Test Method Component	WB/IL-1 ¹	cryo WB/IL-1	WB/IL-6	PBMC/IL-6 ¹	MM6/IL-6
	Other reagents	<ul style="list-style-type: none"> Validated IL-1β ELISA kit 0.9% NaCl (pyrogen-free) Water for injection Endotoxin standard 	<ul style="list-style-type: none"> Validated IL-1β ELISA kit 0.9% NaCl (pyrogen-free) Water for injection Liquid nitrogen Endotoxin standard 	<ul style="list-style-type: none"> Validated IL-6 ELISA kit 0.9% NaCl (pyrogen-free) Water for injection Endotoxin standard 	<ul style="list-style-type: none"> Validated IL-6 ELISA kit 0.9% NaCl (pyrogen-free) Water for injection Trypan blue Endotoxin standard 	<ul style="list-style-type: none"> Sodium pyruvate (1 mM) Validated IL-6 ELISA kit 0.9% NaCl (pyrogen-free) Water for injection Trypan blue DMSO Endotoxin standard
Dose selection procedures		Interference testing performed to determine the lowest dilution of the test product necessary to achieve an acceptable endotoxin spike recovery (i.e., 50% to 200% recovery) ³	Interference testing performed to determine the lowest dilution of the test product necessary to achieve an acceptable endotoxin spike recovery (i.e., 50% to 200% recovery) ³	Interference testing performed to determine the lowest dilution of the test product necessary to achieve an acceptable endotoxin spike recovery (i.e., 50% to 200% recovery) ³	Interference testing performed to determine the lowest dilution of the test product necessary to achieve an acceptable endotoxin spike recovery (i.e., 50% to 200% recovery) ³	Interference testing performed to determine the lowest dilution of the test product necessary to achieve an acceptable endotoxin spike recovery (i.e., 50% to 200% recovery) ³
Endpoints measured		IL-1 β release via ELISA	IL-1 β release via ELISA	IL-6 release via ELISA	IL-6 release via ELISA	IL-6 release via ELISA
Exposure of the test substance	Pre-test preparation of cells	Collect whole blood, heparinize, and use within 4 hr Plate Method: same collection procedure	<ul style="list-style-type: none"> Collect blood, heparinize, and cryopreserve according to the Konstanz or PEI method Prior to testing, thaw blood at 37°C for 15 min 	Collect whole blood, heparinize, and use within 4 hr	<ul style="list-style-type: none"> Collect whole blood and isolate PBMCs by centrifugation Resuspend PBMCs in RPMI-C (1x10⁶ cells/mL) (use PBMCs within 4 hr of initial blood collection) 	<ul style="list-style-type: none"> Incubate MM6 cells (4x10⁵ cells/mL media) for 24 hr Resuspend cells (2.5x10⁶ cells/mL)³ prior to testing
	Application of the test substance	Tube method: In a microfuge tube mix 1000 μ L 0.9% NaCl+100 μ L sample+100 μ L blood Plate method: In a 96-well plate mix 200 μ L 0.9% NaCl+20 μ L sample+20 μ L blood	Konstanz method: In a 96-well plate mix 200 μ L RPMI+20 μ L sample+20 μ L blood PEI Method: In a 96-well plate mix 180 μ L RPMI + 20 μ L sample+40 μ L blood	In a 96-well plate: Mix 50 μ L standards/samples+100 μ L 0.9% NaCl+50 μ L blood	In a 96-well plate: Mix 50 μ L standards/samples+100 μ L RPMI-C+100 μ L PBMCs	In a 96-well plate: Mix 50 μ L standards/samples+100 μ L RPMI-C+100 μ L cells in suspension
	Duration of exposure	10-24 hr	10-24 hr (NOTE: for Konstanz method, after incubation, freeze wells at -20°C or -80°C and then thaw at 37°C)	16-24 hr	16-24 hr	16-24 hr
	Material used for ELISA	Tube method: centrifuge 2 min @ 10,000g-test supernatant	Blood/RPMI/sample mixture	Blood/saline/sample mixture	Cell supernatant	Cell supernatant

	Test Method Component	WB/IL-1 ¹	cryo WB/IL-1	WB/IL-6	PBMC/IL-6 ¹	MM6/IL-6
		Plate method: mix each well be pipetting and test resuspended mixture				
Known limits of use		Intended for parenteral pharmaceuticals, biological products, and medical devices that have been qualified through interference testing	Intended for parenteral pharmaceuticals, biological products, and medical devices that have been qualified through interference testing	Intended for parenteral pharmaceuticals, biological products, and medical devices that have been qualified through interference testing	Intended for parenteral pharmaceuticals, biological products, and medical devices that have been qualified through interference testing	Intended for parenteral pharmaceuticals, biological products, and medical devices that have been qualified through interference testing
Nature of the response assessed		Pyrogenic substances induce the release of proinflammatory cytokines (e.g., IL-1 β) from monocytoid cells present in human whole blood	Pyrogenic substances induce the release of proinflammatory cytokines (e.g., IL-1 β) from monocytoid cells present in human whole blood	Pyrogenic substances induce the release of proinflammatory cytokines (e.g., IL-6) from monocytoid cells present in human whole blood	Pyrogenic substances induce the release of proinflammatory cytokines (e.g., IL-6) from PBMCs	Pyrogenic substances induce the release of proinflammatory cytokines (e.g., IL-6) from these immortalized monocytoid cells
Appropriate controls	Positive control (PC)	0.5 EU/mL WHO-LPS 94/580 [<i>E. coli</i> 0113:h10:K-] ⁵	0.5 EU/mL WHO-LPS 94/580 [<i>E. coli</i> 0113:h10:K-] ⁵	0.5 EU/mL WHO-LPS 94/580 [<i>E. coli</i> 0113:h10:K-] ⁵	0.5 EU/mL WHO-LPS 94/580 [<i>E. coli</i> 0113:h10:K-] ⁵	0.5 EU/mL WHO-LPS 94/580 [<i>E. coli</i> 0113:h10:K-] ⁵
	Negative control (NC)	0.9% NaCl	0.9% NaCl	0.9% NaCl	0.9% NaCl	0.9% NaCl
	Positive product control (PPC)	Test substance spiked with 0.5 EU/mL	Test substance spiked with 0.5 EU/mL	Test substance spiked with 0.5 EU/mL	Test substance spiked with 0.5 EU/mL	Test substance spiked with 0.5 EU/mL
	Negative product control (NPC)	Test substance spiked with 0.9% NaCl	Test substance spiked with 0.9% NaCl	Test substance spiked with 0.9% NaCl	Test substance spiked with 0.9% NaCl	Test substance spiked with 0.9% NaCl
Acceptable ranges of control responses	Positive controls	<ul style="list-style-type: none"> PC OD 1.6-fold>NC OD PPC OD 1.6-fold>NPC OD PPC OD should be within 50% to 200% of the PC OD 	<ul style="list-style-type: none"> PC OD 1.6-fold>NC OD PPC OD 1.6-fold>NPC OD PPC OD should be within 50% to 200% of the PC OD 	PPC OD should be within 50% to 200% of the PC OD	<ul style="list-style-type: none"> PPC OD should be within 50% to 200% of the PC OD 1 EU/mL standard OD>1000 pg/mL IL-6 standard 	<ul style="list-style-type: none"> PC OD\pm20% of the expected value (i.e., 0.5 EU/mL) PPC OD should be within 50% to 200% of the PC OD
	Negative controls	NC OD \leq 0.100	NC OD \leq 0.100	NC OD<200 pg/mL IL-6 standard	NC OD<500 pg/mL IL-6 standard	• NC OD<0.200
Nature of data to be collected and methods used for data collection		<ul style="list-style-type: none"> IL-1β release from monocytoid cells present in whole blood is quantified with an ELISA for all standards and samples. The endotoxin content of a test substance is calculated by comparing the induced IL-1β release with that induced by the 	<ul style="list-style-type: none"> IL-1β release from monocytoid cells present in whole blood is quantified with an ELISA for all standards and samples. The endotoxin content of a test substance is calculated by comparing the induced IL-1β release with that induced by the 	<ul style="list-style-type: none"> IL-6 release from monocytoid cells present in whole blood is quantified with an ELISA for all standards and samples. The endotoxin content of a test substance is calculated by comparing the induced IL-6 release with that induced by the 	<ul style="list-style-type: none"> IL-6 release from PBMCs is quantified with an ELISA for all standards and samples. The endotoxin content of a test substance is calculated by comparing the induced IL-6 release with that induced by the endotoxin standard curve concentrations 	<ul style="list-style-type: none"> IL-6 release from MM6 cells is quantified with an ELISA for all standards and samples. The endotoxin content of a test substance is calculated by comparing the induced IL-6 release with that induced by the endotoxin standard curve concentrations

	Test Method Component	WB/IL-1 ¹	cryo WB/IL-1	WB/IL-6	PBMC/IL-6 ¹	MM6/IL-6
		endotoxin standard curve concentrations	endotoxin standard curve concentrations	endotoxin standard curve concentrations		
Type of media in which data are stored		Electronic files	Electronic files	Electronic files	Electronic files	Electronic files
Measures of variability		Mean±SD of the OD for each test substance/standard	Mean±SD of the OD for each test substance/standard	Mean±SD of the OD for each test substance/standard	Mean±SD of the OD for each test substance/standard	Mean±SD of the OD for each test substance/standard
Decision criteria and basis for their selection		≥0.5 EU/mL indicates a pyrogenic response based on the pyrogen threshold value as detailed in Section 4.2	≥0.5 EU/mL indicates a pyrogenic response based on the pyrogen threshold value as detailed in Section 4.2	≥0.5 EU/mL indicates a pyrogenic response based on the pyrogen threshold value as detailed in Section 4.2⁶	≥0.5 EU/mL indicates a pyrogenic response based on the pyrogen threshold value as detailed in Section 4.2⁷	≥0.5 EU/mL indicates a pyrogenic response based on the pyrogen threshold value as detailed in Section 4.2

Abbreviations: DMSO = dimethylsulfoxide; EU = endotoxin unit; HIFCS = heat-inactivated fetal calf serum; HSA = human serum albumin; MEM = minimum essential medium; MM6 = MONO MAC-6; NC = negative control; NPC = negative product control; OD = optical density; PBS = phosphate buffered saline; PC = positive control; PPC = positive product control

¹As described in **Section 1.1.1**, a catch-up validation studies were also conducted to evaluate the performance of the WB/IL-1 test method using 96-well plates, and the PBMC/IL-6 test method when using cryopreserved PBMCs. The plating procedure (WB/IL-1) and the cryopreservation procedure (PBMC/IL-6) are the only differences in the test method protocols (see **Appendix A**).

²Medium should be qualified for testing by a valid LAL test (e.g., USP28NF23<85>) indicating that the endotoxin contamination is <0.06 IU/mL); FBS concentration for MM6 cells varies based on whether it is for maintenance/propagation (10%) or assay (2%) conditions.

³Dilution of the test material should not exceed the maximum valid dilution (MVD), where MVD = (endotoxin limit concentration)/(detection limit of the assay)

⁴Cell numbers represent viable cells based on trypan blue exclusion

⁵Or another endotoxin calibrated against this standard

⁶Specifies that each substance must "pass" (i.e., be non-pyrogenic) in blood from three different donors.

⁷Includes a sequential decision strategy in which 3 to 4 donors are tested per substance. 1) If all donors show negative - product is non-pyrogenic; 2) If ≥2 donors show a positive - product is pyrogenic;

3) If only one donor shows a positive, an additional 3 to 4 donors are tested and if no more than one donor is positive (out of 6 to 8 donors) - product is non-pyrogenic; otherwise, product is pyrogenic.

952
953
954
955
956
957
958
959
960
961
962
963

964 2.2.1 Methods Used to Analyze the Data, Including Methods to Analyze for Interference
965 With the Assay

966 Once a substance has been tested in the requisite number of donor samples (see **Section**
967 **2.2.2**), the resulting sample test medium (as indicated in **Table 2-1**) is assayed in
968 quadruplicate in the relevant cytokine ELISA. Outliers are identified using Dixon's test
969 ($p>0.05$), and are excluded from the calculations of endotoxin content. Endotoxin standard
970 curves are included in each assay, from which the endotoxin content of each replicate is
971 estimated using a 4-parameter logistic model.

972 As indicated in **Table 2-1**, mean optical density (OD) readings are calculated for the positive
973 and negative control samples, as well as for the relevant positive and negative product
974 controls. The acceptable range of the positive product control (50% to 200% of the positive
975 control response) defines the threshold for interference with the test system. If the positive
976 product control response falls outside of this range, the samples are then assayed at the
977 lowest dilution that does not cause interference.

978 2.2.2 Decision Criteria and the Basis for the Prediction Model Used to Identify a
979 Pyrogenic Substance

980 As described in **Section 4.2**, historical RPT data were used to establish a threshold pyrogen
981 dose (i.e., the endotoxin dose at which fever was induced in 50% of the rabbits), which was
982 determined to be 5 EU/kg. Based on the largest allowable volume for injection in rabbits (10
983 mL/kg), the limit of detection that the *in vitro* pyrogen tests must meet was defined as 0.5
984 EU/mL. Accordingly, the prediction model for each test method was established based on
985 this limit of detection (i.e., a substance is considered pyrogenic if the mean response is
986 greater than or equal to the 0.5 EU/mL standard).

987 For three of the test methods, results from multiple donors (cryo WB/IL-1 [n=5], WB/IL-6
988 [n=3], and PBMC/IL-6 [n=3 to 4]) are required to determine the pyrogenicity of a test
989 substance. In contrast, a single donor sample is used for the WB/IL-1 test method, as is a
990 single cell sample for the MM6/IL-6 test method. As outlined in **Table 2-2**, unlike the cryo
991 WB/IL-1 test method, the WB/IL-6 and PBMC/IL-6 test methods employ a decision strategy
992 that takes into account the individual responses of each donor sample.

993 **Table 2-2 Prediction Model Used for *In Vitro* Pyrogenicity Test Methods**

Test Method	No. Donors	No. Positive	No. Negative	Decision
PBMC/IL-6	4 ¹	4	0	Pyrogenic
		3	1	Pyrogenic
		2	2	Pyrogenic
		1	3	Non-pyrogenic
	3 ¹	0	4	Non-pyrogenic
		3	0	Pyrogenic
		2	1	Pyrogenic
		1	2	Non-pyrogenic
WB/IL-6	3	0	3	Non-pyrogenic
		3	3	Pyrogenic
		2	1	Pyrogenic
		1	2	Non-pyrogenic
cryo WB/IL-1	5 (pooled) ²	0	3	Non-pyrogenic
		1	0	Pyrogenic
WB/IL-1	1	0	1	Non-pyrogenic
		1	0	Pyrogenic
MM6/IL-6	NA ³	0	1	Non-pyrogenic
		1	0	Pyrogenic

994 Abbreviations: cryo = Cryopreserved; IL-1 = Interleukin-1; IL-6 = Interleukin-6; MM6 = Mono Mac 6; NA = Not
 995 applicable; No. = Number; WB = Whole blood

996 ¹Samples are collected from four donors for the PBMC/IL-6 test method, one donor sample may be excluded based on
 997 quality criteria, in which case the prediction model may be applied to results from three donors.

998 ²Samples are collect from five donors for the cryo WB/IL-1 test method and pooled prior to cryopreservation

999 ³Not applicable, because source material is obtained from an immortalized cell line.

1000

1001 2.2.3 Information and Data to be Included in the Study Report and Availability of1002 Standard Forms for Data Collection and Submission

1003 The test report should include the following information, if relevant to the conduct of the
 1004 study:

1005 *Test Substances and Control Substances*

- 1006 • Name, and type (e.g., pharmaceutical, biological product, medical device
 1007 eluate, etc.) of test product
- 1008 • Purity and composition of the test substance or preparation
- 1009 • Physicochemical properties such as physical state, volatility, pH, stability,
 1010 chemical class, water solubility relevant to the conduct of the study
- 1011 • Treatment of the test/control substances prior to testing, if applicable (e.g.,
 1012 vortexing, sonication, warming; solvent used)
- 1013 • Stability, if known

1014 *Justification of the Specific Protocol(s) Used*

1015 *Test Method Integrity*

- 1016 • The procedure used to ensure the integrity (i.e., accuracy and reliability) of the
- 1017 test method over time
- 1018 • If the test method employs proprietary components, documentation of the
- 1019 procedure used to ensure their integrity from “lot-to-lot” and over time

1020 *Criteria for an Acceptable Test*

- 1021 • Acceptable concurrent positive control ranges based on historical data from
- 1022 the testing laboratory (which should be included in the report)
- 1023 • Acceptable negative control data, including historical control ranges from the
- 1024 testing laboratory (which should be included in the report)

1025 *Test Conditions*

- 1026 • Cell system used; donor information, if relevant
- 1027 • Calibration information for the equipment used for measuring cytokine release
- 1028 (e.g., spectrophotometer)
- 1029 • Details of test procedure used
- 1030 • Description of modifications of the test procedure made by the testing
- 1031 laboratory for the substance being tested
- 1032 • Reference to the laboratory’s historical data for the cell system and protocol
- 1033 • Description of data and quality assurance evaluation criteria used

1034 *Results*

- 1035 • Tabulation of data from individual test samples

1036 *Description of Other Effects Observed*

1037 *Discussion of the Results*

1038 *Conclusion*

1039 *A Good Laboratory Practice (GLP) Quality Assurance Statement*

1040 • This statement addresses all GLP inspections and audits made during the
1041 study, and the dates the results were reported to the Study Director. This
1042 statement also serves to confirm that the final report reflects the raw data.

1043 Reporting requirements for GLP-compliant studies are provided in the relevant guidelines
1044 (e.g., OECD 1998; EPA 2003a, 2003b; FDA 2003).

1045 **2.3 Basis for Selection of the Test Method Systems**

1046 One of the difficulties associated with the currently required pyrogenicity test methods (i.e.,
1047 RPT and BET) is that both require extrapolation of the response from a non-human system to
1048 the human. In contrast, and as discussed in **Section 1.1.1**, all five of these test methods
1049 employ human cells in an attempt to mimic the human fever response *in vitro*. Because these
1050 test methods are conducted using cells of human origin, it is proposed that they will better
1051 reflect the human physiological response than current, non-human methods, and thus more
1052 effectively predict human adverse effects.

1053 The whole blood test methods (i.e., WB/IL-1, cryo WB/IL-1, WB/IL-6) offer the
1054 convenience of performing the assay directly on a human blood sample, with minimal pretest
1055 preparation. The cryo WB/IL-1 test method was developed to offer the convenience of an
1056 increased time interval between the time of blood collection and the time a test is initiated
1057 (since the fresh blood methods require testing within four hours of collection), as well as
1058 increased standardization through the pooling of five donor samples to produce a larger
1059 sample bank of cells to use in the test. The MM6/IL-6 test method provides increased
1060 standardization by using an immortalized cell line that may be maintained in the laboratory
1061 indefinitely, and transferred among laboratories. Finally, the PBMC/IL-6 test method was
1062 developed in an attempt to pyrogen detection sensitivity by using the monocyte fraction of
1063 whole blood, which is considered to be the most sensitive human blood cell type to the
1064 presence of endotoxin.

1065 **2.4 Proprietary Components**

1066 Data from the test methods that use the IL-6 endpoint (i.e., WB/IL-6, PBMC/IL-6, MM6/IL-
1067 6) were obtained using a Novartis-developed IL-6 ELISA assay. The monoclonal anti-IL-6
1068 antibody used in this ELISA was developed using an in-house IL-6 clone at Novartis (i.e.,

1069 clone 6). Unless the Novartis assay is made publicly available, alternative, commercially
1070 available IL-6 ELISA assays will be necessary, and would have to be individually validated
1071 for use in these procedures.

1072 The original MM6 cell line used in the MM6/IL-6 test method is maintained by Prof. H.
1073 Ziegler-Heitbrock at the University of Munich. According to Section 2.4 of the ECVAM
1074 MM6/IL-6 BRD (see **Appendix A**), these cells are available for research purposes by request
1075 to Prof. Ziegler-Heitbrock directly, or the German Collection of Microorganisms and Cell
1076 Cultures in Germany. Once obtained, the testing laboratory may establish a working cell-
1077 bank. However, the BRD also states that conditions for licensing of the MM6 cells are to be
1078 negotiated with Prof. Ziegler-Heitbrock directly.

1079 According to the US Patent and Trademark Office (USPTO), patents are held for "Test for
1080 determining pyrogenic effect of a material" (U.S. 5,891,728, Apr 6, 1999), and "Pyrogenicity
1081 test for use with automated immunoassay systems" (U.S. 6,696,261 B2, Feb 24, 2004). These
1082 patents cover the WB/IL-1 and WB/IL-6 test methods, respectively. In addition, and related
1083 to the cryo WB/IL-1 test method, there is a patent pending for "Test procedure with
1084 biological system - Preparations containing deep-frozen blood are used for determining blood
1085 response" (USPTO 436518000).

1086 There are several measures in the study validity criteria that may be used to verify the
1087 integrity of proprietary components. As outlined in **Table 2-1**, an endotoxin standard curve is
1088 established for each assay, which is in turn used to define the endotoxin activity of the test
1089 substances. In addition, positive and negative controls, along with positive and negative
1090 product controls, are used for interference testing, and serve as internal controls for each
1091 assay.

1092 **2.5 Number of Replicates**

1093 2.5.1 Number of Donors

1094 There is no rationale provided for the number of donors included for each test method. As
1095 described in **Section 2.2.2**, samples from multiple donors are for three of the test methods.
1096 The cryo WB/IL-1 test method uses blood from five different donors and the WB/IL-6 and
1097 PBMC/IL-6 test methods use blood from at least three donors, which are tested individually.

1098 In contrast, a single donor sample is used for the WB/IL-1 test method, as is a single cell
1099 culture for the MM6/IL-6 test method.

1100 2.5.2 Number of Assay Replicates

1101 Once each substance has been tested in the requisite number of donor samples (see **Section**
1102 **2.2.2**), the resulting sample test medium is assayed in quadruplicate in the relevant cytokine
1103 ELISA. As indicated in **Section 2.2.1**, Dixon's test is used to detect outliers among the
1104 replicates. Section 2.5 of the ECVAM BRDs states that four replicates were chosen as it is
1105 considered the minimum number for inclusion in Dixon's test.

1106 **2.6 Modifications to the Test Method Protocols Based on ECVAM Validation** 1107 **Study Results**

1108 In the MM6/IL-6 test method, prevalidation studies demonstrated that pre-incubation of the
1109 cells at a defined initial concentration of 2×10^7 cells/50 mL RPMI-C for 24 hr greatly
1110 improved test method performance. Therefore, this modification was included in the
1111 validation study, and subsequently carried forward to the recommended MM6/IL-6 test
1112 method protocol.

1113 For the PBMC/IL-6 test method, a single blood donor was initially used as a source of
1114 PBMCs. However, the use of PBMCs from four separate donors (assayed individually) was
1115 shown to reduce variability, and this modification was carried forward in the recommended
1116 PBMC/IL-6 test method protocol.

1117 No modifications were made to the WB/IL-1, cryo WB/IL-1, and WB/IL-6 test method
1118 protocols as a result of the prevalidation or validation testing experiences.

1119 **2.7 Differences with Comparable Validated Test Methods with Established** 1120 **Performance Standards**

1121 The differences between the *in vitro* pyrogenicity test methods and the currently accepted
1122 pyrogenicity test methods (i.e., BET and RPT) are described in **Sections 1.2.3** and **1.3.2**.

1123 **3.0 SUBSTANCES USED FOR THE VALIDATION OF *IN VITRO***
1124 **PYROGENICITY TEST METHODS**

1125 **3.1 Rationale for the Substances or Products Selected for Testing**

1126 *In vitro* pyrogenicity test method validation studies should evaluate an adequate sample of
1127 substances and products of the types that are intended to be tested with these methods. The
1128 list of test substances selected for inclusion in the ECVAM validation studies consists solely
1129 of marketed parenteral pharmaceuticals that have been labeled as free from detectable
1130 pyrogens. No specific rationale was provided for the selection of these test substances. For
1131 evaluating test method performance, each test substance was spiked with a Gram-negative
1132 endotoxin standard (WHO-LPS 94/580 [*E. coli* O113:H10:K-]). Endotoxin was selected as a
1133 “model” pyrogen for inclusion based on its availability in a standardized form and because of
1134 the known ability of these cells to respond to endotoxin-based pyrogens.. No non-endotoxin-
1135 based pyrogenic substances are presently available in a standardized form.

1136 Each sample contained the appropriate endotoxin spike concentration when tested at its
1137 Maximum Valid Dilution (MVD). The MVD takes into account the endotoxin limit
1138 concentration (ELC) and the detection limit of the particular test method. The U.S. and
1139 European Pharmacopoeias assign ELCs for drugs based on their specific administered dose,
1140 route of administration, and dosing regimen. Based on the selected threshold pyrogen dose of
1141 0.5 EU/mL (see **Section 4.0**), and the decision criteria used in the validation studies to
1142 identify a pyrogenic response (≥ 0.5 EU/mL, see **Section 5.0**), a concentration of 0.5 EU/mL
1143 was used as the detection limit for the *in vitro* test methods when calculating the MVDs for
1144 each of the test substances.

1145 **3.2 Number of Substances**

1146 A total of 13 substances were included in the performance analysis of each of the five *in vitro*
1147 test methods. Ten substances, each spiked with five concentrations of endotoxin, were used
1148 to evaluate accuracy. For the evaluation of reproducibility, three substances were each spiked
1149 with four concentrations of endotoxin.

1150

1151 3.3 Identification and Description of Substances Tested

1152 As indicated in Section 3.1, the test substances selected for use in the validation studies were
 1153 marketed parenteral pharmaceuticals. Table 3-1 lists the 10 test substances used to evaluate
 1154 accuracy, and Table 3-2 lists the three test substances used to evaluate reproducibility.

1155 **Table 3-1 Test Substances (Parenteral Drugs) Used in the Validation Studies for**
 1156 **Determining Test Method Accuracy¹**

Test Substance ²	Source	Active Ingredient	Indication	MVD (-fold)
Beloc®	Astra Zeneca	Metoprolol tartrate	Heart dysfunction	140
Binotal®	Aventis	Ampicillin	Antibiotic	140
Ethanol 13% (w/w)	B. Braun	Ethanol	Diluent	35
Fenistil®	Novartis	Dimetindenmaleat	Antiallergic	175
Glucose 5% (w/v)	Eifel	Glucose	Nutrition	70
MCP®	Hexal	Metoclopramid	Antiemetic	350
Orasthin®	Aventis	Oxytocin	Initiation of delivery	700
Sostril®	GSK	Ranitidine	Antiacidic	140
Drug A - 0.9% NaCl	-	0.9% NaCl	-	35
Drug B - 0.9% NaCl	-	0.9% NaCl	-	70

1157 Abbreviations: MVD = Maximum valid dilution; GSK = GlaxoSmithKline; NaCl = Sodium chloride;
 1158 w/w = Weight/weight; w/v = Weight/volume

1159 ¹Each substance was tested in all five *in vitro* pyrogenicity test methods.

1160 ²Each test substance was spiked with 0, 0.25, 0.5, 0.5, or 1.0 EU/mL of endotoxin (WHO-LPS 94/580

1161 [*E. coli* O113:H10:K-]). Each sample contained the appropriate spike concentration when tested at its Maximum Valid
 1162 Dilution (MVD).

1163

1164 **Table 3-2 Test Substances (Parenteral Drugs) Used in the Validation Studies for**
 1165 **Determining Test Method Reproducibility¹**

Test Substance ²	Source	Agent	Indication
Gelafundin®	Braun Melsungen	Gelatin	Transfusion
Jonosteril®	Fresenius	Electrolytes	Infusion
Haemate®	Aventis	Factor VIII	Hemophilia

1166 ¹Each substance was tested in all five *in vitro* pyrogenicity test methods.

1167 ²Each test substance was spiked with 0, 0, 0.5, or 1.0 EU/mL of endotoxin (WHO-LPS 94/580 [*E. coli* O113:H10:K-]). Each
 1168 sample contained the appropriate spike concentration when tested at its Maximum Valid Dilution (MVD).

1169

1170 3.4 Sample Coding Procedure

1171 According to the ECVAM BRDs (Section 3.4), the ten test substances and the five spike
 1172 concentrations used for the evaluation of accuracy were blinded to the testing laboratories.

1173 For the reproducibility analyses, although the four spike concentrations were blinded to the
1174 participating laboratories, the identities of the three test substances were not.

1175 **3.5 Rationale for the Selection of the Recommended Reference Substances**

1176 Reference substances are used to assess the accuracy and reliability of a proposed,
1177 mechanistically and functionally similar test method and are a representative subset of those
1178 used to demonstrate the reliability and the accuracy of the validated reference test method (in
1179 this case, the RPT). These substances should:

- 1180 • represent the range of responses that the validated test method is capable of
1181 measuring or predicting
- 1182 • have produced consistent results in the validated test method
- 1183 • produce responses that reflect the accuracy of the validated test method
- 1184 • have well-defined chemical structures and/or compositions
- 1185 • be readily available
- 1186 • not be associated with excessive hazard or prohibitive disposal costs

1187 Pyrogenicity test methods provide a unique situation with regard to reference substances for
1188 use in validation studies because the only available standardized reference pyrogen is the
1189 international reference standard endotoxin (WHO-LPS 94/580 [*E. coli* 0113:h10:K-]). This
1190 reference standard was used to spike each of the test substances used in the ECVAM
1191 validation studies (see **Tables 3-2** and **3-3**). This reference standard is also used as a positive
1192 control and for qualifying the test methods during interference testing, and is also used when
1193 performing the BET.

1194 As described in **Section 4.0**, the response of the reference test method (i.e., RPT) to
1195 endotoxin is well documented. For this reason, the threshold pyrogen dose used for
1196 establishing the decision criteria for the *in vitro* test methods was based on historical RPT
1197 data for rabbits treated with endotoxin in the laboratory that supplied the RPT data.

1198

1198

1199

1200

1201

1202

1203

1204

1205

1206

1207

1208

1209

1210

1211 *[This Page Intentionally Left Blank]*

1212

1213

1214 **4.0** ***IN VIVO* REFERENCE DATA FOR THE ASSESSMENT OF TEST**
1215 **METHOD ACCURACY**

1216 **4.1** **Description of the Protocol Used to Generate *In Vivo* Data**

1217 4.1.1 The Rabbit Pyrogen Test (RPT)

1218 The RPT protocols most widely accepted by regulatory agencies are outlined in the USP,
1219 (USP23 NF28<151>), the U.S. Code of Federal Regulations (21CFR610.13), the European
1220 Pharmacopoeia (EP, EP5.0 2.6.8), and the Japanese Pharmacopoeia (JP, JP XIV), as
1221 summarized in **Table 4-1**. The RPT involves measuring the temperature increase in rabbits
1222 following an intravenous injection (via the ear vein) of a test substance in a dose not to
1223 exceed 10 mL/kg injected within a period of not more than 10 minutes. Initially, three rabbits
1224 are injected and the increase (or decrease) in temperature relative to the baseline value is
1225 measured at 30-min intervals for up to three hours. The resulting data are used to calculate an
1226 overall temperature increase by summing the results from all three animals, which is then
1227 used to assign a label of pyrogenic or non-pyrogenic.

1228 **Table 4-1 Test Guidelines for the Rabbit Pyrogen Test**

RPT Protocol Component	Reference			
	21CFR610.13	EP5.0 2.6.8	JP XIV	USP28 NF23 <151>
Number of rabbits	3 or 8 ¹	3, 6, 9, or 12 ¹	3 or 8 ¹	3 or 8 ¹
Rabbit species/strain	Not specified	Not specified	Not specified	Not specified
Exclusion criteria for rabbits during the initial selection of rabbits	<ul style="list-style-type: none"> Used in a negative pyrogen test in the preceding 2 days Used in a pyrogen test in which its temperature rose $\geq 0.6^{\circ}\text{C}$ in the preceding 2 weeks 	<ul style="list-style-type: none"> Weight < 1.5 kg Decreased weight in the preceding week Used in a negative pyrogen test in the preceding 3 days Used in a positive pyrogen test in the preceding 3 weeks 	<ul style="list-style-type: none"> Weight < 1.5 kg Decreased weight in the preceding week Previously used in a positive pyrogen test Rabbits from negative pyrogen tests may be reused only when a "as long a resting period as possible is taken" 	<ul style="list-style-type: none"> Used in a negative pyrogen test in the preceding 2 days Used in a pyrogen test in which its temperature rose $\geq 0.6^{\circ}\text{C}$ in the preceding 2 weeks
Testing room conditions	20 to 23°C	Within 3°C of the housing quarters (temperature not specified)	20 to 27°C and constant humidity	20 to 23°C
Food/water during test	Food withheld during the test, but water available at all times	Food withheld overnight and until end of the test. Water withheld during the test.	Food withheld beginning several hrs. prior to first temperature recording and until the end of the test.	Food withheld during the test period, but water available at all times
Depth of temperature probe in rectum	Not less than 7.5 cm	Approximately 5 cm	60-90 mm	Not less than 7.5 cm
Preliminary test	≤ 7 days prior to main test, perform all procedures used for the main test except the injection.	<ul style="list-style-type: none"> 1-3 days prior to main test, treat test animals with an injection of warmed (38.5°C) pyrogen-free saline Record temperature at 90 min prior to injection and every 30 min thereafter up to 3 hr. Exclude any rabbits with an increase of $> 0.6^{\circ}\text{C}$ 	Not specified	≤ 7 days prior to main test, perform all procedures used for the main test except the injection.
Baseline temperature	<ul style="list-style-type: none"> Record temperature ≤ 30 min prior to injection For any group of rabbits, use only if baseline temperatures do not vary $> 1^{\circ}\text{C}$ among rabbits Exclude rabbits with baseline temperature $> 39.8^{\circ}\text{C}$ 	<ul style="list-style-type: none"> Mean of two temperature recordings at 40 minutes and 10 minutes prior to injection Exclude rabbits if variation $> 0.2^{\circ}\text{C}$ between measurements noted Exclude rabbits with initial temperature $> 39.8^{\circ}\text{C}$ or $< 38.0^{\circ}\text{C}$ 	<ul style="list-style-type: none"> Record temperature three times at one-hr intervals prior to injection Assuming no appreciable variability among recordings, use the last recording as the baseline value. Exclude animals if 2nd and 3rd temperature measurements exceed 39.8°C 	<ul style="list-style-type: none"> Record temperature ≤ 30 min prior to injection For any group of rabbits, use only if baseline temperatures do not vary $> 1^{\circ}\text{C}$ among rabbits Exclude rabbits with baseline $> 39.8^{\circ}\text{C}$
Injection volume	≥ 3 mL/kg <i>BUT</i> ≤ 10 mL/kg	≥ 0.5 mL/kg <i>BUT</i> ≤ 10 mL/kg	10 mL/kg, unless otherwise specified	≤ 10 mL/kg

Injection time	≤10 min	≤4 min, unless otherwise indicated	Not specified, but injection should occur within 15 min of the third pretest temperature recording	≤10 min
Injection site	Marginal ear vein	Marginal ear vein	Marginal ear vein	Marginal ear vein
Pre-warming of test material	37°C±2°C	38.5°C	37°C	37°C±2°C
Temperature recording intervals after injection	30 min intervals for 1 to 3 hr	≤30 min intervals for 3 hr	1 hr intervals for 3 hr	30 min intervals for 1 to 3 hr

Abbreviations: CFR = U.S. Code of Federal Regulations; EP = European Pharmacopoeia; JP = Japanese Pharmacopoeia; RPT = Rabbit pyrogen test; USP = United States Pharmacopoeia

¹Each test is initially conducted with three animals and additional animals are tested to resolve equivocal results in the first three animals

1229
1230
1231

1232 4.1.2 Current In Vivo Pyrogenicity Test Method Protocols

1233 As indicated in **Table 4-1**, U.S. and international regulatory agencies have tailored the RPT
1234 test method protocol to suit their specific needs and goals in protecting human health. The
1235 current test method protocols (i.e., 21CFR610.13; EP5.0 2.6.8; JP XIV; USP28NF23<151>)
1236 recommend using healthy, adult rabbits with no specific breed/strain requirements. Rabbits
1237 are to be adequately acclimated to their surroundings and housed in an environment free from
1238 excessive external stimuli. Each rabbit is conditioned prior to the test with a sham test that
1239 includes all of the procedural steps except the injection. Reuse of test rabbits is permitted
1240 only after an appropriate withdrawal period has been completed (2-3 days for a negative
1241 result, 2-3 weeks for a positive result).

1242 The test is conducted in a room that is designated solely for pyrogen testing, in which the
1243 temperature is within 3°C of the uniform temperature of the housing room (i.e., 20°C±3°C).
1244 Food is withheld during the test, but access to water is continuous. The baseline temperature,
1245 which is used to calculate the increase in temperature during the test, is measured 30-40
1246 minutes prior to injection of the test substance. In each group of rabbits tested, the variation
1247 in baseline temperature among the rabbits should not vary more than 1°C, and rabbits with an
1248 initial temperature greater than 39.8°C are excluded from testing.

1249 The test substance is pre-warmed to approximately 37°C and injected (≤ 10 mL/kg) into the
1250 marginal ear vein, completing each injection within 10 minutes. The rectal temperature is
1251 recorded at 30-minute intervals for up to three hours after the injection. The decision criteria
1252 outlined in **Table 4-2** are then used to determine a pyrogenic response. As shown in **Table 4-**
1253 **2**, these decision criteria by which labels of pyrogenic or non-pyrogenic are assigned vary
1254 among the USP, CFR, EP, and JP test guidelines.

1255

1255 **Table 4-2 Decision Criteria for Determining a Pyrogenic Response in the Rabbit**
 1256 **Pyrogen Test**

RPT Protocol	No. Rabbits	Product passes if:	Product fails if:
USP28 NF23<151>	3	0/3 rabbits show an increase of $\geq 0.5^{\circ}\text{C}$	NA ¹
	5 ¹	$\leq 3/8$ rabbits show an increase of $\geq 0.5^{\circ}\text{C}$ <u>AND</u> the summed responses ≤ 3.3 .	$> 3/8$ rabbits show an increase of $\geq 0.5^{\circ}\text{C}$ <u>AND/OR</u> the sum of all responses > 3.3
21CFR610.13	3	0/3 rabbits show an increase of $\geq 0.5^{\circ}\text{C}$	NA ¹
	5 ¹	$\leq 3/8$ rabbits show an increase of $\geq 0.6^{\circ}\text{C}$ <u>AND</u> the summed responses ≤ 3.7 .	$> 3/8$ rabbits show an increase of $\geq 0.6^{\circ}\text{C}$ <u>AND/OR</u> the summed responses > 3.7
EP5.0 2.6.8	3	Summed responses $\leq 1.15^{\circ}\text{C}$	Summed responses $> 2.65^{\circ}$
	6 ²	Summed responses $\leq 2.80^{\circ}\text{C}$	Summed responses $> 4.30^{\circ}$
	9 ²	Summed responses $\leq 4.45^{\circ}\text{C}$	Summed responses $> 5.95^{\circ}$
	12	Summed responses $\leq 6.60^{\circ}\text{C}$	Summed responses $> 6.60^{\circ}$
JP XIV	3	3/3 rabbits show an increase of $< 0.6^{\circ}\text{C}$ <u>AND</u> the summed responses $\leq 1.4^{\circ}\text{C}$	$\geq 2/3$ rabbits show an increase $\geq 0.6^{\circ}\text{C}$
	5 ³	$\geq 4/5$ rabbits show an increase $< 0.6^{\circ}\text{C}$	$\geq 2/5$ rabbits show an increase $\geq 0.6^{\circ}\text{C}$

1257 CFR = U.S. Code of Federal Regulations; EP = European Pharmacopoeia; JP = Japanese Pharmacopoeia; USP = United
 1258 States Pharmacopoeia; NA = Not applicable

1259 ¹If $\geq 1/3$ rabbits show an increase of $\geq 0.5^{\circ}\text{C}$, continue test with an additional five rabbits.

1260 ²Three additional animals are tested when the summed responses falls in between the previous range.

1261 ³Five additional animals are tested when neither criteria is met, and results are based on these five animals only.

1262

1263 4.2 Reference Data Used to Assess *In Vitro* Test Method Accuracy

1264 There were no direct comparisons using the same test substances in the proposed *in vitro* test
 1265 methods and the RPT. Rather, historical data from 171 rabbits tested with endotoxin (0, 5,
 1266 10, 15, 20 EU/kg in 1 mL/kg) were obtained from a single laboratory. Two different sources
 1267 of endotoxin have been used: *E. coli* EC5; and *E. coli* EC6 that were reported to be identical
 1268 to the WHO standard used in the validation studies). These historical data were used to
 1269 establish a threshold pyrogen dose (i.e., the endotoxin dose at which fever was induced in
 1270 50% of the rabbits), which was determined to be 5 EU/kg. Based on the largest allowable
 1271 volume for injection in rabbits (10 mL/kg), the limit of detection that new pyrogen tests must
 1272 meet was defined as 0.5 EU/mL.

1273 **4.3 Availability of Original Records for the *In Vivo* Reference Data**

1274 **Section 4.1** of the ECVAM BRDs indicates that the Paul-Ehrlich Institute, Germany,
1275 provided the historical RPT data.

1276 **4.4 *In Vivo* Data Quality**

1277 Ideally, all data supporting the validity of a test method should be obtained and reported from
1278 studies conducted in accordance with GLP guidelines (OECD 1998; EPA 2003a, 2003b;
1279 FDA 2003). These guidelines provide an internationally standardized approach for the
1280 reporting requirements of studies designed for regulatory submissions, internal audits of
1281 laboratory records and data summaries, the archive of study data and records, and
1282 information about the test protocol and laboratory personnel, to provide assurances regarding
1283 the integrity, reliability, and accountability of the study. According to Section 4.4 of the
1284 ECVAM BRDs (with the exception of the WB/IL-6 BRD, which states that this section is
1285 "not applicable"), all RPT procedures were GLP-compliant.

1286 **4.5 Availability and Use of Toxicity Information from the Species of Interest**

1287 As stated in **Section 1.2.1**, the major regulatory requirement for pyrogenicity testing is for
1288 end-product release of human and animal parenteral drugs, medical devices, and human
1289 biological products. Results from such testing are used to limit to an acceptable level the
1290 risks of febrile reactions to the injection and/or implantation of the product. Therefore, for
1291 protection of both human and veterinary health, it is vital that the test method employed
1292 provide an accurate estimation of the potential for a pyrogenic reaction in humans and other
1293 species.

1294 The estimation of the potential for a pyrogenic reaction has traditionally been obtained in
1295 rabbits using the RPT, and more recently in horseshoe crabs with the BET (for Gram-
1296 negative endotoxins only). A number of studies have compared febrile responses between
1297 rabbits and humans and showed that rabbits tended to be equal to, or more sensitive than
1298 humans with respect to their threshold pyrogenic responses, but higher doses were more
1299 pyrogenic in humans (Co Tui and Schrifft 1942; Westphal 1956; Keene et al. 1961). Greisman
1300 and Hornick (1969) compared three purified endotoxin preparations in mature New Zealand
1301 white rabbits and in male volunteers and showed that the pyrogenic response induced was

1302 similar in both species. Based on these studies, the rabbit is considered to be predictive of the
1303 human response (and may often over predict the response).

1304 **4.6 Information on the Accuracy and Reliability of the *In Vivo* Test Method**

1305 Hoffmann et al. (2005) modeled the sensitivity and specificity of the RPT. They reviewed
1306 historical data from 171 rabbits challenged with endotoxin in a single laboratory in order to
1307 establish a threshold pyrogen dose (i.e., the endotoxin dose at which fever was induced in
1308 50% of the rabbits). A threshold value of 0.5 EU/mL was defined by regression analysis of
1309 the data. The performance characteristics of the RPT (i.e., sensitivity and specificity) were
1310 then determined using a 2 x 2 contingency table, and incorporating the parameters obtained
1311 from the regression analysis. The authors considered the prevalence of the endotoxin spikes
1312 included in the ECVAM accuracy evaluations in the validation studies (i.e., 0 EU/mL: 20%;
1313 0.25 EU/mL: 20%; 0.5 EU/mL: 40%; 1.0 EU/mL: 20%) and applied the threshold pyrogen
1314 dose of 0.5 EU/mL to calculate theoretical values for sensitivity (58%) and specificity (83%)
1315 of the RPT.

1316 The accuracy and reliability of the RPT for endotoxin testing has been considered adequate
1317 for U.S. and international regulatory needs for many years. Since its inclusion in the USP in
1318 1941, the RPT had been extensively used, and is the preferred method for detection of
1319 pyrogenicity for product development, because of the inability of the BET to detect non-
1320 endotoxin pyrogens.

1321 **5.0 TEST METHOD DATA AND RESULTS**

1322 **5.1 Test Method Protocol**

1323 The standard operating procedures (SOPs) used during the ECVAM validation studies are
1324 included in **Appendix A**. As described in **Section 2.1**, there are many similarities among the
1325 protocols for each of the *in vitro* pyrogenicity test methods, with very few notable differences
1326 other than the type of cells used (i.e., whole blood cells, PBMCs, monocytoid cell line) and
1327 the proinflammatory cytokine assayed (i.e., IL-1 β or IL-6). These similarities and differences
1328 are outlined in **Table 2-1**. An internationally accepted endotoxin standard (i.e., WHO-LPS
1329 94/580 [*E. coli* 0113:h10:K-]) was used to spike samples of saline or marketed parenteral
1330 pharmaceuticals. The same pharmaceuticals were used to create the spiked samples for all
1331 five test methods (see **Table 3-1** and **3-2**). These samples were included in a series of studies
1332 designed to determine the relevance and reliability of each of the *in vitro* test methods.

1333 **5.2 Availability of Copies of Original Data Used to Evaluate Test Method**
1334 **Performance**

1335 ECVAM provided raw data from the validation studies in an electronic format (Excel[®]
1336 spreadsheets) that consisted of OD₄₅₀ measurements for all replicates included in each of the
1337 validation studies.

1338 NICEATM attempted to obtain additional *in vitro* and/or *in vivo* pyrogenicity test method
1339 data. A *Federal Register (FR)* notice (Vol. 70, No. 241, pp. 74833-4, December 16, 2005;
1340 available at <http://iccvam.niehs.nih.gov/methods/pyrogen.htm>) was published requesting
1341 original *in vitro* pyrogenicity test method data and reference data from the currently used
1342 pyrogenicity test methods (i.e., RPT and/or BET). In addition, the *FR* notice was sent directly
1343 to more than 100 interested stakeholders internationally. Despite these efforts, no additional
1344 data were submitted.

1345 **5.3 Description of the Statistical Approaches Used to Evaluate the Resulting Data**

1346 Details of the statistical approaches used to evaluate the accuracy and reliability of each of
1347 the five test methods are included in Section 5.3 of each ECVAM BRD. Briefly, as indicated
1348 in **Section 3.2**, ten substances (each spiked with five concentrations of endotoxin) were

1349 tested in each test method to evaluate accuracy, while three substances (each spiked with four
1350 concentrations of endotoxin) were used to evaluate test method reproducibility. Varying
1351 concentrations of endotoxin-spiked saline were tested for the analysis of intralaboratory
1352 repeatability.

1353 The evaluation of intralaboratory repeatability included a coefficient of variation (CV)
1354 analysis of the OD₄₅₀ measurements for the replicates of each endotoxin concentration.
1355 Boxplots were also generated to demonstrate variability among OD₄₅₀ measurements for each
1356 concentration. Similar analyses were conducted for the three substances used to assess intra-
1357 and interlaboratory reproducibility.

1358 The reproducibility analysis procedures used incorporated the decision criteria that were
1359 developed to differentiate between pyrogenic and non-pyrogenic materials (using a threshold
1360 value of 0.5 EU/mL). In all reproducibility analyses, a single run consisted of each of the
1361 substances (as described above and in **Section 3.2**) assayed in quadruplicate. Acceptability
1362 criteria for each run included a CV analysis to remove highly variable samples from the
1363 analyses. This criterion ranged from CV<0.25 to CV<0.45, depending on the test method
1364 being considered. For the measurement of intralaboratory reproducibility, pair-wise
1365 comparisons between the runs were determined and the correlations between runs expressed
1366 as a percentage of agreement between two individual laboratories. It should be noted that this
1367 analysis takes into account the agreement of the resulting pyrogenicity call (i.e., pyrogenic or
1368 non-pyrogenic), but does not consider whether the call is correct. All three possible
1369 combinations were compared (i.e., run 1 vs. run 2; run 1 vs. run 3; run 2 vs. run 3) and a
1370 mean value calculated, which was intended to provide an overall proportion of inter-run
1371 agreement.

1372 Similar analyses were conducted for an assessment of interlaboratory reproducibility, in
1373 which pairwise comparisons between laboratories were determined and the correlations were
1374 expressed as a percentage of agreement. This analysis included each run from each
1375 laboratory (n=3 per laboratory) and all possible interlaboratory combinations were compared.
1376 Similar to the intralaboratory analysis, this analysis takes the resulting pyrogenicity call from
1377 each run in each laboratory into consideration, but does not consider whether the call is
1378 correct. **Section 7.0** provides additional details and the resulting data from these analyses.

1379 For the accuracy analysis, 2 x 2 contingency tables were constructed using the decision
1380 criteria defined in **Table 4-2** to assign a pyrogenicity call. Each run for each sample from
1381 each laboratory was considered independently. Accordingly, the *in vitro* call was compared
1382 to the "true status" (based on the known endotoxin spike concentration) of the sample. The
1383 resulting accuracy statistics were calculated based on the overall database for each test
1384 method. Similar to the reproducibility analyses, acceptability criteria for each run included a
1385 CV analysis to remove highly variable samples from the analyses, for which a range of
1386 CV<0.25 to CV<0.45 was used, depending on the test method being considered. **Section 6.0**
1387 provides additional details and the resulting data from these analyses.

1388 **5.4 Summary of Results**

1389 Graphical representations of the repeatability and reproducibility analyses are provided in
1390 Section 5.2 of each ECVAM BRD (see **Appendix A**). The tabulated results from which the
1391 intra- and interlaboratory reproducibility analyses and accuracy analyses can be conducted
1392 are provided in Section 5.4 of the ECVAM BRDs. The tables in that section include the test
1393 substance name, the endotoxin spike concentration, the pyrogenicity call for each *in vitro*
1394 run, and the "true status" of each test substance.

1395 **5.5 Use of Coded Chemicals and Compliance with GLP Guidelines**

1396 Ideally, all data supporting the validity of a test method should be obtained using coded
1397 chemicals and reported in accordance with GLP guidelines (OECD 1998; EPA 2003a,
1398 2003b; FDA 2003). As described in **Section 8.1**, all studies were carried out in accordance
1399 with GLP guidelines based on the information available in the ECVAM BRDs (**Appendix**
1400 **A**), **Section 3.4** indicates that the ten test substances and the five spike concentrations used
1401 for the accuracy evaluation were blinded to the testing laboratories. However, although the
1402 four spike concentrations were blinded to the participating laboratories for the reproducibility
1403 studies, the identity of the three test substances was not.

1404 **5.6 Lot-to-Lot Consistency of Test Substances**

1405 Lot-to-lot consistency of test substances is evaluated to ensure that the same substance, with
1406 the same physicochemical properties, is being evaluated over the duration of the study. In
1407 these studies, the test substances used were from released clinical lots of parenteral

1408 pharmaceuticals, which would imply that they had been subjected to rigorous chemical
1409 manufacturing control analyses to verify that the compositions are consistent. In addition, the
1410 international standard for Gram-negative endotoxin, WHO-LPS 94/580 (*E. coli*
1411 O113:H10:K-), is used as the spike solution, which provides a measure of consistency for the
1412 positive control substance, and the spike substance.

1413 **5.7 Availability of Data for External Audit**

1414 As described in **Section 8.4**, all records are stored and archived by the participating
1415 laboratories and are available for inspection.

1416 **6.0 RELEVANCE OF THE *IN VITRO* PYROGENICITY TEST METHODS**

1417 **6.1 Accuracy of the *In Vitro* Pyrogenicity Test Methods**

1418 A critical component of an ICCVAM evaluation of the validation status of a method is an
1419 assessment of its relevance. The measure of relevance used in this evaluation is the
1420 performance of the new test in identifying pyrogens as compared to the performance of the
1421 current reference method (ICCVAM 2003). This aspect of assay performance is typically
1422 evaluated by calculating:

- 1423 • Concordance (also referred to as Accuracy): the proportion of correct
1424 outcomes (positive and negative) of a test method
- 1425 • Sensitivity: the proportion of true positive substances that are correctly
1426 classified as positive
- 1427 • Specificity: the proportion of true negative substances that are correctly
1428 classified as negative
- 1429 • Positive predictivity: the proportion of correct positive responses among
1430 substances testing positive
- 1431 • Negative predictivity: the proportion of correct negative responses among
1432 substances testing negative
- 1433 • False positive rate: the proportion of true negative substances that are falsely
1434 identified as positive
- 1435 • False negative rate: the proportion of true positive substances that are falsely
1436 identified as negative.

1437 The ability of the *in vitro* pyrogenicity test methods to correctly identify the presence of
1438 Gram-negative endotoxin was evaluated using parenteral pharmaceuticals spiked endotoxin
1439 (WHO-LPS 94/580 [*E. coli* O113:H10:K-]). As described in **Section 3.2**, ten substances (see
1440 **Table 3-1**) spiked with five concentrations of endotoxin were used for the evaluation. The
1441 individual spike concentrations in each substance were tested once, using each test method,
1442 in three different laboratories, providing a total of 150 runs (i.e., 10 substances x 5 spike
1443 solutions x 3 laboratories = 150) for evaluation in each test method. The quality criteria

1444 outlined in **Table 2-1** were used to identify run outliers that were subsequently excluded
1445 from the evaluation resulting in fewer than 150 runs included in the evaluation.

1446 As described in **Section 4.2**, no RPT assays were conducted in parallel with the *in vitro*
1447 pyrogenicity test methods during the ECVAM validation studies. Instead, historical RPT data
1448 from rabbits tested with endotoxin were used to establish a threshold pyrogen dose (i.e., the
1449 endotoxin dose at which fever was induced in 50% of the rabbits), which was subsequently
1450 used to establish the limit of detection (i.e., 0.5 EU/mL) that the *in vitro* test methods being
1451 validated must meet. Accordingly, the *in vitro* call was compared to the "true status" (based
1452 on the known endotoxin spike concentration) of the sample. The resulting calls were used to
1453 construct 2x2 contingency tables, which were used to calculate the resulting test performance
1454 values.

1455 6.1.1 Relevance of the Cryo WB/IL-1 Test Method

1456 Of the 150 available runs for the cryo WB/IL-1 test method, 10 showed excessive variability
1457 among the four replicates (i.e., CV >45%), resulting in their exclusion from the analysis. An
1458 additional 20 runs did not qualify according to one or more of the criteria outlined in **Table**
1459 **2-1**. Therefore a total of 120 runs were used in the performance analysis which showed that
1460 the cryo WB/IL-1 test method has a concordance of 92% (110/120), a sensitivity of 97%
1461 (75/77), a specificity of 81% (35/43), a false negative rate of 3% (2/77), and a false positive
1462 rate of 19% (8/43) (see **Table 6-1**).

1463

1463 **Table 6-1 Performance Analysis for Five *In Vitro* Pyrogenicity Test Methods¹**

Test Method	Concordance ²	Sensitivity	Specificity	False Negative Rate	False Positive Rate
cryo WB/IL-1	91.7% (110/120)	97.4% (75/77)	81.4% (35/43)	2.6% (2/77)	18.6% (8/43)
MM6/IL-6	93.2% (138/148)	95.5% (85/89)	89.8% (53/59)	4.5% (4/89)	10.2% (6/59)
PBMC/IL-6	93.3% (140/150)	92.2% (83/90)	95.0% (57/60)	7.8% (7/90)	5.0% (3/60)
PBMC/IL-6 (cryo) ³	91.9% (136/148)	88.8% (79/89)	96.6% (57/59)	11.2% (10/89)	3.4% (2/59)
WB/IL-6	91.9% (136/148)	88.8% (79/89)	96.6% (57/59)	11.2% (10/89)	3.4% (2/59)
WB/IL-1	81.0% (119/147)	72.7% (64/88)	93.2% (55/59)	27.3% (24/88)	6.8% (4/59)
WB/IL-1 (96-well plate method) ⁴	92.8% (129/139)	98.8% (83/84)	83.6% (46/55)	1.2% (1/84)	16.4% (9/55)

1464 Abbreviations: cryo = Cryopreserved; IL-1 = Interleukin-1; IL-6 = Interleukin -6; MM6 = Mono Mac 6; PBMC = Peripheral
 1465 blood mononuclear cells; WB = Whole blood

1466 ¹Based on results of 10 parenteral drugs tested in each of three different laboratories; samples of each drug were tested with
 1467 or without being spiked with a Gram-negative endotoxin standard (0, 0.25, 0.5, 0.5, or 1.0 EU/mL).

1468 ²Percentage (Number of correct runs/total number of runs)

1469 ³A modification of the PBMC/IL-6 test method using cryopreserved PBMCs.

1470 ⁴A modification of the WB/IL-1 test method using 96-well plates instead of tubes for the test substance incubation.

1471

1472 6.1.2 Relevance of the MM6/IL-6 Test Method

1473 Of the 150 available runs for the MM6/IL-6 test method, two showed excessive variability
 1474 among the four replicates (i.e., CV >25%), resulting in their exclusion from the analysis. No
 1475 runs were excluded based on the criteria outlined in **Table 2-1**. Therefore a total of 148 runs
 1476 were used in the performance analysis. Based on this analysis, the MM6/IL-6 test method has
 1477 a concordance of 93% (138/148), a sensitivity of 96% (85/89), a specificity of 90% (53/59), a
 1478 false negative rate of 4% (4/89), and a false positive rate of 10% (6/59) (see **Table 6-1**).

1479 6.1.3 Relevance of the PBMC/IL-6 Test Method

1480 None of the 150 available runs for the PBMC/IL-6 test method showed excessive variability
 1481 (i.e., CV >40%) and all runs met the criteria outlined in **Table 2-1**. Therefore all 150 runs
 1482 were included in the performance analysis. Based on this analysis, the PBMC/IL-6 test
 1483 method has a concordance of 93% (140/150), a sensitivity of 92% (83/90), a specificity of
 1484 95% (57/60), a false negative rate of 8% (7/90), and a false positive rate of 5% (3/60) (see
 1485 **Table 6-1**).

1486 6.1.3.1 *Relevance of the PBMC/IL-6 Method When Using Cryo PBMCs*

1487 As indicated in **Table 2-1**, the PBMC/IL-6 test method protocol was also conducted using a
1488 modified protocol that included cryo PBMCs. None of the 150 available runs for this
1489 modification of the PBMC/IL-6 test method showed excessive variability (i.e., CV >40%)
1490 and all runs met the criteria outlined in **Table 2-1**. Therefore all runs were included in a
1491 performance analysis. Based on this analysis, the PBMC/IL-6 test method, when using cryo
1492 PBMCs, has a concordance of 87% (130/150), a sensitivity of 93% (84/90), a specificity of
1493 77% (46/60), a false negative rate of 7% (6/90), and a false positive rate of 23% (14/60). The
1494 high false positive rate can be attributed to a large number of false positives (50% [10/20]) in
1495 one of the three laboratories (the false positive rate in the remaining two laboratories is 10%).

1496 6.1.4 Relevance of the WB/IL-6 Test Method

1497 None of the 150 available runs for the WB/IL-6 test method showed excessive variability
1498 (i.e., CV >45%) and all runs met the criteria outlined in **Table 2-1**. However, two samples
1499 were mishandled by one of the testing laboratories, and thus the two associated runs were
1500 excluded from the analysis. As a result, 148 runs were included in the performance analysis
1501 for the detection of Gram-negative endotoxin. Based on this analysis, the WB/IL-6 test
1502 method has a concordance of 92% (136/148), a sensitivity of 89% (79/89), a specificity of
1503 97% (57/59), a false negative rate of 11% (10/89), and a false positive rate of 3% (2/59) (see
1504 **Table 6-1**).

1505 6.1.5 Relevance of the WB/IL-1 Test Method

1506 Of the 150 available runs for the WB/IL-1 test method, three showed excessive variability
1507 among the four replicates (i.e., CV >45%), resulting in their exclusion from the analysis. No
1508 runs were excluded based on the criteria outlined in **Table 2-1**. Therefore a total of 147 runs
1509 were used in the performance analysis. Based on this analysis, the WB/IL-1 test method has a
1510 concordance of 81% (119/147), a sensitivity of 73% (64/88), a specificity of 93% (55/59), a
1511 false negative rate of 27% (24/88), and a false positive rate of 7% (4/59) (see **Table 6-1**).

1512 6.1.5.1 *Relevance of the WB/IL-1 Test Method When Using 96-Well Plates*

1513 As indicated in **Table 2-1**, the WB/IL-1 test method protocol was also conducted using a
1514 modified protocol that used 96-well plates instead of individual tubes. Of the 150 available

1515 runs for this modification of the WB/IL-1 test method, 11 showed excessive variability (i.e.,
1516 CV >45%). No runs were excluded based on the criteria outlined in **Table 2-1**. Therefore, a
1517 total of 139 runs were included in a performance analysis. Based on this analysis, the WB/IL-
1518 1 test method, when using 96-well plates, has a concordance of 93% (129/139), a sensitivity
1519 of 99% (83/84), a specificity of 84% (46/55), a false negative rate of 1% (1/84), and a false
1520 positive rate of 16% (9/55).

1521 **6.2 Summary of the Performance Statistics for *In Vitro* Pyrogenicity Test**

1522 **Methods**

1523 The performance of five *in vitro* pyrogenicity test methods for the detection of Gram-
1524 negative endotoxin (based on 10 parenteral pharmaceuticals, each spiked with five
1525 concentrations of endotoxin) was evaluated. As outlined in **Table 6-1**, this analysis indicated
1526 that concordance among the test methods ranged from 81% to 93%, sensitivity ranged from
1527 89% to 97%, specificity ranged from 81% to 97%, false negative rates ranged from 3% to
1528 27%, and false positive rates ranged from 3% to 19%¹⁰.

1529 **6.2.1 Discordant Results**

1530 It was not possible to make a direct comparison between the RPT and *in vitro* pyrogenicity
1531 test results without the availability of parallel testing data (i.e., same test substance tested
1532 using the *in vitro* and *in vivo* methods). Therefore, *in vitro* results that are discordant from the
1533 RPT could not be identified with these studies. Discordant results reflect a failure of the *in*
1534 *vitro* test method(s) to identify Gram-negative endotoxin spiked into a test substance at the
1535 threshold concentration (0.5 EU/mL) established based on historical data from the RPT (see
1536 **Section 4.2**).

1537 **6.2.2 Strengths and Limitations of *In Vitro* Pyrogenicity Test Methods**

1538 The limitations of these test methods have not been fully explored and identified. For this
1539 reason, pre-testing product specific validation will be necessary to establish if a particular test
1540 substance/material is appropriate for evaluation using these *in vitro* test methods. A

¹⁰ Including the cryopreservation modification in the PBMC/IL-6 test method protocol resulted in a false positive rate of 23% and a false negative rate of 7%. Including the 96-well plate modification in the WB/IL-1 test method resulted in a false positive rate of 16% and a false negative rate of 1%.

1541 recognized limitation of the *in vitro* methods is the lack of data to determine their responses
1542 to, and suitability for, non-endotoxin pyrogens that are known to be detected by the RPT.
1543 However, an advantage to these *in vitro* test methods is that they are derived from human
1544 tissues, and thus avoid potential uncertainty associated with cross-species extrapolation.

1545

1546

1547

1548

1549

1550

1551

1552

1553

1554

1555

1556

1557

1558

1559 **7.0 RELIABILITY OF THE *IN VITRO* PYROGENICITY TEST METHODS**

1560 An assessment of test method reliability (intralaboratory repeatability and intra- and
1561 interlaboratory reproducibility) is an essential element of any evaluation of the performance
1562 of an alternative test method (ICCVAM 2003). Repeatability refers to the closeness of
1563 agreement among test results obtained within a single laboratory when the procedure is
1564 performed on the same substance under identical conditions within a given time period
1565 (ICCVAM 1997, 2003). Intralaboratory reproducibility refers to the determination of the
1566 extent to which qualified personnel within the same laboratory can replicate results using a
1567 specific test protocol at different times. Interlaboratory reproducibility refers to the
1568 determination of the extent to which different laboratories can replicate results using the
1569 same protocol and test chemicals, and indicates the extent to which a test method can be
1570 transferred successfully among laboratories. A reliability assessment includes a quantitative
1571 and/or qualitative analysis of repeatability and intra- and interlaboratory reproducibility. In
1572 addition, measures of central tendency and variation are summarized for historical control
1573 data (negative, vehicle, positive), where applicable.

1574 An evaluation of intralaboratory repeatability and reproducibility could be conducted because
1575 *in vitro* pyrogenicity test data were available from replicate wells within individual
1576 experiments, and from replicate experiments within the individual laboratories. In addition,
1577 comparable data were available from each of the three laboratories that performed the
1578 validation studies, which allowed an evaluation of interlaboratory reproducibility.

1579 **7.1 Selection Rationale for the Substances Used to Evaluate the Reliability of *In***
1580 ***Vitro* Pyrogenicity Test Methods**

1581 The quality of a reliability evaluation depends on the extent to which the substances tested
1582 adequately represent the range of physicochemical characteristics and response levels that the
1583 test method should be capable of evaluating.

1584 The rationale for selecting the substances used in the validation studies was discussed in
1585 **Section 3.0**. In brief, substances that were used in the ECVAM validation studies were
1586 marketed parenteral pharmaceuticals, labeled as free from detectable pyrogens (i.e., they had
1587 passed the RPT or BET test). No rationale for the selection of these specific test substances

1588 was provided. Each test substance was spiked with a series of concentrations of Gram-
1589 negative endotoxin standard (WHO-LPS 94/580 [*E. coli* O113:H10:K-]). Endotoxin was
1590 selected as a “model” pyrogen for inclusion based on its availability in a standardized form.
1591 No other pyrogenic substances are presently available in a standardized form.

1592 7.2 Analysis of Intralaboratory Repeatability and Reproducibility

1593 Intralaboratory repeatability analyses were performed using the OD values obtained for each
1594 test with each spiked sample. All analyses of intra- and interlaboratory reproducibility were
1595 performed on the classifications of pyrogenic or non-pyrogenic, rather than on the absolute
1596 OD values generated in each run. Analyses of intralaboratory reliability typically include a
1597 coefficient of variation (CV) analysis, which is a statistical measure of the deviation of a
1598 variable from its mean (e.g., Holzhütter et al. 1996).

1599 7.2.1 Intralaboratory Repeatability

1600 In the ECVAM validation study, intralaboratory repeatability of each test method was
1601 evaluated by testing saline and various endotoxin spikes (0.06 to 0.5 EU/mL) in saline and
1602 evaluating the closeness of agreement among optical density readings for cytokine
1603 measurements at each concentration. Each experiment was conducted up to three times for
1604 each test method. Up to 20 replicates per concentration were tested and results indicated that
1605 variability in OD measurements increased with increasing endotoxin concentration, but the
1606 variability was not so great to interfere with distinguishing the 0.5 EU/mL spike
1607 concentration (i.e., the threshold for pyrogenicity) from the lower concentrations. **Table 7-1**
1608 details the study designs for each of these evaluations. At least four different study designs
1609 were employed for each test method with the exception of the cryo WB/IL-1. In the ECVAM
1610 cryo WB/IL-1 BRD (see **Appendix A**) Appendix D indicates that because intralaboratory
1611 reliability was extensively evaluated in the WB/IL-1, only a subset (n=2) of these studies was
1612 conducted as a part of a "catch-up validation" study. Based on the "acceptable"
1613 intralaboratory performance in this subset of studies, additional studies were not considered
1614 necessary.

1615

1616 **Table 7-1 Intralaboratory Repeatability Assessed with Saline Spiked with WHO-**
 1617 **LPS 94/580**

Experiment	Study Design	Test Method				
		MM6/IL-6	PBMC/IL-6	WB/IL-1	WB/IL-6	cryo WB/IL-1 ¹
1A	Endotoxin concentration (EU/mL)	0, 0.25, 0.5	0, 0.25, 0.5	0, 0.5	0, 0.5	0, 0.5
	N (per spike)	20	20	32	20	32
	Repetitions of experiment	1	1	1	1	1
1B	Endotoxin concentration (EU/mL)	0, 0.063, 0.125, 0.25, 0.5	0, 0.063, 0.125, 0.25, 0.5	0, 0.063, 0.125, 0.25, 0.5	0, 0.063, 0.125, 0.25, 0.5	0, 0.063, 0.125, 0.25, 0.5
	N (per spike)	12	12	12	10	12
	Repetitions of experiment	1	1	1	1	1
2A	Endotoxin concentration (EU/mL)	0, 0.25, 0.5	0, 0.5	0, 0.5	0, 0.25, 0.5	ND
	N (per spike)	20	8	12	8	ND
	Repetitions of experiment	3	3	3	3	ND
2B	Endotoxin concentration (EU/mL)	0, 0.25, 0.5	0, 0.063, 0.125, 0.25, 0.5	0, 0.25, 0.5	0, 0.5	ND
	N (per spike)	20	8	8	5	ND
	Repetitions of experiment	3	3	3	8	ND
2C	Endotoxin concentration (EU/mL)	ND	0, 0.125, 0.25, 0.5	0, 0.5	ND	ND
	N (per spike)	ND	8	5	ND	ND
	Repetitions of experiment	ND	8	8	ND	ND

1618 Abbreviations: cryo = Cryopreserved; EU = Endotoxin unit; IL-1 = Interleukin-1; IL-6 = Interleukin-6; LPS =
 1619 Lipopolysaccharide; MM6 = Mono Mac 6; N = number of replicates; ND = Not done; PBMC = Peripheral
 1620 blood mononuclear cells; WB = Whole blood; WHO = World Health Organization

1621 ¹Assessed using saline spiked with WHO-LPS 94/580 endotoxin.

1622 ²The cryo WB/IL-1 test method was included in a catch-up validation study to assess intralaboratory reliability
 1623 in a subset of experiments (n=2).

1624

1625 7.2.2 Intralaboratory Reproducibility

1626 Intralaboratory reproducibility was evaluated using three marketed pharmaceuticals spiked
 1627 with various concentrations of endotoxin (see **Table 3-2**). Three identical, independent runs
 1628 conducted in each of the three testing laboratories, with the exception of the cryo WB/IL-1

1629 test method¹¹. All three possible combinations were compared (i.e., run 1 vs. run 2; run 1 vs.
1630 run 3; run 2 vs. run 3) and a mean value calculated, intended to provide an overall proportion
1631 of inter-run agreement. In all reproducibility analyses, a single run consisted of each of the
1632 substances assayed in quadruplicate. Acceptability criteria for each run included a CV
1633 analysis to remove highly variable responses from the analyses. The criterion used to identify
1634 outliers ranged from CV <0.25 to CV <0.45, depending on the method being considered, and
1635 was arbitrarily set based on results using saline spiked with endotoxin. For example, for the
1636 MM6/IL-6 test method, the CV for any single spike concentration was ≤ 0.12 , and therefore,
1637 the outlier criterion was set at 0.25.

1638 Agreement between different runs was determined for each substance in three laboratories.
1639 As shown in **Table 7-2**, the agreement across three runs in an individual lab ranged from
1640 75% to 100%.

¹¹ The ECVAM cryo WB/IL-1 test method BRD states that there was no direct assessment of intralaboratory reproducibility because such an evaluation was performed in the WB IL-1 (fresh blood) test method, and the authors assume that variability is not affected by the change to cryopreserved blood assayed in 96-well plates.

1641 **Table 7-2 Intralaboratory Reproducibility of *In Vitro* Pyrogenicity Test Methods**

Run Comparison ¹	WB/IL-1			Cryo WB/IL-1			WB/IL-6			PBMC/IL-6			MM6/IL-6		
	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3	Lab 1	Lab 2	Lab 3
1 vs 2	92% (11/12)	100% (8/8)	100% (12/12)	ND ³	ND	ND	75% (9/12)	92% (11/12)	100% (12/12)	92% (11/12)	100% (12/12)	100% (12/12)	100% (12/12)	92% (11/12)	100% (12/12)
1 vs 3	83% (10/12)	88% (7/8)	92% (11/12)	ND	ND	ND	100% (12/12)	92% (11/12)	100% (12/12)	100% (12/12)	100% (12/12)	92% (11/12)	100% (12/12)	92% (11/12)	92% (11/12)
2 vs 3	92% (11/12)	NA ⁴	92% (11/12)	ND	ND	ND	75% (9/12)	92% (11/12)	100% (12/12)	92% (11/12)	100% (12/12)	92% (11/12)	100% (12/12)	100% (12/12)	92% (11/12)
Mean	89%	-	95%	ND	ND	ND	83%	92%	100%	95%	100%	95%	100%	95%	95%
Agreement ² across 3 runs	83%	-	92%	ND	ND	ND	75%	92%	100%	92%	100%	94%	100%	92%	92%

1642 Abbreviations: Cryo = Cryopreserved; IL-1 = Interleukin-1; IL-6 = Interleukin-6; MM6 = Mono Mac 6; NA = Not assessed; ND = Not done; WB = Whole blood

1643 ¹Comparison between 3 individual runs within each laboratory

1644 ²All possible combinations of runs among the 3 laboratories were compared

1645 ³Not done. The cryo WB/IL-1 BRD states that an assessment of intralaboratory reproducibility was performed using the WB IL-1 (fresh blood) test method, and it was assumed
1646 that intralaboratory variability would not be affected by the change to cryopreserved blood assayed in 96-well plates.

1647 ⁴Not assessed due to lack of sufficient data. The sensitivity criteria were not met for 1/3 substance in run 2, and 1/3 substance in run 3.

1648 7.2.3 Interlaboratory Reproducibility

1649 Interlaboratory reproducibility was evaluated in two different studies. In both studies, each
 1650 run from one laboratory was compared with all runs of another laboratory. The proportions
 1651 of similarly classified samples provide a measure of reproducibility. In the first study, in a
 1652 similar manner as the evaluation of intralaboratory reproducibility, the interlaboratory
 1653 reproducibility was evaluated using results from three marketed pharmaceuticals spiked with
 1654 endotoxin and tested three times in each of the three laboratories. As shown in **Table 7-3**, the
 1655 agreement across three laboratories for each test method (where three runs per laboratory
 1656 were conducted) ranged from 58% to 86%, depending on the test method considered. In
 1657 comparison, the agreement across three laboratories for the cryo WB/IL-1 test method, for
 1658 which only one run per laboratory was conducted, was 92%.

1659

1660 **Table 7-3 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity Test Methods**

Lab Comparison ¹	Agreement Between Laboratories ¹				
	WB/IL-1	Cryo WB/IL-1	WB/IL-6	PBMC/IL-6	MM6/IL-6
1 vs 2	92% (77/84) ²	92% (11/12) ³	72% (78/108)	81% (87/108)	97% (105/108)
1 vs 3	77% (83/108)	92% (11/12) ³	75% (81/108)	86% (93/108)	89% (96/108)
2 vs 3	68% (57/84) ²	92% (11/12) ³	97% (105/108)	89% (96/108)	86% (93/108)
Mean	79%	92%	81%	85%	90%
Agreement across 3 labs ⁴	58% (167/288) ²	92% (11/12) ³	72% (234/324)	78% (252/324)	86% (279/324)

1661 Abbreviations: Cryo = Cryopreserved; IL-1 = Interleukin-1; IL-6 = Interleukin-6; MM6 = Mono Mac 6; WB = Whole blood

1662 ¹Data from three substances (see **Table 3-2**) spiked with endotoxin (WHO-LPS 94/580 [*E. coli* O113:H10:K-]) at 0, 0, 0.5
 1663 and 1.0 EU/mL tested three times in three different laboratories, with the exception of cryo WB/IL-1 (only the preliminary
 1664 run from each laboratory used for analysis)1665 ²Some of the runs did not meet the assay acceptance criteria and therefore were excluded from the analysis.1666 ³For the cryo WB/IL-1 test method, each substance tested only once in each laboratory.1667 ⁴All possible combinations of runs among the 3 laboratories were compared (with the exception of cryo WB/IL-1, which
 1668 was only tested once in each laboratory, resulting in only one possible combination per substance).
 1669

1670 In the second study, reproducibility was evaluated with the same ten substances used for
 1671 evaluating accuracy. In this study, each of the substances was spiked with five concentrations
 1672 of endotoxin and tested once in each of three laboratories. As shown in **Table 7-4**, the
 1673 agreement across three laboratories for each test method ranged from 57% to 88%,
 1674 depending on the test method considered. The levels, and order of agreement among

1675 laboratories was the same for both studies; the WB/IL-1 test method showed the least
1676 agreement (57-58%) and the cryo WB/IL-1 test method showed the most (88-92%)¹².

1677

1678 **Table 7-4 Interlaboratory Reproducibility of *In Vitro* Pyrogenicity Test Methods**

Lab Comparison ¹	Agreement Between Laboratories ¹				
	WB/IL-1 ²	Cryo WB/IL-1	WB/IL-6	PBMC/IL-6 ³	MM6/IL-6
1 vs 2	73% (35/48)	84% (38/45)	85% (41/48)	84% (42/50)	90% (45/50)
1 vs 3	82% (40/49)	88% (21/24)	85% (41/48)	86% (43/50)	90% (43/48)
2 vs 3	70% (33/47)	100% (25/25)	88% (44/50)	90% (45/50)	83% (40/48)
Mean	75%	91%	86%	87%	88%
Agreement across 3 labs	57% (27/47)	88% (21/24)	79% (38/48)	80% (40/50)	81% (39/48)

1679 Abbreviations: Cryo = Cryopreserved; IL-1 = Interleukin-1; IL-6 = Interleukin-6; MM6 = Mono Mac 6; WB = Whole blood
1680 ¹Data from 10 substances spiked with endotoxin (WHO-LPS 94/580 [*E. coli* O113:H10:K-]) at 0, 0.25, 0.5, 0.5, and 1.0
1681 EU/mL tested once in three different laboratories

1682 ²Interlaboratory reproducibility was also evaluated for the modified WB/IL-1 test method protocol (i.e., 96-well plates, see
1683 **Table 2-1**) in which the reproducibility between laboratories ranged from 83% to 92% (mean=89%) and the agreement
1684 across three labs was 83%.

1685 ³Interlaboratory reproducibility was also evaluated for the modified PBMC/IL-6 test method protocol (i.e., cryo PBMCs, see
1686 **Table 2-1**) in which the reproducibility between laboratories ranged from 76% to 96% (mean=84%) and the agreement
1687 across three labs was 76%.

1688

1689 7.3 Historical Positive and Negative Control Data

1690 No historical control data were provided for any of the five *in vitro* pyrogenicity test
1691 methods. However, the intralaboratory repeatability analysis described in **Section 7.2.1**
1692 included repeat testing of both spiked (0.5 EU/mL endotoxin) and unspiked saline, and the
1693 accumulated positive and negative control values, respectively for each of the methods. As a
1694 result, the database that was accumulated during the ECVAM validation studies provides an
1695 indication of the range and variability in responses for the positive and negative controls.

¹² Interlaboratory reproducibility was also evaluated for the modified PBMC/IL-6 test protocol (i.e., cryo PBMCs, see **Table 2-1**) in which the reproducibility between laboratories ranged from 76% to 96% (mean=84%) and the agreement across three labs was 76%.

1696 **8.0 TEST METHOD DATA QUALITY**

1697 **8.1 Adherence to National and International GLP Guidelines**

1698 Ideally, all data supporting the validity of a test method should be obtained and reported in
1699 accordance with GLP guidelines, which are nationally and internationally recognized rules
1700 designed to produce high-quality laboratory records. GLPs provide a standardized approach
1701 to report and archive laboratory data and records, and information about the test protocol to
1702 insure the integrity, reliability, and accountability of a study (OECD 1998; U.S. EPA 2003a,
1703 2003b; FDA 2003).

1704 The ECVAM validation studies for the five proposed methods were carried out in accordance
1705 with GLP guidelines and all deviations from GLP compliance are noted in the BRDs.

1706 Although no direct statement of GLP compliance was provided for the historical RPT data
1707 provided in the ECVAM BRDs, the studies were conducted at the Paul Ehrlich Institute
1708 (PEI), which is a German Federal Agency for Sera and Vaccines that provides regional
1709 support for German regulatory authorities, is qualified for granting marketing approval of
1710 certain marketed biological products (e.g., sera, vaccines, test allergens), and is a WHO
1711 collaborating center for quality assurance of blood products and *in vitro* diagnostics
1712 (<http://www.pei.de>).

1713 **8.2 Data Quality Audits**

1714 Formal assessments of data quality, such as a quality assurance (QA) audit, generally involve
1715 a systematic and critical comparison of the data provided in a study report with the laboratory
1716 records generated for the study. No attempt was made to formally audit the quality of the
1717 data presented in the five ECVAM BRDs. However, as indicated in **Section 5.2**, the raw data
1718 from the validation studies are available from the participating laboratories for a quality
1719 analysis.

1720 **8.3 Impact of Deviations from GLP Guidelines**

1721 The impact of the deviations from the GLP guidelines, as reported in the ECVAM BRDs,
1722 was not evaluated.

1723

1724 **8.4 Availability of Laboratory Notebooks or Other Records**

1725 All records are stored and archived by the participating laboratories and are available for
1726 inspection.

1727 **8.5 Need for Data Quality**

1728 Data quality is a critical component of the validation process. To ensure data quality,
1729 ICCVAM recommends that all data generated during the validation of a method be available,
1730 along with the detailed protocol(s) under which the data were produced. Original data should
1731 be available for examination, as should supporting documentation such as laboratory
1732 notebooks. Ideally, the data should adhere to GLP guidelines (ICCVAM 1997). Data
1733 protocols for the validation studies summarized here are available from ECVAM (see
1734 **Appendix A**), and the data from the individual laboratories are available for inspection, as
1735 indicated in **Section 8.4**.

1736

1737

1738 **9.0 OTHER SCIENTIFIC REPORTS AND REVIEWS**

1739 **9.1 Summaries of *In Vitro* Pyrogenicity Test Methods and Data from Published**
1740 **and Unpublished Studies**

1741 Data from the *in vitro* pyrogenicity methods proposed in the ECVAM submission are
1742 provided in **Appendix A** and were used in the performance analyses described in **Section 6.0**
1743 and **Section 7.0**. A *FR* notice (Vol. 70, No. 241, pp. 74833-74834, Dec 16, 2005) was
1744 published requesting the submission of data from the RPT, the BET, and from *in vitro*
1745 pyrogenicity testing with the five test methods described in this BRD. No data were received
1746 in response to this notice.

1747 NICEATM conducted a prescreen evaluation of the ECVAM BRDs to verify that the
1748 information contained fulfilled the requirements prescribed in the ICCVAM submission
1749 guidelines (ICCVAM 2003). Based on this evaluation, the PWG requested additional data for
1750 the proposed methods to support many of the claims in the ECVAM BRD (e.g., the ability to
1751 detect both endotoxin and non-endotoxin pyrogens). In response to this request, ECVAM
1752 provided supplemental information that included unpublished data from *in vitro* pyrogenicity
1753 tests as an attempt to address these issues (see **Appendix B**).

1754 As indicated in **Section 1.5**, NICEATM conducted an online literature search for relevant
1755 information on the proposed test methods using multiple internet databases (i.e., PubMed,
1756 SCOPUS, TOXLINE, Web of Science). This search revealed nine additional scientific
1757 publications that contained data from *in vitro* pyrogenicity product testing. Some of these
1758 studies contained comparisons of the results obtained in the *in vitro* method(s) with those
1759 obtained in the RPT and/or BET (see **Tables 9-1 to 9-8**). These studies were not included in
1760 previous sections of the BRD because they used a different method or protocol, or because
1761 they lacked sufficient information for an evaluation of accuracy and reliability (e.g., an
1762 adequate validation study design was not included, a standardized reference pyrogen was not
1763 used). Summaries of these published studies and available data from the *in vitro* pyrogenicity
1764 methods are presented below.

1765 9.1.1 Andrade et al. (2003)

1766 The authors evaluated the utility of human PBMCs and diluted WB for *in vitro* pyrogenicity
1767 tests and compared the responses to those obtained in the BET and RPT for the same diverse
1768 sampling of parenteral pharmaceuticals and biological products (see **Tables 9-1** and **9-2**).

1769 Interference testing of each substance was performed with spikes of the international
1770 endotoxin standard (i.e., WHO-LPS 94/580). These studies established an endotoxin
1771 detection limit of 0.06 EU/mL for both *in vitro* assays, and the results were consistent with
1772 those from the BET and RPT. The authors concluded that both the PBMC and WB methods
1773 were comparable to the BET and the RPT in their ability to detect and quantify the presence
1774 of endotoxin. In addition, the WB test method was able to detect concentration-dependent IL-
1775 6 release on exposure of WB to non-endotoxin pyrogens and pyrogens from Gram-positive
1776 organisms (i.e., *Candida albicans* and *Staphylococcus aureus*).

1777

1777

1778 **Table 9-1 Results of Pyrogenicity Testing of Pharmaceutical/Biological Products in**
 1779 **the Human PBMC Assay, the BET, and the RPT¹**

Product	Number of Batches ²	PBMC (EU/mL)	BET (EU/mL)	RPT
Ampicillin - 1000 mg/5 ml A	1	<6	<0.06	Pass
Ampicillin - 1000 mg/5 ml A	1	<6	<0.06	Pass
Gentamycin - 80 mg/2 mL	2	<3	<0.06	Pass
Oxacillin - 500 mg/5 mL	2	<3	<0.06	Pass
Enoxaparin - 100 mg/mL	3	<1.2	<0.06	Pass
Insulin 100 - U/mL	2	<3	<0.06	Pass
Tenoxicam - 40mg/2 mL	1	<6	<0.06	Pass
Metoclopramide - 10 mg/2 mL	4	<3	<0.06	Pass
Calcium folinate - 50 mg/5 mL	1	<2.4	<0.06	Pass
Ranitidine - 25 mg/mL	2	<6	1.2-2.4	Pass
Pantoprazol - 40 mg/10 mL	1	<3	<0.06	Pass
Human serum albumin - 20%	1	<4.8	0.48-0.96	Pass
Erythropoietin - 4000 IU/vial A	1	<1.2	0.48-0.96	Pass
Erythropoietin - 2000 IU/vial B	1	112 ± 10 ¹	491-983	Fail
Erythropoietin - 4000 IU/vial C	1	<1.2	<0.06	Pass
recG-CSF - 200 µg/vial A	3	<0.6	<0.06	Pass
Saline solution - 0.9% A	1	<0.3	<0.06	Pass

1780 Abbreviations: BET = Bacterial Endotoxin Test; CSF = Colony Stimulating Factor; PBMC = Peripheral blood mononuclear
 1781 cells; rec = Recombinant; RPT = Rabbit pyrogen test; SD = Standard deviation

1782 ¹From Andrade et al. (2003)

1783 ²Batch results were combined; PBMC and BET study values represent a mean±SD value or consensus detection limits (n=3
 1784 donors; 4 replicates from each donor).

1785

1786

1786
1787
1788**Table 9-2 Results of Pyrogenicity Testing of Pharmaceutical/Biological Products by the Human WB Culture Assay, the BET, and the RPT¹**

Product	Number of Batches ²	WB Culture (EU/mL)	BET (EU/mL)	RPT
Dipyron - 500 mg/mL	3	<24	<0.06	Pass
Amikacin - 500 mg/2 mL	2	<12	<0.06	Pass
Ampicillin - 1000 mg/5 ml A	1	<6	<0.06	Pass
Ampicillin - 1000 mg/5 ml A	1	<6	<0.06	Pass
Gentamycin - 80 mg/2 mL	2	<6	<0.06	Pass
Oxacillin - 500 mg/5 mL	2	<6	<0.06	Pass
Vancomycin - 500 mg/5 mL	2	<6	<0.06	Pass
Enoxaparin - 100 mg/mL	3	<0.6	<0.06	Pass
Heparin - 5000 IU/mL	2	<0.6	<0.06	Pass
Insulin - 100 U/mL	3	<6	<0.06	Pass
Ketoprofen - 100 mg/2mL	1	<6	<0.06	Pass
Diclofenac - 75 mg/3 mL	1	<12	<0.06	Pass
Tenoxicam - 40 mg/2 mL	2	<6	<0.06	Pass
Metoclopramide - 10 mg/2 mL	3	<3	<0.06	Pass
Cytarabine - 100 mg/5mL	1	<1.2	<0.06	Pass
Calcium folinate - 50 mg/5 mL	1	<0.6	<0.06	Pass
Ranitidine - 25 mg/mL	1	<6	1.2-2.4	Pass
Pantoprazol - 40 mg/10 mL	1	<6	<0.06	Pass
Furosemide - 10 mg/mL	2	<0.6	<0.06	Pass
rec-hGH - 4 IU/vial A	2	<0.2	<0.06	Pass
rec-hGH - 4 IU/vial B	1	12.4±2.5 ¹	15.84-31.68	Pass
Human serum albumin - 20%	1	<2.4	0.48-0.96	Pass
Erythropoietin - 4000 IU/vial A	1	0.76	0.48-0.96	Pass
Erythropoietin - 2000 IU/vial B	1	141±2.8 ¹	491-983	Fail
Erythropoietin - 4000 IU/vial C	1	<0.6	<0.06	Pass
recG-CSF - 300 µg/vial	3	<0.6	<0.06	Pass
Saline solution 0.9% A	2	<0.3	<0.06	Pass
Saline solution 0.9% B	1	44.8±5 ¹	48-96	Fail
Glucose - 0.5%	1	2054±95 ¹	1920-3840	Fail
Vitamin K - 10 mg/mL	2	<6	<0.06	Pass

1789
1790
1791
1792
1793
1794
1795

Abbreviations: BET = Bacterial Endotoxin Test; CSF = Colony Stimulating Factor ; hGH = Human growth hormone; rec = recombinant; RPT = Rabbit pyrogen test; SD = Standard deviation; WB = Whole Blood

¹From Andrade et al. (2003)

²Batch results were combined; PBMC and BET study values represent a mean±SD value or consensus detection limits (n= 3 donors; 4 replicates from each donor).

1795

1796 9.1.2 Bleeker et al. (1994)

1797 This study measured IL-6 release from PBMCs as an indicator of pyrogenicity for *in vitro*
1798 safety testing of hemoglobin (Hb) solutions. The authors demonstrated that pure,
1799 polymerized Hb produced under aseptic conditions did not induce or inhibit IL-6 production,
1800 whereas production under non-aseptic conditions led to IL-6 release, which was also seen
1801 with the BET. Based on these findings, the authors concluded that IL-6 release from isolated
1802 PBMCs provides a sensitive indicator of endotoxin contamination in Hb solutions. The
1803 observed detection limit for endotoxin in Hb solutions (below 0.4 EU/mL), led the authors to
1804 suggest that this test method would be more sensitive to the presence of endotoxin than the
1805 RPT.

1806 9.1.3 Carlin and Viitanen (2003)

1807 Using WB and MM6-based *in vitro* pyrogenicity methods, this study evaluated the pyrogenic
1808 potential of a multivalent vaccine, Infanrix[®] (GlaxoSmithKline) that contains protein and
1809 polysaccharide components from both Gram-positive and Gram-negative bacteria. The five
1810 Infanrix[®] vaccines studied (e.g., Infanrix[®], Infanrix[®] Hep B, Infanrix[®] polio, Infanrix[®] hexa,
1811 and Infanrix[®] polio Hib) contain Gram-positive bacterial components that are potentially
1812 pyrogenic but not detectable in the BET. IL-6 production in the WB/IL-6 test method varied
1813 among the seven donor blood samples in response to each of the five vaccines. Some donor
1814 samples produced a weak or no IL-6 release and others produced a large release (**Table 9-3**).
1815 However, IL-6 production from any single donor was similar to all vaccines when tested at
1816 various times. The variability in the magnitude of response to each vaccine among donors
1817 and the consistency of the response of any single donor was also seen when IL-1 β was used
1818 as a marker. IL-6 release from WB was also examined following exposure to three
1819 concentrations of endotoxin standard (0.2, 2, and 20 pg/1.2 mL). All donor WB samples
1820 released IL-6 in a concentration responsive manner.

1821 The IL-6 release from MM6 cells (**Table 9-4**) exposed to the five Infanrix[®] vaccines was
1822 measured using an ELISA assay and compared to the responses induced by three
1823 concentrations of endotoxin standard (0.2, 2, and 20 pg/1.2 mL) in three separate
1824 experiments. The MM6 cells produced minimal responses to the vaccines when compared to

1825 WB, but released significant amounts of IL-6 in response to high concentrations of
1826 endotoxin. However, IL-6 induction by two different endotoxin standards in MM6 cells was
1827 strongly attenuated (>80% inhibition) when either of two vaccines (Infanrix[®] and Infanrix[®]
1828 Hep-B) was present (data not included in **Table 9-4**). Based on these studies, the authors
1829 suggested that a BET or RPT result might not correlate with the human fever response one
1830 might expect in humans immunized with such vaccines, because the production of
1831 proinflammatory cytokines may be compromised by various components in the vaccine
1832 product, and because Gram-positive components in the vaccines would not be detected in the
1833 BET.

1834 **Table 9-3 IL-6 Production from WB after Exposure to Endotoxin or Five Infanrix® Vaccines^{1,2}**

Experiment (Blood Donor)	Endotoxin		Vaccine (Absorbance in ELISA; n=4 ³)					
	Endotoxin (pg/1.2 mL)	Absorbance (ELISA)	Dilution (µL vaccine/1.2 mL)	Infanrix®	Infanrix® Hep- B	Infanrix® Hexa	Infanrix® Polio	Infanrix® Polio Hib
1	0.2	0.47	0.03	0.945	1.052	1.069	0.869	1.082
	2	0.971	0.3	1.826	2.055	2.014	1.832	1.919
	20	1.116	3	2.826	2.587	2.638	2.609	2.2
2	0.2	0.001	0.03	0.149	0.256	0.231	NT	0.284
	2	0.127	0.3	0.869	0.847	1.095	NT	0.933
	20	0.764	3	1.998	1.986	2.187	NT	1.685
3	0.2	-0.007	0.03	0.005	0.037	0.009	0.007	0.208
	2	0.09	0.3	0.275	0.457	0.282	0.321	0.261
	20	0.811	3	0.941	1.057	0.795	1.284	1.325
4	0.2	0.006	0.03	0.056	0.053	0.028	0.088	0.104
	2	0.043	0.3	0.165	0.312	0.44	0.309	0.533
	20	0.458	3	1.229	1.489	1.476	1.181	1.242
5	0.2	0.043	0.03	-	0.071	-0.003	-0.003	0.011
	2	0.024	0.3	0.007	0.014	0.004	0.03	0.05
	20	0.435	3	0.042	0.164	0.008	0.08	0.12
6	0.2	0.013	0.03	-0.009	-0.018	-0.01	-0.022	0.012
	2	0.022	0.3	-0.007	-0.008	0.005	-0.019	-0.007
	20	0.569	3	0.132	0.411	0.042	0.132	0.188
7	0.2	0.036	0.03	-0.012	-0.012	-0.01	-0.014	0.07
	2	0.014	0.3	-0.01	-0.01	-0.012	-0.011	-0.013
	20	0.436	3	0.183	0.274	0.045	0.183	0.525

1835 Abbreviations: ELISA = Enzyme-linked immunosorbent assay; Hep = Hepatitis; IL-6 = Interleukin-6; NT = Not tested; WB = Whole blood

1836 ¹From Carlin and Viitanen (2003)

1837 ²WB was challenged with endotoxin standard or vaccine in pyrogen-free water to provide the final concentration and incubated overnight at 37°C.

1838 ³Duplicate samples were run in two separate experiments.

1839

1840

1841 **Table 9-4 IL-6 Production by MM6 Cells after Exposure to Endotoxin or Five Infanrix[®] Vaccines^{1,2}**

MM6 Batch	Endotoxin			Vaccine (Absorbance in ELISA; 250,000 MM6 cells); n=4 ³				
	Endotoxin (pg/1.2 mL)	Absorbance in IL-6 ELISA	Dilution (µL vaccine/1.2 mL)	Infanrix [®]	Infanrix [®] Hep-B	Infanrix [®] Hexa	Infanrix [®] Polio	Infanrix [®] Polio Hib
1	0.2	-0.001	0.3	0.013	0.014	0.001	0.002	-0.001
	2	0.026	3	0.078	0.158	0.06	0.105	0.07
	20	0.383	30	0.054	0.052	0.053	0.106	0.089
2	0.2	-0.001	0.3	0.004	0.01	0.001	0.003	0.004
	2	0.025	3	0.033	0.062	0.019	0.037	0.032
	20	0.4	30	0.013	0.012	0.018	0.038	0.038
3	0.2	-0.009	0.3	-0.012	-0.017	-0.021	-0.014	-0.019
	2	0.03	3	0.019	0.05	0.01	0.043	0.026
	20	0.192	30	-0.018	-0.012	-0.007	0	0.005

1842 Abbreviations: ELISA = Enzyme-linked immunosorbent assay; IL-6 = Interleukin-6; MM6 = Mono Mac 6

1843 ¹From Carlin and Viitanen (2003)1844 ²MM6 cells were stimulated with endotoxin standard or vaccine in pyrogen-free water to provide the final concentration and incubated overnight at 37°C.1845 ³n = Duplicate samples were run in two separate experiments.

1846

1847

1848 9.1.4 Carlin and Viitanen (2005)

1849 This study provides support for the findings from a previous study (Carlin and Viitanen
1850 (2003) in which the authors demonstrated IL-6 release by a WB method in response to
1851 pyrogenic or spiked multivalent vaccine preparations that were inactive in the BET. It also
1852 confirms that IL-6 was released from WB of some, but not all donors. The present study
1853 demonstrates that IL-6 release in susceptible donors was caused by toxoids from Gram-
1854 positive diphtheria, and to a lesser extent, from tetanus bacterial components of the vaccines.
1855 The WB donors were studied for two years and their responses to the individual vaccines,
1856 whether responsive or non-responsive, were consistent. The responses of these donors to
1857 Gram-negative endotoxin or LTA from Gram-positive bacteria were consistent and
1858 confirmed the findings of Fennrich et al. (1999) with respect to the consistency of responses
1859 among several hundred blood donors to endotoxin. The authors concluded that individual
1860 donor-specific differences in IL-6 release from WB exposed to the multivalent vaccines
1861 resulted from toxoids present in the diphtheria or tetanus component, and noted that these
1862 donor-specific responses to the vaccines were not observed in the BET.

1863 9.1.5 Daneshian et al. (2006)

1864 This study describes the development of a modification to the WB/IL-1 method termed
1865 AWIPT (Adsorb, Wash, *In Vitro* Pyrogen Test). The authors indicate that this modification is
1866 intended to increase sensitivity to the presence of endotoxin contamination by isolating
1867 endotoxin from WB. To accomplish this, the sample containing endotoxin (naturally
1868 occurring or spiked) is treated with HSA covalently linked to macroporous acrylic beads. The
1869 HSA-treated beads bind the endotoxin, which is subsequently eluted from the beads. The
1870 WB/IL-1 test method is therefore performed using a slightly modified protocol in which the
1871 diluted WB is incubated overnight with the sample in the bead suspension.

1872 The results showed that HSA-coated beads bind endotoxin in a concentration-dependent
1873 manner (when spiked with 0, 25, 50, and 100 pg/mL LPS), but little or none was bound to
1874 unmodified beads. The test showed a detection limit of 25 pg/mL LPS (i.e., 0.25 EU/mL),
1875 which is less sensitive than the BET (3 pg/mL) and slightly lower than the RPT (50 pg/mL).
1876 IL-1 β secretion in response to either LPS or LTA was generally higher using the AWIPT
1877 procedure, but the concentrations of LPS or LTA needed to induce a response were similar;

1878 thus the sensitivity of this test modification was comparable to that of the unmodified
1879 WB/IL-1 test method (*Note: Additional data received from ECVAM suggests that*
1880 *modification of the AWIPT is capable of increasing the sensitivity for detection of endotoxin*
1881 *from the current 0.25 EU/mL to 0.0001 EU/mL [see **Appendix B**]).*

1882 Daneshian et al. studied the kinetics of cytokine release from WB in response to a challenge
1883 with 2 pg/mL endotoxin. IL-1 β release in the AWIPT-treated samples lagged slightly behind
1884 that of the standard WB/IL-1 test in the 0 to 8 hr time period, whereas more IL-1 β was
1885 produced in the AWIPT-treated samples in the 10 to 30 hr time period. Some
1886 immunomodulatory or toxic cancer drug samples tested in the WB/IL-1 method interfered
1887 with the WB/IL-1 assay and required a higher dilution (1/10 to 1/100) to detect IL-1 β .
1888 Detection of endotoxin spiked into these test samples (measured as IL-1 β release) generally
1889 occurred at lower dilutions in AWIPT than in the WB/IL-1 test method, suggesting that the
1890 interfering substances were removed by the procedure. For example, five dilutions (ranging
1891 from 1/3 to 1/316) of liposomal daunorubicin were spiked with 25 pg/mL of endotoxin and
1892 detection of IL-1 β was compared between the two methods. This cytokine was not detectable
1893 in the WB/IL-1 method (< 30% of the IL-1 β released by endotoxin) at any drug dilution,
1894 whereas in the AWIPT, IL-1 β was detected at drug dilutions of 1/32, 1/100, and 1/316
1895 (>78% of the IL-1 β released by endotoxin).

1896 The authors concluded that the inclusion of endotoxin adsorption and washing steps in the
1897 WB/IL-1 method (i.e., the AWIPT) to remove potentially interfering substances improved
1898 the detection of pyrogenic contaminants in immunomodulatory and toxic cancer drug
1899 samples. They suggest that the AWIPT method offers an improvement for safety testing of
1900 products administered to patients, and for batch control in pharmaceutical processing.

1901 9.1.6 Eperon et al. (1996, 1997)

1902 Eperon and colleagues developed an *in vitro* test system for measuring pyrogenic substances
1903 using two clones derived from MM6 cells (H.W.L. Ziegler-Heitbrock, University of Munich)
1904 and one from a THP-1 cell line (European Collection of Animal Cell Cultures (Porton Down,
1905 Salisbury, United Kingdom). These clones are reported to be phenotypically stable over time
1906 with respect to their superior responsiveness to endotoxin than the parent cell lines.

1907 Endotoxin content was measured by the release of TNF- α using an immunoassay. These

1908 clones demonstrate high LPS sensitivity when non-pyrogenic fetal calf serum is used in the
1909 assay as a serum supplement. Enhanced expression of the cell-surface endotoxin receptor
1910 CD14 was obtained by pretreatment of the cells for two days with calcitrol. Purified
1911 endotoxin (i.e., LPS; smooth strain and rough mutant), other cellular components from
1912 Gram-negative or Gram-positive bacteria, and Mycobacteria were tested. The MM6 clones
1913 responded to these pyrogenic products in an order of potency of detection equivalent to that
1914 found in the RPT and similar to that observed in the BET (i.e., Gram-negative endotoxin >
1915 Gram-positive material > non-endotoxin pyrogens). The response of the THP-1 clone was
1916 similar to that of the MM6 clones, except that the THP-1 clone did not respond to
1917 diphosphoryl lipid A, a structural component of LPS.

1918 Pyrogenicity testing of a panel of stable blood products including albumin and
1919 Immunoglobulin G (IgG) for parenteral use produced similar results in the RPT and MM6 or
1920 THP-1 clones when tested as received (i.e., free of detectable pyrogens) or positive results
1921 when spiked with 20 EU/mL endotoxin (**Table 9-5**) with a few exceptions. For example, in
1922 the cell-based test, there was one borderline, but significant positive result in an unspiked
1923 sample (i.e., false positive relative to the RPT). In the BET, 4 of 13 (31%) unspiked samples
1924 tested positive (i.e., false positive). The results suggest that the cell-based assays may
1925 produce fewer false positives than the BET.

1926 When ten bacterial and viral vaccine preparations were evaluated, the monocytoid cell-based
1927 test method (e.g., combined results from two experiment with each cell line) correlated well
1928 with the RPT (positive or negative for endotoxin) with the exception of one preparation that
1929 produced nearly 10-fold less TNF- α than the other samples, and was near the limit of
1930 detection, but not statistically significant different from the negative control (**Table 9-6**). The
1931 authors suggest that these cloned monocytoid cell-based test methods are valid *in vitro*
1932 alternatives for detection of endotoxin in commercial preparations, and produce results
1933 comparable to the RPT and BET.

1934

1934

1935 **Table 9-5 Pyrogenic Activity of Blood Preparations for Parenteral Use¹**

Preparation	Endotoxin Spike	RPT ²	BET ^{3,4}	Cell Test ^{5,6}
IgG for i.v. use	20 EU/mL	+	+	+
	20 EU/mL	+	+	+
	-	-	-	-
	-	-	-	-
	-	-	-	-
	-	-	-	-
	-	-	-	-
	-	-	-	-
Albumin	20 EU/mL	+	+	+
	20 EU/mL	+	+	+
	-	-	-	-
	-	-	-	-
	-	-	-	-
	-	-	-	-
<i>Test threshold⁷</i>		$\Delta t = 1.5^\circ\text{C}$	300 pg/mL LPS	50 pg/mL TNF

1936 Abbreviations: BET = Bacterial Endotoxin Test; ELISA = Enzyme-linked immunosorbent assay; EU =
 1937 Endotoxin Units; IgG = Immunoglobulin G; i.v. = Intravenous; LPS = lipopolysaccharide; RPT = Rabbit
 1938 pyrogen test; TNF = Tumor necrosis factor

1939 ¹From Eperon et al. (1997)

1940 ²n=3

1941 ³n=2

1942 ⁴Haemachem BET assay (St. Louis)

1943 ⁵n=4 [Note: cell type not specified, although author's claim that either MM6 or THP-1 are equally capable of
 1944 endotoxin detection]

1945 ⁶TNF induction was determined using a commercial TNF ELISA.

1946 ⁷RPT threshold was obtained from the European Pharmacopoeia; the threshold for the BET and cell-based test
 1947 methods was considered to be equal to 2 standard deviations from the mean of a set of negative samples.

1948

1949

1949 **Table 9-6 Pyrogenic Activity of Vaccine Preparations¹**

Vaccine Preparation ²	Batch ³	Cell Test ⁴	Pyrogenicity
IgG for i.m. use	A-1	n.d. ⁵	-
	A-2	n.d.	-
Bacterial vaccines	B-1	10.8±0.3	+
	C-1	6.0±3.6	+
	D-1	1.4±1.8	-
Viral vaccines	E-1	n.d.	-
	E-2	n.d.	-
	F-1	n.d.	-
	F-2	n.d.	-
	G-1	21.2±3.2	+

1950 Abbreviations: IgG = Immunoglobulin G; i.m. = Intramuscular; n.d. = Non-detectable

1951 ¹From Eperon et al. (1997)

1952 ²Vaccine solutions were tested at 1/20 (v/v)

1953 ³Letters refer to distinct types of vaccine preparations; numbers to different lots

1954 ⁴TNF- α production in ng/mL \pm S.E.M (n=3) [Note: cell type not specified, although author's claim that either

1955 *MM6* or *THP-1* are equally capable of endotoxin detection]

1956 ⁵No measurable quantity of cytokine was detected.

1957

1958

1959 9.1.7 Pool et al. (1998)

1960 This study describes a WB assay for the detection of pyrogens in blood products. IL-6 release
 1961 from WB in response to endotoxin is used to define a pyrogenic response. This assay was
 1962 highly responsive to *Escherichia coli* endotoxin (i.e., the limit of detection of endotoxin was
 1963 1.25 EU/mL), and also responded to whole bacteria (*E. coli* and *Bacillus subtilis*
 1964 [*B.subtilis*]). There was considerable variation in IL-6 levels released from WB between
 1965 donors following exposure to endotoxin, but each donor response was always linear. The
 1966 pyrogenicity of production batches of human serum albumin (HSA), fibronectin (Fn), and
 1967 stabilized human serum (SHS) solutions were evaluated using the WB method and compared
 1968 to the BET and RPT. Spike recovery in batches of these samples varied between 90 and
 1969 116% for *E. coli* endotoxin, 74 to 111% for *B. subtilis*, and 61 to 99% for *E. coli* and the
 1970 products tested did not interfere with the IL-6 assay system. Good correlations were found
 1971 among the WB, BET, and RPT assay results (**Table 9-7**). Of 22 products tested, the WB
 1972 assay and the RPT were in agreement (i.e., pass or fail) for all tests, while one sample was
 1973 classified as negative in the BET, but positive in both the WB method and the RPT. The
 1974 detection limit for endotoxin by the WB method was 1.25 EU/mL, which is lower than the
 1975 established pyrogen cut-off level (as stated in the European Pharmacopoeia) for the products

1976 under investigation (i.e., 2 EU/mL for HSA and SHS; 4.5 EU/mL for Fn). The authors
 1977 concluded that the WB assay was able to detect both Gram-negative and a Gram-positive
 1978 pyrogens and exhibited greater sensitivity to endotoxin than the RPT.
 1979

1980 **Table 9.7 Comparison of the WB test, BET, and the RPT for Detecting the**
 1981 **Pyrogenicity of Production Batches of Biological Products¹**

Product	Batch	WB (EU/mL) ²	BET	RPT
Fibronectin - 0.5 mg/mL	Fn3195	<0.05	Pass	Pass
	Fn3296	<0.05	Pass	Pass
	Fn3596	1.28	Pass	Pass
Human serum albumin - 200 mg/mL	B274	29.4	Fail	Fail
	B291	<0.05	Pass	Pass
	B293	<0.05	Pass	Pass
	B294	<0.05	Pass	Pass
	B295	<0.05	Pass	Pass
	B296	<0.05	Pass	Pass
	B297	<0.05	Pass	Pass
	B298	1	Pass	Pass
	B299	1.1	Pass	Pass
	B300S	1	Pass	Pass
	B301	<0.05	Pass	Pass
	B302	>20	Pass ³	Fail
Stabilized human serum - 50 mg/mL	SS349	0.7	Pass	Pass
	SS350	<0.05	Pass	Pass
	SS351	<0.05	Pass	Pass
	SS352	0.5	Pass	Pass
	SS353	<0.05	Pass	Pass
	SS354	0.6	Pass	Pass
	SS355	0.5	Pass	Pass

1982 Abbreviations: BET = Bacterial Endotoxin Test; ELISA = Enzyme-linked immunosorbent assay; EU = Endotoxin Units; IL-
 1983 6 = Interleukin-6; WB = Whole blood; RPT = Rabbit pyrogen test

1984 ¹From Pool et al. (1998)

1985 ²Result based on IL-6 secretion in human WB using an ELISA calibrated to an *E. coli* endotoxin standard (Kabi
 1986 Diagnostica).

1987 ³False negative relative to the RPT response
 1988
 1989

1990

1990 9.1.8 Taktak et al. (1991)

1991 This paper summarizes the development of an *in vitro* pyrogenicity test method based on IL-
1992 6 release from MM6 cells. A detectable level of IL-6 was released in response to 2.5 pg/mL
1993 of endotoxin, yielding a level of sensitivity of 25 pg/mL when testing 5% HSA at a 1/10
1994 dilution for the presence of endotoxin. Three batches of a therapeutic HSA that caused fever
1995 in humans were positive in the MM6/IL-6 method, whereas the same substances were
1996 negative in the BET and the RPT (**Table 9-8**). As in the BET, the samples required a 1/10
1997 dilution to remove interfering substances, and had assay sensitivity equal to that of the BET
1998 (25 pg/mL) and 40-fold greater than the RPT (1000 pg/mL). The authors suggest that the
1999 MM6/IL-6 method represents an important alternative to the existing pyrogen tests and may
2000 be a more appropriate end product test for the detection of pyrogens in parenteral products,
2001 such as HSA, that cannot be detected in the BET.

2002

2002

2003 **Table 9-8 Results of Pyrogenicity Testing of Batches of Therapeutic HSA Using the**
 2004 **MM6/IL-6, BET, and RPT¹**

Batch of HSA	Endotoxin Quantitation by IL-6 Release (pg/mL)	Endotoxin Quantitation by the BET		RPT Result
		IU/mL	pg/mL ⁶	
1 ²	97±2.3 ^{3,4}	1.0-2.0	140-280	Pass
2 ²	30±2.8 ⁴	2.4-3.2	336-448	Pass
3 ²	31±2.3 ⁴	0.5-0.75	70-105	Pass
4	<25 ⁵	<0.24	<34	Pass
5	<25 ⁵	3.6-4.8 ⁷	504-762	Pass
6	<25 ⁵	<0.26	<36	Pass

2005 Abbreviations: BET = Bacterial Endotoxin Test; HSA = Human serum albumin; IL-6 = Interleukin-6; IU =
 2006 International units; pg = Picograms; RPT = Rabbit pyrogen test

2007 ¹From Taktak et al. (1991)

2008 ²Batch of HSA used that caused fever in humans.

2009 ³Mean±S.E.M.

2010 ⁴Values are significantly different from subthreshold concentrations of endotoxin (<2.5 pg/mL; p<0.001).

2011 ⁵Values below the detection limit of the test system (25 EU/mL); preparations of HSA were tested at a dilution
 2012 of 1/10 and 2.5 pg/mL endotoxin was the lowest concentration of endotoxin tested that evoked a significant
 2013 release of IL-6.

2014 ⁶1.0 IU=0.14 ng for preparation used.

2015 ⁷False positive relative to RPT.

2016

2017

2017

2018 **9.2 Conclusions from Scientific Literature Based on Independent Peer-Reviewed**
2019 **Reports and/or Reviews**

2020 Hartung et al. (2001) provided a summary report of an ECVAM-sponsored workshop to
2021 review the current status of pyrogen testing, to review the capabilities of new pyrogen tests,
2022 and to provide recommendations for their continued development. The need for alternatives
2023 to the RPT and the BET was discussed, and their respective limitations were highlighted. The
2024 workshop compared the utility of the various methods (i.e., *in vitro* pyrogenicity test
2025 methods, BET, RPT) for testing a variety of pyrogenic materials. Workshop conclusions
2026 indicated a need for alternative test methods to address the limitations of the BET and RPT,
2027 but stressed the need for appropriate validation of any new method.

2028 An additional nine reports describing studies of cell-based *in vitro* pyrogenicity methods
2029 were obtained from the literature search described in **Section 9.1** and representative of the
2030 alternative *in vitro* pyrogenicity methods considered in the ECVAM workshop. Although
2031 these reports did not include data on test substances that could be used in the performance
2032 analysis in **Section 6.0** and **Section 7.0**, they did evaluate the use of the *in vitro* pyrogenicity
2033 test methods for sensitivity to endotoxin (i.e., endotoxin detection limit), specificity of the
2034 response to endotoxin and/or non-endotoxin pyrogens (i.e., spectrum and relative potency of
2035 various pyrogens detected), and/or the impact of interfering substances. A summary of each
2036 study is presented below.

2037 9.2.1 De Groote et al. (1992)

2038 The authors measured the release of various cytokines (IL-1 β , IL-6, TNF- α , IL-2, IFN- γ ,
2039 and granulocyte-macrophage colony stimulating factor [GM-CSF]) in response to endotoxin
2040 or phytohaemagglutinin (PHA) stimulation of WB and PBMC cultures. Endotoxin stimulated
2041 IL-1 β , TNF- α , and IL-6 release, while PHA stimulated IL-2, IFN- γ , and GM-CSF release.
2042 There was a significant correlation between production of the three endotoxin-induced
2043 cytokines and the number of monocytes in the challenged culture, suggesting that monocytes
2044 are the major source of these cytokines: the other cytokines did not correlate with any of the
2045 cell types. The data also suggested that WB produced less variable levels of cytokines than

2046 PBMC on exposure to endotoxin. Consistent results were obtained with the WB test using
2047 more than 50 different blood donors. The authors suggest that WB is a more appropriate
2048 choice for studying cytokine production *in vitro* and its modulation by exogenous or
2049 endogenous factors, because natural cell-to-cell interactions are preserved, immune
2050 mediators are available, and cytokine levels obtained with PBMC were more variable.

2051 9.2.2 Fennrich et al. (1999)

2052 Fennrich and colleagues compared a commercially available human WB/IL-1 pyrogenicity
2053 assay (PyroCheck[®] from DPC Biermann, Bad Nauheim distributed by Millenia, UK) to the
2054 BET and RPT. There was a concentration-dependent IL-1 β release in WB that was incubated
2055 with nitrocellulose filters containing live *E. coli* bacteria and *E. coli* killed by heat or by
2056 antibiotics. The authors also tested air conditioning filters from a veterinary sheep facility
2057 and identified filters to be contaminated with bacteria that were later confirmed by
2058 microbiological tests (the data and the identity of the organisms identified were not
2059 presented). The authors compared the PyroCheck[®], BET, and the RPT (**Table 9-9**) and
2060 concluded that PyroCheck[®] is a simple, accurate test that detects a wider range of pyrogens
2061 than the BET.

2062

2062

2063 **Table 9-9 Comparison of the Application Spectra of the RPT, the BET, and the**
 2064 **Human WB Assay (PyroCheck®)¹**

Test		Applications		
		PyroCheck®	BET	RPT
Pyrogens	Gram-negative	+	+	+
	Gram-positive	+	-	+
	Fungi	+	-	+
Product pyrogenicity	Biologicals	+	-	+
	Pharmaceuticals	+	+	+
	Medical devices	+	+ ²	-
	Air quality	+ ²	+ ²	-
	Blood products	+	-	-

2065 Abbreviations: BET = Bacterial Endotoxin Test; RPT = Rabbit pyrogen test; WB = Whole blood

2066 ¹From Fennrich et al. (1999)

2067 ²Based on preliminary data

2068

2069

2070 9.2.3 Hansen and Christensen (1990)

2071 This study compared the results from PBMC exposed to endotoxin or ultraviolet light-killed

2072 *S. aureus* as an index of pyrogenicity, and then compared these results to the BET and the

2073 RPT. The authors used human PBMC obtained from heparinized peripheral blood and

2074 measured IL-1-like material in culture supernatants by evaluating co-mitogenic activity on

2075 PHA-stimulated murine thymocytes (measured in units of IL-1 β where 1 unit is defined as

2076 the concentration that gives 50% of the maximal incorporation of ³H-thymidine in the

2077 thymocyte assay). The endpoint is referred to as an IL-1-like material because other

2078 cytokines such as IL-2, IL-6, and TNF α may also stimulate the proliferative response of the

2079 thymocytes. When exposed to endotoxin, PBMC secreted cytokines in a concentration-

2080 dependent manner that provided a limit of detection of 200 pg/mL of endotoxin. In

2081 comparison, the BET can normally detect 10 to 100 pg/mL of endotoxin, while the RPT can

2082 detect 500 pg/mL. Therefore, the PBMC procedure had a level of detection of endotoxin 2.5-

2083 fold lower than that of the RPT and 2-fold higher than the BET. The PBMCs also responded

2084 with greater sensitivity to the Gram-positive pyrogen *S. aureus* (10⁵ cells/mL), which was not

2085 detected in the BET (10⁹ cells/mL). Based on these results, the authors proposed that the

2086 PBMC test be used as an alternative *in vitro* test to the BET and RPT.

2087 9.2.4 Hartung and Wendel (1996)

2088 The authors stimulated human WB with various inflammatory agents to release endogenous
2089 cytokines (IL-1 β , TNF α) and inflammatory mediators (PGE₂) as an *in vitro* method for the
2090 detection of pyrogenic materials. Cytokines were released in a concentration-dependent
2091 manner following exposure to endotoxin or LTA. Heat-killed Gram-positive bacteria (*S.*
2092 *aureus*) or components of these organisms (i.e., muropeptides, LTA, enterotoxins,
2093 streptolysin O) and plant mitogens such as phorbol myristate acetate (PMA) and PHA also
2094 produced a cytokine response. Higher concentrations (three orders of magnitude) of the
2095 Gram-positive pyrogens were needed to elicit a response as compared to Gram-negative
2096 pyrogenic material.

2097 Studies to determine the variability among the responses of different donor WB samples
2098 were also performed. Only two of the 18 donor samples released IL-1 β in response to 1
2099 pg/mL of endotoxin, but all responded to 10 pg/mL endotoxin. The release of IL-1 β from the
2100 WB samples of 45 individual donors exposed to 100 ng/mL of endotoxin was also consistent.
2101 Based on these results, the authors suggested using the WB/IL-1 test method as an *in vitro*
2102 alternative to the RPT.

2103 9.2.5 Moesby et al. (1999)

2104 Moesby and colleagues compared pyrogenicity testing using MM6 cells, isolated PBMC, and
2105 the BET. LPS and ultraviolet light (UV)-killed Gram-negative *Staphylococcus typhimurium*
2106 (*S. typhimurium*) or Gram-positive *S. aureus* produced concentration-dependent increases in
2107 IL-6 production in MM6 or PBMC cultures. PBMC, but not MM6 cells, were able to
2108 differentiate UV-irradiated yeast (*C. albicans*) and mold (*Aspergillus niger*) pyrogens, as
2109 evidenced by statistically significant increases in IL-6 production. The BET can detect Gram-
2110 negative endotoxin, but not Gram-positive endotoxin or LTA (the pyrogenic component of
2111 Gram-positive bacteria), and it may weakly detect yeast or viral pyrogens that the MM6
2112 assay could not detect. Therefore, the authors suggest that pyrogen testing using MM6 cells
2113 would be a useful supplement to the BET for the detection of both Gram-negative and Gram-
2114 positive bacteria.

2115 9.2.6 Nakagawa et al. (2002)

2116 Nakagawa and colleagues describe an *in vitro* pyrogen test system based on pro-
2117 inflammatory cytokine release from a sub-clone of MM6 cells (i.e., MM6-CA8) and compare
2118 this response to a human WB culture system and the RPT. Similar to MM6 cells, MM6-CA8
2119 were developed for superior reactivity to both endotoxin and peptidoglycan (PG). The MM6-
2120 CA8 cells release IL-6, TNF- α , and IL-1, but in greater quantities than MM6 cells in the
2121 range of 1 to 1000 pg/mL of endotoxin (up to 4-fold greater) or to 1 to 1000 ng/mL PG (up to
2122 10-fold greater) compared to MM6 cells. The range of responses of human WB to the various
2123 pyrogens was similar to that of the MM6-CA8 cells. The relative potencies of the various
2124 pyrogens in the RPT were similar to those of the cytokine-induction potencies in the WB and
2125 MM6-CA8 methods, except for polyinosinic:polycytidylic acid (Poly [I:C]), which was
2126 reported to be 10,000-fold more potent as a pyrogen injected in rabbits when compared to
2127 humans. The authors conclude that these results suggest MM6-CA8 cells can detect a variety
2128 of pyrogens using IL-6 as the marker, and that these responses are highly relevant to the
2129 prediction of human fever reactions.

2130 9.2.7 Pool et al. (1999)

2131 This article describes a method to differentiate between endotoxin and non-endotoxin
2132 pyrogens when testing HSA solutions in a WB culture assay. Detection limits for four Gram-
2133 positive (*Bacillus stearothermophilus* [*B. stearothermophilus*], *B. subtilis*, *Micrococcus*
2134 *luteus*, and *S. aureus*) and four Gram-negative bacteria (*E. coli*, *Klebsiella pneumoniae*,
2135 *Proteus vulgaris*, and *Pseudomonas aeruginosa*) were expressed as the number of whole
2136 bacteria required to produce a pyrogenic response equal to that of 1.25 EU/mL endotoxin. *B.*
2137 *stearothermophilus* and *E. coli* produced concentration-dependent increases in IL-6
2138 production. The cationic antibiotic Polymyxin B, which inhibits the binding of endotoxin to
2139 the CD14 receptor, produced concentration-dependent inhibition of IL-6 release following
2140 exposure to 10 EU/mL endotoxin in the WB assay at concentrations up to 1 EU/mL and
2141 completely inhibited IL-6 release at concentrations above 2 EU/mL. In contrast, Polymyxin
2142 B had no effect on IL-6 release following exposure to *B. subtilis*. These data suggest that
2143 Polymyxin B may be useful for differentiating endotoxin and non-endotoxin pyrogenic
2144 contaminants. The data also suggests that binding of endotoxin to Polymyxin B (e.g., by
2145 linkage to an affinity column) may be used in the depyrogenation process.

2146 9.2.8 Poole et al. (2003)

2147 This paper describes a rapid single-plate *in vitro* test for the presence of pyrogenic substances
2148 based on monocyte activation. The assay uses polyclonal antibodies to IL-6 or TNF α
2149 cytokines, coated and stabilized onto 96-well plates. Monocytoid cells (e.g., PBMC, MM6 or
2150 THP-1 cells), endotoxin standard (LPS), test sample, and a second biotinylated antibody
2151 specific for the cytokine (e.g., either IL-6 or TNF- α) are incubated for 2 to 4 hr in the
2152 antibody-coated wells. An ELISA for one of the cytokines is then performed on the washed
2153 plate. IL-6 is preferred and provides a limit of detection of 0.015 EU/mL with PBMC, 0.05
2154 EU/mL in MM6 cells, and 0.03 EU/mL with diluted WB. The amount of TNF- α released in
2155 WB in response to endotoxin was approximately 50 to 70% lower than IL-6, but was released
2156 earlier (i.e., 2 hr vs. 4 hr). The amount of IL-6 released on exposure to endotoxin tended to be
2157 greater in this single plate test when compared to the traditional two plate test (i.e., in which
2158 the supernatant from one plate is transferred to a second plate for the ELISA) using PBMCs,
2159 MM6 cells, THP-1 cells, or WB. The authors report that this single plate assay using IL-6
2160 release as the endpoint can be completed in 5 hr, and that this time could be reduced to 3 hr
2161 using TNF α as the endpoint (because it is released earlier from the cells). The authors also
2162 suggest that this single plate test method is readily adaptable to high throughput assays.

2163 9.2.9 Schindler et al. (2004)

2164 The authors optimized conditions for use of cryopreserved human WB in pyrogen testing to
2165 obviate the need for fresh WB. The release of IL-1 β from fresh and cryo WB collected from
2166 five donors was used as the measure of endotoxin presence. Challenge with 0.5 or 1.0
2167 EU/mL endotoxin resulted in IL-1 β release in bloods from all donors, although kinetic
2168 studies suggested that IL-1 β release was delayed one hour in the cryopreserved samples.
2169 Cryopreservation did not appear to alter the spectrum of detectable pyrogens or immune
2170 stimuli when results were compared to that of fresh WB, and no cytokine release was
2171 measured in materials that fresh WB did not respond to. Seven clinical-grade (i.e., endotoxin-
2172 free) parenteral products spiked with 0.5 EU/mL of endotoxin revealed that there was less
2173 interference in cryo WB than in fresh WB based on lower minimal interference dilutions that
2174 were always at or below the maximal valid dilution (MVD) for each product. The data
2175 showed that a broad variety of drugs could be tested for pyrogenic contaminants using cryo

2176 WB while maintaining the endotoxin limit concentration (ELC) established in the various
2177 Pharmacopoeias.
2178
2179
2180

2180

2181

2182

2183

2184

2185

2186

2187

2188

2189

2190

2191

2192

2193 ***[This Page Intentionally Left Blank]***

2194

2195

2196 **10.0 ANIMAL WELFARE CONSIDERATIONS (REFINEMENT, REDUCTION,**
2197 **AND REPLACEMENT)**

2198 **10.1 How the Five *In Vitro* Test Methods Will Refine, Reduce, or Replace Animal**
2199 **Use**

2200 ICCVAM promotes the scientific validation and regulatory acceptance of new methods that
2201 refine, reduce, or replace animal use where scientifically feasible. Refinement, Reduction,
2202 and Replacement are known as the three "Rs" of animal protection. These principles of
2203 humane treatment of laboratory animals are described as:

- 2204 • refining experimental procedures such that animal suffering is minimized;
- 2205 • reducing animal use through improved science and experimental design; and
- 2206 • replacing animal models with non-animal procedures (e.g., *in vitro*
2207 technologies), where possible (Russell and Burch 1959)

2208 The *in vitro* pyrogenicity test methods address each of these animal welfare considerations.
2209 The assays use monocytoïd cells of human origin, obtained either from whole blood
2210 donations or from an immortalized cell line. The currently accepted pyrogenicity test
2211 methods require the use of either rabbits (RPT) or horseshoe crab haemolymph (BET).
2212 Because isolated cells are treated in these human cell assays, treatment-related pain and
2213 suffering are avoided in live animals. The capability of these five *in vitro* assays to detect
2214 Gram-negative endotoxin, suggests that they may reduce or replace the use of rabbits and/or
2215 horseshoe crabs for pyrogen testing. However, the RPT will detect classes of pyrogens that
2216 have not been examined/validated in the human cell tests, and therefore may still need to be
2217 used in certain circumstances.

2218 **10.2 Requirement for the Use of Animals**

2219 10.2.1 Rationale for the Use of Animals

2220 Non-human animal species are not used for these *in vitro* tests. As indicated above, the *in*
2221 *vitro* pyrogenicity methods use monocytoïd cells of human origin obtained from either an
2222 immortalized cell line (MM6/IL-6) or whole blood donations (WB/IL-1, WB/IL-6, cryo
2223 WB/IL-1, and PBMC/IL-6) within a short time before the test is to be performed. Therefore,

2224 human volunteers are required for four of the five methods. Standard phlebotomy techniques
2225 are used to obtain the blood samples. The only pain and distress experienced by the donor
2226 would be associated with the collection procedure (i.e., needle stick), which is commonplace
2227 in medical procedures. The use of an experienced phlebotomist perform the collection should
2228 minimize pain and distress. As indicated in **Section 2.2.2**, blood samples from up to five
2229 donors are required for a single assay, depending on the method under consideration.

2230

2231

2232 **11.0 PRACTICAL CONSIDERATIONS**

2233 Several issues are taken into account when assessing the practicality of using an *in vitro* test
2234 method in place of an *in vivo* test method. In addition to reliability and accuracy evaluations,
2235 assessments of the laboratory equipment and supplies needed to carry-out the *in vitro* test
2236 method, level of personnel training, labor costs, and the time required to complete the test
2237 method relative to the *in vivo* test method are necessary. The time, personnel cost, and effort
2238 required to conduct the proposed test method(s) must be considered to be reasonable when
2239 compared to the *in vivo* test method it is intended to replace.

2240 **11.1 Transferability of the *In Vitro* Pyrogenicity Test Methods**

2241 Test method transferability addresses the ability of a method to be accurately and reliably
2242 performed by multiple laboratories (ICCVAM 2003), including those experienced in the
2243 particular type of procedure, and otherwise competent laboratories with less or no experience
2244 in the particular procedure. The degree of transferability of a test method can be evaluated by
2245 its interlaboratory reproducibility. ECVAM measured the transferability (i.e., interlaboratory
2246 reproducibility) of each assay among experienced laboratories. The results obtained, and
2247 presented in **Tables 7-3** and **7-4**, provide an estimate of the minimum variability to be
2248 expected; interlaboratory variability is anticipated to be greater (i.e., lower transferability)
2249 among laboratories having less experience with the assays.

2250 **11.1.1 Facilities and Major Fixed Equipment**

2251 A standard laboratory facility for sterile tissue culture is necessary for performing the *in vitro*
2252 pyrogenicity methods described here. The major equipment necessary to conduct the tests are
2253 readily available and include, a laminar flow hood, tissue culture incubator, water bath, and
2254 spectrophotometric microplate reader.

2255 In contrast, the RPT requires a facility that meets applicable State and Federal regulations for
2256 the care and housing of laboratory animals. The primary expense for equipping a facility to
2257 conduct the RPT would be the acquisition of an adequate animal room and associated
2258 housing (e.g., cages, bedding, food, water, etc.) for boarding animals during the study, and
2259 specifically trained animal care support personnel.

2260

2261 11.1.2 General Availability of Other Necessary Equipment and Supplies.

2262 The equipment and supplies necessary to conduct the *in vitro* pyrogenicity test methods (e.g.,
2263 micropipetters, sterile tissue culture vessels, disposable plastic ware, assay reagents) are
2264 readily available in most scientific laboratories, or can be obtained from any of several
2265 scientific laboratory equipment vendors.

2266 The RPT requires fewer general laboratory supplies. Those that are needed are readily
2267 available in most toxicity testing laboratories, or could be readily obtained from any of a
2268 number of scientific laboratory equipment vendors.

2269 **11.2 Personnel Training Considerations**

2270 Training considerations are defined as the level of instruction needed for personnel to
2271 conduct the test method accurately and reliably (ICCVAM 2003). Evaluation of the levels of
2272 training and expertise needed to conduct the test method, as well as the training requirements
2273 needed to insure that personnel are competent in the test procedures, are discussed below.

2274 11.2.1 Required Training and Expertise Needed to Conduct the *In Vitro* Pyrogenicity Test
2275 Methods

2276 Laboratory personnel require training in the relevant enzyme immunoassay protocols and the
2277 aseptic techniques associated with mammalian tissue culture. The quality criteria associated
2278 with each *in vitro* test method may be used to ensure that personnel are competent in the
2279 performance of the various procedures. When a technician has mastered all aspects of the
2280 protocol, and can independently conduct the assay such that the quality criteria have been
2281 met, the individual is considered to have demonstrated proficiency in the assay.

2282 The RPT requires training in the care and handling of laboratory animals, and the collection
2283 of accurate rectal temperature measurements at the appropriate time intervals from each
2284 rabbit. The laboratory personnel must be adequately trained to maintain the animals, and to
2285 accurately and consistently record the proper body temperature. It is not known what, if any,
2286 proficiency requirements are in place for the RPT.

2287 11.3 Cost Considerations

2288 In addition to the major fixed equipment, there are three additional factors that contribute to
2289 the cost of the proposed *in vitro* methods: 1) cost of the monocytoid cell line (i.e., Mono Mac
2290 6); 2) cost of the reagents for the ELISA procedure; and 3) labor costs for laboratory
2291 personnel.

2292 Because the proposed *in vitro* test methods are relatively more labor-intensive than the RPT,
2293 it is estimated that the cost of any these *in vitro* methods would be more than that of the BET
2294 or the RPT. However, because these test methods are amenable to high throughput screening
2295 in a properly equipped laboratory, these increased costs could be considerably reduced.

2296 Another consideration is the need for a laboratory animal veterinarian, the housing and care
2297 of the rabbits before, during, and after the test, and the cost of replacement rabbits, where
2298 necessary. These costs could offset the one-time costs of the tissue culture equipment and
2299 microplate reader.

2300 11.4 Time Considerations

2301 The *in vitro* pyrogenicity methods require two working days for completion. On the first day,
2302 blood is drawn and the test materials are prepared and incubated with the monocytoid cells.
2303 On the second day, the cytokine release from the cells is determined by immunoassay. The
2304 BET and RPT can both be completed within one working day. However, depending on the
2305 specific protocol employed, the RPT could require additional testing in up to 12 animals,
2306 which would extend the time to completion.

2307

2308

2308

2309

2310

2311

2312

2313

2314

2315

2316

2317

2318

2319

2320

2321 ***[This Page Intentionally Left Blank]***

2322

2323

2324 **12.0 REFERENCES**

- 2325 Andrade SS, Silveira RL, Schmidt CA, Brum Junior L, Dalmora SL. 2003. Comparative
2326 evaluation of the human whole blood and human peripheral blood monocyte tests for
2327 pyrogens. *Int J Pharmaceutics* 265:115-124.
- 2328 Beutler B, Rietschel ET. 2003. Innate immune sensing and its roots: the story of endotoxin.
2329 *Nat Rev Immunol* 3(2):169-176.
- 2330 Bleeker WK, de Groot EM, den Boer PJ, BiesselsPTM, Asrden LA, Bakker JC. 1994.
2331 Measurement of interleukin-6 production by monocytes for *in vitro* safety testing of
2332 hemoglobin solutions. *Art Cells Blood Subs Immob Biotech* 22(3):835-840.
- 2333 Carlin G, Viitanen E. 2003. In vitro pyrogenicity of a multivalent vaccine; infanix.
2334 *PHARMEUROPA* 15(3):418-423.
- 2335 Carlin G, Viitanen E. 2005. In vitro pyrogenicity of the diphtheria, tetanus and acellular
2336 pertussis components of a trivalent vaccine. *Vaccine* 23(28):3709-3715.
- 2337 CFR. 2005. Code of Federal Regulations. Title 21, Volume 7. Revised as of April 1, 2005.
2338 21CFR610.13.
- 2339 Cooper JF, Levin J, Wagner HN, Jr. 1971. Quantitative comparison of in vitro and in vivo
2340 methods for the detection of endotoxin. *J Lab Clin Med* 78(1):138-148.
- 2341 Co Tui, Shrift, MH. 1942. Alternative test for pyrogens in infusion fluids. *Proc Soc Exp Biol*
2342 *Med* 49:320-330.
- 2343 Daneshian M, Guenther A, Wendel A, Hartung T, von Aulock S. 2006. *In vitro* pyrogen test
2344 for toxin or immunomodulatory drugs. *J Immunol Methods* 213(1-2):169-175.
- 2345 De Groote D, Zangerle PF, Gevaert Y, Fassotte MF, BeguinY, Noizat-Pirenne, Pirene J,
2346 Gathy R, Lopex M, Dehart I, Igo D, Baudrihaye M, Delacroix D, Franchimont. 1992.
2347 *Cytokine* 4(3):239-248.
- 2348 Dinarello CA, O'Conner JV, LoPreste G, Swift RL. 1984. Human leukocytic pyrogen test for
2349 detection of pyrogenic material in growth hormone produced by recombinant *Escherichia*
2350 *Coli*. *J. Clin. Microbiol.* 20:323-329.
- 2351 Dinarello CA. 1999. Cytokines as endogenous pyrogens. *Journal of Infectious Diseases*
2352 179(Suppl. 2).
- 2353 Duff GW, Atkins E. 1982. The detection of endotoxin by in vitro production of endogenous
2354 pyrogen: comparison with limulus amebocyte lysate gelation. *J Immunol Methods* 52(3):323-
2355 331.
- 2356 EP. 2005. Biological Tests, Supplement 5.6, 2.6.8 Pyrogens. In: European Pharmacopoeia.
2357 Strasbourg, France:Council of Europe.
- 2358 EPA. 2003a. Good Laboratory Practice Standards. Toxic Substances Control Act. 40 CFR
2359 792.
- 2360 EPA. 2003b. Good Laboratory Practice Standards. Federal Insecticide, Fungicide, and
2361 Rodenticide Act. 40 CFR 160.

- 2362 Eperon S, Jungi TW. 1996. The use of human monocytoid lines as indicators of endotoxin. J
2363 Immunol Methods 194:121-129.
- 2364 Eperon S, de Groote DD, Werner-Felmayer G, Jungi T. 1977. Human monocytoid cell lines
2365 as indicators of endotoxin: comparison with rabbit pyrogen and *Limulus* amoebocyte lysate
2366 assay. J Immunol Methods 207:135-145.
- 2367 FDA. 2003. Good Laboratory Practices for Nonclinical Laboratory Studies. 21 CFR 58.
- 2368 Fennrich S, Wendel A, Hartung T. 1999. New applications of the human whole blood
2369 pyrogen assay (PyroCheck). ALTEX 16(3):146-149.
- 2370 Greisman SE, Hornick RB. 1969. Comparative pyrogenic reactivity of rabbit and man to
2371 bacterial endotoxin. Proc Soc Exp Biol Med 131(4):1154-1158.
- 2372 Hansen E Wind, Christensen JD. 1990. Comparison of cultured human mononuclear cells,
2373 *Limulus* amoebocyte lysate and rabbits in the detection of pyrogens. J Clin Pharm Therap
2374 15:425-433.
- 2375 Hartung T, Wendel A. 1996. Detection of pyrogens using human whole blood. In Vitro
2376 Toxicol 9(4):353-359.
- 2377 Hartung T, Aaberge I, Berthold S, Carlin G, Charton E, Coecke S, et al. 2001. Novel pyrogen
2378 tests based on the human fever reaction. The report and recommendations of ECVAM
2379 Workshop 43. European Centre for the Validation of Alternative Methods. European Centre
2380 for the Validation of Alternative Methods. Altern Lab Anim 29(2):99-123.
- 2381 Hoffmann S, Peterbauer A, Schindler S, Fennrich S, Poole S, Mistry Y, et al. 2005.
2382 International validation of novel pyrogen tests based on human monocytoid cells. J Immunol
2383 Methods 298(1-2):161-173.
- 2384 Holzhütter HG, Archer G, Dami N, Lovell DP, Saltelli A, Sjostrom M. 1996.
2385 Recommendations of the application of biostatistical methods during the development and
2386 validation of alternative toxicological methods. ECVAM Biostatistics Task Force Report 1.
2387 ATLA, 24:511-530.
- 2388 ICCVAM. 1997. Validation and Regulatory Acceptance of Toxicological Test Methods: A
2389 Report of the ad hoc Interagency Coordinating Committee on the Validation of Alternative
2390 Methods. NIH Publication No.: 97-3981. Research Triangle Park:National Toxicology
2391 Program.
- 2392 ICCVAM. 2003. ICCVAM Guidelines for the Nomination and Submission of New, Revised,
2393 and Alternative Test Methods. NIH Publication No: 03-4508. Research Triangle
2394 Park:National Toxicology Program.
- 2395 ICCVAM Authorization Act of 2000. 2000. 42 U.S.C. § 2851-2, 2851-5. Available:
2396 <http://iccvam.niehs.nih.gov/about/overview.htm> [accessed 15 September 2006].
- 2397 JP. 2001. Bacterial Endotoxins Test. In: The Japanese Pharmacopoeia 14th ed. Part I.
2398 Shibuya, Shiubya-ku, Tokyo, Japan:The Society of Japanese Pharmacopoeia.
- 2399 Keene WR, Silberman HR, Landy M. 1961. Observations on the pyrogenic response to and
2400 its application to the bioassay of endotoxin. J Clin Invest 40(21):295-301.

- 2401 Moesby L, Jensen S, Hansen EW, Christensen JD. 1999. A comparative study of Mono Mac
2402 6 cells, isolated mononuclear cells and Limulus amoebocyte lysate assay in pyrogen testing.
2403 Int J Pharmaceutics 191:141-149.
- 2404 Moesby L, Hansen EW, Christensen JD. 2000. Endotoxin testing of proteins for parenteral
2405 administration using the Mono Mac 6 assay. Journal of Clinical Pharmacy and Therapeutics
2406 25(4):283.
- 2407 Moesby L, Hansen EW, Christensen JD, Høyer CH, Juhl GL, Olsen HB. 2005. Dry and
2408 moist heat sterilization cannot inactivate pyrogenicity of Gram positive microorganisms. Eur
2409 J Pharm Sci 26:318-323.
- 2410 Nakagawa Y, Maeda H, Murai T. 2002. Evaluation of the *in vitro* pyrogen test system based
2411 on proinflammatory cytokine release from human monocytes: Comparison with a human
2412 whole blood culture test system and with the rabbit pyrogen test. Clin Diag Lab Immunol
2413 9(3):588-597.
- 2414 Ochiai M, Kataoka M, Toyozumi H, Yamamoto A, Horiuchi Y. 2001. A quantitative in vitro
2415 assay for detecting biological activity of endotoxin using rabbit peripheral blood. Jpn J Infect
2416 Dis 54(4):131-136.
- 2417 OECD. 1998. OECD Series on Principles of Good Laboratory Practice and Compliance
2418 Monitoring Number 1: OECD principles on Good Laboratory Practice (as revised in 1997).
2419 ENV/MC/CHEM(98)17. OECD, Paris.
- 2420 Pool EJ, Johaar G, James S, Petersen I, Bouic P. 1998. The detection of pyrogens in blood
2421 products using an *ex vivo* whole blood culture assay. J Immunoassay 10(2&3):95-111.
- 2422 Pool EJ, Johaar G, James S, Petersen I, Bouic P. 1999. Differentiation between endotoxin
2423 and non-endotoxin pyrogens in human albumin solutions using an *ex vivo* whole blood
2424 culture assay. J Immunoassay 20(1&2):79-89.
- 2425 Poole S, Mistry Y, Ball C, Gaines Das RE, Opie LP, Tucker G, et al. 2003. A rapid 'one-
2426 plate' in vitro test for pyrogens. J Immunol Methods 274(1-2):209-220.
- 2427 Poole S, Thorpe R, Meager A, Gearing AJ. 1988. Assay of pyrogenic contamination in
2428 pharmaceuticals by cytokine release from monocytes. Dev Biol Stand 69:121-123.
- 2429 Russell WMS, Burch RL. 1959. The principles of humane experimental technique. Potters
2430 Bar, Herts, England:Universities Federation for Animal Welfare.
- 2431 Schindler S, Asmus S, von Aulock S, Wendel A, Hartung T, Fennrich S. 2004.
2432 Cryopreservation of human whole blood for pyrogenicity testing. J Immunol Methods
2433 294:89-100.
- 2434 Taktak YS, Selkirk S, Bristow AF, Carpenter A, Ball C, Rafferty B, Poole S. 1991. Assay of
2435 pyrogens by Interleukin-6 release from monocytic cell lines. J Pharm Pharmacol 43:578-582
- 2436 Tsuchiya S, Yamabe M, Yamaguchi Y, Kobayashi Y, Konno T, Tada K. 1980. Establishment
2437 and characterization of a human acute monocytic leukemia cell line (THP-1). Int J Cancer
2438 26(2):171-176.
- 2439 USP. 1999. The U.S. Pharmacopeia. USP23-NF28<151>. Ed. The U.S. Pharmacopeial
2440 Convention. Rockville, MD:The U.S. Pharmacopeial Convention.

- 2441 USP. 2005. The U.S. Pharmacopeia. USP28-NF23<85>. Ed. The U.S. Pharmacopeial
2442 Convention. Rockville, MD:The U.S. Pharmacopeial Convention.
- 2443 Westphal O. 1956. In: Polysaccharides in Biology (GF Springer, ed.). Trans. 2nd Conf. New
2444 York:Josiah Macy Jr. Foundation. 115
- 2445 Westphal O, Jann K. 1965. Bacterial Lipopolysaccharides: extraction with phenol-water and
2446 further applications of the procedure. In: Methods in carbohydrate chemistry (RL Whistler
2447 and ML Wolfrom, eds) Vol. 5. New York:Academic Press Inc, 83-91.
- 2448 Ziegler-Heitbrock HW, Thiel E, Futterer A, Herzog V, Wirtz A, Riethmuller G. 1988.
2449 Establishment of a human cell line (Mono Mac 6) with characteristics of mature monocytes.
2450 Int J Cancer 41(3):456-461.
- 2451
- 2452
- 2453

2454 **13.0 GLOSSARY¹³**

2455 **Amebocytes:** The blood cells of the horseshoe crab (*Limulus polyphemus* or *Tachypleus*
2456 *tridentatus*) that contain the active components of the reagent used in the bacterial endotoxin
2457 test (i.e., amebocyte lysate).

2458 **Assay²:** The experimental system used. Often used interchangeably with "test" and "test
2459 method."

2460 **Bacterial Endotoxin Test (BET)³:** A test used to quantify endotoxins of Gram-negative
2461 bacterial origin using amebocyte lysate from the horseshoe crab (*Limulus polyphemus* or
2462 *Tachypleus tridentatus*). There are two types of techniques for this test: the gel-clot
2463 techniques, which are based on gel formation, and the photometric techniques. The
2464 photometric techniques include the turbidimetric technique, which is based on the
2465 development of turbidity after cleavage of an endogenous substrate, and a chromogenic
2466 method, which is based on the development of color after cleavage of a synthetic peptide-
2467 chromogen complex.

2468 **Coded substances:** Substances labeled by code rather than name so that they can be tested
2469 and evaluated without knowledge of their identity or anticipation of test results. Coded
2470 substances are used to avoid intentional or unintentional bias when evaluating laboratory or
2471 test method performance.

2472 **Coefficient of variation:** A statistical representation of the precision of a test. It is expressed
2473 as a percentage and is calculated as follows:

2474

$$2475 \left(\frac{\textit{standard deviation}}{\textit{mean}} \right) \times 100\%$$

2476

¹³ The definitions in this Glossary are restricted to their uses with respect to the Rabbit Pyrogen Test, the *in vitro* pyrogenicity test methods included in this BRD, and the Limulus Amebocyte Lysate (LAL) assay.

²From ICCVAM (2003)

³From USP (2005)

2477 **Concordance²**: The proportion of all substances tested that are correctly classified as
2478 positive or negative. It is a measure of test method performance and one aspect of
2479 “relevance.” The term is often used interchangeably with “accuracy” (see also “two-by-two”
2480 table). Concordance is highly dependent on the prevalence of positives in the population
2481 being examined.

2482 **Endogenous pyrogens**: Various cytokines including interleukins (e.g., IL-1 α , IL-1 β), tumor
2483 necrosis factor (i.e., TNF- α , TNF- β), and interferon (IFN- γ) released from leukocytes in
2484 response to external stimuli (e.g., endotoxin) capable of causing an increase in body
2485 temperature above the normal level.

2486 **Endotoxin Limit Concentration (ELC)**: The concentration at which endotoxin is
2487 considered to be pyrogenic. It is expressed as the ratio of the threshold pyrogen dose (K) and
2488 the rabbit pyrogen test dose or the maximum human dose administered on a weight (kg) basis
2489 in 1 hr (M) defined as K/M. The ELC varies based on M.

2490 • The FDA ELC for non-intrathecal medical devices is 0.5 EU/mL

2491 • The FDA ELC for intrathecal medical devices is 0.06 EU/mL

2492 **Endpoint²**: The biological or chemical process, response, or effect assessed by a test method.

2493 **False negative²**: A substance incorrectly identified as negative by a test method.

2494 **False negative rate²**: The proportion of all positive substances falsely identified by a test
2495 method as negative (see “two-by-two” table). It is one indicator of test method accuracy.

2496 **False positive²**: A substance incorrectly identified as positive by a test method.

2497 **False positive rate²**: The proportion of all negative substances that are falsely identified by
2498 a test method as positive (see “two-by-two” table). It is one indicator of test method
2499 accuracy.

2500 **Fever**: Elevation of body temperature above the normal level.

2501 **Good Laboratory Practices (GLP)²**: Regulations promulgated by the U.S. Food and Drug
2502 Administration and the U.S. Environmental Protection Agency, and principles and
2503 procedures adopted by the Organization for Economic Cooperation and Development and

2504 Japanese authorities that describe record keeping and quality assurance procedures for
2505 laboratory records that will be the basis for data submissions to national regulatory agencies.

2506 **Hazard²**: The potential for an adverse health or ecological effect. A hazard potential results
2507 only if an exposure occurs that leads to the possibility of an adverse effect being manifested.

2508 **Interlaboratory reproducibility²**: A measure of whether different qualified laboratories
2509 using the same protocol and test substances can produce qualitatively and quantitatively
2510 similar results. Interlaboratory reproducibility is determined during the prevalidation and
2511 validation processes and indicates the extent to which a test method can be transferred
2512 successfully among laboratories.

2513 **Intralaboratory repeatability²**: The closeness of agreement between test results obtained
2514 within a single laboratory when the procedure is performed on the same substance under
2515 identical conditions within a given time period.

2516 **Intralaboratory reproducibility²**: The first stage of validation; a determination of whether
2517 qualified people within the same laboratory can successfully replicate results using a specific
2518 test protocol at different times.

2519 **In Vitro**: In glass. Refers to assays that are carried out in an artificial system (e.g., in a test
2520 tube or Petri-dish) and typically use single-cell organisms, cultured cells, cell-free extracts, or
2521 purified cellular components.

2522 **In Vivo**: In the living organism. Refers to assays performed in multi-cellular organisms.

2523 **Lipopolysaccharide (LPS)**: A complex of lipid and carbohydrate (endotoxin) released from
2524 the cell walls of Gram-negative organisms that is pyrogenic and capable of producing septic
2525 shock.

2526 **Lipoteichoic acid (LTA)**: A polyol phosphate polymer bearing a strong negative charge that
2527 is covalently linked to the peptidoglycan in Gram-positive bacteria. It is strongly antigenic,
2528 but is generally absent in Gram-negative bacteria and therefore is considered to be the
2529 primary pyrogenic component of Gram-positive bacteria.

2530 **Minimum Valid Concentration (MVC):** The concentration of a product when it is diluted
2531 to the maximum valid dilution (MVD) expressed as $\lambda M/K$, where:

2532 • λ = The sensitivity of the LAL reagent used expressed as EU/mL. The value
2533 varies with the method used. For the gel-clot method it is the labeled LAL
2534 sensitivity (EU/mL). For the chromogenic, turbidometric, or kinetic-
2535 turbidometric LAL test methods it is the lowest point used in the standard
2536 curve.

2537 • M = The maximum human dose for pyrogenicity administered on a weight
2538 basis (kg) in 1 hr, or the rabbit pyrogen test dose (whichever is larger). It is
2539 one of the variables used to define the Endotoxin Limit Concentration (ELC)
2540 defined as the ratio of K/M , where K is the threshold pyrogen dose in rabbits
2541 or humans.

2542 • K = See **threshold pyrogen dose**.

2543 **Maximum Valid Dilution (MVD):** When a USP Endotoxin Limit Concentration (ELC) is
2544 defined, the MVD is the ratio of the product of the ELC and the product potency to the LAL
2545 reagent sensitivity (λ) expressed as $([ELC \times \text{Product Potency}]/\lambda)$. If there is no official USP
2546 ELC defined, then the MVD is the ratio of the Product Potency/Minimum Valid
2547 Concentration (MVC).

2548 **Monocytoid cells:** Cells obtained from peripheral blood or grown in culture that
2549 phenotypically resemble monocytes or macrophages.

2550 **Negative control:** An untreated sample containing all components of a test system, except
2551 the test substance solvent, which is replaced with a known non-reactive material, such as
2552 water. This sample is processed with test substance-treated samples and other control
2553 samples to determine whether the solvent interacts with the test system.

2554 **Negative predictivity²:** The proportion of correct negative responses among substances
2555 testing negative by a test method (see “two-by-two” table). It is one indicator of test method
2556 accuracy. Negative predictivity is a function of the sensitivity of the test method and the
2557 prevalence of negatives among the substances tested.

- 2558 **Parenteral:** Introduction into the body by some means other than through the
2559 gastrointestinal tract; referring particularly to intravenous (i.v.), intramuscular (i.m.),
2560 subcutaneous (s.c.), or intrathecal (i.t.) injection.
- 2561 **Performance²:** The accuracy and reliability characteristics of a test method (see “accuracy,
2562 reliability”).
- 2563 **pH:** A measure of the acidity or alkalinity of a solution. A pH of 7.0 is neutral; higher pHs
2564 are alkaline, lower pHs are acidic.
- 2565 **Positive control:** A sample containing all components of a test system and treated with a
2566 substance known to induce a positive response, which is processed with the test substance-
2567 treated and other control samples to demonstrate the sensitivity of each experiment and to
2568 allow for an assessment of variability in the conduct of the assay over time.
- 2569 **Positive predictivity²:** The proportion of correct positive responses among substances
2570 testing positive by a test method (see “two-by-two” table). It is one indicator of test method
2571 accuracy. Positive predictivity is a function of the sensitivity of the test method and the
2572 prevalence of positives among the substances tested.
- 2573 **Prevalence²:** The proportion of positives in the population of substances tested (see “two-by-
2574 two” table).
- 2575 **Protocol²:** The precise, step-by-step description of a test method, including a listing of all
2576 necessary reagents, criteria and procedures for evaluation of the test data.
- 2577 **Pyrogen:** A substance that causes a rise in body temperature above normal or that produces
2578 a fever. Gram-negative, Gram-positive, and acid-fast bacteria, molds, viruses, and yeast and
2579 some of their cellular constituents are pyrogenic.
- 2580 **Quality assurance²:** A management process by which adherence to laboratory testing
2581 standards, requirements, and record keeping procedures is assessed independently by
2582 individuals other than those performing the testing.
- 2583 **Rabbit Pyrogen Test (RPT)³:** A test designed to limit to an acceptable level the risks of
2584 febrile reaction in the patient to the administration, by injection, or the product concerned.

2585 The test involves measuring the rise in temperature of rabbits following the intravenous
2586 injection of a test solution.

2587 **Reduction alternative²:** A new or modified test method that reduces the number of animals
2588 required.

2589 **Reference test method²:** The accepted *in vivo* test method used for regulatory purposes to
2590 evaluate the potential of a test substance to be hazardous to the species of interest.

2591 **Refinement alternative²:** A new or modified test method that refines procedures to lessen
2592 or eliminate pain or distress in animals or enhances animal well-being.

2593 **Relevance²:** The extent to which a test method correctly predicts or measures the biological
2594 effect of interest in humans or another species of interest. Relevance incorporates
2595 consideration of the “accuracy” or “concordance” of a test method.

2596 **Reliability²:** A measure of the degree to which a test method can be performed reproducibly
2597 within and among laboratories over time. It is assessed by calculating intra- and
2598 interlaboratory reproducibility and intralaboratory repeatability.

2599 **Replacement alternative²:** A new or modified test method that replaces animals with non-
2600 animal systems or one animal species with a phylogenetically lower one (e.g., a mammal with
2601 an invertebrate).

2602 **Reproducibility²:** The consistency of individual test results obtained in a single laboratory
2603 (intralaboratory reproducibility) or in different laboratories (interlaboratory reproducibility)
2604 using the same protocol and test substances (see intra- and interlaboratory reproducibility).

2605 **Sensitivity²:** The proportion of all positive substances that are classified correctly as
2606 positive in a test method. It is a measure of test method accuracy (see “two-by-two” table).

2607 **Specificity²:** The proportion of all negative substances that are classified correctly as
2608 negative in a test method. It is a measure of test method accuracy (see “two-by-two” table).

2609 **Test²:** The experimental system used; used interchangeably with “test method” and “assay.”

2610 **Test method²:** A process or procedure used to obtain information on the characteristics of a
2611 substance or agent. Toxicological test methods generate information regarding the ability of a
2612 substance or agent to produce a specified biological effect under specified conditions. Used
2613 interchangeably with “test” and “assay.” See also “validated test method” and “reference
2614 test.”

2615 **Test method component:** Structural, functional, and procedural elements of a test method
2616 that are used to develop the test method protocol. These components include unique
2617 characteristics of the test method, critical procedural details, and quality control measures.

2618 **Threshold pyrogen dose:** The dose level at which a product is considered to be pyrogenic or
2619 non-pyrogenic. It is one of the variables (K) used to calculate the Endotoxin Limit
2620 Concentration (ELC) defined as K/M , where M is the rabbit pyrogen test dose or the
2621 maximum human dose administered in 1 hr (whichever is larger).

- 2622 • The threshold pyrogen dose for non-intrathecal use in rabbits and humans is
2623 5.0 EU/kg
- 2624 • The threshold pyrogen dose for intrathecal use in rabbits and humans is 0.2
2625 EU/kg

2626 **Tiered testing:** A testing strategy where all existing information on a test substance is
2627 reviewed, in a specified order, prior to *in vivo* testing. If the irritancy potential of a test
2628 substance can be assigned, based on the existing information, no additional testing is
2629 required. If the irritancy potential of a test substance cannot be assigned, based on the
2630 existing information, a step-wise animal testing procedure is performed until an unequivocal
2631 classification can be made.

2632 **Transferability²:** The ability of a test method or procedure to be accurately and reliably
2633 performed in different, competent laboratories.

2634 **Two-by-two table²:** The two-by-two table can be used for calculating accuracy (concordance)
2635 $([a+d]/[a+b+c+d])$, negative predictivity $(d/[c+d])$, positive predictivity $(a/[a+b])$, prevalence
2636 $([a+c]/[a+b+c+d])$, sensitivity $(a/[a+c])$, specificity $(d/[b+d])$, false positive rate $(b/[b+d])$,
2637 and false negative rate $(c/[a+c])$.

		NEW TEST OUTCOME		
		Positive	Negative	Total
Reference Test Outcome	Positive	a	c	a + c
	Negative	b	d	b + d
	Total	a + b	c + d	a + b + c + d

2638

2639 **Validated test method²:** An accepted test method for which validation studies have been
 2640 completed to determine the relevance and reliability of this method for a specific proposed
 2641 use.

2642 **Validation²:** The process by which the reliability and relevance of a procedure are
 2643 established for a specific purpose.

2644 **Weight of evidence (process):** The strengths and weaknesses of a collection of information
 2645 are used as the basis for a conclusion that may not be evident from the individual data.

2646

2647

2648

2649

Appendix A

2650

2651 **Individual BRDs Submitted by ECVAM on Five *In Vitro* Pyrogenicity Test Methods**

2652

2653

2654 **The Human Whole Blood (WB)/IL-1 *In Vitro* Pyrogen Test:**

2655 **Application of Cryopreserved Human WB..... TAB A1**

2656 **An Alternative *In Vitro* Pyrogenicity Test Using the Monocytoid**

2657 **Cell Line Mono Mac 6 (MM6)/IL-6..... TAB A2**

2658 **The Human Peripheral Blood Mononuclear Cell (PBMC)/IL-6**

2659 ***In Vitro* Pyrogen Test TAB A3**

2660 **The Human WB/IL-1 *In Vitro* Pyrogen Test..... TAB A4**

2661 **The Human WB/IL-6 *In Vitro* Pyrogen Test..... TAB A5**

2662

2663

2664

2665

2665

2666

2667

2668

2669

2670

2671

2672

2673

2674

2675

2676

2677

2678 ***[This Page Intentionally Left Blank]***

2679

2680

2680

2681

2682

2683

2684

Appendix B

2685

2686

ECVAM Response to ICCVAM Questions

2687

2688 **ECVAM Information and Additional Unpublished Data.....TAB B**

2689

2690

2691

2691

2692

2693

2694

2695

2696

2697

2698

2699

2700

2701

2702

2703

2704 ***[This Page Intentionally Left Blank]***

2705

2706

2707
2708
2709
2710
2711
2712
2713
2714
2715
2716
2717
2718

Appendix C

Supplemental Information

Key References TAB C1
Guidelines for Pyrogenicity Testing TAB C2

2718

2719

2720

2721

2722

2723

2724

2725

2726

2727

2728

2729

2730

2731

[This Page Intentionally Left Blank]

2732

Appendix C1

2733

2734

Key References

2735

2736 **List of key references** (see **Section 12.0** for a complete list of references):

2737

2738 Bleeker WK, de Groot EM, den Boer PJ, BiesselsPTM, Asrden LA, Bakker JC. 1994.

2739 Measurement of interleukin-6 production by monocytes for *in vitro* safety testing of

2740 hemoglobin solutions. *Art Cells Blood Subs Immob Biotech* 22(3):835-840.

2741

2742 Dinarello CA, O'Connor JV, LoPreste G, Swift RL. 1984. Human leukocyte pyrogen test for

2743 detection of pyrogenic material in growth hormone produced by recombinant *Escherichia*

2744 *coli*. *J Clin Microbiol* 20: 323-329.

2745

2746 Hartung T, Wendel A. 1996. Detection of pyrogens using human whole blood. *In Vitro*

2747 *Toxicol* 9:353-359.

2748

2749 Hartung T, Aaberge I, Berthold S, Carlin G, Charton E, Coecke S, et al. 2001. Novel pyrogen

2750 tests based on the human fever reaction. The report and recommendations of ECVAM

2751 Workshop 43. European Centre for the Validation of Alternative Methods. European Centre

2752 for the Validation of Alternative Methods. *Altern Lab Anim* 29(2):99-123.

2753

2754 Hoffmann, S., Luederitz-Puechel, U., Montag-Lessing, T., Hartung, T. 2005. Optimisation of

2755 pyrogen testing in parenterals according to different pharmacopoeias by probabilistic

2756 modeling. *J of Endotoxin Res* 11(1):26-31.

2757

2758 Hoffmann S, Peterbauer A, Schindler S, Fennrich S, Poole S, Mistry Y, Montag-Lessing T,

2759 Spreitzer I, Löschner B, van Aalderen M, Bos R, Gommer M, Nibbeling R, Werner-

2760 Felmayer G, Loitzl P, Jungi T, Brcic M, Brügger P, Frey E, Bowe G, Casado J, Coecke S, de

2761 Lange J, Mogster B, Næss LM, Aaberge IS, Wendel A, Hartung T. 2005. International

2762 validation of novel pyrogen tests based on human monocytoid cells. *J Immunol Methods*

2763 298:161-173.

2764 Poole S, Mistry Y, Ball C, Gaines Das RE, Opie LP, Tucker G, et al. 2003. A rapid 'one-

2765 plate' *in vitro* test for pyrogens. *J Immunol Methods* 274(1-2):209-220.

2766

2767 Schindler S, Asmus S, von Aulock S, Wendel A, Hartung T, Fennrich S. 2004.

2768 Cyropreservation of human whole blood for pyrogencity testing. *J Immunol Methods*

2769 294:89-100.

2770

2771 Taktak YS, Selkirk S, Bristow AF, Carpenter A, Ball C, Rafferty B, Poole S. 1991. Assay of

2772 pyrogens by Interleukin-6 release from monocytic cell lines. *J Pharm Pharmacol* 43:578-582.

2773

2773

2774

2775

2776

2777

2778

2779

2780

2781

2782

2783

2784

2785

2786 ***[This Page Intentionally Left Blank]***

2787

2788

2789

APPENDIX C2

2790

2791

2792

Guidelines for Pyrogenicity Testing

2793

2794

2795 **List of included guidelines:**

2796

2797 United States Pharmacopeia

2798 Bacterial Endotoxins Test

2799 Pyrogen Test

2800 Biologics

2801

2802 European Pharmacopoeia

2803 Pyrogens

2804 Bacterial Endotoxins

2805 Parenteral Preparations

2806

2806

2807

2808

2809

2810

2811

2812

2813

2814

2815

2816

2817

2818

2819 *[This Page Intentionally Left Blank]*

2820

2821