

APPENDIX I
OHIO WATER MICROBIOLOGY LABORATORY
ANALYSIS OF *CLOSTRIDIUM PERFRINGENS* IN ENVIRONMENTAL WATER SAMPLES

1. Prepare mCP basal agar as follows (or order Oxoid CM0992):
 - 900 mL reagent water
 - 30.0 g tryptose
 - 20.0 g yeast extract
 - 5.0 g sucrose
 - 1.0 g L-cysteine hydrochloride
 - 0.1 g MgSO₄·7H₂O
 - 0.04 g bromcresol purple
 - 15.0 g Bacto agar
2. pH adjust agar to 7.6 (with 1 N NaOH)
3. Label each batch of media with an assigned number and record in the media QC log book.
4. Dispense 100-mL aliquots into dilution bottles and autoclave for 15 minutes at 121°C. Store bottles in the refrigerator for up to 6 months.
5. To prepare plates, melt the agar using a beaker with water on a hot plate or by placing in the autoclave for a 5-minute cycle. After the agar is tempered and before pouring plates, add the following ingredients to each 100-mL bottle of agar:
 - 0.04 g D-cycloserine
 - 0.0025 g polymyxin-B sulfate
 - 0.2 mL 4.5% FeCl₃·6H₂O solution
(filter-sterilize 0.45 g FeCl₃ in 10 mL reagent water, store at 4°C for up to 2 months)
 - 2.0 mL 0.5% Phenolphthalein diphosphate solution*
(filter-sterilize 0.5 g phenolphthalein diphosphate in 100 mL reagent water, store at 4°C)
 - 8.0 mL 0.075% indoxyl β-D glucoside solution**
(dissolve 0.006 g indoxyl β-D glucoside in 8.0 mL sterile reagent water)

*Alternatively, 0.01 g phenolphthalein diphosphate can be added directly without making the solution.
**Alternatively, 0.006 g indoxyl β-D glucoside can be added directly without making the solution.

To make a large batch of agar, combine the 100-mL bottles into a sterile flask, calculate the amounts of additives for the total volume of agar, and add each additive in one batch.
6. Label each sub-batch with an assigned number and record in media logbook.
7. Store plates inverted in the refrigerator for up to 1 month.

3/12/2007

M:\QC manual\Appendix I - Clostridium perfringens.doc

8. Analyze samples by membrane filtration using appropriate volumes. For turbid surface water, use 3-, 10-, and 30-mL volumes. For clean surface water, use 10-, 30-, and 50-mL volumes. Record media batch and sub-batch numbers on the Results Worksheet.

Note: To quantify the concentration of spores only, heat treat the sample at 60°C for 15 min prior to filtration.

9. Run a filter blank using buffer water before filtering each sample.
10. Plate a positive control (10 and 30 mL of 10⁻² dilution raw sewage) for every 20th sample or when a new analyst is processing or reading plates.
11. Incubate plates anaerobically in a GasPak bag or chamber at 44.5°C for 24 hours for the analysis of vegetative bacteria or at 41°C for 24 hours for the analysis of spores.
12. In the fume hood using protective gloves and forceps, place membrane filters with straw-yellow colonies onto cellulose pads saturated with NH₄OH. Wait for approximately 15 seconds, then examine for positive colonies.
13. Magenta colonies that are approximately 1 to 2 mm in diameter are enumerated as *C. perfringens*. New analysts may compare size and color of colonies to those in the positive control.
14. Calculate the number of colony-forming units per 100-mL (CFU/100 mL) sample using the following equation:

$$CFU/100mL = \frac{Colony\ count * 100}{Volume\ plated(mL)}$$