



# Laboratory Methods for the Diagnosis of Epidemic Dysentery and Cholera

Centers for Disease Control and Prevention  
Atlanta, Georgia 1999



# **Laboratory Methods for the Diagnosis of Epidemic Dysentery and Cholera**

**Centers for Disease Control and Prevention  
Atlanta, Georgia 1999**

This manual was prepared by the National Center for Infectious Diseases (NCID), Centers for Disease Control and Prevention (CDC), Atlanta, Georgia, USA, in cooperation with the World Health Organization Regional Office for Africa, (WHO/AFRO) Harare, Zimbabwe.

Jeffrey P. Koplan, M.D., M.P.H., Director, CDC

James M. Hughes, M.D., Director, NCID, CDC

Mitchell L. Cohen, M.D., Director, Division of Bacterial and Mycotic Diseases, NCID, CDC

Ebrahim Malek Samba, M.B.,B.S., Regional Director, WHO/AFRO

Antoine Bonaventure Kabore, M.D., M.P.H., Director Division for Prevention and Control of Communicable Diseases, WHO/AFRO

***The following CDC staff members prepared this report:***

Cheryl A. Bopp, M.S.

Allen A. Ries, M.D., M.P.H.

Joy G. Wells, M.S.

***Production:***

J. Kevin Burlison, Graphics

James D. Gathany, Photography

Lynne McIntyre, M.A.L.S., Editor

**Cover:** From top, *Escherichia coli* O157:H7 on sorbitol MacConkey agar, *Vibrio cholerae* O1 on TCBS agar, and *Shigella flexneri* on xylose lysine desoxycholate agar.

# Acknowledgments

Funding for the development of this manual was provided by the U.S. Agency for International Development, Bureau for Africa, Office of Sustainable Development.

This manual was developed as a result of a joint effort by the World Health Organization Regional Office for Africa, WHO Headquarters, and the Centers for Disease Control and Prevention as part of the activities of the WHO Global Task Force on Cholera Control. In particular, the staff of the project for Improving Preparedness and Response to Cholera and Other Epidemic Diarrhoeal Diseases in Southern Africa have worked closely with many laboratorians and epidemiologists in southern Africa to develop an integrated approach to the laboratory diagnosis of cholera and dysentery upon which this manual is based.

We also appreciate the valuable assistance of Ms. Katherine Greene, Dr. Eric Mintz, Ms. Nancy Puhr, Dr. Nancy Strockbine, Dr. Robert Tauxe, and Dr. Fred Tenover, Centers for Disease Control and Prevention, Atlanta, Georgia, USA; Dr. Lianne Kuppens, World Health Organization, Geneva, Switzerland; Dr. Elizabeth Mason, World Health Organization, Harare, Zimbabwe; and Ms. Catherine Mundy, Liverpool School of Tropical Medicine, Liverpool, UK.

# Introduction

Cholera and dysentery have afflicted humankind for centuries. The epidemics they cause have affected the outcome of wars and the fates of countries. In much of the world, epidemic cholera and dysentery are uncommon, but during the past decade these two diseases have re-emerged as causes of significant morbidity and mortality in many developing countries.

Only a few pathogens cause epidemic diarrhea, although there are many that cause sporadic diarrhea. In developing countries, two etiologic agents are responsible for most epidemic diarrhea: toxigenic *Vibrio cholerae* serogroup O1, which causes watery diarrhea, and *Shigella dysenteriae* serotype 1, which causes bloody diarrhea. Recently, two additional organisms have emerged to cause epidemic diarrhea, *Vibrio cholerae* serogroup O139, which causes watery diarrhea, and *Escherichia coli* serotype O157:H7, which causes bloody diarrhea. The latter is a common agent of diarrhea only in developed countries.

This manual focuses on the epidemiology of these four organisms and the laboratory methods used to identify them and to test their susceptibility to antimicrobial agents in the epidemic setting. The laboratory techniques and study methodology described provide accurate and useful information for the control of epidemics using a minimum of resources. The manual emphasizes coordination of the activities of the microbiologist and the epidemiologist in order to obtain information that can be generalized to develop effective treatment policies for these epidemic diarrheal diseases. It encourages focused studies to determine the organisms causing epidemics and their antimicrobial susceptibility patterns rather than relying on random information that may not accurately represent a situation.

Often the countries that face the challenge of responding to an epidemic are those with the least resources. Therefore, the microbiology laboratory must use its resources wisely in order to have the greatest impact on reducing morbidity and mortality during an epidemic. There may be several ways to reach the end result of identifying the organism causing the outbreak or the epidemic. Often, however, a small added benefit requires a much larger expenditure of materials and time. In this manual this problem is addressed specifically. The procedures described are not new; most have been used for a number of years. However, these procedures were specifically selected for testing specimens from outbreaks rather than for general use in a clinical microbiology laboratory. The selected procedures minimize the materials needed by the laboratory while deriving the most useful information.

# Table of Contents

## Acknowledgments

## Introduction

<b>Chapter 1.</b>	<b>The Public Health Role of Clinical Laboratories</b> . . . . .	<b>1</b>
	A. Epidemic Diarrhea . . . . .	1
	B. Public Health Role of the Laboratory . . . . .	2
<b>Chapter 2.</b>	<b>Collection and Transport of Fecal Specimens</b> . . . . .	<b>7</b>
	A. Collection of Stool . . . . .	7
	B. Preparing Specimens for Shipment . . . . .	10
<b>Chapter 3.</b>	<b>Epidemiology of Dysentery Caused by <i>Shigella</i></b> . . . . .	<b>13</b>
	A. Epidemiology of <i>Shigella</i> . . . . .	13
	B. Clinical Manifestations . . . . .	14
	C. Treatment . . . . .	14
<b>Chapter 4.</b>	<b>Isolation and Identification of <i>Shigella</i></b> . . . . .	<b>17</b>
	A. Isolation Methods . . . . .	17
	B. Biochemical Screening Tests . . . . .	20
	C. Serologic Identification of <i>Shigella</i> . . . . .	26
	D. Media for Isolation and Identification of <i>Shigella</i> . . . . .	28
<b>Chapter 5.</b>	<b>Etiology and Epidemiology of Cholera</b> . . . . .	<b>37</b>
	A. Historical Background . . . . .	37
	B. Clinical Manifestations . . . . .	38
	C. Treatment . . . . .	39
	D. Epidemiology . . . . .	39
	E. Cholera Vaccine . . . . .	40
<b>Chapter 6.</b>	<b>Isolation and Identification of <i>Vibrio cholerae</i></b> <b>Serogroups O1 and O139</b> . . . . .	<b>41</b>
	A. Isolation Methods . . . . .	41
	B. Serologic Identification of <i>V. cholerae</i> O1 and O139 . . . . .	49
	C. Media and Reagents for <i>V. cholerae</i> . . . . .	51
<b>Chapter 7.</b>	<b>Epidemiology of <i>Escherichia coli</i> Serotype O157:H7</b> . . . . .	<b>55</b>
<b>Chapter 8.</b>	<b>Isolation and Identification of <i>Escherichia coli</i></b> <b>Serotype O157:H7</b> . . . . .	<b>57</b>
	A. Isolation and Identification Methods . . . . .	57
	B. Preparation and Quality Control of Sorbitol-MacConkey Agar . . . . .	60
<b>Chapter 9.</b>	<b>Antimicrobial Susceptibility Testing</b> <b>(Agar Disk Diffusion Method)</b> . . . . .	<b>61</b>
	A. Considerations for Antimicrobial Susceptibility Testing . . . . .	61
	B. Procedure for Agar Disk Diffusion . . . . .	61
	C. Special Considerations for Susceptibility Testing of <i>Vibrio cholerae</i> . . . . .	71
	D. Preparation and Quality Control of Media and Reagents . . . . .	71

<b>Chapter 10.</b>	<b>Storage of Isolates</b> . . . . .	<b>75</b>
	A. Short-Term Storage . . . . .	75
	B. Long-Term Storage . . . . .	75
<b>Chapter 11.</b>	<b>Quality Control of Media and Reagents</b> . . . . .	<b>77</b>
	A. Quality Control of Media . . . . .	77
	B. Quality Control of Reagents . . . . .	78
	C. Advantages of Centralized Acquisition of Media and Reagents . . . . .	79
<b>Chapter 12.</b>	<b>Standard Safety Practices in the Microbiology Laboratory</b> . . . . .	<b>81</b>
	A. Standard Microbiological Safety Practices . . . . .	81
	B. Special Practices . . . . .	83
	C. Protective Clothing and Equipment . . . . .	84
<b>Chapter 13.</b>	<b>Packing and Shipping of Clinical Specimens and Etiologic Agents</b> . . . . .	<b>87</b>
	A. Preparation for Transport of Infectious Specimens and Cultures . . . . .	87
	B. Transport and Shipment of Cultures and Specimens . . . . .	87
<b>Annex A:</b>	<b>Diagnostic Supplies Needed for 1 Year for Laboratory Confirmation of Outbreaks and for Laboratory-Based Surveillance for <i>Vibrio cholerae</i> O1/O139 Antimicrobial Susceptibility</b> . . . . .	<b>93</b>
<b>Annex B.</b>	<b>Supplies Needed for Laboratory Identification of <i>Shigella dysenteriae</i> 1 During an Outbreak</b> . . . . .	<b>95</b>
<b>Annex C.</b>	<b>Guidelines for Establishing a Public Health Laboratory Network for Cholera Control</b> . . . . .	<b>97</b>
<b>Annex D.</b>	<b>International Reference Laboratories</b> . . . . .	<b>101</b>
<b>Annex E.</b>	<b>Designing a Survey to Examine Antimicrobial Susceptibility of Organisms Causing Epidemic Diarrhea</b> . . . . .	<b>103</b>
<b>Annex F.</b>	<b>Stool Specimen Data Sheet for Epidemic Diarrhea</b> . . . . .	<b>105</b>
<b>Annex G.</b>	<b>Most Frequently Encountered Reactions in Screening Biochemicals</b> . . . . .	<b>106</b>
<b>Annex H.</b>	<b>Diagnostic Laboratory Supplies for Isolation and Presumptive Identification of <i>Escherichia coli</i> O157:H7 During an Outbreak (Sufficient for 100 Specimens)</b> . . . . .	<b>107</b>