

Module 2: Safety Precautions for Tuberculosis Microscopy

Learning Objectives

At the end of this module, you will be able to

- Explain airborne transmission of TB
- Describe the risks when collecting sputum
- Describe personal health and safety practices
- Describe why there should be 3 areas in the TB laboratory
- Describe methods for the disposal of contaminated material
- Describe chemical safety precautions in the laboratory

Content Overview

- **Transmission of TB bacilli**
- **Proper collection of sputum**
- **Laboratory Arrangement**
- **Safe work practices in the TB microscopy laboratory**
- **Safe disposal of infectious waste**
- **Chemical Safety**

Importance of Laboratory Safety

- Prevents laboratory-acquired infection
- Specialized equipment may aid good laboratory practice but does NOT replace it

Transmission of TB Bacilli

- ***Mycobacterium tuberculosis* is almost always transmitted by patients with active pulmonary disease**
 - TB patient expels bacilli in small droplets of respiratory secretions
 - Secretions quickly evaporate leaving “droplet nuclei” less than 5 μm in diameter
 - **Droplet nuclei of this size containing 1–3 bacilli can remain suspended indefinitely in the air**
 - **Following inhalation, droplet nuclei are able to reach deep into the lungs to produce infection**

Aerosol Formation: Spread of droplets




Coughing

Singing

Sneezing

Talking

Relative Risk from Exposure to Infectious TB Case



Safe Specimen Collection

- **If a patient presents to the laboratory, and is coughing, ask the patient to cover his or her mouth**
- **Collect specimens outside where**
 - **air movement will dilute droplet nuclei**
 - **sunlight will rapidly inactivate TB bacilli**
- **Stand clear of patients when specimens are collected**

Never Stand in Front of the Patient During Collection



Safety Practices: Airflow

- Establish airflow in working areas that will direct potentially infectious particles away from personnel
 - air is exhausted into a remote area
 - extraction fans can be useful
 - **laboratory coats will not help against infection with TB**

Personal Protective Equipment



- **Masks**
- **Gloves**
- **Lab Coats**



Personnel Protective Equipment: Masks



- Surgical masks do not filter out infectious droplet nuclei



Personnel Protective Equipment: Gloves

- Not necessary to perform sputum microscopy
- Lack of availability does NOT mean that sputum smears cannot be prepared!
- Wearing gloves can give technicians a false sense of safety
- Do not re-use gloves
- Remove gloves before using or operating equipment to avoid contamination (e.g., microscope or telephone)
- Never wear gloves outside the laboratory



Personnel Protective Equipment: Laboratory Coats

- A lack of laboratory coats does NOT mean that sputum microscopy cannot be performed
- If laboratory coats are worn:
 - Leave at worksite for organization to clean
 - Tie at the back, not the front
 - Use appropriate size
 - DO NOT wear outside of the laboratory

Laboratory Design: Working Area

- Record keeping and storage [Clean]
- Smear preparation and staining [Dirty]
- Performing microscopy [Clean]

Laboratory Design: Bench Area

- Record keeping and storage [Clean]
 - area for entering data into the register and for storing slides
- Smear preparation and staining [Dirty]
 - should be well lit area
 - near an open window to ensure adequate ventilation during smear preparation
 - sink with running water also required

Laboratory Design: Microscopy Bench

- Microscopy area [Clean]
 - Use a flat bench or stable table for microscope
- Place microscope in area where lighting is subdued, preferably
- If no electricity is available:
 - Use daylight as the light source
 - Place the microscope directly in front of a window

Biological Safety Cabinets

- Not required for performing sputum smear microscopy
- Are only necessary for cultures and drug susceptibility testing (DST)
 - large amount of organism being handled
- BSC are very expensive to purchase and maintain
 - Require yearly maintenance
 - Require filter replacement



Waste Disposal

- Discard specimens by one of the following methods:
 - burning
 - burying
 - autoclaving


**All materials used
should be considered
contaminated!**

Chemical Safety

- Alcohols are flammable; avoid flame
- Phenol is toxic
 - Avoid direct contact with the skin or mucus membranes
 - Reduce exposure to fumes, work in ventilated area
- Acids are corrosive
 - Use personal protection equipment while handling acids
 - Work in ventilated area
 - Avoid direct contact with the skin, mucus membranes, clothes and paper

Chemical Safety: Handling Acids

**ALWAYS ADD
ACID TO WATER**



**NEVER ADD
WATER TO ACID**



Summary

- How is TB transmitted from person to person?
- What precautions must be taken when collecting sputum specimens?
- What universal precaution must you take when handling specimens?
- Why do lab coats, gloves and surgical masks offer little protection?
- What precautions must you take when handling acids?