

United States Department of Agriculture Food Safety and Inspection Service Office of Public Health Science Laboratory QA/QC Division 950 College Station Road Athens, GA 30605

### **Laboratory Guidebook Notice of Change**

Chapter new, **revised**, or archived: MLG 5.04

Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat

Products

Effective Date: 1/28/08

Description and purpose of change(s):

FSIS validated the use of a replacement enrichment broth, modified Tryptone Soya or Trypticase Soy Broth with novobiocin plus casaminoacids (mTSB+n) to replace modified EC broth with novobiocin (mEC+n). Based on validation study results, the time and temperature range of the incubation was modified. Additional changes to this method include a revised limit of detection, instructions for preparing beef trim/ground beef component subsamples and environmental sponge samples for analysis, and an expanded description of how FSIS confirms and reports *E. coli* O157:H7.

The methods described in this guidebook are for use by the FSIS laboratories. FSIS does not specifically endorse any of the mentioned test products and acknowledges that equivalent products may be available for laboratory use.

| MLG 5.04  | Page 1 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

#### Procedure Outline

- 5.1 Introduction
  - 5.1.1 General
  - 5.1.2 Limits of Detection
- 5.2 Safety
- 5.3 Quality Control Practices
- 5.4 Equipment, Materials, Media, Reagents and Test Kits
  - 5.4.1 Equipment
  - 5.4.2 Media, Reagents and Cultures
  - 5.4.3 Test Kits
- 5.5 Detection Procedure
- 5.6 Isolation Procedure
- 5.7 Identification and Confirmation
- 5.8 Storage of Cultures
- 5.9 Selected References

#### 5.1 Introduction

#### 5.1.1 General

The following method is used for the analysis of raw and ready-to-eat meat products, and laboratory environmental sponge samples for *Escherichia coli* O157:H7. The method is based on enrichment in a selective broth medium, application of a rapid screening test, immunomagnetic separation (IMS) in paramagnetic columns, and plating on a highly selective medium.

Unless otherwise stated all measurements cited in this method have a tolerance of  $\pm 2\%$ .

#### **5.1.2** Limits of Detection

This test method has been shown to detect 0.23 colony forming units (cfu)/g in a 25 g sample of 75/25 (lean/fat) ground beef.

### 5.2 Safety

E. coli O157:H7 is a human pathogen with a low infectious dose. (Ingestion of 100 cells can cause disease.) The use of gloves and eye protection is mandatory and all work surfaces must be disinfected prior to and immediately after use. Laboratory personnel must abide by CDC guidelines for manipulating Biosafety Class II pathogens. A Class II laminar flow biosafety cabinet is recommended for activities with potential for producing aerosols of pathogens. All available Material Safety Data Sheets (MSDS) should be obtained from the manufacturer for

| MLG 5.04  | Page 2 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

the media, chemicals, reagents and microorganisms used in the analysis. The personnel who will handle the materials should read all MSDS sheets.

### **5.3** Quality Control Practices

- a. Rainbow<sup>®</sup> Agar plates have a shelf life of 2 weeks.
- b. All media and E-Buffer must be pre-warmed to 18-35°C prior to use.
- c. The fluorescent strain of *E. coli* O157:H7 is used in this procedure to monitor for cross contamination. The protocol for the use of the fluorescent strain of *E. coli* O157:H7 as a positive control follows:

Wild-type strains of *E. coli* O157:H7 transformed with pGFP produce a green fluorescent protein. As a result of this transformation, fluorescent strains of *E. coli* O157:H7 possess the unique property of expressing bright green fluorescence visible in the dark when illuminated by long-wave UV light. This property, which sets them apart from typical *E. coli* O157:H7, makes them useful positive controls for analyses of meat samples for *E. coli* O157:H7. At different steps in the procedure, both test samples and (fluorescent) positive controls can be tested for the bright green fluorescence as a Quality Control measure to make sure that positive sample isolates actually came from the test sample and not from accidental contamination by the positive control cultures.

Results of studies done at the FSIS Beltsville Microbial Pathogens Laboratory showed that these fluorescent cultures can be subjected to *E. coli* O157:H7 isolation and identification procedures without losing their fluorescent properties. These strains retain their fluorescent properties when grown in SOB media with added ampicillin (SOB + A). These cultures must be transferred every 7 days to fresh SOB + A media, according to the protocol outlined below. The fluorescent colonies are ready to be used as positive controls on day 3 of the following protocol, and for the next 6 consecutive days without losing their fluorescent properties. If these cultures are not needed on a continuous basis, they can be stored at refrigeration temperatures on SOB + A agar plates in ziplock bags or sealed with parafilm for 1 month and then transferred, or started up again 2 days before needed. *Strict adherence* to the protocol described below is essential, in order to ensure that the fluorescent strains do not lose their ability to express green fluorescence.

i. Test the fluorescent *E. coli* O157:H7 strain (FSIS culture # EC 465-97 or the currently designated control strain) on SOB + A agar plate for fluorescence by illuminating colonies under long-wave UV light in the dark.

| MLG 5.04  | Page 3 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

- ii. Select only fluorescent colonies and inoculate into 10 ml of SOB + A broth in a tube. Incubate at  $35 \pm 2$ °C overnight.
- iii. Streak the culture from the SOB + A broth onto a SOB + A agar plate. Incubate at  $35 \pm 2$ °C overnight.
- iv. Examine colonies on the plate for fluorescence. The fluorescent colonies are ready to be inoculated into modified TSB broth + novobiocin (mTSB+n) at this stage. These cultures on SOB + A agar plates can be stored refrigerated and be used as positive controls for 6 more days. Incubate the inoculated mTSB+n positive control culture at  $42 \pm 1^{\circ}$ C overnight, along with the test samples.
- v. Continue analysis per Sections 5.5-5.7 and test the Blood Agar Plates of the fluorescent positive controls and any positive sample cultures for fluorescence.

### 5.4 Equipment, Materials, Media, Reagents and Test Kits

#### 5.4.1 Equipment

- a. Balance, sensitivity of 0.1 g
- b. Stomacher<sup>™</sup> 400 or 3500 with appropriate sizes of sterile Stomacher bags, with or without mesh. (Tekmar Co., Cincinnati, Ohio), or equivalent bag mixer and bags
- c. Incubators, static  $42 \pm 1$  °C,  $35 \pm 2$  °C, and  $22 \pm 2$  °C
- d. Micropipettors to deliver 15-1000 µl with sterile disposable filtered micropipet tips
- e. Mechanical Pipettor with 1.0 ml, 5.0 ml, 10.0 ml sterile pipettes
- f. Inoculating loops, "hockey sticks" or spreaders, and needles
- g. UV light (long-wave, e.g. VWR # 36553-124, or equivalent)
- h. Filter unit, 0.2 μm, nylon, sterile
- i. Infrared thermometer
- j. LabQuake<sup>®</sup> Agitator (or equivalent) with clips to hold microcentrifuge tubes
- k. Sterile disposable 12 x 75 mm polypropylene tubes (e.g. Fisher # 14-956-1B, or equivalent)
- 1. Microcentrifuge and sterile 1.5 ml microcentrifuge tubes
- m. Sterile 50 ml conical tubes (e.g. Falcon® # 2070, or equivalent) or sterile bottles
- n. Sterile 40 µm Cell Strainer (Falcon® # 2340, or equivalent)
- o. MACS<sup>®</sup> Large Cell Separation Columns (Miltenyi Biotec # 422-02, or equivalent)
- p. OctoMACS<sup>®</sup> Separation Magnet (Miltenyi Biotec # 421-09, or equivalent)

| MLG 5.04  | Page 4 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

- q. Multistand to support OctoMACS® Separation Magnet (Miltenyi Biotec # 423-03, or equivalent)
- r. Tray, autoclavable, approximately 130 mm x 83 mm (e.g. VWR # 62663-222, or equivalent) for use with the OctoMACS®
- s. VITEK® system

### 5.4.2 Media, Reagents and Cultures

- a. Modified TSB broth with novobiocin plus casaminoacids (mTSB+n) Oxoid, or equivalent
- b. Rainbow<sup>®</sup> Agar O157 (Biolog Inc., Hayward California, 94545) containing 10 mg/L novobiocin plus 0.8 mg/L potassium tellurite, or equivalent selective medium
- c. Tryptic soy agar with 5% sheep blood
- d. SOB + A Medium
- e. E Buffer, approximately 7 ml per sample (Buffered Peptone Water, Bovine Albumin Sigma # A7906-500G, or equivalent and Tween-20<sup>®</sup>, or equivalent)
- f. Disinfectant (Lysol® I. C., 2.0%, or equivalent)
- g. Dynal® # 710.04 anti-E. coli O157 antibody-coated paramagnetic beads (Dynal Inc., Lake Success, NY 11042), or equivalent
- h. E. coli O157:H7 strain 465-97 (positive control used throughout method)
- i. E. coli ATCC strain 25922 (negative control for bead capture and screen tests)

### **5.4.3** Test Kits

- a. *E. coli* O157:H7 latex agglutination test kit (RIM® *E. coli* O157:H7 Latex Test Kit, REMEL, 12076 Santa Fe Drive, Lenexa, KS 66215), or equivalent
- b. Biochemical test kit and system, GNI and GNI Plus cards (VITEK® system, bioMerieux Vitek, Inc., 595 Anglum Drive, Hazelwood, MO 63042-2395), or equivalent
- c. Shiga toxin test kit [Premier® EHEC, cat. # 608096 (Meridian Diagnostics, Inc., 3471 River Hills Dr., Cincinnati, OH, 45244), or equivalent

#### **5.5** Detection Procedure

a. Sample Preparation

Note: Disinfect the intact sample package(s) surface prior to opening.

| MLG 5.04  | Page 5 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

- i. Raw ground beef microbiological testing programs.
  - Randomly collect five  $65 \pm 2$  g sub-samples (total of  $325 \pm 10$  g) that are representative of the entire sample. Place each  $65 \pm 2$  g sub-sample in a sterile Strainer Stomacher<sup>TM</sup> bag. Add 585 ml mTSB+n broth and pummel for two minutes in a Stomacher<sup>TM</sup>.
- ii. Cooked meat patties and semi-dry and dry fermented sausages. Randomly prepare five 65 ± 2 g sub-samples (total of 325 ± 10 g) that are representative of the entire sample. When appropriate, sample representative portions from both the outer surface (shell) and inner section (core) of RTE products, especially semi-dry and dry fermented sausages. Place each 65 ± 2 g sub-sample in a sterile Strainer Stomacher<sup>™</sup> bag. Add 585 ml mTSB+n broth and pummel for two minutes in a Stomacher<sup>™</sup>.
- iii. Outbreak-related samples. Randomly collect thirteen  $25 \pm 1$  g sub-samples (total of  $325 \pm 13$  g) that are representative of the entire sample. When possible avoid the exterior surfaces of non-intact samples. Place each  $25 \pm 1$  g sub-sample in a sterile Strainer Stomacher bag and add 225 ml of mTSB+n broth. Pummel for 2 minutes in a Stomacher.
- iv. <u>Beef trim/trim components</u>. Prepare a representative sample of the product by taking approximately equal portions from a variety of locations. For non-comminuted components, FSIS inspectors submit samples that generously represent the surface of trimmings. For intact chunk beef or cuts, favor sampling of the surface to mimic submitted beef trim product and to provide the best opportunity for detection.
- v. <u>Environmental sponge samples</u>. Add  $50 \pm 5$  ml of mTSB+n broth to each bagged sponge sample and stomach 2 minutes.
- vi. <u>Most Probable Numbers (MPN) Determination.</u> Follow MPN instructions given in the specific program protocol or see MLG Appendix 2, Most Probable Number Procedure and Tables.
- b. Incubate all bags (static) with their contents for 15 to 22 h at  $42 \pm 1^{\circ}$ C. Include a positive, negative, and uninoculated medium control for each group of samples tested. Use the fluorescent *E. coli* O157:H7 strain (FSIS culture # EC 465-97) as a positive control and *E. coli* ATCC strain 25922 as the negative control.

| MLG 5.04  | Page 6 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

- c. From the enrichment cultures in the Stomacher  $^{\text{\tiny TM}}$  bags, perform the screening test for *E. coli* O157:H7 following the manufacturer's instructions. The enrichment culture may be analyzed immediately upon removal from the incubator without waiting for tempering to room temperature. To prevent clogging the pipette tip, be sure to collect the appropriate size sample from the enrichment culture outside the inner strainer bag.
- d. Samples negative by the screening test can be reported as negative for *E. coli* O157:H7 and discarded.
- e. Samples positive by the screening test should be reported as potential positives. Begin isolation procedures from the enrichment culture in the Stomacher  $^{\text{\tiny TM}}$  bag.

#### **5.6** Isolation Procedure

Note: Steps a.-l. may be performed in a sequence that is convenient to the laboratory personnel.

- a. Prepare E Buffer by mixing 0.5 g Bovine Albumin and 50 µl Tween-20<sup>®</sup> into 100 ml Buffered Peptone Water (BPW). Filter sterilize (0.2 µm) and store at 2-8°C.
- b. Remove Rainbow<sup>®</sup> Agar plates from 2-8°C storage, allowing 3 plates for each screen-positive culture and each control. Be sure that plates have no visible surface moisture at the time of use. If necessary, dry plates (e.g. for up to 30 minutes in a laminar flow hood with the lids removed) prior to use. Dried plates that are not used should be labeled "dried", placed in bags and returned to 2-8°C.
- c. Remove a bottle of E Buffer from 2-8°C storage. Decant 7 ml of E Buffer for each culture and each control into a sterile tube or bottle and allow it to warm to at least 18°C. (Return the stock E Buffer to 2-8°C.)
- d. For each positive control, negative control and screen-positive culture to be analyzed, order and label 50 ml conical centrifuge tubes so that the positive control is first, followed by the negative control, then all cultures. Maintain this order for subsequent steps.
- e. For each positive control, negative control, and screen-positive culture, label two sterile 1.5 ml microcentrifuge tubes (for step g and step s), one 50 ml conical centrifuge tube (for step h.) and two 12 x 75 mm capped tubes (one for step p.). For each pair of 12 x 75 mm tubes, label one tube and add 0.9 ml E Buffer (for step q.).
- f. Prepare the Dynal #710.04 E. coli O157:H7 immunomagnetic bead suspension

| MLG 5.04  | Page 7 of 13 |                    |  |
|---|--------------|--------------------|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |              |                    |  |
| Revision: 04 Replaces: MLG 5.03   |              | Effective: 1/28/08 |  |

by following Table 1 below. Be sure to include the positive and negative controls in the total number of cultures. Use the bead suspension immediately (step g), or hold at 2-8°C. Return the stock vial of Dynal #710.04 *E. coli* O157:H7 immunomagnetic beads to 2-8°C.

- g. Vortex the bead solution briefly (2-3 seconds), then add 50  $\mu$ l to a labeled microcentrifuge tube (from step e), one for each control and screen-positive culture. Use immediately or hold these tubes at 2-8°C.
- h. Place a 40  $\mu$ m Cell Strainer on a labeled 50-ml conical centrifuge tube (from step e.). Pipet 5  $\pm$  1 ml of each control and enrichment culture into the respective Cell Strainer and collect at least 1.0 ml of filtrate.
- i. Do not proceed with more than the number of tubes that the OctoMacs® magnet(s) will hold. Transfer 1.0 ml of a filtrate (step h.) to the corresponding microcentrifuge tube containing the immunomagnetic bead suspension (step g.) and place in the clips of the LabQuake® tube agitator. Rotate the tubes for 10-15 min at 18-30°C.
- j. Attach the OctoMACS® Magnet to the Multistand.
- k. Position a tray on the base of the Multistand so that it will collect the filtrate passing through the columns. Add approximately 300 ml of 2% Lysol<sup>®</sup> I. C. (or equivalent) disinfectant to cover the bottom of the tray.
- 1. Label and place the appropriate number of Large Cell Separation columns on the OctoMACS® Magnet. Insert columns from the front making sure the column tips do not touch any surfaces. Leave the plungers in the bags at this time to maintain sterility.
- m. Transfer at least 0.5ml E Buffer to the top of each column and let the buffer run through.
- n. Resuspend, then transfer each culture and control from step i. to its corresponding column.
- o. After a culture or control has drained through, wash the column by applying 1.0 ml of E Buffer to each column and allow to drain. Repeat 3 more times for a total of 4 washes.

| MLG 5.04  | Page 8 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

- p. After the last wash has drained, remove the column from the OctoMACS<sup>®</sup> Magnet and insert the tip into an empty labeled 12 x 75 mm tube (from step e.). Apply 1.0 ml of E Buffer to the column, and using the plunger supplied with the column, *immediately* flush out the beads into the tube. Use a smooth, steady motion to avoid splattering. Cap the tubes. Repeat this for each column. If the OctoMACS<sup>®</sup> magnet is to be used for a second set of cultures, it must be decontaminated as described in step u, below. Repeat steps j.-s. for the additional cultures.
- q. Vortex the tubes from step p. briefly to resuspend the beads. Make a 1:10 dilution of each treated bead suspension by adding 0.1 ml of the bead suspension to a 12 x 75 mm labeled tube containing 0.9 ml E Buffer (from step e.).
- r. Vortex briefly to maintain beads in suspension and plate 0.1 ml from each tube (from step p. and step q.) onto a labeled Rainbow<sup>®</sup> Agar plate. Use a hockey stick or spreader to spread plate the beads, being careful not to spread the beads against the edge of the plate.
- s. Vortex the tubes containing undiluted beads (from step p.) and transfer to a labeled microfuge tube (from step e.) and centrifuge at least one minute using a bench-top microcentrifuge to concentrate the beads. Withdraw and discard the supernatant without disturbing the beads. Add 0.1 ml of E Buffer to the beads, resuspend the beads and transfer the beads to a labeled Rainbow<sup>®</sup> Agar plate. Spread plate the beads as described in step r.
- t. As soon as there is no visible moisture on the agar surface, invert plates and incubate for 24-26 h at  $35 \pm 2$ °C.
- u. Decontaminate the OctoMACS<sup>®</sup> Magnet by applying 2% Lysol<sup>®</sup> I. C. (or equivalent) disinfectant directly to the surface. After approximately ten minutes, rinse with deionized or tap water. Allow the unit to air-dry or use absorbent paper towels to dry the unit.

| MLG 5.04  | Page 9 of 13       |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

Table 1.

| # of Cultures  | ul of Beads*   | ul of E-Buffer | # of Cultures | ul of Beads* | ul of E-Buffer |
|----------------|--|----------------|---------------|--------------|----------------|
| 1              | 15   | 135            | 26            | 145          | 1305           |
| 2              | 20   | 180            | 27            | 150          | 1350           |
| 3              | 25   | 225            | 28            | 155          | 1395           |
| 4              | 30   | 270            | 29            | 160          | 1440           |
| 5              | 35   | 315            | 30            | 165          | 1485           |
| 6              | 40   | 360            | 31            | 175          | 1575           |
| 7              | 45   | 405            | 32            | 180          | 1620           |
| 8              | 50   | 450            | 33            | 185          | 1665           |
| 9              | 55   | 495            | 34            | 190          | 1710           |
| 10             | 60   | 540            | 35            | 195          | 1755           |
| 11             | 65   | 585            | 36            | 200          | 1800           |
| 12             | 70   | 630            | 37            | 205          | 1845           |
| 13             | 75   | 675            | 38            | 210          | 1890           |
| 14             | 80   | 720            | 39            | 215          | 1935           |
| 15             | 85   | 765            | 40            | 220          | 1980           |
| 16             | 90   | 810            | 41            | 230          | 2070           |
| 17             | 95   | 855            | 42            | 235          | 2115           |
| 18             | 100  | 900            | 43            | 240          | 2160           |
| 19             | 105  | 945            | 44            | 245          | 2205           |
| 20             | 110  | 990            | 45            | 250          | 2250           |
| 21             | 120  | 1080           | 46            | 255          | 2295           |
| 22             | 125  | 1125           | 47            | 260          | 2340           |
| 23             | 130  | 1170           | 48            | 265          | 2385           |
| 24             | 135  | 1215           | 49            | 270          | 2430           |
| 25             | 140  | 1260           | 50            | 275          | 2475           |
| * Dynal® anti- | * Dynal® anti-E. coli O157:H7 antibody-coated paramagnetic beads (vortex briefly before use) |                |               |              |                |

Dynal anti-E. coli O15/:H/ antibody-coated paramagnetic beads (vortex briefly before use)

| MLG 5.04  | Page 10 of 13      |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

#### 5.7 Identification and Confirmation

a. After incubation, *E. coli* O157:H7 colonies have black or gray coloration on Rainbow<sup>®</sup> Agar. When *E. coli* O157:H7 colonies are surrounded by pink or magenta colonies, they may have a bluish hue. Mark colonies typical of *E. coli* O157:H7 and perform latex agglutination assays for O157, following manufacturer's instructions. Streak all latex positive colonies, up to a total of five from each sample (one per subsample, if possible) onto Blood Agar plates. Incubate Blood Agar plates for 16-24 h at  $35 \pm 2^{\circ}$ C.

Note: If no typical colonies are present, hold the original Rainbow<sup>®</sup> plates at 20-24°C for an additional 6-24 h then re-examine for typical colonies.

b. After incubation, examine the Blood Agar plates for purity under visible light, and evidence of cross contamination with the positive control by using long wave UV light. Only the positive control culture, *E. coli* O157:H7 strain 465-97, should fluoresce. If the Blood Agar plates appear pure and uncontaminated, perform the following confirmatory tests:

#### i. Biochemical confirmation.

Inoculate VITEK-GNI or GNI Plus cards (if using VITEK® system), or an equivalent biochemical identification testing system. The cytochrome oxidase and gram stain tests are optional.

#### ii. O157 and H7 confirmation.

To confirm the absence or presence of O157 and H7 antigens, use an *E. coli* O157:H7 latex test agglutination kit (RIM<sup>®</sup> *E. coli* O157:H7 Latex Test Kit, or equivalent). Use growth from the Blood Agar plate (from step b). For inconclusive results, genetic testing (e.g. PCR) may be necessary.

#### iii. Shiga toxin/toxin genes confirmation.

The presence of Shiga toxin(s) in a culture isolate is confirmed by the use of a toxin assay, e.g., Meridian Premier<sup>®</sup> EHEC Kit, or equivalent. When Shiga toxin(s) is (are) not demonstrated, detection of one or more toxin genes by PCR is used for confirmation. NOTE: The positive control culture, *E. coli* O157:H7 (FSIS culture # EC 465-97), is toxin-negative.

| MLG 5.04  | Page 11 of 13      |  |  |
|---|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |  |  |
| Revision: 04  | Effective: 1/28/08 |  |  |

- c. FSIS uses the following definitions for reporting *E. coli* O157:H7:
  - <u>Potential Positive</u> a sample that causes a positive reaction with the screen test.
  - <u>Presumptive Positive</u> a sample that has typical colonies, observed on Rainbow Agar, and reacts specifically with O157 antiserum.
  - <u>Confirmed Positive</u> a biochemically identified *Escherichia coli* isolate that is serologically or genetically determined to be "O157" that meets at least one of the following criteria:
    - 1) Positive for Shiga toxin (ST) production
    - 2) Positive for Shiga toxin gene(s) (stx)
    - 3) Genetically determined to be "H7"

If an FSIS laboratory's confirmatory test results are insufficient to allow identification, then the isolate is transferred to the Outbreaks Section of the Eastern Laboratory Microbiology Branch (OSEL), or current FSIS reference laboratory, for further testing prior to reporting.

For example, if a presumptive positive isolate is tested in a FSIS laboratory for Shiga toxin and it is Shiga toxin negative, then the isolate is transferred to OSEL for further testing.

A Shiga toxin gene PCR test is performed in OSEL. If the isolate fails to express Shiga toxin (toxin negative result), but has the genes necessary to express Shiga toxin, then it is considered by the agency to be *E. coli* O157:H7 confirmed positive.

If the isolate was H7 negative and found to be Shiga toxin and gene negative (by EHEC test and genetic test for Shiga toxin genes) additional PCR test(s) for H7 gene(s) are performed. If the H7 PCR test is positive, the isolate is considered *E. coli* O157:H7 confirmed positive.

If the isolate is *E. coli* O157 presumptive positive, but additional tests show it to be H7 negative (by latex agglutination and PCR) and Shiga toxin negative (by EHEC test and genetic test for Shiga toxin genes), then the isolate would be reported as *E. coli* O157:H7 negative.

| MLG 5.04  |                    | Page 12 of 13      |  |  |
|---|--------------------|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |                    |  |  |
| Revision: 04  | Replaces: MLG 5.03 | Effective: 1/28/08 |  |  |

### **5.8** Storage of Cultures

For storage requirements of the fluorescent *E. coli* O157:H7 strain (FSIS culture # EC 465-97 or the currently designated control strain), refer to Section 5.3.c. of this chapter.

Store other "working" E. coli stock cultures on nutrient agar slants. Transfer stocks monthly onto duplicate nutrient agar slants, incubate overnight at  $35 \pm 2$ °C, and then store them at 2-8°C. Use one of the slants as the working culture. Use the other slant for sub-culturing to reduce the opportunity for contamination. Cultures may be subcultured up to 5 times. After this period the culture must be re-confirmed biochemically or a new culture initiated.

For long term storage freeze cultures using cryo-beads i.e. Cryostor<sup>TM</sup> or lyophilize.

#### **5.9** Selected References

Deng, M.Y., P.M. Fratamico, M.A. Jensen, and G. Tice. 1997. DNA amplification using tableted PCR reagents for identification of *Escherichia coli* O157:H7 isolated from foods. Journal of Rapid Methods and Automation in Microbiology. 5:61-74.

Ewing, W. H. 1986. Edwards and Ewing's Identification of Enterobacteriaceae, 4th Edition. Elsevier Science Publishing Co., Inc., New York.

Feng P., S. D. Weagant, and M. A. Grant. FDA Bacteriological Analytical Manual *Online*, Chapter 4: Enumeration of *Escherichia coli* and the Coliform Bacteria.

Fratamico, P. M., M. Y. Deng, T. P. Strobaugh, and S. A. Palumbo. 1997. Construction and characterization of *Escherichia coli* O157:H7 strains expressing firefly luciferase and green fluorescent protein and their use in survival studies. J. Food Prot. **60**:1167-1173.

Harrison, B. and D. Warburton. 1997. Identification of *Escherichia coli* verotoxins by the Meridian Premier EHEC kit<sup>®</sup>. Laboratory procedure MFLP-93 *In* The Compendium of Analytical Methods, Vol. 3. Health Protection Branch, Health Canada, Ottawa, Canada.

Hitchins, A. D., P. A. Hartman, and E. C. D. Todd. 1992. Coliforms-*Escherichia coli* and its toxins, p. 325-369. *In* C. Vanderzant and D. F. Splittstoesser (ed.), Compendium of Methods for the Microbiological Examination of Foods. 3rd Edition. Amer. Publ. Hlth. Assoc., Washington, D.C. 20005.

| MLG 5.04  |                    | Page 13 of 13      |  |  |
|---|--------------------|--------------------|--|--|
| Title: Detection, Isolation and Identification of Escherichia coli O157:H7 from Meat Products |                    |                    |  |  |
| Revision: 04  | Replaces: MLG 5.03 | Effective: 1/28/08 |  |  |

Park, C. H., K. M. Gates, N. M. Vandl, and D. L. Hixon. 1996. Isolation of Shiga-like toxin producing *Escherichia coli* (O157 and non-O157) in a community hospital. Diagn. Microbiol. Infect. Dis. **26**:69-72.

Richmond, J.Y. and R.W. McKinney (ed.). 2007. Biosafety in Microbiological and Biomedical Laboratories, 5<sup>th</sup> ed. U.S. Government Printing Office, Washington, D.C.

Taormina, P. J., M. Rocelle, S. Clavero, and L. R. Beuchat. 1998. Comparison of selective agar media and enrichment broths for recovering heat-stressed *Escherichia coli* O157:H7 from ground beef. Food Microbiol. **15**:631-638.

Weagant, S. D., J. L. Bryant, and K. G. Jinneman. 1995. An improved rapid technique for isolation of *Escherichia coli* from foods. J. Food Prot. **58**:7-12.