

Chapter 8

Isolation and Identification of *Escherichia coli* Serotype O157:H7

Isolation and identification of *Escherichia coli* serotype O157:H7 can be greatly enhanced when optimal laboratory media and techniques are employed. The methods presented here are intended to be economical and to offer laboratorians some flexibility in choice of protocol and media. Laboratories that do not have sufficient resources to adopt the methods described below should consider sending specimens or isolates to other laboratory facilities that routinely perform these procedures. Laboratory supplies required for diagnosis of *E. coli* O157:H7 are listed in Annex H.

A. Isolation and Identification Methods

E. coli O157:H7 rapidly ferments lactose and is indistinguishable from most other *E. coli* serotypes on traditional lactose-containing media. However, unlike approximately 80% of other *E. coli*, nearly all isolates of *E. coli* O157:H7 ferment D-sorbitol slowly, or not at all. Sorbitol-MacConkey (SMAC) agar was developed to take advantage of this characteristic by substituting the carbohydrate sorbitol for lactose in MacConkey agar, and it is the medium of choice for isolation of *E. coli* O157:H7. Sorbitol-negative colonies will appear colorless on SMAC (Figure 8-1).




Figure 8-1. *E. coli* O157:H7 colonies are colorless on SMAC. Non-O157 *E. coli* colonies are pink.

Isolation and Identification of Escherichia coli Serotype O157:H7

Enrichment for *E. coli* O157:H7 is not usually necessary for isolation of the organism from acutely ill patients.

Figure 8-2 illustrates the procedure for recovery of *E. coli* O157:H7 from fecal specimens. SMAC is inoculated as described in Chapter 4 (Figure 4-2). Incubate 18 to 24 hours at 35° to 37°C. After 18 to 24 hours' incubation, the amount and type of growth (e.g., sorbitol-positive or sorbitol-negative) on SMAC should be recorded on data sheets for each patient specimen (Figure 8-3). Colonies suspicious of *E. coli* O157:H7 will appear colorless and about 2 to 3 mm in diameter (Figure 8-1).

Test sorbitol-negative colonies selected from SMAC with *E. coli* O157 anti-serum or latex reagents (O157 antibody-coated latex and control latex) according to the procedures recommended by the manufacturer. Suspected colonies may be tested with antisera directly from the SMAC plate or subcultured to a nonselective medium (HIA, for example) and tested the next day. This provides more growth on which to perform the agglutination assay (however, some manufacturers of O157 latex reagents recommend testing only colonies taken directly from the plate). If colonies are tested directly from the plate, a colony that is positive in O157 latex reagent should also be transferred to another medium for subsequent testing. Once one colony from a plate has been identified as O157-positive, no further colonies from the same plate need to be tested.




Figure 8-2. Procedure for recovery of *E. coli* O157:H7 from fecal specimens

Escherichia coli O157:H7 Worksheet

SPECIMEN NUMBER	MEDIA	SORBITOL -	SORBITOL +	COLONY	SEROLOGY		BIOCHEMICAL ID <i>E. coli</i> YES/NO	PRESUMPTIVE IDENTIFICATION
					O157 LATEX	CONTROL LATEX		
	SMAC			SM1				
				SM2				
				SM3				
	SMAC			SM1				
				SM2				
				SM3				
	SMAC			SM1				
				SM2				
				SM3				
	SMAC			SM1				
				SM2				
				SM3				
	SMAC			SM1				
				SM2				
				SM3				
	SMAC			SM1				
				SM2				
				SM3				

Figure 8-3. Escherichia coli O157:H7 worksheet

