

UNITED STATES ARMY MEDICAL RESEARCH INSTITUTE OF INFECTIOUS DISEASES



Specimen Collection and Submission Manual

Special Pathogens Laboratory

Diagnostic Systems Division

Fort Detrick, Maryland 21702

May 2012



This manual is designed to provide detailed instructions for submission of samples that will be analyzed in the SPL. Tests that are not listed may require special preparation or advance notice. Please notify the SPL in advance for tests not listed in this manual. Please call the laboratory manager to obtain additional information.

Table of Contents

Introduction.....3

SPL Contact Information.....4

Technical Guidelines and General Information.....5–7

List of Tests.....8-14


Diagnostic Specimen Requirements.....15-23

Environmental samples.....23

Acronyms and definitions.....24-25

INTRODUCTION

1. This manual has been prepared to assist the Special Pathogens Laboratory (SPL) customer base. It contains information about procedures performed in the SPL. Included are information about specimen types and special instructions about particular tests and procedures. Changes to this handbook will be published as needed in the form of revised pages and distributed through normal channels.
2. Your comments, suggestions, and assistance in regards to the improvement of this manual or of services in general will be appreciated and can be addressed to me, my managerial and supervisory staff, or any other member of the SPL.

A handwritten signature in black ink, appearing to read 'William Dorman', written over a horizontal line.

William Dorman

Manager, Special Pathogens Laboratory

Disclaimer:

Opinions, interpretations, conclusions, and recommendations are those of the authors and are not necessarily endorsed by the US Army or the Department of Defense. The mention of specific commercial entities does not imply endorsement by the US Army or the Department of Defense.

SPECIAL PATHOGENS LABORATORY CONTACT LISTING

Commercial Area Code (301), DSN 343

USAMRIID EMERGENCY HOTLINE	(888) 872-7443
USMRIID Security (24hrs), SPL after hours notification	619-2257
USAMRIID Commander	619-2833
Division Chief	619-4721
Director, Special Pathogens Laboratory	619-4738
Manager, Special Pathogens Laboratory	619-3318
FAX	619-2492

Normal hours of operation are Monday through Friday 0730 to 1630. Support during non-duty hours is available by prior arrangement.

TECHNICAL GUIDELINES AND GENERAL INFORMATION

Acceptable Specimen Volumes

Test entries list minimum acceptable specimen volumes. The minimum volume is defined as the absolute minimum needed to run a validated test algorithm. If there is insufficient specimen volume for testing, there may be delays, and the request may be referred to management for testing approval using modified procedures. If modified procedures are required to complete testing, the client will be notified before testing.

Accreditation/Licensure

For clinical microbiology samples, the SPL is CLIP (Clinical Laboratory Improvement Program) certified. SPL also maintains current licenses, permits, and registrations required by federal regulations.

For environmental samples, the SPL adheres to ISO 17025 standards as a sample testing laboratory.

For additional information or copies of certificates, please contact the SPL.

General Bacterial Isolation Requirements

SPECIMEN	CONTAINER OF TRANSPORT DEVICE	PROCEDURE
Throat	Swab	Swab tonsillar areas, posterior pharynx, and areas of inflammation, ulceration, or exudation. Send properly labeled swab along with request to the laboratory.
Nasopharynx	Flexible calcium alginate swab	Insert swab until it reaches posterior nasopharynx; rub mucosa gently; then withdraw and return to sterile container. Send to the laboratory.
Urine	Sterile screw-capped container or sterile urine container	Have patient clean penile or labial area carefully; collect midstream sample in sterile container cover container immediately; send to the laboratory promptly; refrigerate if storing longer than one (1) hour.
Superficial wounds and abscesses	Swab	Clean the site around the wound, Vigorously swab area with Culturette; place swab in swab transport media, send to the lab as soon as possible. If anaerobe is suspected, a third swab should be collected and placed in an anaerobic transport system.
Rectal / Oral	Swab	Specimen should be collected into a sterile container; stool must not be retrieved from the toilet or contaminated with urine; specimen must be sent to the lab within 1 hour or refrigerated no longer than 24 hours if transit is delayed. (<i>Salmonella</i> , <i>Shigella</i> and <i>Campylobacter</i> may not be recovered from refrigerated specimens).
Blood	Blood culture bottles (aerobic and/or anaerobic) or isolator tubes	Prepare patient's skin for venipuncture by swabbing with alcohol prep then iodine prep; allow to dry 1 minute; do not touch skin after preparation; prepare tops of blood culture bottles in the same manner as patient's skin; use syringe to collect 20 ml of blood; change needles before inoculating 10 ml of blood into each bottle. Send to the laboratory within 24 hours.

Crisis Contingency Plan

SPL maintains a corporate contingency plan for crisis recovery and business continuation. The purpose of this plan is to ensure prompt recovery of SPLs' critical business functions in the event of a crisis affecting any aspect of our continued patient care service. In the event of a local, regional, or national crisis that adversely affects timely delivery of specimens to SPL facilities, SPL will expeditiously initiate specific client-notification procedures to provide clients with necessary information and instruction on prearranged transportation and testing alternatives.

Health Insurance Portability and Accountability Act (HIPAA)

SPL is committed to complying with privacy and security standards promulgated in the Health Insurance Portability and Accountability Act (HIPAA). MEDCOM has implemented policies, processes, and procedures designed to ensure compliance with the standards. Compliance is continuously monitored for effectiveness. Workforce training is completed annually.

SPL complies with security standards by ensuring that systems, policies, and procedures meet or exceed all required and addressable implementation specifications. Internet and interface connectivity is encrypted and/or password protected, and electronic access is limited to authorized entities. Breaches of protected health information (PHI) or other confidential business information are reviewed and reported to the Department of Health and Human Services (DHHS) as necessary.

Inappropriate Submissions

All specimens should be collected, labeled, transported, and processed according to this procedure. Review the appropriate container type, volume, and special handling requirements needed for analysis before the specimen is collected. Due to the nature of typical specimens submitted to the SPL, less than optimal specimens may be tested upon consultation with SPL management. If any of the guidelines for these processes are not met, the specimen may be rejected or the test may be canceled. SPL will contact the client for resolution. The following list represents some possible causes for specimen rejection or test cancellation:

- Inappropriate specimen type
- Insufficient volume for analysis
- Improperly labeled specimen
- Inappropriate specimen container
- Improper specimen transport
- Specimen has leaked in transit
- Specimen has been submitted in incorrect or expired transport media
- Incomplete test request
- Test order without a specimen
- Specimen without a test order
- No specimen type provided
- No source provided*
- Compromised specimen (e.g., hemolysis, lipemic, or clotted specimens)

* The source of specimen, when appropriate, must be included on the paper or electronic request. The source of specimen is required for all infectious disease testing, including PCR tests.

Laboratory Result Reporting

SPL communicates laboratory results to clients by several means, including printed reports, email, and verbal results. Clients may request a phone notification or fax report by written notification when tests are ordered.

Preliminary results may be offered for infectious diseases and other tests in which a final report follows. Final results are generated at the completion of testing and may contain updated information from the preliminary result. When critical results are obtained, results are called to the physician or requesting lab.

Critical Values/Normal Values

SPL's testing laboratories operate Monday through Friday 8am-4:30pm. Critical results are reported as soon as testing has been completed and a critical value has been identified. SPL reports critical values immediately to the contact(s) provided by our client(s) in accordance with the Laboratory Accreditation Program Inspection Checklists.

The normal value for most infectious agent detection is "not-detectable." A confirmed detectable/identified agent will be reported immediately to the client as well as federal agencies; as appropriate.

Referral Testing

One of SPL's service goals is to support clients by providing comprehensive service for all laboratory testing. To accomplish this goal, SPL may on occasion, select additional resources to perform additional tests not performed at SPL. Primary referral resources are typically intrinsic to USAMRIID and share the same aspects of service, quality, reliability, and turnaround time.

Every effort is made to test specimens at SPL, although the referral services may not fall under accredited services. Those tests not performed by SPL are clearly identified in the sample report.

Shipment Address

Samples to be shipped by commercial carrier should be addressed to:

USAMRIID
Attn: SPL MCMR-UID-F
1425 Porter Street
Fort Detrick MD 21702-5011

It is highly recommended that SPL personnel or the USAMRIID Command be notified before shipment.

Samples should be clearly marked for special testing and preferably the chain of custody procedures maintained. Each sample should be placed in a water-tight receptacle made of glass, plastic, or metal with a screw cap closure. It is recommended that the screw cap be reinforced with adhesive tape. Liquid specimens should be placed individually in a second plastic vial or zip-top bag to prevent leakage. Absorbent material sufficient to absorb the entire content of the primary receptacle is placed between the primary and secondary packaging; for liquid specimens placed in a plastic vial, the absorbent material should be placed between the plastic vial and another secondary packaging material. The secondary packaging should be of material that prevents leakage.

For transportation, the samples should be packaged in an International Air Transportation Association (IATA) or Department of Transportation 49, Code of Federal Regulation 173 approved container, accordingly.

Additional guidance on packing and shipping infectious substances can be found through American Society for Microbiology:

<http://www.asm.org/images/pdf/Clinical/pack-ship-7-15-2011.pdf>

LIST OF TESTS

Arbovirus - Detection in clinical samples

Methodology: culture

Tests May Include: West Nile virus (WNV), eastern equine encephalitis virus, Venezuelan equine encephalitis virus, yellow fever

Specimen: see specimen requirements

Shipping: ship at refrigeration temperature.

Turnaround: within 10 days of specimen receipt

Arbovirus - PCR panel

Methodology: PCR

Tests May Include: West Nile virus (WNV), eastern equine encephalitis virus, Venezuelan equine encephalitis virus, yellow fever

Specimen: see specimen requirements

Shipping: ship at refrigeration temperature

Turnaround: 2-3 days

Bacillus anthracis - Detection in clinical samples

Methodology: culture, biochemical testing, JBAIDS and Laboratory Response Network protocols

Tests May Include: biochemical, molecular, and immunological methods

Specimen: see specimen requirements

Shipping: ship at refrigeration temperature.

Turnaround: 2-3 days

Note: The JBAIDS PCR kit is FDA cleared for the identification of *Bacillus anthracis* from whole blood, blood culture samples and isolated organism grown on agar plates.

Bacillus anthracis - Identification/confirmation of referred isolate

Methodology: biochemical testing, PCR and Laboratory Response Network protocols

Tests May Include: biochemical, molecular, and immunological methods

Specimen: actively growing pure culture on suitable medium

Shipping: ship at room temperature.

Turnaround: within 1-2 days of specimen receipt

Note: This organism has been designated as a select agent (Select Agent Regulation, 42 CFR, 73, Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.

Bacteria, Aerobic - Identification

Methodology: biochemical testing, Typing or grouping if appropriate, Sequencing of 16S ribosomal DNA if indicated

Specimen: actively growing pure culture on suitable medium

Shipping: ship at room temperature.

Turnaround: not available

Bacteria, anaerobic- identification

See: *Clostridium botulinum*, *Clostridium perfringens*

Comments: anaerobic bacteria other than *Clostridium* spp. not identified at SPL

Brucella spp - detection in clinical samples

Methodology: culture and Laboratory Response Network protocols

Tests May Include: biochemical, molecular, and immunological methods

Specimen: see specimen requirements

Shipping: ship at refrigeration temperature.

Turnaround: negative results available after 7-21 days of incubation

Brucella spp - identification/confirmation of referred isolate

Methodology: biochemical testing and laboratory response network protocols

Tests May Include: biochemical, molecular, and immunological methods

Specimen: actively growing pure culture on suitable medium

Shipping: ship at room temperature.

Turnaround: within 1 week of specimen receipt

Note: *B. abortus*, *melitensis*, and *suis* have been designated as select agents (Select Agent Regulation, 42 CFR, 73, Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.

Burkholderia mallei - detection in clinical samples

Methodology: culture, biochemical testing and Laboratory Response Network protocols

Tests May Include: biochemical, molecular, and immunological methods

Specimen: see specimen requirements

Shipping: ship at refrigeration temperature.

Turnaround: 4-7 days

Burkholderia mallei - identification/confirmation of referred isolate**Methodology:** biochemical testing, PCR and Laboratory Response Network protocols**Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** actively growing culture on suitable medium, see specimen requirements**Shipping:** ship at room temperature.**Turnaround:** 4 working days**Note:** This organism has been designated as a select agent (Select Agent Regulation, 42 CFR, 73, Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.***Burkholderia pseudomallei*** - detection in clinical samples**Methodology:** culture, biochemical testing, PCR and Laboratory Response Network protocols**Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** see specimen requirements**Shipping:** ship at refrigeration temperature.**Turnaround:** 4-7 days***Burkholderia pseudomallei*** - identification/confirmation of referred isolate**Methodology:** biochemical testing, PCR and Laboratory Response Network protocols**Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** actively growing culture on suitable medium, see specimen requirements**Shipping:** sShip at room temperature.**Turnaround:** 4 working days**Note:** This organism has been designated as a Select Agent (Select Agent Regulation, 42 CFR, 73, Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.***Clostridium botulinum*** - detection in clinical and environmental samples and food**Methodology:** culture**Specimen:** see specimen requirements**Shipping:** see specimen requirements**Turnaround:** within 10 working days of specimen receipt**Note:** This organism has been designated as a select agent (Select Agent Regulation, 42 CFR, 73, Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.

Clostridium botulinum - identification/confirmation of referred isolate

Methodology: biochemical testing, toxin testing

Specimen: actively growing pure culture on suitable medium

Shipping: ship at room temperature.

Turnaround: not available

***Clostridium botulinum* toxin** - detection in clinical and environmental samples and food

Methodology: toxin neutralization assay, mouse bio-assay

Specimen: see specimen requirements

Shipping: see specimen requirements

Turnaround: within 7 working days of specimen receipt

Clostridium perfringens - detection in stool or implicated food

Methodology: culture

Specimen: stool in ParaPak C&S, modified Carey-Blair or equivalent - fill to line (approximately 5 ml)
Implicated food - minimum of 10 g in original container or transferred to sterile container using sterile instruments.

Shipping: ship stool at room temperature, food at refrigeration temperature.

Turnaround: within 5 working days of specimen receipt

Clostridium perfringens - identification/confirmation of referred isolate

Methodology: biochemical testing, toxin testing

Specimen: actively growing pure culture on suitable medium

Shipping: ship at room temperature.

Turnaround: Within 5 working days of specimen receipt

Coxiella burnetii - detection of DNA in clinical samples

Methodology: Laboratory Response Network Protocols

Tests May Include: PCR or immunoassays

Specimen: tissue or bone marrow (100 mg)

Whole EDTA blood or serum (0.5 ml)

Nasopharyngeal or throat swab, dry or in transport medium

Sputum, bronchial/tracheal washings (0.5 ml)

Lesion exudates

Shipping: ship at refrigeration temperature.

Turnaround: 2 days

Francisella tularensis - detection in clinical samples**Methodology:** culture and Laboratory Response Network protocols**Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** see specimen requirements**Shipping:** ship at refrigeration temperature.**Turnaround:** negative results available after 5-7 days of incubation**Note:** The JBAIDS PCR kit is FDA cleared for the identification of *Bacillus anthracis* from whole blood, blood culture samples and isolated organism grown on agar plates..***Francisella tularensis*** - identification/confirmation of referred isolate**Methodology:** culture and Laboratory Response Network protocols**Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** actively growing pure culture on suitable medium**Shipping:** ship at room or refrigeration temperature.**Turnaround:** 3-5 days**Note:** This organism has been designated as a select agent (Select Agent Regulation, 42 CFR, 73, Interim Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.**Influenza virus** - detection of RNA in clinical samples**Methodology:** real-time PCR**Specimen of choice:** nasopharyngeal swab is the preferred specimen; other acceptable specimens include throat swab; nasal aspirate; combined nasal swab with oropharyngeal swab**Shipping:** ship at refrigeration temperature.**Turnaround:** 1 week**Ricin toxin** - detection in non-clinical samples**Methodology:** immunological and molecular methods**Specimen:** food, soil, water, environmental surface wipe**Shipping:** ship at refrigeration temperature.**Turnaround:** 1-2 days**Staphylococcal enterotoxin A & B** - detection in non-clinical samples**Methodology:** immunological and molecular methods**Specimen:** food, soil, water, environmental surface wipe**Shipping:** sShip at refrigeration temperature.**Turnaround:** 1-2 days

Vaccinia virus - detection in clinical samples**Methodology:** culture**Specimen:** roof of lesion in a sterile container swab of lesion, dry or in transport medium. Contact lab for details. Vesicular fluid sample in viral transport medium.**Shipping:** ship at refrigeration temperature.**Comments:** this test is recommended for patients exhibiting acute, generalized, vesicular or pustular rash illness. For details on evaluating patients see the CDC poster "Evaluating Patients for Smallpox" (<http://www.bt.cdc.gov/agent/smallpox/diagnosis/index.asp>)**Shipping:** ship at refrigeration temperature.**Turnaround:** 1 – 3 weeks**Vaccinia virus** - detection of DNA in clinical samples**Methodology:** PCR (Laboratory Response Network protocols)**Specimen:** roof of lesion in a sterile container swab of lesion, dry or in transport medium. Contact lab for details. Touch-prep (slide) of vesicular fluid.**Shipping:** ship at refrigeration temperature.**Turnaround:** 1-2 days**Comments:** This test is recommended for patients exhibiting acute, generalized, vesicular or pustular rash illness. For details on evaluating patients see the CDC poster "Evaluating Patients for Smallpox" (<http://www.bt.cdc.gov/agent/smallpox/diagnosis/index.asp>)**Comments:** This test is for research use only.**Varicella zoster virus** - detection of DNA in clinical samples**Methodology:** PCR (Laboratory Response Network protocols)**Specimen:** roof of lesion in a sterile container swab of lesion, dry or in transport medium. Contact lab for details. Touch-prep (slide) of vesicular fluid.**Shipping:** ship at refrigeration temperature.**Turnaround:** 1-2 days**Comments:** These tests are recommended for patients exhibiting acute, generalized, vesicular or pustular rash illness. For details on evaluating patients see the CDC poster "Evaluating Patients for Smallpox" (<http://www.bt.cdc.gov/agent/smallpox/diagnosis/index.asp>)**Comments:** This test is for research use only.**Variola virus****Restrictions:** If smallpox is suspected, call SPL for instructions.**Viral hemorrhagic fevers** – detection in clinical samples**Methodology:** culture and serology**Tests May Include:** Ebola, Marburg, CCHF, RVF, Hantavirus, Lassa, Junin, Machupo**Specimen:** see specimen requirements**Shipping:** ship at refrigeration temperature.**Turnaround:** Usually 10 to 21 days

Viral hemorrhagic fevers – PCR**Methodology:** PCR**Tests May Include:** Ebola, Marburg, CCHF, RVF, Hantavirus, Lassa, Junin, Machupo**Specimen:** see specimen requirements**Shipping:** ship at refrigeration temperature.**Turnaround:** 1-2 days***Yersinia pestis* - detection in clinical samples****Methodology:** culture, JBAIDS and Laboratory Response Network protocols**Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** see specimen Requirements**Shipping:** ship at refrigeration temperature.**Turnaround:** 5 days**Note:** The JBAIDS PCR kit is FDA cleared for the identification of *Yersinia pestis* from whole blood, sputum, blood culture samples and isolated organism grown on agar plates.***Yersinia pestis* - identification/confirmation of referred isolate****Tests May Include:** biochemical, molecular, and immunological methods**Specimen:** actively growing pure culture on suitable medium**Shipping:** ship at room or refrigeration temperature.**Turnaround:** 2-3 days**Note:** This organism has been designated as a select agent (Select Agent Regulation, 42 CFR, 73, Interim Final Rule). Special handling criteria apply. Please contact the laboratory for special instructions.

Diagnostic Specimen Requirements for *Bacillus anthracis*

Type of Infection	Specimen type	Minimum Volume	Collection Comments
Cutaneous anthrax	Vesicle Swab	2 swabs	Vesicle should be unroofed and 2 sterile, dry swabs should be soaked in the vesicular fluid.
	Vesicle Aspirate	1 ml	An aspirate of the fluid is also an appropriate specimen.
	Eschar Swab	2 swabs	Roll swabs beneath the edge of the eschar without removing it.
	Fresh/Frozen Tissue	1 punch biopsy	
Gastro-intestinal anthrax	Stool	5 g	If unable to obtain stool, obtain rectal swab by inserting swab 1 inch beyond anal sphincter.
	Rectal Swab	1 swab	
Inhalation anthrax	Nasal Swab	1 swab	For epidemiologic purposes only. Useful only within 24 hours of exposure.
	Sputum	1 ml	If patient has a productive cough, this is the specimen of choice in the early course of the disease.
	Tracheal Aspirates, Bronchoalveolar Wash, etc.	1 ml	
Meningitis	CSF	1 ml	Centrifuge ≥ 1 ml of fluid
Blood	3.2% Sodium Citrate Blood (Blue Top)	1 ml	For molecular testing.
	Serum/Plasma	2 ml	
	Blood Culture	Refer to manufacturer's recommendation	Collect appropriate blood volume and number of sets per submitting lab's protocol. Collect prior to antibiotic use if possible. Most likely to be positive in later stages of disease.
	Blood in Sodium Polyanethol Sulfonate (SPS)	8 ml	Blood collected in SPS may be inoculated into blood culture bottles upon receipt at SPL.
Other	Pleural Fluid	1 ml	
Environmental	Swab	1 swab	

Diagnostic Specimen Requirements for *Brucella* spp.

Specimen Type	Minimum Volume	Collection Comments
Blood culture	Refer to manufacturer's recommendation	Collect appropriate blood volume and number of sets per submitting lab's protocol. Collect prior to antibiotic use if possible. Multiple specimens increase possibility of obtaining a positive culture.
Blood in sodium polyanethol sulfonate (SPS)	8 ml	Blood collected with SPS may be inoculated into blood culture bottles upon receipt at SPL.
Bone marrow blood culture in bottle fluid	1 ml	Collect appropriate bone marrow volume per manufacturer's recommendation. Some blood culture systems are appropriate for bone marrow.
Bone marrow in sodium SPS	1 ml	Bone marrow collected in SPS may be inoculated into culture bottles upon receipt at SPL.
Abscess material	1 ml	Collect as needed based on clinical presentation. Appropriate postmortem specimen.
Lymph node, liver/spleen biopsy	1-5 g	Collect as needed based on clinical presentation. Appropriate postmortem specimen.
Synovial fluid, CSF, other body fluids	1 ml	Collect as needed based on clinical presentation. Appropriate postmortem specimen.
Whole blood	1 ml	3.2% sodium citrate blood (blue top)
Nasal swab	1 swab	For epidemiologic purposes only. Useful only within 24 hours of exposure.
Serum: Acute and Convalescent	2 ml	Acute-phase specimen should be collected ASAP after onset of disease. Convalescent-phase specimen should be collected >14 days after the acute specimen.
Environmental	1 swab	

Diagnostic Specimen Requirements for *Burkholderia* spp.

Specimen Type	Minimum Volume	Collection Comments
Abscess material, tissues	1 ml	Collect tissues and fluids rather than swabs, when possible. Collect as needed based on clinical presentation. Appropriate postmortem specimen.
CSF, other body fluids	1 ml	Collect as needed based on clinical presentation. Appropriate postmortem specimen.
Sputum	1 ml	
Skin swab	1 swab	
Urine	1 ml	Collect a midstream clean-catch or a catheterized specimen.
Blood culture	Refer to manufacturer's recommendation	Collect appropriate blood volume and number of sets per submitting lab's protocol. Collect prior to antibiotic use if possible.
Blood in sodium polyanethol sulfonate (SPS)	8 ml	Blood collected in SPS may be inoculated into blood culture bottles upon receipt at SPL.
Throat or Nasal swab	1 swab	For epidemiologic purposes only. Useful only within 24 hours of exposure.
Whole blood	1 ml	3.2% sodium citrate blood (Blue top)
Serum	1 ml	Collect in serum separator tube (SST™) or red top tube. For molecular testing.
Environmental	1 swab	

Collection and Transport of Diagnostic Samples for Botulism Testing

Suspected Food borne Botulism

Acceptable Specimens	Required Volume/Comments	Toxin Assay (T) or Culture (C) Performed
Serum (priority sample type)	5 ml (less results in incomplete testing)	T
Gastric contents	20 ml, anaerobic transport system	T, C
Vomitus	20 ml, anaerobic transport system	T, C
Stool	25-50 g (walnut-size) collected before anti-toxin treatment.	T, C
Sterile water enema	Collect using a minimal amount of water, before anti-toxin treatment.	T, C
Implicated consumed food – commercial or home-prepared	Leave foods in their original containers, if possible, or transfer to sterile, leak-proof containers. Empty containers with remnants of food are acceptable.	T, C
Unopened home-prepared food from the batch consumed by the patient	Leave foods in their original containers.	T, C
Unopened commercial products	Products are referred immediately to the FDA.	T, C

Suspected Infant Botulism

Acceptable Specimens	Required Volume/Comments	Toxin Assay (T) or Culture (C) Performed
Stool (priority sample type)	As above.	T, C
Sterile water enema	As above.	T, C
Serum	3 ml (0.5 ml will allow screening but incomplete testing)	T
Rectal Swab	More useful for culture than toxin detection.	C
Potential Sources	Include honey, opened formula (Unopened commercial products are referred to the FDA.), other foods/liquids fed to the infant. Environmental sampling is discouraged.	C

Suspected Wound Botulism

Acceptable Specimens	Required Volume/Comments	Toxin Assay (T) or Culture (C) Performed
Serum (priority sample type)	5 ml (less results in incomplete testing)	T
Wound swab	Anaerobic transport system	C
Tissue or exudate	Anaerobic transport system	C
Stool	To rule out food borne botulism. Collect as above	T,C

Suspected Intentional Toxin Release

Acceptable Specimens	Required Volume/Comments	Toxin Assay (T) or Culture (C) Performed
Clinical material	Serum, stool, sterile water enema as above.	T
Food	As above.	T
Environmental swabs	Send in individual clean, dry containers	T

Shipping Requirements (Botulism)

Notify the laboratory in advance.

- Collect and transport clinical samples in sterile, leak-proof containers.
- Leave foods in their original containers, if possible, or place in sterile leak-proof, unbreakable containers. Place each container in a separate sealed plastic bag to prevent cross-contamination during shipping. Label completely.
- Ship by the most rapid means available.
- Store and ship specimens in anaerobic transport systems at room temperature.
- Store and ship all other specimens at 4°C.
- Freezing should be avoided as it decreases recovery of *C. botulinum* and may decrease toxin activity. However, if a delay of more than several days cannot be avoided, freeze samples for storage and ship frozen.

Diagnostic Specimen Requirements for *Francisella tularensis*

Type of Infection	Specimen type	Minimum Volume	Collection Comments
Pulmonary	Sputum, throat swab, tracheal aspirates, bronchoalveolar wash, etc.	1 ml	
	Nasal swab	1 swab	For epidemiologic purposes only. Useful only within 24 hours of exposure.
Ulceroglandular	Ulcer scraping, biopsy, or swab (eye)	1 g 1 swab	Specimen from advancing edge of the lesion not central necrotic area, which is usually secondarily infected.
Glandular	Lymph node aspirate, tissue	1ml 1-5g	
Septicemia	Blood culture	Refer to manufacturer's recommendation	Collect appropriate blood volume and number of sets per submitting lab's protocol. Collect before antibiotic use if possible. Most likely to be positive in later stage of disease.
	Blood in sodium polyanethol sulfonate (SPS)	8 ml	Blood collected in SPS may be inoculated into blood culture bottles upon receipt at SPL.
Meningitis	CSF	1 ml	Centrifuge ≥ 1 ml of fluid
Misc/Other	Whole blood	1 ml	3.2% sodium citrate blood (Blue top)
	Serum/plasma	2 ml	Acute-phase specimen should be collected ASAP after onset of disease. Convalescent-phase specimen should be collected >14 days after the acute specimen.
Postmortem	Lymph, lung, liver, spleen tissue, bone marrow, CSF	1-5g 1 ml	
Environmental	Swab	1 swab	

Diagnostic Specimen Requirements for *Yersinia pestis*

Type of Infection	Specimen type	Minimum Volume	Collection Comments
Bubonic Plague	Lymph node (bubo) aspirate	2 ml	
Septicemic Plague	Blood culture	Refer to manufacturer's recommendation	A series of three venipuncture specimens taken 15-30 minutes apart is most effective. Collect before antibiotic use if possible.
	Blood in sodium polyanethol sulfonate (SPS)	8 ml	Blood collected in SPS may be inoculated into blood culture bottles upon receipt at SPL.
Pneumonic Plague	Sputum	1 ml	"Bloody" sputum is a hallmark of this disease.
	Tracheal aspirates, bronchoalveolar wash, etc.	1 ml	Bronchial or tracheal aspirates are the specimens of choice.
	Nasal/Throat swab	1 swab	For epidemiologic purposes only. Useful only within 24 hours of exposure.
Misc/Other	Whole blood	1 ml	3.2% sodium citrate blood (Blue top)
	Serum/plasma	2 ml	Acute-phase specimen should be collected ASAP after onset of disease. Convalescent-phase specimen should be collected >14 days after the acute specimen.
Postmortem	Lymph and lung tissue, bone-marrow	1-5 g 1ml	

Diagnostic Specimen Requirements for Viruses

Type of Infection	Specimen type	Minimum Volume	Collection Comments
Equine Encephalitis VEE, EEE and WEE	0-24 hr Nasal swabs, induced respiratory secretions	1 swab or 1ml	Place swab or fluid in viral transport media.
	24-72 h Serum; throat swabs, CSF	1 swab or 1ml	If submitting for culture place in viral transport media.
	> 6 days Serum; tissue	1ml 1-5g	Serum for IgM; Pathology samples plus brain.
Viral Hemorrhagic Fevers	0-24 hr Nasal swabs, induced respiratory secretions	1 swab or 1ml	Place swab or fluid in viral transport media.
	24-72 h Serum	1ml	If submitting for culture place in viral transport media.
	> 6 days Serum; tissue	1ml 1-5g	Serum for IgM; Pathology samples plus adrenal gland.
Environmental	Swab	1 swab	

Environmental samples can be collected to determine the nature of a bio-aerosol either during, shortly after, or considerably after an event. Obviously, the sooner that the environmental sample is taken, in conjunction with early postexposure clinical samples, can help to identify the agent.

Samples taken well after an event may allow identification of the agent used. While this information would likely be too late for useful prophylactic treatment, when combined with other information, may be used in the prosecution of war crimes or other criminal proceeding. What, where, when, how, etc of the sample collection should be documented both in writing and with pictures. The types of samples taken can be extremely variable. Some of the possible samples are:

- Aerosol collections in buffer solutions
- Soil
- Swabs (Dacron or macro-foam; cotton has been shown to be PCR inhibitory)
- Dry powders
- Container of unknown substance
- Vegetation
- Food / Water
- Dead animals
- Human remains

What is collected will depend on the situation. At a minimum, anything that appears to be contaminated can be sampled with swabs or with absorbent paper or cloth. Samples should ideally be double bagged in Ziploc bags (the outside of the inner bag decontaminated with dilute bleach before placing in the second bag) labeled with time and place of collection along with any other pertinent data. As always, please contact the SPL with any questions regarding sample collection or submission.

Acronyms/Definitions

The following terms and acronyms, associated with the Special Pathogens Laboratory (SPL), are used throughout this manual. These are defined here to provide clarity of meaning.

Chain of Custody (CoC) Chain of custody refers to the act or acts undertaken by the SPL to guard and trace the flow of samples and documents from the time of receipt through permanent archiving or destruction. The SPL accomplishes chain of custody through a series of standard operating procedures.

Clinical Laboratory Improvement Program (CLIP) Program derived by the Department of Defense that allows military owned clinical laboratories to meet the Clinical Laboratory Improvement Amendments (CLIA) of 1988 standards, with certain exceptions to meet military operational requirements.

Health Insurance Portability and Accountability Act of 1996 (HIPAA) The HIPAA Privacy Rule provides federal protections for personal health information held by covered entities and gives patients an array of rights with respect to that information. At the same time, the privacy rule is balanced so that it permits the disclosure of personal health information needed for patient care and other important purposes.

ISO 17025 The ISO 17025 is one of several international standards developed by the International Organization for Standardization based in Geneva, Switzerland. The ISO 17025 establishes a laboratory's technical competence by requiring that a laboratory have a quality system in place and with such quality system in place, requires criteria specific to the go of ensuring valid test data.

Joint Biological Agent Identification and Diagnostic System (JBAIDS) JBAIDS is the United States DoD standard platform used to identify biological agents in a dual purpose role: for diagnostic applications in a clinical setting and for environmental and food sample confirmatory testing. The ruggedized JBAIDS is an open platform that analyzes 32 samples and is deployed in field hospitals, mobile analytical labs, shipboard medical labs, food and water safety test centers, research labs, and other mobile scenarios.

Laboratory Response Network (LRN) The LRN is a collaborative effort within the US federal government involving the Association of Public Health Laboratories and the Centers for Disease Control and Prevention (CDC). Most state public health laboratories participate as reference laboratories of the LRN. These facilities support hundreds of sentinel laboratories in local hospitals throughout the United States and can provide sophisticated confirmatory diagnosis and typing of biological agents that may be used in a bioterrorist attack or other bio-agent incident. The LRN was established in 1999.

Polymerase Chain Reaction (PCR) The PCR is a scientific technique in molecular biology to amplify a single or a few copies of a piece of DNA across several orders of magnitude, generating thousands to millions of copies of a particular DNA sequence.

Protected Health Information (PHI) The privacy rule protects all "*individually identifiable health information*" held or transmitted by a covered entity or its business associate, in any form or media, whether electronic, paper, or oral. The privacy rule calls this information "protected health information."

Sodium polyanethol sulfonate (SPS) SPS tubes are used for blood culture specimen collection in microbiology. Eight gentle tube inversions will prevent the blood from clotting. The blood can remain in the SPS tube before it has to be transferred to a blood culture bottle.

Special Pathogens Laboratory (SPL) Special Pathogens Laboratory receives and analyzes clinical, environmental and biological material for the presence of biological threat agents and disease causing agents. The laboratory works within the chain of command of the Diagnostic Systems Division at USAMRIID.

United States Army Medical Research Institute of Infectious Diseases (USAMRIID) USAMRIID is a basic and applied research facility located at Fort Detrick, Maryland, whose mission is to support military readiness through defensive medical research and development.

USAMRIID

SPECIAL PATHOGENS LABORATORY

*“The **Specimen Collection and Submission Manual** is a comprehensive reference guide of testing services, shipping information and submission requirements. Efforts have been made to be concise and current with all information, but as testing procedures and regulations change, some information may not be current. Many of these tests are only available with prior approval from SPL. It is impossible to address all situations in this guide, so please contact us with any questions.”*

Ft. Detrick

Maryland

