

QUALITY ASSURANCE PROJECT PLAN (QAPjP) and QA Report for Pacific 2001

Prepared by: Cristian M. Mihele
Date: June 28, 2001

2.	Principal Investigator	3
3.	Team Members	3
4.	Measurement Program.....	3
5.	Measurement Species and Units	3
6.	Representative Size Range (if PM)	3
7.	Measurement Platform (surface, airborne).....	3
8.	Measurement Sites (surface only).....	3
9.	Measurement Objective(s)	3
10.	Measurement Details.....	3
10.1.	Field Measurements	3
10.1.1.	Measurement Principle	3
10.1.2.	Instrumentation (Manufacturer/Model).....	3
10.1.3.	Flow System	3
10.1.4.	Inlet Height Above Ground (if surface)	4
10.1.5.	Nominal Flow Rate.....	4
10.1.6.	Flow Measurement/Control.....	4
10.1.7.	Flow Temperature and Pressure	4
10.1.8.	Sampling Times/Period/Frequency.....	4
10.1.9.	Sampling Methods	4
10.1.10.	Filter Type/Coating Type/Reagent Type	4
10.1.11.	Planned Changes to Instruments or Methods During Study.....	4
10.2.	Laboratory Measurements (If Applicable).....	4
10.2.1.	Laboratory Name and Address	4
10.2.2.	Analytical Method(s)	4
10.2.3.	Sample Extraction or Work-up.....	4
10.2.4.	Analytical Detection Limits.....	4
11.	Quality Assurance/Quality Control	5
11.1.	Field Quality Assurance/Quality Control	5
11.1.1.	Traceability.....	5
11.1.2.	Calibration	5
11.1.3.	Zeros and spans	5
11.1.4.	Blanks	5
11.1.5.	Field Quality Control procedures	5
11.1.6.	Precision determination	5
11.1.7.	Comparison with other measurements	5
11.1.8.	Inspections and Audits.....	5
11.2.	Laboratory Quality Assurance/Quality Control.....	5
11.2.1.	Traceability.....	5
11.2.2.	Calibration procedures.....	5
11.2.3.	Blanks	5
11.2.4.	Other lab QC.....	6
11.2.5.	Precision determination	6
11.2.6.	Comparison with other methods	6
11.2.7.	Audits	6
12.	Data Management and Quality Control	6
12.1.	Raw Data Recording.....	6
12.2.	Final Data Reporting.....	6
12.3.	Data Quality Control and Validation.....	6

QAPjP / QA Report

12.4.	Validity Flags.....	6
12.5.	Below Method Detection Limit Values	6
12.6.	Derived Parameters	6
12.7.	Explanation of Zero or Negative Data.....	7
13.	Data Quality Objectives (Pre-Study)	7
13.1.	Accuracy	7
13.2.	Precision	7
13.3.	Comparability	7
13.4.	Representativeness	7
13.5.	Completeness	7
13.6.	Other Quality Information.....	7
14.	Significant Changes to Site, Instruments or Methods During Study	8
15.	Post-study Data Quality Indicators (DQIs)	8
15.1.1.	Accuracy	8
15.1.2.	Precision	8
15.1.3.	Comparability	8
15.1.4.	Representativeness	8
15.1.5.	Completeness	8
15.2.	Blank correction (describe whether done and method used):	8
15.3.	Other Quality Information.....	8
16.	References:	8

1. Principal Investigator

Cristian M. Mihele, Meteorological Service of Canada, 4905 Dufferin Street,
Toronto, Ontario, M3H 5T4

2. Team Members

Donald R. Hastie, York University, 4700 Keele Street, Toronto, Ontario, M3J 1P3
Douglas A. Lane, Meteorological Service of Canada, 4905 Dufferin Street, Toronto,
Ontario, M3H 5T4

3. Measurement Program

Measurements of RO_x radicals by chemical amplification

4. Measurement Species and Units

RO_x=HO₂+RO₂+OH+RO where R is an organic radical (e.g. CH₃)

Units: pptv

5. Representative Size Range (if PM)

N/A

6. Measurement Platform (surface, airborne)

The sampling will take place 3 m above ground level on trailer

7. Measurement Sites (surface only)

Langley

8. Measurement Objective(s)

9. Measurement Details

9.1. Field Measurements

9.1.1. Measurement Principle

RO_x radicals: chemical amplification

9.1.2. Instrumentation (Manufacturer/Model)

The air sample is doped with high concentrations of NO and CO so that the radicals enter a chemical amplification cycle. As a result a large number (about 150) of NO₂ molecules are produced for each radical that enters the system. This number is called the chain length or amplification factor. The amount of NO₂ produced is the measured with a luminol detector. The chain length is determined using a radical calibration source based on water photolysis. The system has been built in-house and developed at York University.

9.1.3. Flow System

The air sample is draw and mixed with NO and CO in a ¼" teflon FEP reactor, passed through the luminol detector and then exhausted thorough a CO to CO₂ converter for safety reasons.

9.1.4. Inlet Height Above Ground (if surface)

3 m above ground level on top of trailer

9.1.5. Nominal Flow Rate

Flow through the Luminol detector 1.6 L min⁻¹

Reagents flows: NO 16 cm³ min⁻¹

CO 64 cm³ min⁻¹

9.1.6. Flow Measurement/Control

Mass flow controllers are used to control the sampling flow and the flow of reagents.

9.1.7. Flow Temperature and Pressure

Ambient temperature and pressure

9.1.8. Sampling Times/Period/Frequency

1 minute values averaged to 10 minute values

9.1.9. Sampling Methods

The chemical amplifier samples continuously and the raw data is acquired and stored by a PC.

9.1.10. Filter Type/Coating Type/Reagent Type

9.1.11. Planned Changes to Instruments or Methods During Study

The chemical amplifier is adversely affected by the ambient relative humidity. Changes are under way to diminish this effect. If these changes are not successful, a correction based on ambient relative humidity will be applied to the final radical data.

9.2. Laboratory Measurements (If Applicable)

9.2.1. Laboratory Name and Address

N/A

9.2.2. Analytical Method(s)

N/A

9.2.3. Sample Extraction or Work-up

N/A

9.2.4. Analytical Detection Limits

Reported detection limits: 2-6 pptv.

The detection limit for this instrument is determined by the variation of ambient ozone. It will be determined based on the radical data collected from 2am to 5am, when the radicals are not expected to be present.

10. Quality Assurance/Quality Control

10.1. Field Quality Assurance/Quality Control

10.1.1. Traceability

Flow rates will be referenced to an MSC MKS primary flow standard.
The luminol detector will be calibrated against a traceable NO₂ permeation source.

10.1.2. Calibration

The luminol detector will be calibrated against a traceable NO₂ permeation source once every 3 days or each time the luminol solution will be changed.
The chain length will be determined in the field once every 3 days with the radical calibration source based on water photolysis. If the weather is not favorable for sampling (rain) additional calibrations will be performed.

10.1.3. Zeros and spans

N/A

10.1.4. Blanks

N/A

10.1.5. Field Quality Control procedures

The luminol solution for the NO₂ detector will be prepared on site; the luminol detector will be operated in the linear regime; the raw data will be examined daily to check for possible malfunctions.

10.1.6. Precision determination

10.1.7. Comparison with other measurements

Since the number of other radical detectors in use is relatively small (1 to 2 in North America and 4 to 6 in Europe), comparisons with other similar measurements are difficult because of the budget constraints. The research groups involved in radical measurements by chemical amplification participated in a laboratory intercomparison in 1997 (PRICE) and results show that this good consistency between these instruments.

10.1.8. Inspections and Audits

10.2. Laboratory Quality Assurance/Quality Control

10.2.1. Traceability

10.2.2. Calibration procedures

10.2.3. Blanks

10.2.4. Other lab QC

10.2.5. Precision determination

10.2.6. Comparison with other methods

10.2.7. Audits

11. Data Management and Quality Control

11.1. Raw Data Recording

The output from the luminol detector is collected with a 16 bit data acquisition board from Data Translation via Labtech. The 1-second NO₂ raw data is averaged over 1 minute and, along with the chain length, is used to compute the radical concentrations.

11.2. Final Data Reporting

The final data will be reported as 10 min average.

11.3. Data Quality Control and Validation

11.4. Validity Flags

The NARSTO flags will be used for the radical data

- V0 Valid value
- V1 Valid value but comprised wholly or partially of below-MDL data
- V2 Valid estimated value
- V3 Valid interpolated value
- V4 Valid value despite failing to meet some QC or statistical criteria
- V5 Valid value but qualified because of possible contamination (e.g.,
pollution source, laboratory contamination source)
- V6 Valid value but qualified due to non-standard sampling conditions (e.g.,
instrument malfunction, sample handling)
- M1 Missing value because no value is available
- M2 Missing value because invalidated by data originator
- H1 Historical data that have not been assessed or validated

11.5. Below Method Detection Limit Values

Reported detection limits: 2-6 pptv.

The detection limit for this instrument is determined by the variation of ambient ozone. It will be determined based on the radical data collected from 2am to 5am, when the radical are not expected to be present.

11.6. Derived Parameters

11.7. Explanation of Zero or Negative Data

12. Data Quality Objectives (Pre-Study)

12.1. Accuracy

See section 11.1.7

12.2. Precision

See section 11.1.7

12.3. Comparability

See section 11.1.7

12.4. Representativeness

12.5. Completeness

RO_x radicals completeness objective =70%.

12.6. Other Quality Information

End of Pre-Study QAPjP

Start of Post-Study QA Report

13. Significant Changes to Site, Instruments or Methods During Study

14. Post-study Data Quality Indicators (DQIs)

14.1.1. Accuracy

14.1.2. Precision

14.1.3. Comparability

14.1.4. Representativeness

14.1.5. Completeness

14.2. Blank correction (describe whether done and method used):

14.3. Other Quality Information

15. References: