

<b>X17B2 SAFETY CHECKLIST</b>			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
<b>Rev: B</b>	<b>Effective: 5/24/07</b>	<b>Page 1 of 1</b>	Number: LS-SCL-0038
Reviewed by: Thomas McDonald		Reviewed by:	<b>Approved by: J. Aloï</b>

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

- 1. Valid Padlock Index and log, no locks open
- 2. X17B2 hutch – threshold door piece (painted yellow) bolted in place under each of 2 doors going into hutch.
- 3. X17B2 hutch – lead shielding mat against threshold door piece in each of 2 doors going into hutch.
- 4. X17B2 hutch interior – water cooled backstop/ Bremsstrahlung shield #4 downstream end of hutch (see photo).
  - Banded and in the beam path (**for X17B2 operation**)

OR

  - Banded and removed from beam path (**for X17B3 operation**)
- 5. X17B1 Safety Checklist completed.

OPCO / Beamline Staff: \_\_\_\_\_ Date: \_\_\_\_\_