

<b>X17B3 SAFETY CHECKLIST</b>			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
<b>Rev: B</b>	<b>Effective: 5/24/07</b>	<b>Page 1 of 1</b>	Number: LS-SCL-0039
Reviewed by: Thomas McDonald		Reviewed by:	<b>Approved by: J. Aloï</b>

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

- 1. X17B1 Safety Checklist complete and in place.
- 2. X17B2 Safety Checklist complete and in place.
- 3. X17B3 hutch – threshold door piece (painted yellow) bolted in place under each of 2 doors going into hutch.
- 4. X17B3 hutch – lead shielding mats against threshold door piece in each of 2 doors going into hutch.
- 5. X17B3 hutch interior – water cooled backstop/ Bremsstrahlung shield #5 (see photo) in place and banded at downstream end of hutch.

OPCO / Beamline Staff: \_\_\_\_\_ Date: \_\_\_\_\_