

X3A SAFETY CHECKLIST			BROOKHAVEN NATIONAL LABORATORY NATIONAL SYNCHROTRON LIGHT SOURCE
Rev: B	Effective: 5/24/2007	Page 1 of 1	Number: LS-SCL-0089
Reviewed by: Thomas McDonald		Reviewed by:	Approved by: J. Aloï

Original signatures on file.

The only official copy of this file is the one on-line in the NSLS Quality Assurance website. Before using a printed copy, verify that it is the most current version by checking the document effective date on the NSLS QA website.

1. Valid Padlock Index and log, no padlocks open.

THE FOLLOWING MUST BE CHECKED STARTING FROM THE SHIELD WALL

- 2. Bremsstrahlung shield (BS-1) in place and secure under EZ-1.
- 3. Exclusion Zone (EZ-1) in place as per photo.
- 4. Exclusion Zone (EZ-2) in place and secure.
- 5. Bremsstrahlung shield (BS-2) in place and secure as per photo.
- 6. Exclusion Zone (EZ-3) in place and secure.
- 7. Bremsstrahlung shield (BS-3) in place and secure as per photo.
- 8. Monochromator view port covered with leaded glass.
- 9. Bremsstrahlung shield (BS-4) in place and secure as per photo
- 10. Mirror tank view ports (2) covered with leaded glass.
- 11. Lead shielding in place around X3B beampipe inside X3A hutch (upstream wall).
- 12. Lead shielding in place around X3B beampipe inside X3A hutch (downstream wall).

OPCO / Beamline Staff: _____ Date: _____