



SOUTHERN RESEARCH
INSTITUTE

Affiliated with The University of Alabama at Birmingham



Spore Lot Consistency: Points to Ponder

Pat Fellows
NIAID Aerosol Challenge Workshop
3 December 2003



Considerations

- Challenge Material
- Standardized Protocols/Procedures
- Media
- Characterizations
- Purity
- Viability
- Storage Conditions



Challenge Material

- Will NIAID provide seed stock of challenge material or will they provide challenge material “ready to use”
- Protocols for harvest and purification
- Protocols to characterize lots
- Qualification standards



Media

- Agar vs. Broth
 - agar harvest more time consuming; total percent sporulation may be lower
 - differences in virulence
 - LD₅₀ values higher when spores grown on agar

	LD ₅₀ -blood agar	LD ₅₀ -L&D
Vollum 1B	2501	307
Ames	488	175
Texas-2	740	519

Ivins & Fellows, unpublished data




Characterizations

- Genetic Characterizations
 - genotypic analysis
- Phenotypic Characterizations
 - pX01 toxin gene expression
 - Protein gel analysis to identify PA, LF, EF
 - Biological activity - rat lethality; rabbit/guinea pig skin test
 - pX02 capsule gene expression
 - mucoid phenotype




Sterne - non-encapsulated



Ames - encapsulated



Mixture of encapsulated and non-encapsulated colonies - not Ames, however, this is indicative of loss of pX02 and therefore loss in virulence




India Ink Staining identifies the presence of capsule



Purity

- Colony morphology
 - free of contaminants
- Hypaque purification
 - density gradient that separates spores from vegetative cells & debris
- Microscopic evaluation
 - percent refractility
 - Malachite Green staining





SOUTHERN RESEARCH
INSTITUTE

Affiliated with The University of Alabama at Birmingham




Photo courtesy of Joel Bozue



SOUTHERN RESEARCH
INSTITUTE

Affiliated with The University of Alabama at Birmingham





Photo courtesy of Joel Bozue



Photos courtesy of Sue Welkos



Viability

- Viable counts
 - periodic plate counts
- LD₅₀ Determinations
 - initial determination
 - periodic
 - yearly
 - lot to lot bridging



Storage Conditions

- Aliquots of known concentration can be stored at -70°C with glycerol
- Batches can be stored in 1% phenol at 4°C
 - highly concentrated; not recommended for “working stock” concentrations lower than 10^7
 - type of container important
 - glass and polycarbonate appropriate storage



Challenge

- Spores must be prepared at the appropriate concentration required for challenge as well as taking into account “spray efficiency factor”
- Are spores heat shocked just prior to challenge
 - synchronizes spores for germination
 - get rid of any vegetative cells



Conclusions

Protocols and standards should be developed that consider:

- Procedures for spore harvest and purification
- Genotypic & Phenotypic characteristics
- Viability and LD₅₀ determinations
- Storage conditions
- Procedures for challenge