

Total Coliform	Total Coliform Fermentation Technique	9221 A, B, C
	Total Coliform Membrane Filter Technique	9222 A, B, C
	ONPG-MUG Test	9223
Fecal Coliforms	Fecal Coliform Procedure	9221 E
	Fecal Coliform Filter Procedure	9222 D
Heterotrophic bacteria	Pour Plate Method	9215 B
Turbidity	Nephelometric Method	2130 B

Alternative testing methods for disinfectant residuals listed at 40 CFR 141.74(a)(2)

Residual	Methodology	SM 21 st Edition ¹
Free Chlorine	Amperometric Titration	4500-Cl D
	DPD Ferrous Titrimetric	4500-Cl F
	DPD Colorimetric	4500-Cl G
	Syringaldazine (FACTS)	4500-Cl H
	Amperometric Titration	4500-Cl D
Total Chlorine	Amperometric Titration (Low level measurement)	4500-Cl E
	DPD Ferrous Titrimetric	4500-Cl F
	DPD Colorimetric	4500-Cl G
	Iodometric Electrode	4500-Cl I
	Amperometric Titration	4500-Cl O ₂ C
Chlorine Dioxide	Amperometric Titration	4500-Cl O ₂ E
	Indigo Method	4500-O ₃ B

Alternative testing methods for contaminants listed at 40 CFR 141.131(b)(1)

Contaminant	Methodology	SM 21 st Edition ¹
HAA5	LLE (diazomethane)/GC/ECD	6251 B