

Method 1622/1623 Spiking Suspension Enumeration Form: *Cryptosporidium Giardia* (accessibility version of form - text only)

Laboratory name:

Laboratory ICR ID (if applicable):

Internal spiking suspension enumeration ID number (note on Line 7 of bench sheet)

Source of stock *Cryptosporidium* oocyst/ *Giardia* cyst suspension

Stock *Cryptosporidium* oocyst/ *Giardia* cyst suspension lot number

Date that stock *Cryptosporidium* oocyst/ *Giardia* cyst suspension was received

Date of this spiking suspension enumeration

Analyst performing enumeration

header 1: Well number (or hemacytometer chamber number if not using well-slide enumeration)

1, 2, 3, 4, 5, 6, 7, 8, 9, 10

header 2: Oocysts/cysts counted

Time enumeration completed (samples must be spiked within 24 hours of enumeration)

Mean of the 10 counts

Standard deviation of the 10 counts

Relative standard deviation

Volume used per well (in microliters )

Comments:

---

July 16, 1998, Revision - Draft