# Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)



Bruce G. Weniger, MD, MPH Assistant Chief for Vaccine Development Immunization Safety Branch National Immunization Program Centers for Disease Control and Prevention

Conference: Innovative Administration Systems for Vaccines Rockville, Maryland, USA, 18-19 December 2003



## Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)

- Historical usage, efficacy and reactogenicity
- Public health need for high-speed devices for mass vaccination
- Safety concerns for multi-use-nozzle jet injectors (MUNJIs)
- New-generation disposable-cartridge jet injectors (DCJIs)
- Universal standards for DCJI cartridges



# What is Jet Injection?

- Needle-free delivery of liquid through skin or mucous membrane
- Mechanism
  - Moving piston pressurizes fluid in the dose chamber <\_\_\_</p>
    - >1,000 p.s.i. (>7,000 kPa)
    - Piston driven by release of compressed metal spring or compressed gas
    - Power sources: foot or hand pump/lever, gas cylinder, or electrical motor
  - Fluid ejected through *orifice* (diameter +/- 0.15 mm) in the *nozzle* 
    - High pressure over small surface area
    - Penetrates epidermis or mucous membrane
- Deposition to targeted tissues
  - Intramuscular
  - Subcutaneous
  - Intradermal
  - Penetration dependent on various factors
    - skin thickness, orifice size, pressure, speed, etc.
- Dose volumes often variable
  - Human devices: 0.1 to 1.0 mL
  - Veterinary devices: <1.0 to  $\ge 2$  mL



Further information: www.cdc.gov/nip/dev/jetinject.htm



# Historical origins of jet injection - I

- Technology invented in France, 1860s
  - Various names used for technique:
    - l'aqua-puncture
    - l'hydropuncture
    - *douche filiforme* [threadlike spray]
  - H. Galante et C<sup>ie</sup> device
    - Manufactured in Paris
    - Clinical studies: Dr. L. Sales-Girons, Pierrefonds-les-Bains spa
    - Treated "uncontrolled neuralgia" and other conditions
    - Injected plain water and other liquids

Citation:

www.cdc.gov/nip/dev/jetinject.htm#bibliography

X. M. BÉCLARD présente à l'Académie, de la part de M. Galante, un petit instrument pour pratiquer l'hydropuncture.

#### CORRESPONDANCE MANUSCRITE.

323

Un jet capillaire de liquide, lorsqu'il est projeté sur un point de la peau et avec la force de vingt-cinq atmosphères au moins, la perce, et le liquide est introduit dans les mailles des tissus organisés.

Il s'agissait de faire un petit instrument qui produisît ce jet et cette force; c'est ce qu'a fait M. le docteur Sales-Girons avec une scringue fabriquée par M. Galante.



Le liquide est contenu dans une petite ampoule de verre; en tirant le piston il entre dans le corps de la pompe, dont le diamètre n'excède pas 3 millimètres de calibre; en poussant avec le simple effet de la main, il sort un filet capillaire avec la pression voulue de 25 à 30 atmosphères.

Appliqué sur le point désigné de la peau, l'opération de l'hydro-puncture se pratique avec la plus grande facilité. Ce procédé est mis à profit aujourd'hui contre les névralgies rebelles, avec de l'eau pure ou des liquides médicamenteux.

#### {Béclard/Galante | 866}

# Historical origins of jet injection - II

- Modern era: "rediscovered" in USA
- Clinicians noted deposition of hydraulic fluid in worker's hand from a pinpoint leak in pressurized fluid line
- First commercial device:
  - Hypospray<sup>®</sup> (R.P. Scherer Corp., Detroit)
    - Introduced late 1940s
    - Targeted to childhood diabetes to overcome needle phobia
    - Self-administration
    - Generally single-patient use





# Historical origins of jet injection - III

• *Multi-use-nozzle jet injectors* (MUNJIs) introduced 1950s



- Developed under U.S. military contract for mass vaccination of recruits
- 600-1,000 injections/hour
- Billions of vaccine doses administered worldwide
  - INF, MEA, MEN<sub>ps</sub>, POL<sub>IPV</sub>, SMA, YEL, inter alia
- Ped-O-Jet<sup>®</sup> brand: earliest and most widely used MUNJI worldwide



Plate 17.5. Jet injectors were used to administer smallpox vaccine in the left arm and measles vaccine in the right arm. Sometimes the injections were given simultaneously. Aluminium foil covers the vaccine vials to prevent the inactivation of virus by exposure to light.

# Medicaments Delivered by Jet Injectors

- Anesthetics, local and systemic, sedation
  - e.g., lidocaine, morphine, midazolam
- Antibiotics
  - e.g., penicillin, streptomycin
- Corticosteroids
- Heparin
- Hormones
  - e.g., insulin, erythropoetin, chorionic gonadotropin, growth hormone
- Immunostimulating drugs
  - e.g., interferon
- Tuberculin
- Vaccines
- Vitamins



# Immunogenicity of Jet-injected Vaccines

Generally equivalent to or better than needle/syringe

- Effect of antigen-presenting dendritic (langerhans) cells in skin ?
- Live vaccines
  - Bacille Calmette-Guérin
  - Measles
  - Mumps
  - Measles-Mumps-Rubella
  - Measles-Smallpox (vaccinia)
  - Rubella
  - Smallpox (vaccinia)
  - Yellow fever

- Inactivated vaccines
  - Botulism
  - Cholera
  - Diphtheria tetanus pertussis
  - Hepatitis A
  - Hepatitis B
  - Influenza
  - Meningococcus A, C
  - Plague
  - Polio
  - Tetanus
  - Tularemia Typhoid





Source: Jackson, et al. Vaccine 19:4703-4709, 2001.

# Local Adverse Effects of Jet-injected Vaccines

- Slightly more reactogenic than needle-syringe (N-S)
  - Immediate local reactions
    - erythema
    - hematoma
  - Delayed local reactions
    - soreness
    - erythema
    - induration/edema
    - ecchymosis
  - Inactivated vaccines usually more reactogenic than live vaccines
- Pain generally less than or similar to N-S
  - Results vary by study
  - Patients report "no pain", "less pain", and "equal pain to N/S"
- Other local adverse events
  - Bleeding at injection site more common than with N-S
  - Laceration if movement occurs during injection
  - Traumatic injury (e.g., neuropathy) reported



# Advantages and Disadvantages of Jet Injection

## Advantages

- Administers existing, off-the-shelf vaccines
- Long history of efficacious use with many vaccine types
- Avoids need and delay in reformulating and licensing new vaccines intended for other needle-free methods
- Both end-user filling *and* manufacturer prefilling possible
- Eliminates needlestick dangers and sharps-waste burden
- Potential high rates of vaccination (>600 patients / hour)
- Disadvantages
  - Increased pain for irritating medications
    - e.g., adjuvanted vaccine
  - Potential laceration injury from improper technique
  - Potential blood cross-contamination via MUNJI devices
  - Requires reconstitution of dried vaccine into liquid



Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)

- Historical usage, efficacy and reactogenicity
- Public health need for high-speed devices for mass vaccination
- Safety concerns for multi-use-nozzle jet injectors (MUNJIs)
- New-generation disposable-cartridge jet injectors
  (DCJIs)
- Universal standards for DCJI cartridges



# Public Health Need for High-speed Devices for Mass Vaccination

# Global disease eradication

- Smallpox eradicated (1975) with MUNJIs and bifurcated needles
- Polio almost eradicated with (high-speed) oral polio vaccine
  - Next target for eradication: measles
    - Will use strategy of National Immunization Days, as with polio
    - Parenteral vaccine requires injections
    - Constraints of limited time and skilled manpower for using needle-syringes

# Pandemics and regional/local epidemics

- Influenza
- Meningococcal meningitis
- Yellow fever

# Rapid response to bioterrorism

- Parenteral vaccines: anthrax, plague, glanders, tularemia, etc.
- But not smallpox (bifurcated needles simpler and just as fast)



# Reported Measles Cases, Brazil, 1990-1996\*



Data courtesy: Glacus de Souza Brito, Ministério da Saúde, São Paulo



# Vaccine Doses Distributed in the U.S., 1962-2000



"Swine" Influenza Vaccine (A/New Jersey/8/76 - Hsw1Nsw1) 1976-1977 Mass Campaigns with Multi-use-nozzle Jet Injectors (MUNJIs)





12 Doses 0.5 mil. each D.5 mil. eac

## Published Russian Studies of Potential Biowarfare Agent Vaccination with High-speed Jet Injectors

#### ANTHRAX vaccine

Saltykov 1971, Burgasov 1973, Cherkasski 1973, Pilipenko 1974 (+ plague, tularemia), Burgasov 1976, Loktev 1980 (+ plague, tularemia)

#### BOTULISM vaccine

Agafonov 1974 (+ typhoid, tetanus)

#### PLAGUE vaccine

Agafonov 1972, Pilipenko 1974 (+ anthrax, tularemia), Agafonov 1978 (oral plague + JI typhoid, oral smallpox), Lebedinskii 1979, Loktev 1980 (+ anthrax, tularemia), Gapochko 1981 (+ smallpox), Gapochko 1992 (+ typhus, smallpox)

#### SMALLPOX vaccine

Agafonov 1969, Lev 1972, Marennikova 1972, Agafonov 1973 (+ tularemia), Balabaniuk 1974, Agafonov 1977, Agafonov 1978 (oral smallpox + typhoid, oral plague), Gapochko 1981 (+ plague), Gapochko 1992 (+ plague, typhus)

#### TYPHUS vaccine

Gapochko 1992 (+ plague, smallpox)



Прививка оспенной вакцины с помощью безыгольного инъектора.

#### TULAREMIA vaccine

Agafonov 1973 (+ smallpox), Pilipenko 1974 (+ anthrax, plague), Uglovoi 1975, Loktev 1980 (+ anthrax, plague)

#### VENEZUELAN EQUINE ENCEPHALITIS vaccine

Lukin 1977

www.cdc.gov/nip/dev/jetinject.htm#bibliography

## Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)

- Historical usage, efficacy and reactogenicity
- Public health need for high-speed devices for mass vaccination
- Safety concerns for multi-use-nozzle jet injectors (MUNJIs)
- New-generation disposable-cartridge jet injectors
  (DCJIs)
- Universal standards for DCJI cartridges



Types of Multi-use-nozzle Jet Injectors (MUNJIs) Most fill from multi-dose vial attached "on-tool"







Hypospray Professional (L) & K3 (R) models











BIP-4 BIP-100









#### CENTERS FOR DISEASE CONTROL



#### June 13, 1986 / Vol. 35 / No. 23

- 373 Hepatitis B Associated with Jet Gun Injection — California
- 376 Recommendations for Providing Dialysis Treatment to Patients Infected with HTLV-III/LAV
- 383 Gastroenteritis Outbreaks on Two Caribbean Cruise Ships
- 384 Rapid Nutritional Status Evaluation During Drought Conditions — Republic of Niger
- 387 Salmonellosis Outbreaks Associated with Commercial Frozen Pasta — Massachusetts, New Jersey, New York

CDC

#### Epidemiologic Notes and Reports

#### Hepatitis B Associated with Jet Gun Injection — California

In March 1985, during routine investigation of hepatitis B (HB) case reports, an epidemiologist at the Long Beach (California) Department of Public Health noted that three HB patients had each received injections at the same weight-reduction clinic (clinic A) before disease onset. When review of previous case records and questioning of newly reported HB patients identified five additional HB cases among clinic attendees, the California Department of Health Services joined in the investigation of the clinic on July 1, 1985.

Clinic A belonged to a chain of 29 weight-reduction clinics located throughout southern California. Attendees at the clinics typically received a series of daily parenteral injections of human chorionic gonadotropin (HCG). Injections were usually given by jet injectors (Med-E-Jet Corp, Cleveland, Ohio), although some attendees received injections with single-use disposable needles and syringes. A standard regimen consisted of 30 injections; however, individuals varied considerably in duration of treatment and number of injections received.

The investigation focused on a cohort of 341 persons who attended clinic A during the first 6 months of 1985. Clinical history, review of risk factors for acquiring hepatitis B virus



CDC

FIGURE 1. Cases of acute clinical hepatitis B in a weight-reduction clinic, by month of onset — California, 1984-1985



Source: MMWR 1986; 35: 373-376







Effect of Acetone Swabbing in Laboratory Study of HBsAg Contamination of Jet Injector Devices

## Hepatitis B surface antigen (HBsAg) positives Number / total sampled (%)

	Ped-O-Jet®		<u>Med-E-Jet®</u>		
Site Tested	<u>No Swabbing</u>	<u>Swabbing</u>	<u>No Swabbing</u>	<u>Swabbing</u>	
Next ejectate	19/50 (38%)	A 3/50 (6%)	40/50 (80%)	B 29/45 (64%)	
Nozzle exterior	8/10 (80%)	A 3/10 (30%)	9/10 (90%)	B 7/9 (78%)	
Nozzle interior	0/10 (0%)	0/10 (0%)	8/10 (80%)	B 5/9 (56%)	
Nozzle tip	not applicable	not applicable	9/10 (90%)	B 6/9 (67%)	

<sup>A</sup> Ped-O-Jet - swab vs. no swab: p < 0.05, Fishers Exact Test <sup>B</sup> Med-E-Jet - swab vs. no swab: not significant

Canter et al 1990; unpublished data, CDC, 1986-1987, courtesy Walter W. Bond



# Jet Injectors Study, Brazil, 1989-1991 © Occult Blood By Dipstick\* In Vaccine Samples

	No. of	Occult blood	
Location	vaccinees	detected	<b>Percent</b>
Sao Paulo	1 1 0 0		
and Recife	1193	2	0.2
Amazon	30	2	6.6
Sao Paulo	1662	24	1.4
TOTAL	2885	28	1.0

\* sensitivity 0.002-0.01 uL BLOOD/0.5 mL DOSE

de Souza Brito, et al. Int'l Conf AIDS, Yokohama, 1996 abst. PC0132



# Jet Injection of Liquid Vaccines: Latest Policies on Multi-use-nozzle Jet Injectors

## US Department of Defense - 1997

"Accordingly, we have sent a recommendation through the DoD Medical Products Quality Control System (MPQCS) that the <u>use</u> of these products, regardless of manufacturer, <u>be discontinued</u> until assurances of their safety are received." [Defense Logistics Agency, DPSC-M, 9 December 1997]

## WHO - 1998

"...until safe needle-free injectors are identified through independent safety testing, only needles and syringes should be used for immunization, according to current WHO policy 'Safety of Injections WHO/EPI/LHIS/97.04 and 'Safety of Injections in Immunization Programmes' WHO/EPI/LHIS/96.05." [Note for the Record, 25 Mar 1998]

## • CDC (ACIP) - 2002

"...the <u>use</u> of existing multiple-use-nozzle JIs <u>should be limited</u> .... [and] considered when the theoretical risk for bloodborne disease transmission is outweighed by the benefits of rapid vaccination with limited manpower in responding to serious disease threats (e.g., pandemic influenza or bioterrorism event), and by any competing risks of iatrogenic or occupational infections resulting from conventional needles and syringes." [MMWR 2002;51(RR-2):12-13]



# In Vivo MUNJI Safety Evaluation Methodology

## Animal models

calves, piglets, human HBsAg+ carriers

## Study locations:

- Public Health Laboratory Serv., London
- Univ. Florida Sch. Vet. Med., Gainesville
- São Paulo state, Brazil
- Procedure
  - Use MUNJI to inject subject with buffered saline
  - Next ejectate(s) into specimen container
    - Represents what next vaccinee would receive



- Repeat, to collect ~100 ejectates (fewer in human subjects)
  - Sterilize injector head between human subjects
- Analyze ejectate for extremely small quantities of blood
  - Direct and indirect ELISA
  - Serum albumin is target antigen as surrogate for blood
  - Plus HBsAg+ in human study



# A model to assess the infection potential of jet injectors used in mass immunisation

P.N. Hoffman<sup>a,\*</sup>, R.A. Abuknesha<sup>b</sup>, N.J. Andrews<sup>c</sup>, D. Samuel<sup>d</sup>, J.S. Lloyd<sup>e,1</sup>

<sup>a</sup> Laboratory of Hospital Infection, Central Public Health Laboratory, 61 Colindale Avenue, London NW9 5HT, UK <sup>b</sup> Division of Life Sciences, King's College London, University of London, London SE1 8WA, UK <sup>c</sup> PHLS Statistics Unit, 61 Colindale Avenue, London NW9 5EQ, UK <sup>d</sup> Enteric and Respiratory Virus Laboratory, Central Public Health Laboratory, London NW9 5HT, UK <sup>c</sup> GPV, WHO, 1211 Geneva 27, Switzerland

Received 16 October 2000; received in revised form 8 March 2001; accepted 13 March 2001

#### Abstract

Jet injectors are needleless injectors that penetrate skin with high-pressure fluid. They have potential advantages over needles and syringes in mass immunisation programs, but concerns over their capacity to transfer blood-borne viruses have been a barrier to acceptance. Hepatitis B infection can transmit in 10 pl of blood; detection of such low volumes presents severe difficulties to such assessments. A model to assess jet injector safety was developed using injection of an inert buffer into calves and assaying the next injector discharge, representing the next dose of vaccine, for blood using a highly sensitive ELISA. Four injectors were tested: two with reusable heads and direct skin contact, one with single-use injector heads and one where the injector head discharged at a distance from the skin. All injectors tested transmitted significant (over 10 pl) volumes of blood; the volumes and frequency of contamination varied with injector. The source of the contamination was consistent with contamination by efflux of injected fluid and blood from the pressurised pocket in tissue that is formed during injection. This insight should inform the design of safe jet injectors. © 2001 Elsevier Science Ltd. All rights reserved.



# In vivo Model for Evaluating Safety of Mult-use-nozzle Jet Injectors

Blood quantity (pL) detected in next discharge (? >10 unsafe ?)

Injector	Injections Sampled	0-9.9	10-49.9	<u>&gt;</u> 50
"A"	114	75	20	19
		(66%)	(18%)	(17%)
"B"	48	2	5	41
		(4%)	(10%)	(85%)
"C"	119	69	26	24
		(58%)	(22%)	(20%)

Source: Hoffman, et al. *Vaccine* 2001;19:4020-27.



# In vivo Model for Evaluating Safety of Multi-use-nozzle Jet Injectors

Ejectates with blood >10pL / 0.5mL

Injector	No. of Animal Injections	Pre- Inject	Post 1 <sup>st</sup>	Post 2 <sup>nd</sup>	Post 3 <sup>rd</sup>	Post 4 <sup>th</sup>	Post 5 <sup>th</sup>	Post 6 <sup>th</sup>
"D"	20	7	13	8	1	1	2	1
	(100%)	(35%)	(65%)	(40%)	(5%)	(5%)	(10%)	(5%)

Source: Hoffman, et al. *Vaccine* 2001;19:4020-27.



In vivo Model for Evaluating Safety of Multi-use-nozzle Jet Injectors

## Results

- 5 devices, 3 models, 3 sites
  - bovine published
  - porcine, human unpublished
- All produced ≥1 "unclean" ejectate(s) following up to ~100 in vivo injections

# Caveats and limitations

To be discussed in later presentation from WHO



## Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)

- ✓ Historical usage, efficacy and reactogenicity
- Public health need for high-speed devices for mass vaccination
- Safety concerns for multi-use-nozzle jet injectors (MUNJIs)
- New-generation disposable-cartridge jet injectors (DCJIs)
- Universal standards for DCJI cartridges



Slow-speed, Disposable Cartridge Jet Injectors (DCJIs, "dick-jees") Most filled by end-user from "off-tool" single- or multi-dose vial:



Some investigational devices intended for manufacturer pre-filling: Entire device 1-use disposable (high speed?):







## High-speed, Disposable Cartridge Jet Injectors (DCJIs) Early MEDiVAX<sup>™</sup> prototype by PATH filled from multi-dose vial "on tool":



### Investigational LectraJet<sup>®</sup> by DCI, Inc. filled from multi-dose vial "off tool":



Prototype Manufacturer-prefilled Auto-reconstitution Disposable Cartridges for Jet Injection Developed under CDC SBIR contract p.o. 65819 by Creare, Inc., Hanover, NH, USA (www.creare.com)





Vaccine reconstituted with diluent

After injection, obturator jams to prevent refilling



## Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)

- ✓ Historical usage, efficacy and reactogenicity
- Public health need for high-speed devices for mass vaccination
- Safety concerns for *multi-use-nozzle jet injectors* (MUNJIs)
- New-generation disposable-cartridge jet injectors (DCJIs)
- Universal standards for DCJI cartridges





Manufacturer-pre-filled *Imule*<sup>™</sup> cartridges/vials for *Mini-Imojet*<sup>™</sup> jet injector (Institut Mérieux, later Pasteur Mérieux Connaught, now Aventis Pasteur)



Imule<sup>™</sup> cartridge/vial for use with Mini-Imojet<sup>™</sup> injector (Institut Mérieux, now Aventis Pasteur)



# 10 20 30 40







## Jet Injection of Vaccines: Overview and Challenges for Mass Vaccination with Jet Injectors (JIs)

- Historical usage, efficacy and reactogenicity
- Public health need for highspeed devices for mass vaccination
- Safety concerns for multi-usenozzle jet injectors (MUNJIs)
- New-generation disposablecartridge jet injectors (DCJIs)
- Universal standards for DCJI cartridges





