



Acquisition of digital chest images for pneumoconiosis classification: Methods, procedures, and hardware

Ehsan Samei, PhD

Duke Advanced Imaging Laboratories
Departments of Radiology, Physics,
Biomedical Engineering, and Medical Physics
Duke University Medical Center

Transition to the digital era

Letter



Email

Hardcopy



Softcopy

Analog



Digital

Why digital radiography?

- Four function of film:
 - Image capture
 - Image processing
 - Image display
 - Image storage
- Separation of function enabling optimization of each towards improved image quality

Advantages of digital radiography

- Improved dynamic range
 - Toleration for over-/under-exposure
- Image post-processing
 - Improved visualization
- Digital format
 - Enabling quantification and digital analysis
 - Electronic archival and distribution


Digital radiography in classification of pneumoconiosis?

Notable advantages in
potentially providing accessible,
standardized image data for
visual interpretation or
automated classification


Digital radiography

- Radiography in which the image data are **sampled** into discrete elements in spatial and intensity dimension

Sampling



Pixel = a single spatial/intensity element
(~ atoms in matter)




Common aspects of DR

1. Imaging geometry
2. Detrimental effect of x-ray scatter
3. Analog x-ray capture followed by digitization
4. Digital image
5. Required image pre- and post-processing

Differences between DR systems



1. Detector technologies
2. Inherent image quality attributes
3. Reported exposure indicator
4. Required exposure
5. Image post-processing and appearance

DR technologies




Computed Radiography (CR)	Indirect Flat-panel	Direct Flat-panel	CCD or CMOS-based
PSL phosphor (eg, BaFBr)	CsI or GSO phosphor	a-Se	CsI or GSO phosphor
PSL light-guide	Contact layer	None	Lens or fiber-optic taper
PSL signal digitiz.	Photo diode /TFT arrays	TFT arrays	CCD or CMOS


DR offerings




CR 101


- Invented 1975 by GW Luckey at Kodak
- First commercial unit in 1983 (Fuji FCR 101)
- Most common digital radiographic modality
- Easy retrofit and portable application
- The most economical digital option








ready for use



unexposed
image plate




erasure



fluorescent
lamp

ADC




display

How CR Works?


Photostimulable luminescence (PSL)

1. Phosphor excitation (by x-rays)
2. Formation of latent image
3. Stimulation by light (scanning laser)
4. Emission of light (PSL signal)
5. Light collection (PSL signal)
6. Signal digitization and image formation



CR image quality

- Inherent image quality governed by lateral spread of the laser light, and phosphor thickness




Flat-panel 101


- Technology development facilitated with advancements in TFT array fabrications for flat-panel displays
- First detectors introduced in late 90's
- Quality and speed out-performs CR
- Cost is still the prohibitive utilization factor

Two types of flat-panel detectors


Direct:






Indirect:



Detector structure



Detector structure




How flat-panel detectors work?

1. Pixels are formed by the photodiode/TFT elements coupled with unpixelated capture medium
2. X-ray energy conversion to charge (with intermediary light conversion for indirect detectors)
3. Charge collected in pixel capacitors
4. Pixels read line-by-line after the exposure

Flat-panel image quality

- Inherent image quality governed by lateral spread of light/charge, and phosphor/photo-conductor thickness




CCD/CMOS-based DR 101


- First systems introduced in 90's
- Phosphor-based
- Based on lower-cost Charged Couple Device (CCD) and Complementary metal-oxide-semiconductor (CMOS) light sensors
- Small sensor size requires demagnification
- Quality comparable to CR
- Lower-cost alternative to flat-panel DR

CCD/CMOS-based DR structure

Fiber-optic coupler



Lens coupler




How CCD/CMOS-based DR detectors work?

1. X-ray capture by a phosphor
2. Recording the emitted light by the CCD/CMOS sensor
3. Pixels read line-by-line after the exposure

CCD/CMOS-based DR image quality

- Inherent image quality governed by lateral spread of light, phosphor thickness, and loss of light (“quantum sink”)




Effect of scatter


- An ever present attribute of chest radiography (analog AND digital)
 - Scatter Fraction: 60-90%
 - Scatter Fraction: 40-70% with grid
- Reduces image quality
 - Reduces SNR^2/dose by 2.5-5 times!

Samei et al, *Med Phys*, Sept 2004; Boyce and Samei, *Med Phys*, April 2006.


Full-field DR



Slot-scan DR




DQE




- Full Field System
0.24 mR, 120 kVp
- Slot-Scan System
0.35 mR, 117 kVp
- Slot-Scan System
0.26 mR, 140 kVp

Scatter Fractions



- Full-Field, 120 kVp No Grid
- Full-Field, 120 kVp
- Full-Field, 140 kVp
- Slot-Scan, 117 kVp
- Slot-Scan, 140 kVp


Effective DQE (eDQE)



- Full-Field, 120 kVp No Grid
- Full-Field, 120 kVp
- Full-Field, 140 kVp
- Slot-Scan, 117 kVp
- Slot-Scan, 140 kVp

Demonstration

Full-field, 120 kVp, 0.02 mSv



Slot-scan, 140 kVp, 0.02 mSv




Conclusions

- Digital radiography offers distinct advantages
- Current commercial offerings represent differing technologies with differing image quality attributes
- An initiative involving diverse systems needs to
 - take their similarities and differences into consideration
 - employ controlled unifying conditions

Recommendations

- Standardized image acquisition and processing protocols
- Robust quality control and preventative maintenance programs
- Facility and equipment accreditation programs

Thank you for your attention



DAILabs.duhs.duke.edu