

National Institute of Standards & Technology

Certificate

Standard Reference Material 422C Radioactivity Standard for Liquid Scintillation Counting

Radionuclide Carbon-14-<u>n</u>-hexadecane

Source identification 4222C

Source description 5-ml of solution in NIST borosilicate-

glass ampoule (1)*

Solution composition 14 C-labeled <u>n</u>-hexadecane in inactive

n-hexadecane (2)

Radioactivity concentration 5.402 x 10⁴ Bq g⁻¹

Reference time 1200 EST September 3, 1990

Overall uncertainty 0.81 percent (3)

Photon-emitting impurities None observed (4)

Measuring instrument Liquid-scintillation counter (5)

Half life 5760 ± 50 years ⁽⁶⁾

This standard reference material was prepared in the Center for Radiation Research, Ionizing Radiation Division, Radioactivity Group, Dale D. Hoppes, Group Leader.

Gaithersburg, MD 20899 January, 1991 William P. Reed, Acting Chief Office of Standard Reference Materials

*Notes on back

NOTES

(1) Approximately five milliliters of solution. Ampoule specifications:

body diameter wall thickness barium content lead oxide content other heavy elements $16.5 \pm 0.5 \text{ mm}$ $0.60 \pm 0.04 \text{ mm}$ less than 2.5 percent less than 0.02 percent trace quantities

- The density of <u>n</u>-hexadecane is 0.7709 ± 0.001 g/cm³ at 25.0 °C.
- The overall uncertainty was formed by taking three times the quadratic combination of standard deviations of the mean, or approximations thereof, for the following:

a) liquid-scintillation measurements for 11 vials times 5 runs		0.03 percent
b) standardization utilizing efficiency		_
tracing technique (5)		0.20 percent
c) gravimetric measurements		0.15 percent
d) quenching in the liquid scintillator	•	0.10 percent

Limits of detection for photon-emitting impurities are:

0.02 ys⁻¹ between 100 and 1900 keV.

- The liquid-scintillation counter was standardized using a ³H radioactivity standard, by comparing the theoretical and observed spectra for ¹⁴C and ³H, as described in Coursey, B.M., et al. Appl. Radiat. Isot., <u>37</u>, 405 (1986).
- (6) Mann, W.B., Marlow, W.F., and Hughes, E.E., Int. J. Appl. Radiat. Isot. <u>11</u>, 57 (1961).

For further information please contact J.M. Calhoun, (301) 975-5538, or Dr. B.M. Coursey, (301) 975-5584.