## SUPPLEMENT 2: CHIROPRACTIC SYSTEMS

## FIELD TEST RECORD EDIT CHECKS

When performing this supplementary procedure, the following changes to the edit checks on page AR-11 of the procedures manual apply:

- 1. Verify that a "B" is entered in data item 1.
- 2. Verify that the manufacturer, model number, and serial number for the vertical cassette holder are entered in data items 4 and 5.
- 3. Verify that the certification status of the vertical cassette holder is indicated in the "other" box of data item 6.
- 4. The chosen SIDs (edit check #11) can be greater than 50 inches (125 cm) at data items 52, 57, 62, and 65. Also, an SID of 72 inches is acceptable on the field test record.
- 5. If the indicated SID at data item 9 is 72 inches (182.9 cm), then data items 57, 58, 59, 65, 66, and 67 should be blank.