

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
 PUBLIC HEALTH SERVICE  
 FOOD AND DRUG ADMINISTRATION  
**C-ARM FLUOROSCOPIC AND SPOT-FILM SYSTEMS**  
**FIELD TEST RECORD**

Print Legibly. Use Black  
 Ball Point Pen. Enter One  
 Character per Box. Do Not  
 Write in Shaded Area.

FIELD TEST SERIAL NO: (1-8)

CF

REGIONAL REVIEW (NAME)

(Use Form FDA 2782, Field Test Record Continuation, if more space is needed.)

Card No.

(9-10)

Test Procedure:

1	CF
11	13

2. System Type  M--Mobile  S--Stationary  
 14

Beam Limiting Device Information

3. BLD Manufactured after 5/22/79  Y--yes  N--no  
 15

4. BLD Type  1--Fixed Aperture  2--Stepless  
 16

Component Certification Information

5. Indicate the Status of Each as Follows:

C--Certified

V--Certified with a Variance

N--Not Certified

X--Not Present or not applicable

17 Tub Housing Assembly

18 X-Ray Control

19 Spot Film Device  
(after 4/77)

20 Beam Limiting Device

21 High Voltage Generator

22 Fluoroscopic Imaging  
Assembly

23 Table


24 Image Intensifier (after 4/77)

(before 4/77)

10

Initial Set-up (Fluoroscopic Mode) And Surveyor Protection Test

MDH (Exposure)



6. System Hazardous  Y--yes  N--no  
 25

If yes, describe hazard in remarks and discontinue testing.

7. Warning Label Present  Y--yes  N--no  
 26

Fluoroscopic X-ray Field/Image Receptor Alignment




Image Dimension

X-ray Field Dimension

1/4 8.  27  29 in

12.  30  32 cm

2/1 9.  33  35 in

13.  36  38 cm

3/2 10.  39  41 in

14.  42  44 cm

4/3 11.  45  47 in

15.  48  50 cm

NO MAGNIFICATION

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
 PUBLIC HEALTH SERVICE  
 FOOD AND DRUG ADMINISTRATION  
**C-ARM FLUOROSCOPIC AND SPOT-FILM SYSTEMS**  
**FIELD TEST RECORD**

Print Legibly. Use Black  
 Ball Point Pen. Enter One  
 Character per Box. Do Not  
 Write in Shaded Area.

FIELD TEST SERIAL NO: (1-8)

CF

REGIONAL REVIEW (NAME)

(Use Form FDA 2782, Field Test Record Continuation, if more space is needed.)

10

**Fluoroscopic X-ray Field/Image Receptor Alignment**

	Image Dimension	X-ray Field Dimension		
1/4	16. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> in 51 53	20. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> cm 54 56	24. Visible Area: <input type="text"/> 1--Circular 75 2--Rectangular	
2/1	17. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> in 57 59	21. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> cm 60 62		25. Tube Potential and Current Continuously Indicated During Exposure <input type="text"/> Y--yes N--no 76
3/2	18. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> in 63 65	22. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> cm 66 68		
4/3	19. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> in 69 71	23. <input type="text"/> <input type="text"/> . <input type="text"/> <input type="text"/> cm 72 74		26. Maximum Cumulative Setting for Fluoro Timer is 5 minutes or less <input type="text"/> Y--yes N--no 77

**Primary Protective Barrier/X-ray Field Size Comparison**

27. Diameter of the Image Intensifier Housing   .   cm  
11 13

28. Distance from face of I.I. to base of stand   .   cm  
14 16

**Minimum Fluoro Field Size**

29. X-ray field Image length   .   cm  
17 19

30. X-ray field Image width   .   cm  
20 22

31. Field Shape  1--Circular 2--Rectangular  
23

11

**Entrance Exposure rate: (BLD fully open, no magnification)**

32. Fluoroscopic Technique Factor Control Type  
  
 24  
 M--Manual Only  
 A--Automatic Only  
 B--Both Manual and Automatic

MDH (Exposure rate)

MANUAL	33. <input type="text"/> <input type="text"/> kVp 25 27	34. <input type="text"/> <input type="text"/> mA 28 29	35. <input type="text"/> <input type="text"/> R/min 30 33
	36. Is a high level control present? <input type="text"/> Y--YES N--NO 34		37. Continuous Audible Signal Upon Activation of High Level Control <input type="text"/> Y--YES N--NO X--N/A 35
AUTO	38. <input type="text"/> <input type="text"/> kVp 36 38	39. <input type="text"/> <input type="text"/> mA 39 40	40. <input type="text"/> <input type="text"/> R/min 41 44
	41. Is a high level control present? <input type="text"/> Y--YES N--NO 45		42. Continuous Audible Signal Upon Activation of High Level Control <input type="text"/> Y--YES N--NO X--N/A 46

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
 PUBLIC HEALTH SERVICE  
 FOOD AND DRUG ADMINISTRATION  
**C-ARM FLUOROSCOPIC AND SPOT-FILM SYSTEMS**  
**FIELD TEST RECORD**

(Use Form FDA 2782, Field Test Record Continuation, if more space is needed.)

Print Legibly. Use Black  
 Ball Point Pen. Enter One  
 Character per Box. Do Not  
 Write in Shaded Area.

FIELD TEST SERIAL NO: (1-8)

CF

REGIONAL REVIEW (NAME)

11


**Source -- Image Distance and Minimum Source to Skin Distance**

43. Distance from face of SSD spacer (or from face of BLD if spacer is not present) to top of brass strips    .   cm  
 47 49

44. Outside Separation of Image of Focal Spot Strips    .   cm  
 50 52

**Beam Quality**

47.    kVp  
 53 55



MDH (Pulse Exposure  
 4.5mm Al)

46.    mR @ 4.5 mm Al  
 56 58

47.    sec  
 59 61

48.    mR @ 3.5 mm Al  
 62 64

49.    sec  
 65 67

50.    mR @ 2.5 mm Al  
 68 70

51.    sec  
 71 73

52.    mR @ 1.5 mm Al  
 74 76

53.    sec  
 77 79

54.    mR @ 0.0 mm Al  
 11 13

55.    sec  
 14 16

44. At End of Present Cumulative Time Interval Either  
 Continuous Audible Signal and/or Termination of  
 X-rays  Y--YES  
 17 N--NO

12

**Spot-Film Mode**

57. Type of spot-film device   
 18  
 1. Clip-on cassette holder  
 2. Cut-film changer  
 3. N/A (only a fluorographic camera provided)

58. Spot-film length    cm  
 19 21

59. Spot-film width    cm  
 22 24

60. Technique factors fixed or selectable indicated before exposure  Y=yes N=no  
 25

61. A visible beam-on indication during exposure  Y=yes N=no  
 26

62. An audible indication of exposure termination  Y=yes N=no  
 27

63. Exposure terminated at a preset time or mAs  Y=yes N=no  
 28

64. Means provided for maintaining the fluoro field size during spot-filming  Y=yes N=no X-N/A  
 29

DEPARTMENT OF HEALTH AND HUMAN SERVICES  
 PUBLIC HEALTH SERVICE  
 FOOD AND DRUG ADMINISTRATION  
**C-ARM FLUOROSCOPIC AND SPOT-FILM SYSTEMS**  
**FIELD TEST RECORD**

Print Legibly. Use Black  
 Ball Point Pen. Enter One  
 Character per Box. Do Not  
 Write in Shaded Area.

FIELD TEST SERIAL NO: (1-8)


CF

REGIONAL REVIEW (NAME)

(Use Form FDA 2782, Field Test Record Continuation, if more space is needed.)

12

**Spot-Film Reproducibility and Linearity**



65. MDH Threshold Setting.  $\frac{\quad}{\quad}$  .  $\frac{\quad}{\quad}$  cm  
 0.5-3 phase, 0.2-1 phase  $\frac{0}{\quad}$  .  $\frac{\quad}{30}$

66. Generator Type  $\frac{\quad}{31}$  1-single phase  
 3-three phase

67. Technique Factors  
 $\frac{\quad}{32}$  .  $\frac{\quad}{34}$  kVp  $\frac{\quad}{35}$  .  $\frac{\quad}{37}$  mAs  
 $\frac{\quad}{38}$  .  $\frac{\quad}{40}$  mA  $\frac{\quad}{41}$  .  $\frac{\quad}{43}$  sec.

13

**Reproducibility (no AI in beam)**

68. $\frac{\quad}{44}$ . $\frac{\quad}{48}$ mR	69. $\frac{\quad}{49}$ . $\frac{\quad}{51}$ msec	78. $\frac{\quad}{52}$ . $\frac{\quad}{56}$ mR	79. $\frac{\quad}{57}$ . $\frac{\quad}{59}$ msec
70. $\frac{\quad}{60q}$ . $\frac{\quad}{64}$ mR	71. $\frac{\quad}{65}$ . $\frac{\quad}{67}$ msec	80. $\frac{\quad}{68}$ . $\frac{\quad}{72}$ mR	81. $\frac{\quad}{73}$ . $\frac{\quad}{75}$ msec
72. $\frac{\quad}{11}$ . $\frac{\quad}{15}$ mR	73. $\frac{\quad}{16}$ . $\frac{\quad}{18}$ msec	82. $\frac{\quad}{19}$ . $\frac{\quad}{23}$ mR	83. $\frac{\quad}{24}$ . $\frac{\quad}{26}$ msec
74. $\frac{\quad}{27}$ . $\frac{\quad}{31}$ mR	75. $\frac{\quad}{32}$ . $\frac{\quad}{34}$ msec	84. $\frac{\quad}{32}$ . $\frac{\quad}{34}$ mR	85. $\frac{\quad}{40}$ . $\frac{\quad}{42}$ msec
76. $\frac{\quad}{43}$ . $\frac{\quad}{47}$ mR	77. $\frac{\quad}{48}$ . $\frac{\quad}{50}$ msec	86. $\frac{\quad}{51}$ . $\frac{\quad}{55}$ mR	87. $\frac{\quad}{56}$ . $\frac{\quad}{58}$ msec

**Linearity**

88.  $\frac{\quad}{59}$  .  $\frac{\quad}{62}$  mA

If change in mA causes a kVp shift, readjust kVp (if possible) to value selected at Item

89.  $\frac{\quad}{63}$  .  $\frac{\quad}{67}$  mR

90.  $\frac{\quad}{68}$  .  $\frac{\quad}{72}$  mR

91.  $\frac{\quad}{11}$  .  $\frac{\quad}{15}$  mR

92.  $\frac{\quad}{16}$  .  $\frac{\quad}{20}$  mR

14

**X-ray Field/Spot-film Size Comparison**

93. Distance from the spot-film plane to the SSD spacer (or from face of BLD if spacer is not present)  $\frac{\quad}{21}$  .  $\frac{\quad}{23}$  cm

94. Distance from Image receptor to bottom of test stand  $\frac{\quad}{24}$  .  $\frac{\quad}{26}$  cm

95. Length  $\frac{\quad}{28}$  .  $\frac{\quad}{30}$  cm (For circular field enter diameter twice)

96. Width  $\frac{\quad}{31}$  .  $\frac{\quad}{33}$  cm

REMARKS