

CDC Experience with Laboratory Tests Used to Investigate Incidents of Transfusion-Transmitted Malaria

Marianna Wilson
Division of Parasitic Diseases
Centers for Disease Control and Prevention
Atlanta, GA



The findings and conclusions in this presentation have not been formally disseminated by the Centers for Disease Control and Prevention and should not be construed to represent any agency determination or policy.

DIAGNOSTIC TOOLS AVAILABLE AT CDC

- **Blood film examination**
- **Antibody detection - IFA**
- **DNA detection - PCR**

MALARIA IFA TEST

Antigen: washed cell thick smears of RBCs infected with one species

- *P. falciparum*
- *P. malariae*
- *P. ovale*
- *P. vivax*

MALARIA IFA TEST

- **Sensitivity** $89/92 = 98\%$
 - Based on Pv and Pf infections
- **Specificity** $1/184 = 99.5\%$
 - *Babesia* may cross-react


MALARIA IFA TEST

Human Experimental Infections

- Parasitemia precedes antibody
 - *P. falciparum* 4-6 days
 - *P. vivax* 2-6 days
 - *P. malariae* 4-6 days
- If parasitemia is suppressed, may still develop detectable antibody

MALARIA IFA TEST

Antibody Persistence Post-Rx



Typical Malaria IFA Results

<i>Slide</i>	<i>IFA</i>			
	<i>Diagnosis</i>	<i>P. falcip</i>	<i>P. malariae</i>	<i>P. ovale</i>
<i>P. falciparum</i>	1:16384	1:4096	1:1024	1:64
<i>P. falciparum</i>	1:1024	Neg	Neg	Neg
<i>P. ovale</i>	Neg	Neg	1:256	Neg
<i>P. vivax</i>	Neg	Neg	1:64	1:1024
<i>P. malariae</i>	1:1024	1:1024	1:64	1:64

CDC PCR FOR MALARIA

Two-step nested PCR

- **Genus-specific primers followed by**
- **Species-specific primers**
(Snounou et al, 1993)

CDC PCR FOR MALARIA

Sensitivity & Specificity

100%+


based on smear results

ALGORITHM FOR LABORATORY TESTING OF SUSPECT DONORS

- Screen with serology
- PCR on IFA positives
- Blood film examination on IFA positives

TRANSFUSION TRANSMITTED MALARIA, 1963-1998

How 65 implicated donors were identified



TRANSFUSION TRANSMITTED MALARIA, 1963-1998

In implicated donors, where information was available:

- 98% had positive serology
- 33% had positive blood smear


Laboratory investigation of blood donors, 1996-2005

<i>Year</i>	<i>Species</i>	<i>Donors</i>	<i>IFA +</i>	<i>PCR +</i>
1996	Falcip	7	1	1
1998	Falcip	4	1	1
2002	Malariae	10	1	0
2003	Falcip	2	1	0
2005	Falcip	150	1/100	0


Laboratory investigation of blood donors, 1996-2005

<i>Year</i>	<i>Tx date</i>	<i>Spec date</i>	<i>Slide result</i>	<i>PCR result</i>
1996	11/12/96	3/1997	Pos	Pos
1998	1/15/98	?1/15/98	Neg	Pos
2002	5/1/02	7/11/02	Neg	Neg
2003	3/1/03	4/03	Neg	Neg
2005	6/30/05	8/29/05	Neg	Neg

Detection of Parasites/Parasite Products



Detection of Parasites/Parasite Products



Detection of 10 parasites/unit requires a sensitivity:
-4,000 times better than current PCR
-200,000 times better than microscopy

Selected References

IFA

- Collins et al, Amer J Trop Med Hyg 1964;13:256
- Collins et al, Amer J Trop Med Hyg 1964;13:777
- Sulzer et al, Amer J Trop Med Hyg 1969;18:199
- Wilson et al, Amer J Trop Med Hyg 1970;19:401
- Gleason et al, Amer J Trop Med Hyg 1971;20:10

PCR

- Snounou et al, Mol Biochem Parasitol 1993;61:315
- Johnston et al, J Clin Microbiol 2006;44:1087