

# NKDEP/IFCC Conference to Address Standardization of Urine Albumin/Creatinine Measurement and Reporting

March 27-28, 2007
Washington, DC, USA
Minutes on NKDEP/LWG web site



# **Participants**

Greg Miller, PhD

Virginia Commonwealth University Chair LWG, NKDEP

**Andrew Narva, MD** 

Director, National Kidney Disease Education Program

Elisa Gladstone, MPH

Associate Director, National Kidney Disease Education Program

David Bruns, MD

University of Virginia

David Bunk, PhD

National Institute of Standards and Technology

Gary Curhan, MD, ScD

Brigham and Women's Hospital; Harvard Medical School John Eckfeldt, MD, PhD

University of Minnesota

Jim Fleming, PhD, FACB

Laboratory Corporation of America

Neil Greenberg, PhD

Ortho Clinical Diagnostics

Glen Hortin, MD, PhD

National Institutes of Health, Clinical Center

Yoshihisa Itoh, MD, DrSci

Asahikawa Medical College, Japan

John Lieske, MD

Mayo Clinic

Matthew McQueen, MB, ChB, PhD, FCACB, FRCPC

Hamilton Regional Laboratory Medicine Program, Canada Gary Myers, PhD

Centers for Disease Control and Prevention

Mauro Panteghini, PhD

University of Milan, Italy

Sverre Sandberg, MD, PhD

Haukeland University Hospital and The Norwegian Quality Improvement of Primary Care Laboratories (NOKLUS), Norway

Heinz Schimmel, PhD

Institute for Reference Materials and Measurements of the European Commission, Belgium

David Seccombe, MD, PhD, FRCPC

Canadian External Quality
Assessment Laboratory
(CEQAL), Canada

Jack Zakowski, PhD

Beckman Coulter, Inc



### **Objectives: Urine albumin**

- Frame the issues
- Develop a path forward to improve standardization
- Impact successful implementation of clinical practice guidelines



# **Agenda**

- Review current status of measurement and reporting
- Discussion groups to address specific issues
  - Define the measurand; what molecule(s) should be measured
  - Standardization of sample requirement and reporting
  - Calibration traceability; reference system
- Reports from discussion groups
- Consensus of current status and path forward



#### **Current status of measurement**

#### PT/EQA suggested:

- A range of results for the same sample
  - Influenced by non-commutability
- Urine dipstick results were highly variable
- A range of imprecision; with most acceptable within a method
- A variety of reporting units for albumin and albumin/creatinine ratio



#### **Current status of measurement**

#### Practice surveys suggested:

- No uniformity in sample type (timed vs. first morning vs. random collection)
- Lab recommendations for sample type were not followed
- Variability in lab provided decision limits for albumin and albumin/creatinine ratio
- Physicians react to differences that are smaller than is analytically justified
- Guidelines for the diagnosis of albuminuria were not followed
- Variation regarding treatment once albuminuria has been diagnosed



#### Define the measurand

- Fragmented and dimeric forms of albumin are found in urine
  - Unclear if native or due to post collection handling and storage (including freezing)
  - Some polyclonal Ab methods were relatively insensitive to fragmented forms
- Influence of pH, osmolality, contact with sediment and adsorption, centrifugation, and other sample handling factors is not well understood



## Sample collection issues

- Physiologic effects that influence a result
  - Diurnal, postural, exercise influences
  - Biological variation of normal vs. impaired kidney function
  - First morning void vs. random void vs. stress condition
  - Applicable restrictions, e.g. physiological steady state, non-menstruating, free of concurrent infections, drugs, dietary supplements
- Sample handling that influences a result
  - Non-specific binding to the collection container
  - Degradation and fragmentation during storage and freeze-thaw



# Reporting issues

- Nomenclature
  - Urine albumin excretion rate (UAE) is what is physiologically desired and correlates with CKD and CVD
  - "microalbumin" is confused with a different type/size of albumin rather than detecting a low concentration
  - "urine albumin" is simple and may be preferable
- Urine albumin/creatinine ratio (ACR) is used as a surrogate for UAE
  - A variety of units are used (e.g. mg/g, mg/mol) with different numeric values for decision points
- Time period of collection influences decision point
- Age, gender, race influences decision point



#### **Measurement issues**

- Influence of: labels; competitive vs. non-competitive assays; homogeneous vs. heterogeneous assays; polyclonal vs. monoclonal antibodies
- Influence of epitope for the Ab and albumin fragmented forms
- Calibration traceability to diluted CRM 470 (ERM DA470) serum protein reference material
  - Details of traceability design and measurement implementation appear to influence calibration uniformity
- No urine albumin reference material, nor reference measurement procedure at this time
- No urine creatinine reference material; do have a RMP for urine creatinine listed by JCTLM



# Reference system

# Mayo Clinic LC-MS candidate reference measurement procedure

- Measures N-terminal 24 peptide fragment following digestion
- Internal standards: BSA; 15N-labeled recombinant HSA
- Inter-assay CV is 12-15%
- LOQ is 10 mg/L
- Correlation to immunoassay is better than to HPLC
- Further development and characterization is in progress



# Reference system

#### JCCLS candidate urine albumin reference material

- Purified monomeric human serum albumin (>97.5% pure)
- Characterized by amino acid analysis, SDS-PAGE, MS, HPLC
- Prepared in a buffered aqueous matrix
- Value assignment by traceability to diluted CRM 470 (human serum albumin) using selected routine measurement procedures
- Preliminary evaluation of improvement in uniformity of results among routine methods



- Publish a report describing current status and recommendations for addressing the issues
  - Establish clinical requirements for measurement performance
  - Develop recommendations for nomenclature and reporting
- Define specifications for routine method robustness to range of matrix components in urine samples



- Understand current method performance
  - Conduct a round robin evaluation of routine and higher order methods using a panel of native urine samples
    - Include: candidate RMs, diluted CRM470, representative EQA samples, urine containing modified albumin forms
    - Evaluate adsorption to collection, storage, and sample containers
    - Evaluate effect of centrifugation of sediment



- Develop a reference system for urine albumin
  - Define the measurand
    - What is native in fresh urine
    - Can the N-terminal 24 peptide be used
  - Define a reference material for urine albumin
    - Diluted CRM 470
    - JCCLS candidate material
    - Value assignment
  - Establish a reference measurement procedure
    - Continue development of Mayo LC-IDMS method



- Develop a reference system for urine creatinine
  - Reference measurement procedure exists
    - Only 1 lab listed in JCTLM
  - Reference materials in urine matrix are needed
  - Establish performance requirements for routine method result to be used in a ratio with albumin



- Investigate relationship between albumin/creatinine ratio (ACR) and albumin excretion rate (AER)
  - Develop an algorithm to convert ACR to AER
  - Consider gender, age and race factors
  - Requires standardization to be accomplished
  - Interpretive criteria related to risk assessment