

| Substance | CASRN | Form Tested | Testing Lab | REET I Corneal Thickness Measurement | REET 2 Corneal Thickness Score (1-4) | Overall In Vitro Irritancy Classification |
|-----------|-------|--------------------|-------------|--------------------------------------|--------------------------------------|---|
| HZA-1 | NA | Liquid formulation | Unilever | 32.1 | 2.3 | Severe irritant |
| HZB-1 | NA | Liquid formulation | Unilever | 24.8 | 1.7 | Nonsevere irritant |
| HZC-1 | NA | Liquid formulation | Unilever | 21.2 | 1.0 | Nonsevere irritant |
| HZD-1 | NA | Liquid formulation | Unilever | 20.3 | 1.0 | Nonsevere irritant |
| HZE-1 | NA | Liquid formulation | Unilever | 6.4 | 0.0 | Nonsevere irritant |
| HZF-1 | NA | Liquid formulation | Unilever | 14.3 | 1.0 | Nonsevere irritant |
| HZG-1 | NA | Liquid formulation | Unilever | 29.5 | 2.0 | Severe irritant |
| HZH-1 | NA | Liquid formulation | Unilever | 7.7 | 0.0 | Nonsevere irritant |
| HZI-1 | NA | Liquid formulation | Unilever | 28.6 | 1.7 | Severe irritant |
| HZJ-1 | NA | Liquid formulation | Unilever | 16.3 | 1.0 | Nonsevere irritant |
| HZK-1 | NA | Liquid formulation | Unilever | 36.4 | 2.7 | Severe irritant |
| HZL-1 | NA | Liquid formulation | Unilever | 36.2 | 2.7 | Severe irritant |

Gettings et al. (1996) Data Sorted by Substance

| Substance | CASRN | Form Tested | Testing Lab | REET I Corneal Thickness Measurement | REET 2 Corneal Thickness Score (1-4) | Overall In Vitro Irritancy Classification |
|-----------|-------|--------------------|-------------|--------------------------------------|--------------------------------------|---|
| HZM-1 | NA | Liquid formulation | Unilever | 26.0 | 1.7 | Severe irritant |
| HZN-1 | NA | Liquid formulation | Unilever | 37.7 | 2.7 | Severe irritant |
| HZP-1 | NA | Liquid formulation | Unilever | 25.0 | 1.7 | Severe irritant |
| HZQ-1 | NA | Liquid formulation | Unilever | 7.5 | 0.3 | Nonsevere irritant |
| HZR-1 | NA | Liquid formulation | Unilever | 13.4 | 1.0 | Nonsevere irritant |
| HZS-1 | NA | Liquid formulation | Unilever | 33.3 | 2.3 | Severe irritant |
| HZT-1 | NA | Liquid formulation | Unilever | 2.5 | 0.0 | Nonsevere irritant |
| HZU-1 | NA | Liquid formulation | Unilever | 36.7 | 3.0 | Severe irritant |
| HZV-1 | NA | Liquid formulation | Unilever | 25.6 | 1.3 | Severe irritant |
| HZW-1 | NA | Liquid formulation | Unilever | 23.9 | 1.7 | Nonsevere irritant |
| HZX-1 | NA | Liquid formulation | Unilever | 20.9 | 1.3 | Nonsevere irritant |
| HZY-1 | NA | Liquid formulation | Unilever | 18.6 | 1.0 | Nonsevere irritant |

Gettings et al. (1996) Data Sorted by Substance

| Substance | CASRN | Form Tested | Testing Lab | REET I Corneal Thickness Measurement | REET 2 Corneal Thickness Score (1-4) | Overall In Vitro Irritancy Classification |
|------------------|--------------|--------------------|--------------------|---|---|--|
| HZZ-1 | NA | Liquid formulation | Unilever | 2.8 | 0.0 | Nonsevere irritant |

NA = Not Applicable; REET = Rabbit Enucleated Eye Test