Department of Veterans Affairs	DEPLETED URANIUM (DU) PROGRAM CHECKLIST 24-HOUR URINE URANIUM COLLECTION BALTIMORE VA MARYLAND HEALTH CARE SYSTEM			
	CONSULT URINE INSTR	RUCTIONS (REVI	SED 07/98)	
1. PATIENT NAME 2. SOCIAL				2. SOCIAL SECURITY NUMBER
3. ADDRESS				4. SPECIMEN DATE
5. TELEPHONE NUMBER		6. DATE OF BIF	отн	
3. TELETHONE NOMBER			XIII	
7. REFERRING VA MEDICAL CENTER AND ADDRESS				8. MAIL CODE
9. REFERRING PHYSICIAN				
10. PAGER NUMBER	R 11. TELEPHONE NUMBER			12. FAX NUMBER (To receive report)
TO THE EXTENSE OF	TI TEEL HOWE NO MIDER			12.175 (16 liberty
3. REFERRING FACILTY CONTACT NAME 14. RE		14. REFERRING	NG FACILITY CONTACT PHONE NUMBER	
Call the DU Program at 1-800-815-7533 to obtain the specimen collection kits. Each kit includes shipping materials and eight 500 ml				
containers to be used for the collection of urine over a 24 hour period. The materials in the kits conform to shipping regulations				
#173.199 and FEDEX shipping requirement (IATA Packing Instruction 650).				
The shipper of the diagnostic specimens assumes responsibility for complying with these shipping requirements. Using bottles and				
packing materials supplied by the DU Program, and using them according to directions, assures compliance.				
Schedule patient for 24-hour urine collection.  DATE  A. TIME OF FIRST VOID (D. 1988)	NSCARDED) LIRINE DAY	/ 1	R TIME ∩E EIE	RST VOID URINE DAY 2
A. TIME OF TIME VOID (E	SID (DIOCANDED) GIVINE DAT T		B. TIME OF TRST VOID ORINE DAT 2	
Instruct the patient to urinate directly into the collection container(s), filling each container before proceeding to the next. Uranium sticks to the sides of the container. Therefore, do not transfer urine due to potential loss of analyte. Issue 8 containers to patient to				
ensure full 24-hour collection in approved containers.				
Instruct the patient to collect urine beginning after first morning void of Day 1 and end collection after first morning void on Day 2 (the next day).				
<u>Seal containers as tightly as possible.</u> Bag each urine container with absorbent material. Make sure each plastic bag is sealed tightly.  The sample must be mailed in the package provided.				
TIP: CONTACT YOUR LABORATORY SERVICES SUPERVISOR TO ASSIST IN PACKAGING.				
CENID CDECIMENI VIA EEDEV				
SEND SPECIMEN VIA FEDEX.				
FEDEX TRACKING NUMBER:				
SEND TO:				
VA MARYLAND HEALTH CARE				
PATHOLOGY AND LABORATORY MEDICINE SERVICE (113)				
10 N. GREENE STREET				
BALTIMORE, MARYLAND 2120 ATTN: CHEM. LAB (113) (FOR D				
` ' '	· ·	0000D must be fo	avad to the DII o	effice at 410 605 7042
This completed VA Form 10-9009F and the completed copy of VA Form 10-9009D must be faxed to the DU office at 410-605-7943 before sending the urine sample.				
-	Zinlanka wa si	lasada tabuat		
Another copy of this VA Form 10-9009F, sealed in a Ziploc bag, must be enclosed with the sample for identification purposes.				
Notification of the results can be expected in approximately 90 days.				