

SAFETY FEATURE EVALUATION FORM

SAFETY DENTAL SYRINGES



Date: _____ Department: _____ Occupation: _____

Product: _____ Number of times used: _____

Please **circle** the most appropriate answer for each question. Not applicable (N/A) may be used if the question does not apply to this particular product.

agree.....disagree

1. The safety feature can be activated using a one-handed technique 1 2 3 4 5 N/A
2. The safety feature **does not** obstruct vision of the tip of the sharp and the
intraoral injection site. 1 2 3 4 5 N/A
3. Use of this product requires you to use the safety feature 1 2 3 4 5 N/A
4. This product **does not** require more time to use than a non-safety device 1 2 3 4 5 N/A
5. The safety feature works well with a wide variety of hand sizes 1 2 3 4 5 N/A
6. The device is easy to handle while wearing gloves 1 2 3 4 5 N/A
7. The device is easy to handle when wet 1 2 3 4 5 N/A
8. This device accepts standard anesthetic carpules and does not hinder carpule
changing 1 2 3 4 5 N/A
9. The safety feature **does not** restrict visibility of carpule contents intraorally 1 2 3 4 5 N/A
10. This device accepts standard dental needles of all common lengths and gauges,
and does not interfere with needle changing 1 2 3 4 5 N/A
11. The device provides a better alternative to traditional recapping 1 2 3 4 5 N/A
12. Sterilization of this device is as easy as a standard dental syringe 1 2 3 4 5 N/A
13. For syringes with integral needles only: The needle on this syringe **will not** break
while bending and repositioning in the tissue 1 2 3 4 5 N/A
14. This device is no more difficult to break down after use for sterilization than a
standard dental syringe 1 2 3 4 5 N/A
15. The safety feature operates reliably 1 2 3 4 5 N/A
16. The exposed sharp is permanently blunted or covered after use and prior to
disposal 1 2 3 4 5 N/A
17. There is a clear and unmistakable change (either visible or audible) that occurs
when the safety feature is activated 1 2 3 4 5 N/A
18. The user **does not** need extensive training to operate the product correctly 1 2 3 4 5 N/A
19. The design of the device allows for easy removal of the needle from the syringe 1 2 3 4 5 N/A
20. The design of the device allows for easy removal of the carpule from the syringe 1 2 3 4 5 N/A