Brookhaven National Laboratory – Photon Sciences - National Synchrotron Light Source			
Beamline Hazard Analysis - Beamline U4IR			
LS-BHA-U4IR	Revision: 02	Effective: 9/28/11	Page 1 of 2
The following tasks are included in the beamline safety review. Any operations not listed here, or performed by others than those listed, need additional review by ES&H personnel and/or beamline personnel.			

Prepared By:	G.L. Carr	Approved By:	Lori Stiegler

Tasks / Equipment	Performed by:		Controls	
	Beamline Staff	Users		
Inert compressed gas	Υ	Υ	<ul> <li>Compressed gas training – TQ-COMPGAS1 or NSLS User training</li> <li>Transport cylinders using a cylinder cart</li> <li>Properly secure cylinders</li> <li>Wear safety glasses with side shields while connecting lines</li> <li>Leak check gas lines</li> </ul>	
Solvents	Y	Y	<ul> <li>Limit quantities at beamline to &lt; 250 ml</li> <li>Use PPE - safety glasses, nitrile gloves</li> </ul>	
Pouring cryogens < 5L	Y	Y	<ul> <li>Use goggles</li> <li>Long pants and sleeves - avoid clothing which can trap cryogens (pockets, cuffs)</li> <li>Cryogen or heavy gloves</li> <li>Enclosed shoes</li> <li>Never pour from above chest level</li> </ul>	
Using the LN2 Fill Station or pressurized transfer	Y	Y	<ul> <li>LS-PROC-LN2FILL-MANL training for fill station</li> <li>Use full face shield along with goggles or safety glasses</li> <li>Long pants and sleeves - avoid clothing which can trap cryogens (pockets, cuffs)</li> <li>Cryogen or heavy gloves</li> <li>Enclosed shoes</li> </ul>	
Liquid helium transfer	Y	Y	<ul> <li>Use full face shield along with goggles or safety glasses</li> <li>Long pants and sleeves - avoid clothing which can trap cryogens (pockets, cuffs)</li> <li>Cryogen or heavy gloves</li> <li>Enclosed shoes</li> <li>Use guidance/checklist from PRM 5.1.0</li> </ul>	
Class 2/3a laser	Υ	Υ	Laser registration, warning signs	
Ergonomics concerns –	Υ	Υ	Be aware of body posture  Take force and based as	
awkward postures	.,		Take frequent breaks	
Material Handling/Crane	Υ	N	Training	
Visible light	Y	Y	<ul><li>Covers</li><li>Warning signs</li></ul>	
			Training	

Brookhaven National Laboratory – Photon Sciences - National Synchrotron Light Source			
Beamline Hazard Analysis - Beamline U4IR			
LS-BHA-U4IR	Revision: 02	Effective: 9/28/11	Page 2 of 2
The following tasks are included in the beamline safety review. Any operations not listed here, or performed by others than			
those listed, need additional review by ES&H personnel and/or beamline personnel.			
_			

Magnetic fields > 5G	Y	Y	<ul> <li>Warning signs</li> <li>No workers allowed with ferromagnetic implants</li> </ul>
Cutting/razor blades	Y	Y	<ul> <li>Be aware of hand positioning</li> <li>Use safety knife when possible</li> <li>Consider using cut resistant gloves for frequent cutting tasks</li> </ul>