

GIPSA Performance Verified Mycotoxin Test Kits – Effective 1/31/2013

Company	Test Kit	Part Number	Test Type	Test Kit Range	Test Format	Identified Commodities	Detection Method	Certificate Expiration Date
Beacon Analytical Systems, Inc.	Aflatoxin Plate Kit Quantitative Test	20-0017	Aflatoxin (quantitative)	5 – 100 ppb	Microtiter well plate assay	Corn	Stat Fax Reader Model 303 Plus	1/3/2016
Charm Sciences, Inc.	ROSA® Aflatoxin P/N 10 ppb	LF-AFPN	Aflatoxin (qualitative)	Screens @ 10 ppb	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB)	6/13/2015
Charm Sciences, Inc.	ROSA® Aflatoxin P/N 20 ppb	LF-AFPN	Aflatoxin (qualitative)	Screens @ 20 ppb	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB) and visual	6/13/2015
Charm Sciences, Inc.	ROSA® BEST Aflatoxin P/N	LF-AFPN-BEST	Aflatoxin (qualitative)	Screens @ 10 ppb	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB)	7/28/2014
Charm Sciences, Inc.	ROSA® Aflatoxin BEST P/N	LF-AFPN-BEST	Aflatoxin (qualitative)	Screens @ 20 ppb	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB) and visual	7/28/2014
Charm Sciences, Inc.	ROSA® DON P/N 0.5 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 0.5 ppm	Lateral flow strip	Wheat	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 1.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 1 ppm	Lateral flow strip	Wheat	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 2.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 2 ppm	Lateral flow strip	Wheat	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 5.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 5 ppm	Lateral flow strip	Wheat	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 0.5 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 0.5 ppm	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 1.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 1 ppm	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 2.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 2 ppm	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 5.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 5 ppm	Lateral flow strip	Corn	ROSA®-M Reader (LF-ROSAREADER-M-NB)	5/27/2014
Charm Sciences, Inc.	ROSA® DON P/N 1.0 ppm	LF-DONPN	Deoxynivalenol (qualitative)	Screens @ 1 ppm	Lateral flow strip	Barley	ROSA®-M Reader (LF-ROSAREADER-M-NB)	6/13/2015
Charm Sciences, Inc.	ROSA® FAST Aflatoxin Quantitative Test	LF-AFQ-FAST	Aflatoxin (quantitative)	5 – 100 ppb	Lateral flow strip	Corn, barley, basmati rice, brewers rice, brown rice, corn flour, corn germ meal corn gluten meal, corn grits, corn meal, corn screenings, corn/soy blend, distillers dried grains, distillers dried grains with solubles, hominy, milled rice, millet, oats, popcorn, rough rice, rye, sorghum, soybeans, wheat, and wheat flour	ROSA®-M Reader (LF-ROSAREADER-M-NB)	3/12/2015

GIPSA Performance Verified Mycotoxin Test Kits – Effective 1/31/2013

Company	Test Kit	Part Number	Test Type	Test Kit Range	Test Format	Identified Commodities	Detection Method	Certificate Expiration Date
Charm Sciences, Inc.	ROSA® FAST5 DON Quantitative Test	LF-DONQ-FAST5	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Lateral flow strip	Wheat, corn, barley, distillers dried grains with solubles, malted barley, milled rice, oats, rough rice, sorghum, wheat bran, wheat flour, and wheat middlings	ROSA®-M Reader (LF-ROSAREADER-M-NB)	12/23/2014
Charm Sciences, Inc.	ROSA® FAST5 Fumonisin Quantitative Test	LF-FUMQ-FAST5	Fumonisin (quantitative)	0.5 – 5.0 ppm	Lateral flow strip	Corn, barley, flaking corn grits, millet, oats, rough rice, sorghum, and wheat	ROSA®-M Reader (LF-ROSAREADER-M-NB)	8/2/2015
Charm Sciences, Inc.	ROSA® WET™ Aflatoxin Quantitative Test	LF-AFQ-WET	Aflatoxin (quantitative)	5 – 100 ppb	Lateral flow strip	Corn, barley, corn flour, corn germ meal, corn meal, corn/soy blend, distillers dried grains with solubles, hominy, oats, popcorn, rice bran (defatted), rough rice, sorghum, soybeans, and wheat	ROSA®-M Reader (LF-ROSAREADER-M-NB)	9/12/2015
Charm Sciences, Inc.	ROSA® Ochratoxin Quantitative Test	LF-OCHAQ-G	Ochratoxin A (quantitative)	5 – 100 ppb	Lateral flow strip	Wheat, barley, corn, corn gluten meal, malted barley, oats, rye, sorghum, and soybean meal	ROSA®-M Reader (LF-ROSAREADER-M-NB)	10/1/2015
Charm Sciences, Inc.	ROSA® Zearalenone Quantitative Test	LF-ZEARQ	Zearalenone (quantitative)	100 – 1000 ppb	Lateral flow strip	Corn, barley, distillers dried grains with solubles, flaking corn grits, milled rice, oats, rough rice, sorghum, wheat, and wheat flour	ROSA®-M Reader (LF-ROSAREADER-M-NB)	12/30/2014
EnviroLogix, Inc.	QuickTox™ for QuickScan Aflatoxin	AQ 109 BG	Aflatoxin (quantitative)	5 – 100 ppb	Lateral flow strip	Corn and wheat	EnviroLogix QuickScan™ System	6/29/2013
EnviroLogix, Inc.	QuickTox™ for QuickScan DON	AQ 204 BG AQ 204 BG2 AQ 204 BG3	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Lateral flow strip	Corn, wheat, barley, distillers dried grains with solubles, oats, wheat bran, and wheat flour	EnviroLogix QuickScan™ System	11/8/2013
Neogen Corporation	Agri-Screen® for Aflatoxin	8010	Aflatoxin (qualitative)	Screens @ 20 ppb	Microtiter well plate assay	Corn	Visual	10/19/2014
Neogen Corporation	Agri-Screen® for DON	8310	Deoxynivalenol (qualitative)	Screens @ 1 ppm	Microtiter well plate assay	Wheat	Visual	6/17/2014
Neogen Corporation	Agri-Screen® for DON	8311	Deoxynivalenol (qualitative)	Screens @ 1 ppm	Microtiter well plate assay	Corn	Visual	6/17/2014

GIPSA Performance Verified Mycotoxin Test Kits – Effective 1/31/2013

Company	Test Kit	Part Number	Test Type	Test Kit Range	Test Format	Identified Commodities	Detection Method	Certificate Expiration Date
Neogen Corporation	Reveal® for Aflatoxin (MeOH)	8015	Aflatoxin (qualitative)	Screens @ 20 ppb	Lateral flow strip	Corn	AccuScan III Reader and visual	11/14/2014
Neogen Corporation	Reveal® for Aflatoxin (EtOH)	8015	Aflatoxin (qualitative)	Screens @ 20 ppb	Lateral flow strip	Corn	AccuScan III Reader and visual	11/14/2014
Neogen Corporation	Reveal® for DON	8315	Deoxynivalenol (qualitative)	Screens @ 2 ppm	Lateral flow strip	Wheat	Visual	6/22/2014
Neogen Corporation	Reveal® Q+ for Aflatoxin	8085	Aflatoxin (quantitative)	5 – 100 ppb	Lateral flow strip	Corn, brewers rice, corn flaking grits, corn germ meal, corn gluten meal, corn meal, corn screenings, corn/soy blend, corn starch, cracked corn, distillers dried grains with solubles, popcorn, and sorghum	AccuScan III Reader	3/9/2015
Neogen Corporation	Reveal® Q+ for DON	8385	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Lateral flow strip	Wheat, corn, barley, corn gluten meal, distillers dried grains with solubles, malted barley, oats, and rough rice	AccuScan III Reader	4/13/2015
Neogen Corporation	Reveal® Q+ for Fumonisin	8885	Fumonisin (quantitative)	0.5 – 5.0 ppm	Lateral flow strip	Corn	AccuScan III Reader AccuScan® Pro Reader	9/26/2015
Neogen Corporation	Reveal® Q+ for Zearalenone	8185	Zearalenone (quantitative)	100 – 1000 ppb	Lateral flow strip	Corn, corn gluten meal, distillers dried grains with solubles, rough rice, soy flour, and wheat	AccuScan III Reader	7/2/2015
Neogen Corporation	Veratox® Aflatoxin Quantitative Test	8030	Aflatoxin (quantitative)	5 – 100 ppb	Microtiter well plate assay	Corn, corn bran, corn germ meal, corn gluten meal, corn grits, corn meal, corn/soy blend, cottonseed, distillers dried grains with solubles, lightly pearled barley, malted barley flour, milled rice, millet, popcorn, raw rye, rye flour, sorghum, soybeans, wheat, wheat bran, wheat bran aleurone, wheat flour 2 nd clear, and wheat middlings	Stat Fax Reader Model 321 Plus	3/8/2015

GIPSA Performance Verified Mycotoxin Test Kits – Effective 1/31/2013

Company	Test Kit	Part Number	Test Type	Test Kit Range	Test Format	Identified Commodities	Detection Method	Certificate Expiration Date
Neogen Corporation	Veratox® DON 2/3	8335	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Microtiter well plate assay	Wheat, corn, barley, corn bran, corn germ meal, corn gluten meal, corn grits, lightly pearled barley, malted barley flour, oats, raw rye, rough rice, rye flour, sorghum, wheat bran, wheat bran aleurone, wheat flour, wheat flour 2 nd clear, and wheat middlings	Stat Fax Reader Model 321 Plus	2/10/2015
Neogen Corporation	Veratox® DON 5/5	8331	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Microtiter well plate assay	Wheat, corn, barley, corn bran, corn germ meal, corn gluten meal, corn grits, lightly pearled barley, malted barley flour, oats, rice, raw rye, rye flour, wheat bran, wheat bran aleurone, wheat flour, wheat flour 2 nd clear, and wheat middlings	Stat Fax Reader Model 321 Plus	12/29/2014
R-Biopharm Rhône Ltd.	AflaRhône®	P56 / FL01	Aflatoxin (quantitative)	5 – 100 ppb	Immunoaffinity column	Corn	VICAM Series-4	2/28/2013
R-Biopharm, Inc.	RIDASCREEN® FAST Aflatoxin SC	R9002	Aflatoxin (quantitative)	5 – 100 ppb	Microtiter well plate assay	Corn, brown rice, corn bran, corn flour, corn germ meal, corn gluten meal, corn screenings, corn/soy blend, cracked corn, distillers dried grains, distillers dried grains with solubles, flaking corn grits, milled rice, oats, popcorn, rough rice, rye, sorghum, soybeans, wheat	Stat Fax Reader Model 303 Plus	1/28/2016
R-Biopharm, Inc.	RIDASCREEN® FAST DON	R5901 R5902	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Microtiter well plate assay	Wheat, corn, barley, malted barley, and oats	Stat Fax Reader Model 303 Plus	4/17/2015
R-Biopharm, Inc.	RIDASCREEN® FAST DON SC	R5905	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Microtiter well plate assay	Wheat, corn, barley, malted barley, and oats	Stat Fax Reader Model 303 Plus	7/7/2014
R-Biopharm, Inc.	RIDASCREEN® FAST Fumonisin	R5602	Fumonisin (quantitative)	0.5 – 5.0 ppm	Microtiter well plate assay	Corn	Stat Fax Reader Model 303 Plus	9/21/2015

GIPSA Performance Verified Mycotoxin Test Kits – Effective 1/31/2013

Company	Test Kit	Part Number	Test Type	Test Kit Range	Test Format	Identified Commodities	Detection Method	Certificate Expiration Date
Romer Labs, Inc.	AgraStrip® (EtOH)	COKAS1200	Aflatoxin (qualitative)	Screens @ 10 ppb	Lateral flow strip	Corn	Visual	8/11/2014
Romer Labs, Inc.	AgraStrip® (MeOH)	COKAS1200	Aflatoxin (qualitative)	Screens @ 10 ppb	Lateral flow strip	Corn	Visual	8/11/2014
Romer Labs, Inc.	AgraStrip® (EtOH)	COKAS1000 COKAS1000U	Aflatoxin (qualitative)	Screens @20 ppb	Lateral flow strip	Corn	Visual	8/2/2015
Romer Labs, Inc.	AgraStrip® (MeOH)	COKAS1000 COKAS1000U	Aflatoxin (qualitative)	Screens @20 ppb	Lateral flow strip	Corn	Visual	8/2/2015
Romer Labs, Inc.	FluoroQuant® Afla	COKFA1010	Aflatoxin (quantitative)	5 – 100 ppb	Rapid, fluorometric, solid phase cleanup	Corn, brown rice, corn meal, corn/soy blend, milled rice, popcorn, soybeans, and wheat	Romer FQ-Reader™ (EQFFM3000)	4/27/2015
Romer Labs, Inc.	FluoroQuant® Afla IAC	COKFA4010	Aflatoxin (quantitative)	5 – 100 ppb	Immunoaffinity column	Corn, brown rice, corn flour, corn germ meal, corn gluten meal, corn meal, corn screening, corn/soy blend, distillers dried grains with solubles, flaking corn grits, milled rice, oats, popcorn, rice, rye, sorghum, soybeans, and wheat	Romer FQ-Reader™ (EQFFM3000)	11/22/2014
Thermo Fisher Scientific Diagnostic Division	EZ-Tox DON	600120-DON	Deoxynivalenol (quantitative)	0.5 – 5.0 ppm	Rapid enzyme immunoassay	Wheat, corn, barley, and malted barley	BioTek ELx8081U	4/13/2015
VICAM	AflaTest™	12022	Aflatoxin (quantitative)	5 – 100 ppb	Immunoaffinity column	corn, barley, corn bran, condensed distillers solubles, corn flour, corn gluten feed, corn gluten meal, corn meal, corn screenings, corn soy blend, distillers dried grain, distillers dried grains with solubles, flaking corn grits, popcorn, milled rice, rice bran, rough rice, sorghum, soybeans, and wheat.	VICAM Series-4EX	1/3/2016
VICAM	OchraTest™	13012	Ochratoxin A (quantitative)	5 – 100 ppb	Immunoaffinity column	Wheat	VICAM Series-4EX	9/12/2015