

# Rapid DNA Testing Approaches for Reference Samples

Erica L.R. Butts and Peter M. Vallone

National Institute of Standards and Technology (NIST)  
100 Bureau Drive MS 8314, Gaithersburg, MD 20899-8314

Poster # 30 at:  
23<sup>th</sup> ISHI, Nashville, TN  
October 15-18, 2012

Email: [erica.butts@nist.gov](mailto:erica.butts@nist.gov)  
Phone: 301-975-5107

## Rapid DNA (Services)

Such as Bode Rapid DNA Service™. This is a **lab-based** quick turnaround (< 2 h) approach using techniques/methods familiar to a forensic laboratory (PCR, CE separation, etc). **This is not a field portable instrument**. The profiles are for **investigative leads** (not for uploading into CODIS).

## Rapid DNA (Instruments)

New fully integrated (swab in → profile out) platforms developed by IntegenX, NetBio, ZyGem/Lockheed, and Univ of AZ intended for the STR typing of reference samples.

**Throughput:** 5-8 samples  
**Time:** 90 minutes per 5-8 profiles  
**Cost:** >\$200 per sample

The eventual goal for these platforms is to develop STR profiles in a **non-laboratory setting**. This will rely on proper validation and NDIS approval paths if information will be intended for CODIS comparison.

## Current ways your laboratory can perform Rapid DNA Typing:

**What is Rapid DNA Typing?** Recently, various claims have been made to the term 'Rapid DNA'. It is important to understand the distinctions.

Three main areas exist within Rapid DNA Typing: **Services**, **Techniques**, and **Instruments**

**Rapid DNA Techniques:** genotype results generated in less than 2 hours with standard laboratory equipment and protocols are shown below

**Acknowledgements:**  
This work was funded in part through an interagency agreement between NIST and the FBI Biometrics Center of Excellence. We would also like to thank Life Technologies and Promega for providing the STR kits used for testing and analysis.

**NIST Disclaimer:** Certain commercial equipment, instruments and materials are identified in order to specify experimental procedures as completely as possible. In no case does such identification imply a recommendation or it imply that any of the materials, instruments or equipment identified are necessarily the best available for the purpose. **Points of view are mine** and do not necessarily represent the official position of the National Institute of Standards and Technology.

**References:**  
1. Promega Corporation. (2012). PowerPlex Fusion System Technical Manual. Madison, WI  
2. Promega Corporation. (2012). PunchSolution Kit Technical Manual. Madison, WI  
3. Promega Corporation. (2012). SwabSolution Kit Technical Manual. Madison, WI  
4. Applied Biosystems. (2012). GlobalFiler™ Express Kit Early Evaluation Study Guide. Foster City, CA  
5. ZyGEM. Quick-Start Guide  
6. P.M. Vallone, C.R. Hill, D. Podini, J.M. Butler. Rapid amplification of commercial STR typing kits. Forensic Sci. Int. Genet. Sup. Series 2 (2009) 111-112.  
7. A. Foster, N. Laurin. Development of a fast PCR protocol enabling rapid generation of AmpFISTR Identifier profiles for genotyping of human DNA. Investigative Genetics 2012 3:6.

## Fast PCR

Fast PCR protocols can be adopted into a STR typing workflow

### Buccal Swabs (Cotton)

**Setup: 3 min 35 sec**

**Prep-N-Go Buffer**

Incubate swab at room temperature in 400 µL Prep-n-Go Buffer

**forensicGEM Saliva**

One buccal swab washed with 500 µL DNA-free water and incubated at room temperature for 15 minutes. 20 µL elute added to reaction mix.

75 °C for 15 minutes followed by 95 °C for 5 minutes

2 µL solution added to **Rapid Identifier PCR Setup**

**Setup: 6 min 15 sec**  
**Incubation: 20 min**

**Fast cycling:** The use of robust enzymes and faster thermal cyclers to reduce the amount of time required for PCR amplification.

**Current Rapid Identifier PCR Setup [6]:**  
2.0 µL Identifier Primers  
5.0 µL Takara Perfect Real Time Mix  
0.25 µL Takara SpeedSTAR Polymerase  
0.75 µL Water  
2.0 µL DNA

**Cycling Parameters (28 cycles) [7]:**

95°C	95°C	61°C	72°C
1 min	5 s	15 s	1 min

**Thermal Cyclers: ABI 9700 and Streck Philisa**

**ABI 9700**  
Heating mechanism: Peltier block (AI)  
Tube format: 0.2 mL tubes or plates  
96 reactions per instrument

**Streck Philisa**  
Heating mechanism: Peltier block (AI)  
Tube format: proprietary 50 µL tubes  
8 reactions per instrument

Require the use of gel loading tips to load PCR product into CE setup plate due to tube design

**Prep-N-Go**  
PCR Setup: 5 min 42 sec  
**PCR: 36 Minutes**

**forensicGEM Saliva**  
PCR Setup: 7 min 18 sec  
**PCR: 36 Minutes**

**Prep-N-Go**  
PCR Setup: 2 min 7 sec  
**PCR: 14 Minutes**

**forensicGEM Saliva**  
PCR Setup: 6 min 48 sec  
**PCR: 14 Minutes**


### CE Conditions:

1.2 kV for 7 sec

**CE Run Time:**  
**38 minutes**

**Analysis:**  
GeneMapper ID-X v1.2 with an analytical threshold of 50 RFU

STR genotype results were generated in less than 2 hours with standard laboratory equipment and protocols. The validation of Direct PCR STR kits or Fast PCR protocols will allow a forensic lab to have 'Rapid DNA' capabilities.



## Technique

Rapid DNA Typing

**Sample Type**  
8 samples for each set of parameters

**Pre-PCR Techniques**  
Various techniques suggested to prepare multiple substrates for PCR amplification

**PCR Amplification**  
Fast PCR: Identifier  
Direct PCR: GlobalFiler Express and PowerPlex Fusion

**CE & Analysis**  
Applied Biosystems 3500 Genetic Analyzer

**Results & Final STR Typing Times**  
The total amount of time from sampling to complete profile for 8 samples, to include all liquid handling steps, PCR, CE, sample transfer steps, and data analysis.

## Direct PCR

Direct PCR kits speed up processing by bypassing extraction and can perform a direct amplification in less than 90 minutes. PCR products are separated and detected on a traditional CE system.

### Blood FTA

**Punching:**  
**2 min 21 sec**

### Blood S&S 903



**Punching:**  
**2 min 26 sec**

### Buccal Swabs (Cotton)



**None**

**Prep-N-Go Buffer**  
(for use with GlobalFiler)

3 µL added with PCR setup and one 1.2 mm punch

**Setup: 58 sec**  
**Incubation: 30 min**

**Promega PunchSolution**  
(for use with Fusion)

10 µL PunchSolution Reagent incubated at 70 °C for 30 minutes

**forensicGEM Storage Card (Blood)**  
(for use with GlobalFiler or Fusion)

One 1.2 mm punch washed with 100 µL DNA-free water and incubated at room temperature for 15 minutes. Remove water and add reaction mix to the 1.2 mm punch.

75 °C for 15 minutes followed by 95 °C for 5 minutes

2 µL solution added to **GlobalFiler**  
3 µL solution added to **Fusion**

**Setup: 5 min 28 sec**  
**Incubation: 35 min**

**Prep-N-Go Buffer**  
(for use with GlobalFiler)

Incubate swab at room temperature in 400 µL Prep-n-Go Buffer

**Setup: 1 min 40 sec**  
**Incubation: 30 min**

**Promega SwabSolution**  
(for use with Fusion)

One buccal swab incubated at 70 °C for 30 minutes in 1 mL SwabSolution Reagent

**forensicGEM Saliva**  
(for use with GlobalFiler or Fusion)

One buccal swab washed with 500 µL DNA-free water and incubated at room temperature for 15 minutes. 20 µL elute added to the reaction mix.

75 °C for 15 minutes followed by 95 °C for 5 minutes

2 µL solution added to **GlobalFiler**  
3 µL solution added to **Fusion**

**Setup: 6 min 15 sec**  
**Incubation: 20 min**

**GlobalFiler Express**  
6-dye Technology  
24plex  
3500 Instruments only

**PCR Setup (Punches):**  
6 µL Master Mix  
6 µL Primer Set  
3 µL Prep-n-Go

**Setup: 1 min 34 sec**

**Cycling Conditions:**  
**ABI 9700**  
95 °C for 1 minute

94 °C for 3 seconds  
60 °C for 30 seconds  
**Cycle for 26 cycles**

60 °C for 8 minutes

**GlobalFiler**  
**Cycling Time:**  
**42 minutes**  
**15 seconds**

**PCR Setup (Swabs):**  
6 µL Master Mix  
6 µL Primer Set  
3 µL Prep-n-Go or DNA



**Setup: 2 min 7 sec**

**PowerPlex Fusion**  
5-dye Technology  
24plex  
3100 or 3500 Instruments

**PCR Setup (Punches):**  
5 µL Master Mix  
5 µL Primer Set  
15 µL Water

**Setup: 1 min 33 sec**

**Cycling Conditions:**  
**ABI 9700**  
96 °C for 1 minute

94 °C for 10 seconds  
59 °C for 60 seconds  
72 °C for 30 seconds  
**Cycle for 26 cycles**

60 °C for 20 minutes

**Fusion**  
**Cycling Time:**  
**1 hour**  
**35 minutes**  
**45 seconds**

**PCR Setup (Swabs):**  
5 µL Master Mix  
5 µL Primer Set  
13 µL Water  
2 µL Swab Extract/DNA



**Setup: 2 min 19 sec**

### CE Conditions:

**GlobalFiler:**  
1.2 kV for 16 sec  
1550 sec Run Time  
13.0 kV Run Voltage  
15.0 kV PreRun Voltage

**Fusion:**  
1.2 kV for 24 sec  
1300 sec Run Time  
15.0 kV Run Voltage

**GlobalFiler**  
**CE Run Time:**  
**39 minutes**  
**5 seconds**

**Fusion**  
**CE Run Time:**  
**37 minutes**  
**30 seconds**


**Analysis:**  
GeneMapper ID-X v1.2 with an analytical threshold of 200 RFU for both **GlobalFiler** and **Fusion**

**Analysis of quality single source reference samples:**  
**45 sec**



## GlobalFiler Express

## PowerPlex Fusion



Final processing times reflect one run of 8 samples

Final processing times reflect one run of 8 samples