

### Outline of Topics to Discuss

- Background information on Identifiler Plus and the major differences with Identifiler
- Low Template (LT) DNA samples
  - Challenges and limitations with LT-DNA testing
  - Approaches to genotyping low template DNA
  - NIST LT-DNA data and Peak Height Ratios (PHR) with Identifiler Plus
  - LT-DNA mixture samples using Identifiler Plus
- Direct PCR results with Identifiler Plus
- Identifiler Plus data from the 3500x Genetic Analyzer
- Summary and Conclusions



### Similarities with Identifiler

- Primer sequences and concentrations are the same
- Amplicon sizes are the same (<360 bp)
- Same amount of alleles in allelic ladder
- Same dye set (G5)
- 25 µL reaction volume
- Same species specificity and precision

### Differences with Identifiler

- Master Mix and Primer Mix only – no separate Taq/enzyme to add
- 2 recommended cycling protocols
  - 28 cycles: 1 ng optimum input DNA, full profiles with 125 pg
  - 29 cycles: <500 pg input DNA, extra sensitivity for <125 pg
- Shorter thermal cycling conditions (decreased by ~ 1 hour)
  - One denaturation step (20 sec @ 94°C)
  - One annealing/extending step (3 min @ 59°C)
  - Much shorter final extension (10 min @ 60°C)
- Cleaner baseline and improved heterozygote peak balance
- Optimized to overcome inhibition

### Thermal Cycling Conditions

**Identifiler Plus**







| Initial incubation step | Cycle (28 or 29 cycles <sup>1</sup> ) |                |                 | Final extension | Final hold |
|-------------------------|---------------------------------------|----------------|-----------------|-----------------|------------|
|                         | Denature                              | Anneal/Extend  |                 |                 |            |
| HOLD                    | CYCLE                                 |                |                 | HOLD            | HOLD       |
| 95 °C<br>11 min         | 94 °C<br>20 sec                       | 59 °C<br>3 min | 60 °C<br>10 min | 4 °C<br>∞       |            |

**Identifiler** ↗ This step is eliminated

| Initial incubation step | Denature          | Anneal         | Extend         | Final Extension | Final Step           |
|-------------------------|-------------------|----------------|----------------|-----------------|----------------------|
|                         | CYCLE (28 cycles) |                |                |                 |                      |
| HOLD                    | CYCLE             |                |                | HOLD            | HOLD                 |
| 95 °C<br>11 min         | 94 °C<br>1 min    | 59 °C<br>1 min | 72 °C<br>1 min | 60 °C<br>60 min | 4–25 °C<br>(forever) |


**Identifiler Plus is about 1 hour shorter than Identifiler**




## Low Template (LT) DNA Samples

- Some Definitions of Low Template (LT) DNA**
- Working with **<100-200 pg genomic DNA**
  - Considered to be data below stochastic threshold level where PCR amplification is not as reliable (determined by each laboratory; typically 150-250 RFUs)
  - Enhancing the sensitivity of detection (increasing PCR cycles, PCR product clean-up, increasing CE injection/voltage)
  - Having too few copies of DNA template to ensure reliable PCR amplification (allelic or full locus drop-out)
  - Can often be the minor component of mixture samples consisting of low level DNA template amounts

- Challenges of LT-DNA Testing**  
Gill, P. (2001) *Croatian Med. J.* 42(3): 229-232
- Increased chance for contamination (want a sterile lab environment to reduce staff contamination)
  - Data interpretation is more complicated (due to stochastic variation during PCR amplification):
    - Heterozygote peak imbalance
    - Allele drop-out
    - Allele drop-in
    - Increased stutter products
- LT-DNA profiles should be interpreted with careful guidelines



- Suggestions for Optimal Results with LT-DNA**
- Typically at least 2 – 3 PCR amplifications from the same DNA extract are performed to obtain **consensus profiles**
  - An allele cannot be scored (considered real) unless it is present at least twice in replicate samples
  - Extremely sterile environment is required for PCR setup to avoid contamination from laboratory personnel or other sources






### Signal Enhancement Techniques







- **Additional PCR cycles**
- **More sensitive kits** (Identifiler Plus)
- Microcon cleanup to remove salts that interfere with electrokinetic injection
- Lower PCR volume (concentrates amplicon)
- Increase TaqGold/enzyme concentration
- Longer CE injection times and voltage
  - 10 s @ 3 kV = 30
  - 5 s @ 2 kV = 10


### NIST Example LT-DNA Data with Identifiler Plus

### Experimental Design to Study LT-DNA Issues

- Pristine DNA Samples
  - 2 single-source samples
  - **heterozygous for all loci tested** (permits peak height ratio studies)
- **Low DNA Template Amounts**
  - Dilutions made after DNA quantitation against NIST SRM 2372
  - **100 pg, 30 pg, and 10 pg** (1 ng tested for comparison purposes)
- Replicates
  - **5 separate PCR reactions** for each sample
- STR Multiplex Kits
  - **Identifiler Plus** (half-reactions)
- **Increased Cycle Number**
  - Identifiler Plus (**29 cycles and 32 cycles**; 28 for 1 ng)













## Examination of LT-DNA Mixtures with Identifiler Plus

- ### LT-DNA Mixture Samples
- 2 samples (male and female) were mixed together at 1:3 and 1:5 – 1 ng (1:3 and 1:5) or 100 pg (1:5) or 50 pg (1:3) total DNA
  - 3 person mixture (2 males and female) were mixed together at 1:2:3 – 1 ng or 100 pg total DNA
  - Identifiler Plus (28 and 31 cycles) was tested (half reactions)
  - 5 replicates with 3 extra cycles











### Summary of Data Observed

- The results with pristine full heterozygous samples demonstrate that replicate testing can produce reliable information with single source samples at low levels of DNA when consensus profiles are created.
- With 3 extra cycles, there was better recovery at 10 pg of DNA using Identifiler Plus including less allelic and full locus drop-out. However, there is a greater potential for allele drop-in or high stutter.
- Variability of peak heights in replicates was observed with LT-DNA mixtures using Identifiler Plus.
- More minor contributor peaks were called with 3 extra cycles using Identifiler Plus.



### Examination of direct PCR with Identifiler Plus

### Direct PCR with Identifiler Plus Experimental Design

- Full reactions (25 µL total volume)
- 1.2mm punches
  - Buccal and blood samples
    - 903 and FTA paper
    - Two samples for each condition
- Punch added directly to Master Mix
- Manufacturer thermal cycling protocol
  - 28 cycles
















- ### Summary of Data Observed
- Full profiles were obtained from FTA blood, FTA saliva, and 903 blood punches for both samples using Identifiler Plus with direct PCR.
  - Partial profiles were obtained from 903 saliva samples - this could be due to sampling issues.
  - There were some adenylation issues with FTA blood and saliva samples – this could be remedied with lower injection or less PCR cycles.
  - FTA blood, FTA saliva, and 903 blood profiles had excellent heterozygote and locus-to-locus peak height balance.

### Identifiler Plus data on the 3500x/ Genetic Analyzer (@ AFDIL)



- ### Experimental Design
- Identifiler Plus kit was used with full reactions (25  $\mu$ L total volume)
    - 0.5 ng DNA, 28 cycles
      - GS 500 LIZ size standard
      - GS 600 LIZ v2 size standard (normalization)
  - 23 population samples + 1 Allelic Ladder
  - 3500x/ data compared to 3130x/ data
    - Default injection of 1.2 kV, 24 sec (3500x/)
    - Low injection of 1.2 kV, 10 sec (3500x/)
    - Default injection of 3 kV, 10 sec (3130x/)





### Summary of Data Observed

- The RFU scale for the 3500x/ is different than the 3130x/ (30000 RFU vs 8000 RFU).
- The 3500x/ instrument is more sensitive than the 3130x/ – can adjust the injection time and voltage.
- Identifiler Plus profiles on the 3500x/ are well balanced (inter- and intra-locus and between dye channels).
- The GS 600 v2 size standard is for the normalization of data between different instruments in the lab; the data is comparable to data using the GS 500 size standard.

### Conclusions

- Identifiler Plus is a highly sensitive kit that can result in full profiles down to 30 pg of DNA.
- With 3 extra cycles in LT-DNA mixtures, more minor contributor peaks were called using Identifiler Plus.
- Full profiles were obtained from FTA blood, FTA saliva, and 903 blood punches using Identifiler Plus with direct PCR.
- Identifiler Plus profiles on the 3500x/ are well balanced including good inter- and intra-locus balance as well as between dye channels.


### Acknowledgments

**NIST Funding:** Interagency Agreement 2008-DN-R-121 between the [National Institute of Justice](#) and NIST Office of Law Enforcement Standards

**NIST Disclaimer:** Certain commercial equipment, instruments and materials are identified in order to specify experimental procedures as completely as possible. In no case does such identification imply a recommendation or endorsement by the National Institute of Standards and Technology nor does it imply that any of the materials, instruments or equipment identified are necessarily the best available for the purpose.

**Points of view are mine** and do not necessarily represent the official position or policies of the US Department of Justice or the National Institute of Standards and Technology.

**NIST Team for This Work**



John Butler



Dave Duerwer



Erica Butts

**A special thanks to Applied Biosystems for providing the kits used in this study**

**Contact Info:** [becky.hill@nist.gov](mailto:becky.hill@nist.gov), 301-975-4275