
INSTRUMENT

EAM LINE

SPALLATION NEUTRON SOURCE

EQ-SANS — EXTENDED Q-RANGE SMALL-ANGLE NEUTRON SCATTERING DIFFRACTOMETER

The EQ-SANS Diffractometer is designed for the study of materials across length scales ranging from 0.1 to 100 nm. The high intensity provided by the EQ-SANS enables both



high-throughput experiments and time-resolved experiments facilitated by the pulsed source of SNS. EQ-SANS enables measurements over a wide Q-range at a single instrument configuration, providing improved throughput. The high maximum Q of the instrument allows both large-scale and local structure to be studied by the instrument. The versatility of SANS makes EQ-SANS broadly applicable to a wide range of materials from science and industry.

APPLICATIONS

Life science

- Solution structures of proteins, DNA, and other biological molecules and molecular complexes
- Protein-protein and protein-ligand interactions, kinase regulation
- Protein-membrane interaction
- Materials for drug delivery

Polymer and colloidal systems

- Block copolymers and dendrimers
- Micelles and emulsions
- Polyelectrolytes and ion distribution at solid-liquid interfaces

Materials science

- Simultaneous study of domain and crystalline structures
- Crystallization and precipitation
- Nanoparticles

Earth and environmental sciences

- Pore structure in soils
- Structure of geologic materials

SPECIFICATIONS

14 m
3–4.3 Å
Coupled supercritical hydrogen
~10 ⁷ –10 ⁹ n/ cm ² /s
0.004 Å ⁻¹ < Q < 1.5 Å ⁻¹

LOW-ANGLE DETECTOR	
Sample- to-detector distance	1.3–10 m
Detector size	1 x 1 m
Detector resolution	5.5 x 4.3 mm

Status: Available to users

FOR MORE INFORMATION, CONTACT

Instrument Scientist: William Heller, hellerwt@ornl.gov, 865.241.0093 Instrument Scientist: Chris Stanley, stanleycb@ornl.gov, 865.574.6669 Instrument Scientist: Changwoo Do, doc1@ornl.gov, 865.576.6981

neutrons.ornl.gov/eqsans

