

The following spares and accessories are available from your laboratory dealer.

Description	Catalogue Number
Roller	DSRT1014
Spanner	USRT1C2

For a comprehensive list of parts required by service engineers conducting internal repairs and a service manual, please contact the Technical Service Department of Bibby Sterilin Ltd. quoting both the model and serial number.

Only spare parts supplied or specified by Bibby Sterilin Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

If in any doubt, please contact the Technical Service Department of Bibby Sterilin Ltd. or the point of sale.

Bibby Sterilin Ltd
Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0) 1785 812121
Fax: +44 (0) 1785 813748
E-mail equipmnet@bibby-sterilin.com
www.bibby-sterilin.com

Warranty

Bibby Sterilin Ltd warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of two (2) years. In the event of a justified claim Bibby Sterilin will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by incorrect installation, adaptation, modification, fitting of non approved parts or repair by unauthorised personnel.

Technical Specification

	SRT1	SRT2
No. of rollers	5	10
Speed	33rpm	33rpm
Amplitude	16mm	16mm
Roller size (length x diameter)	340 x 30mm	340 x 30mm
Overall dimensions (WxDxH)	500 x 200 x 90mm	500 x 380 x 90mm
Net weight	5kg	7kg