

Neutron Phasing Methods Applicable to D/H SIR Derivatives

David A. Langs, Hongliang Xu and Herbert A. Hauptman

Hauptman-Woodward Medical Research Institute and Department of Structural Biology,
SUNY at Buffalo, 700 Ellicott Street, Buffalo, NY 14203

Abstract: Twenty years ago Kossiakoff reviewed the then recent advances in neutron diffraction technology as applied to macromolecules and predicted that “Clearly the time is close at hand when neutron diffraction will take its place alongside X-ray diffraction as an equal partner in the structure determination of macromolecules.” The first, and perhaps only successful application to determine a large crystal structure by SIR phasing methods occurred when Roger Koeppe enlisted the help of Beno Schoenborn to phase an orthorhombic form of gramicidin A via H/D exchange. A 5Å map based on 8 determined H/D replacement sites revealed the cylindrical shape of the dimer shown in figure 1, but the helical trace of the peptide chain was not discernable. Koeppe also synthesized a deuterated Val¹ derivative to improve the phases by MIR methods, but crystals of sufficient size could not be obtained. In retrospect, it may have been impossible to determine the Val¹ replacement sites since H/D positions on the same methyl group would be only 1.7Å apart, and not resolved by the experimental resolution. In this regard we have developed an SIR structure solution method that does not require us to determine the replacement sites. Instead, low resolution macromolecular phases can be determined directly from a *a priori* phase invariant triples estimates based on the SIR amplitudes.




Figure 1

Negative Density & Direct Methods: For many years most crystallographers believed that map “positivity” was a necessary condition for direct methods validity even though Hauptman (1976) showed that direct methods could be rigorously extended to the case of unequal atoms, even to the extreme that some could have negative scattering density, such as H in neutron diffraction. With this in mind our direct methods program **SnB** was modified to include the most negative peak positions as hydrogen atoms in structure factor calculations to recompute phases. Tests performed on experimental neutron data for cyclosporin A [C₆₂H₁₁₁N₁₁O₁₂·H₂O, 86% complete to 0.915Å; Knott, Schefer & Schoenborn, 1990] gave a 2% NSnB success rate for this 199 atom structure in spite of the fact that only two of the 9000 triples used had A-values greater than 1.5.

Is Negative Density Harmful? : Error-free neutron data were used to compare the fully deuterated structure with the native hydrogenated one. Similarly, error-free X-ray data were tested for the same purpose. *These results were contrary to our expectations.*

| Cyclosporin | success rate | < δφ > | #atoms |
|---------------|--------------|--------|--------|
| NSnB | | | |
| (H-structure) | 6 % | 30° | 199 |
| (D-structure) | 0.3 % | 50° | 199 |
| X-RAY | | | |
| (SnB) | 1.1 % | 13° | 86 |
| (RANTAN) | 0.1 % | 25° | 86 |

In view of this evidence, one is forced to conclude that *negative density is an asset rather than a hindrance* when it comes to direct methods. Tests indicate that structures as large as lysozyme can be solved by NSnB if data < 0.90Å are available. Larger structures, however, will require derivative SIR data to succeed at much lower resolution.

SIR Direct Methods: One tremendous advantage afforded by neutron diffraction is that an H/D substitution strategy can provide an enormous number of *perfect* isomorphous replacement derivatives. The same can not be said of X-ray data. Since estimates for the three-phase triples invariants for SIR data were formulated 20 years ago, it is a fairly straightforward proposition to incorporate these results into NSnB. Error-free neutron data for a pair of cyclosporin structures, S and T, was tested to this end.



SIR triples A-values were computed from the six E-magnitudes S_h, S_k, S_l, T_h, T_k, T_l according to Hauptman’s 1982 formula. The A_{SIR} values ranged from +10 to -10 with 25,000 triples exceeding the 1.5 value indicated for the native neutron data. These large A_{SIR} values indicate an enormous phasing power, such that phasing may be achieved at much lower resolution. In this regard we began with 0.915Å data, and if successful, proceeded to lower the resolution to determine at which point difficulties ensued. In this way we were able to solve the structure at as low as 2.27Å resolution, *which bodes well for potential macromolecular applications using neutron data from spallation sources.*

The 2.27Å Results: 1000 random NSnB trials were each iterated for 150 cycles and sorted in ascending order on R_{min} as is the usual practice with X-ray data. Density modification [ρ=0 if ρ<0.25σ(ρ); ρ=2σ(ρ) if ρ>2σ(ρ)] replaced “peak picking” for the lower resolution ranges. 13 Solutions are noted among the 75 sets with the lowest R_{min} values (figure 2), but are not clearly separated from the non-solution on this basis.




Figure 2

To better identify the 13 solutions a **sine invariant correlation coefficient** was proposed

$$CC(i,j) = [\langle \sin\Phi_{h,k} \sin\Psi_{h,k} \rangle - \langle \sin\Phi_{h,k} \rangle \langle \sin\Psi_{h,k} \rangle] / \sigma(\sin\Phi_{h,k}) \sigma(\sin\Psi_{h,k})$$

where Φ_{h,k} and Ψ_{h,k} are the triples sine invariants of the ith and jth sets, respectively. Large ICC(i,j) values are expected for equivalent phase sets. The top 75 phase sets in Table 1 will generate 75*74/2 = 2775 CC(i,j) values. When these values are sorted in descending order on ICC the topmost entries (Table 1) will form the basis for identifying the groups of phase sets that are equivalent except for the choice of origin and enantiomorph.

Table 1. Top 75 Pairwise comparisons ranked on CC

| rank | amb1 | amb2 | CC | rank | amb1 | amb2 | CC | rank | amb1 | amb2 | CC |
|------|-------|-------|-------|------|-------|-------|-------|------|-------|-------|-------|
| 1 | 677* | 1169* | 0.224 | 26 | 1281* | 1611* | 0.177 | 51 | 553 | 567 | 0.160 |
| 2 | 33* | 281* | 0.221 | 27 | 1045* | 1281* | 0.177 | 52 | 117* | 1281* | 0.160 |
| 3 | 677* | 1045* | 0.218 | 28 | 553 | 1397 | 0.176 | 53 | 117* | 1169* | 0.160 |
| 4 | 655* | 1169* | 0.217 | 29 | 551 | 715 | 0.175 | 54 | 281* | 1169* | 0.158 |
| 5 | 31 | 961 | 0.212 | 30 | 553 | 1957 | 0.173 | 55 | 33* | 1281* | 0.158 |
| 6 | 117* | 1015* | 0.210 | 31 | 1059 | 1921 | 0.171 | 56 | 1711 | 1807 | 0.157 |
| 7 | 827 | 1305 | 0.209 | 32 | 65 | 595 | 0.170 | 57 | 1305 | 1921 | 0.157 |
| 8 | 281* | 677* | 0.209 | 33 | 281* | 1281* | 0.169 | 58 | 589 | 1059 | 0.156 |
| 9 | 371 | 1921 | 0.207 | 34 | 1397 | 1957 | 0.168 | 59 | 371 | 1059 | 0.156 |
| 10 | 677* | 1281* | 0.205 | 35 | 1365 | 1431 | 0.168 | 60 | 33* | 677* | 0.156 |
| 11 | 1169* | 1611* | 0.203 | 36 | 827 | 1921 | 0.168 | 61 | 677* | 1611* | 0.155 |
| 12 | 171 | 1717 | 0.203 | 37 | 33* | 1875* | 0.168 | 62 | 589 | 1345 | 0.155 |
| 13 | 117* | 655* | 0.203 | 38 | 117* | 1611* | 0.167 | 63 | 31 | 803 | 0.155 |
| 14 | 1397 | 1935 | 0.202 | 39 | 1677 | 1711 | 0.166 | 64 | 1305 | 1711 | 0.154 |
| 15 | 281* | 1611* | 0.192 | 40 | 589 | 1921 | 0.165 | 65 | 715 | 1365 | 0.152 |
| 16 | 1717 | 1935 | 0.191 | 41 | 1045* | 1169* | 0.163 | 66 | 595 | 827 | 0.152 |
| 17 | 715 | 1431 | 0.189 | 42 | 1015* | 1281* | 0.162 | 67 | 33* | 1045* | 0.152 |
| 18 | 371 | 1691 | 0.189 | 43 | 677* | 1015* | 0.162 | 68 | 1203 | 1983 | 0.151 |
| 19 | 1169* | 1281* | 0.188 | 44 | 371 | 1305 | 0.162 | 69 | 1169* | 1737* | 0.151 |
| 20 | 371 | 759 | 0.187 | 45 | 281* | 655* | 0.161 | 70 | 1397 | 1717 | 0.150 |
| 21 | 563 | 1711 | 0.186 | 46 | 209 | 1677 | 0.161 | 71 | 589 | 827 | 0.150 |
| 22 | 407 | 1711 | 0.182 | 47 | 181 | 1203 | 0.161 | 72 | 553 | 583 | 0.150 |
| 23 | 33* | 1611* | 0.182 | 48 | 1513 | 1853 | 0.160 | 73 | 551 | 827 | 0.150 |
| 24 | 371 | 827 | 0.179 | 49 | 1345 | 1711 | 0.160 | 74 | 65 | 1921 | 0.150 |
| 25 | 655* | 1045* | 0.178 | 50 | 583 | 1957 | 0.160 | 75 | 1939 | 1999 | 0.149 |

CLUSTER ANALYSIS: Only four equivalent groups of phases can be distilled from the ICCI data indicated in Table 1. The omit-map correlation coefficient and <ρ²> of the electron density are slightly better than R_{min} for identifying the correct solution.

| Group | # of sets | < δφ > | R _{min} | CC(X _h ,F _h) | <ρ ² > |
|-------|-----------|--------|------------------|-------------------------------------|-------------------|
| I | 13 | 35° | 0.423 | 0.581 | 322 |
| II | 5 | 71° | 0.432 | 0.397 | 159 |
| III | 10 | 76° | 0.437 | 0.510 | 167 |
| IV | 41 | 84° | 0.407 | 0.271 | 72 |

When the 13 sets are phase shifted to the same origin and enantiomorph, and then phase averaged, the resulting set has a mean phase error of 38° and this 2.27Å map reveals the peptide backbone and side-chains of the structure shown in figure below.




Figure 3

Larger Macromolecules: The neutron structure of a mutant rubredoxin [Chatake et al. *Acta Cryst D60*, 1364 (2004)] was tested using error-free SIR data at the published resolution of 1.6Å. The fully deuterated crystal was considered to be the native and SIR derivatives were considered by replacing certain deuterated amino acids in the 52 residue sequence. A Pro⁵-mutant replaced the 5 proline residues in the sequence with wild type H-labeled proline. There were a total of 7426 data of which the 5500 largest E-values were phased by the NSnB method. 400 random trials were each run for 300 SnB cycles. CC values were computed for the 400*399/2 = 79800 pairs of phases sets. The 3rd largest CC value identified the phase sets that had the lowest mean phase errors, 68.8 and 76.3°. Additional derivatives would markedly improve

| Amb1 | < δφ > | Amb2 | < δφ > | CC | |
|------|-------------|------|-------------|-----------|-------|
| 561 | 88.6 | 671 | 89.4 | 0.186 | |
| 611 | 88.9 | 661 | 88.4 | 0.172 | |
| 333 | 68.8 | 689 | 76.3 | 0.169 *** | |
| 509 | 88.8 | 671 | 89.4 | 0.142 | |
| 509 | 88.8 | 561 | 88.6 | 0.141 | |
| 203 | 88.1 | 497 | 88.6 | 0.140 | |
| | 37 | 88.3 | 493 | 88.3 | 0.136 |
| 199 | 88.6 | 795 | 89.2 | 0.131 | |
| | 51 | 88.6 | 489 | 88.9 | 0.131 |
| 387 | 89.3 | 537 | 88.8 | 0.128 | |