

Contact: Simon Bystryak Location: Hillsborough, NJ, USA Email: <u>sbystryak@allisystems.com</u> Tel: 973.919.5707

Website: http://www.allisystems.com

Company Profile

Industry Sector: Medical Diagnostics

Company Overview: Allied Innovative Systems (ALLIS) is R&D company with the mission of early stage development and validation of original tecnologies and devices in the field of clinical diagnostics and high throughput drug screening. The company has establised a sound record of bringing new technologies from the level of an academic idea into the stage of functional prototypes ready for commercialization.

Target Market(s): Clinical, research and hospital laboratories worldwide.

Management

Leadership:

Simon Bystryak, President and CEO Natalya Ossina, Chief Strategy Officer Anna Smagulova, Chief Financial Officer

Scientific Advisory Board:

Neil Constantine, Ph.D.: Professor of Biochemistry at Department of Pathology, Medical School, University of Maryland

- Nancy D. Connell, Ph.D.: Professor at Department of Medicine, UMDNJ Center for BioDefense and Center for Emerging Pathogens at New Jersey Medical School
- Armen Sarvazyan, Ph.D., D.Sc.: Chief Scientific Officer of Artann Laboratories, Adjunct Professor of the Surgery Department, Robert Wood Johnson Medical School, UMDNJ, and Honarary Foreign Professor of Physics Faculty of Moscow State University

National Institutes of Health

Larta

National Institutes of Health Commercialization Assistance Program (NIH-CAP)

Key Value Drivers

Technology: AmpliFlux is a platform tecnology for improving the performance of industry standard ELISA-based assays. AmpliFlux is based on the use of a unique photochemical amplification reaction, and consists of two steps: The first step is a conventional ELISA. At the second step of assay, the mixture is irradiated by visible light. Illumination of the samples leads to a drastic increase in the final signal. The secondsep takes just 5 minutes.

Competitive Advantage: A new technology allows one to increase the sensitivity, accuracy, dynamic range, signal-to-noise ratio, save costly reagents and reduce time of analysis. This unique approach has never been used before for improvement of the performance of ELISA-based methods. The outstanding results were achieved for increasing of the assay sensitivity for detection of analytes of different orgin, such as Hepatitis B, HIV, Prostate specific antigen and other physiologically active substances.

Plan & Strategy: Seeking a strategic partner.

*Technology funded by the NIAID and being commericialized under the NIH-CAP.

Product Pipeline



Prototype Products

The company's products include:

- Properly protected photochemical amplification method (PAM) is a product itself.
- AmpliFlux[™] Detection system kits (sets of reagents) for ELISA + PAM
- AmpliFlux[™] Illuminator.

US Patent No. 5,776,703; New patent applications in preparation



