

UNT HEALTH SCIENCE CENTER



Monoclonal Antibodies Against Dopamine Receptor Subtypes D1a and D1b

Learn more!

Robert McClain, PhD
Associate Vice President
robert.mcclain@unthsc.edu
817-735-2618

Research Tool

1997-473 1995-472

Our Inventors

Dr. Robert Luedtke robert.luedtke@unthsc.edu

Publications

"Immunoblot and immunohistochemical comparison of murine monoclonal antibodies specific for the rat D1a and D1b dopamine receptor subtypes" J Neuroimmunol 101(2):170 (1999)

Application

➤ Useful in dopamine receptor ligand studies to differentiate between D1-like dopamine receptor subtypes.

Details

- > SG2-D1a anti D1a dopamine receptor
 - tested with immunoblot and immunohistochemistry techniques
 - o tested in rat
- ➤ SG4-D1b anti D1b dopamine receptor
 - tested with immunoblot and immunohistochemistry techniques
 - o tested in rat

3500 Camp Bowie Blvd Fort Worth, TX 76107 Phone: 817-735-5147 FAX: 817-735-5485 techtransfer@unthsc.edu